dr. ir. Geert M. P. van Kempen dr. ir. Michael van Ginkel dr. ir. Cris L. Luengo Hendriks dr. dr. dipl. phys. Bernd Rieger prof. dr. ir. Lucas J. van Vliet

Quantitative Imaging Group, Department of Applied Sciences, Delft University of Technology

Delft March 20, 2017

Contents

T	inc	lices of functions by subject	4
	1.1	Library Functions	5
		1.1.1 Image Object	5
		1.1.2 Scalar Images	6
		1.1.3 Strings	6
		1.1.4 Arrays	6
		1.1.5 Frameworks	8
		1.1.6 Pixel Tables	8
		1.1.7 Data Structures	9
		1.1.8 Numerical Algorithms	9
		1.1.9 Sorting	9
		1.1.10 Indexing	0
		1.1.11 Memory Management	1
		1.1.12 Support Functions	1
	1.2	File I/O functions	2
		1.2.1 File IO	2
	1.3	Image Processing Functions	4
		1.3.1 Mathematics	4
		1.3.2 Statistics	8
		1.3.3 Manipulation	9
		1.3.4 Interpolation	20
		1.3.5 Painting	20
		1.3.6 Linear Filters	20
		1.3.7 Derivative Filters	21
		1.3.8 Non-Linear Filters	21
		1.3.9 Binary Filters	22
		1.3.10 Mathematical Morphology	22
		1.3.11 Point Operations	23
		1.3.12 Transforms	24
		1.3.13 Distance Transforms	24
	1.4	Application Functions	24

4 Contents

		1.4.1	Smoothing	24
		1.4.2	Sharpening	25
		1.4.3	Line and Edge Detection	25
		1.4.4	Extrema Detection	25
		1.4.5	Object Generation	25
		1.4.6	Noise Generation	26
		1.4.7	Image Restoration	26
		1.4.8	Shift Estimation	27
		1.4.9	Segmentation	27
		1.4.10	Analysis	27
		1.4.11	Measurement	28
		1.4.12	Functions for Microscopy	31
2	Fur	nction	reference	32
_				
3	Ass	sorted	topics	757
	Ass 3.1		topics ary conditions	757 759
		Bound		
	3.1	Bound	ary conditions	759
	3.1 3.2	Bound Compr Deriva	ary conditions	759 760
	3.1 3.2 3.3	Bound Compr Deriva DIPlib	ary conditions	759 760 762
	3.1 3.2 3.3 3.4	Bound Compr Deriva DIPlib Descri	ary conditions	759 760 762 763
	3.1 3.2 3.3 3.4 3.5	Bound Comprise Deriva DIPlib Descrip File fo	ary conditions	759 760 762 763 765
	3.1 3.2 3.3 3.4 3.5 3.6	Bound Compr Deriva DIPlib Descrip File fo	ary conditions	759 760 762 763 765 766
	3.1 3.2 3.3 3.4 3.5 3.6 3.7	Bound Comprise Deriva DIPlib Descrip File for From in	ary conditions	759 760 762 763 765 766
	3.1 3.2 3.3 3.4 3.5 3.6 3.7 3.8 3.9	Bound Comprise Deriva DIPlib Descrip File fo From it General	ary conditions	759 760 762 763 765 766 769 770
	3.1 3.2 3.3 3.4 3.5 3.6 3.7 3.8 3.9 3.10	Bound Comprise Deriva DIPlib Descrip File for From in General General	ary conditions	759 760 762 763 765 766 769 770 771

Chapter 1

Indices of functions by subject

1.1 Library Functions

1.1.1 Image Object

- ChangeDataType Change the data type of an image
- ChangeDimensions Changes the order of the dimensions in an image
- ChangeToOd Make an image zero dimensional
- HasContiguousData Determines whether an image has all data contiguous in memory
- HasNormalStride Determines whether an image has a normal stride
- ImageAssimilate Inherit properties of another image
- ImageCopyProperties Copy the properties of an image
- ImageForge Allocate pixel data for an image
- ImageFree Free an image
- ImageGetData Get the data pointers of a set of images
- ImageGetDataType Read the data type field
- ImageGetDimensionality Read the dimensionality field
- ImageGetDimensions Read the dimensions array
- ImageGetPlane Read the plane number
- ImageGetStride Read the stride array
- ImageGetType Read the type field
- ImageNew Allocate a structure
- ImagesCheck Check properties of several images
- ImagesCheckTwo Check properties of two images
- ImagesCompare Compare properties of several images
- ImagesCompareTwo Compare properties of two images
- ImageSetDataType Set the data type field
- ImageSetDimensions Set the dimensions array
- \bullet ${\tt ImageSetType}$ Set the image type field
- ImagesSeparate Take care of in-place operations
- ImageStrip Restore an image to its initial ("raw") state

1.1.2 Scalar Images

- ConvertDataType Converts the data type of an image
- IsScalar Determines whether an image is a scalar
- ScalarImageNew Allocate a scalar image

1.1.3 Strings

- StringAppend Append a string to another
- StringArrayCopy Copy a string array
- StringArrayFree Array free function
- StringArrayNew Allocate an array of strings
- StringCat Concatenate two strings
- StringCompare Compare two strings
- StringCompareCaseInsensitive Compare two strings without minding case
- StringCopy Copy a String
- StringCrop Crop a string
- StringFree Free a string
- StringNew Allocate a string
- StringReplace Replace the contents of one string with that of another
- UnderscoreSpaces Replace spaces with underscores

1.1.4 Arrays

- ArrayFree Array free function
- ArrayNew Array allocation function
- BooleanArrayCopy Copy an array
- BooleanArrayFind Find value in array
- BooleanArrayFree Array free function
- BooleanArrayNew Array allocation function
- BoundaryArrayFree Array free function
- BoundaryArrayNew Array allocation function
- ComplexArrayCopy Copy an array

- ComplexArrayFind Find value in array
- ComplexArrayFree Array free function
- ComplexArrayNew Array allocation function
- ConvertArray converts the data type of an array
- CoordinateArrayFree Array free function
- CoordinateArrayNew Array allocation function
- DataTypeArrayCopy Copy an array
- DataTypeArrayFind Find value in array
- DataTypeArrayFree Array free function
- DataTypeArrayNew Array allocation function
- FloatArrayCopy Copy an array
- FloatArrayFind Find value in array
- FloatArrayFree Array free function
- FloatArrayNew Array allocation function
- FrameWorkProcessArrayFree Array free function
- FrameWorkProcessArrayNew Array allocation function
- ImageArrayFree Array free function
- ImageArrayNew Array allocation function
- ImageCheckBooleanArray Check a boolean array
- ImageCheckBoundaryArray Check a boundary array
- ImageCheckComplexArray Check a complex array
- ImageCheckFloatArray Check a float array
- ImageCheckIntegerArray Check an integer array
- IntegerArrayCopy Copy an array
- IntegerArrayFind Find value in array
- IntegerArrayFree Array free function
- IntegerArrayNew Array allocation function
- StringArrayCopy Copy a string array
- StringArrayFree Array free function
- StringArrayNew Allocate an array of strings

- VoidPointerArrayCopy Copy an array
- VoidPointerArrayFind Find value in array
- VoidPointerArrayFree Array free function
- VoidPointerArrayNew Array allocation function

1.1.5 Frameworks

- MonadicFrameWork FrameWork for monadic operations
- PixelTableFrameWork FrameWork for PixelTable filters
- ScanFrameWork FrameWork for scanning multiple images
- SeparableFrameWork FrameWork for separable filters
- SingleOutputFrameWork FrameWork for generation functions

1.1.6 Pixel Tables

- BinaryImageToPixelTable Convert a binary image to a pixel table
- GreyValuesInPixelTable Copy greyvalues from image in pixel table
- PixelTableAddRun Add a new run to a pixel table
- PixelTableCreateFilter Create a pixel table from a filter shape
- PixelTableFrameWork FrameWork for PixelTable filters
- PixelTableGetDimensionality Get the dimensionality of a pixel table
- PixelTableGetDimensions Get the dimensions of a pixel table
- PixelTableGetOffsetAndLength Converts the pixel table's runs
- PixelTableGetOrigin Get the origin of the pixel table
- PixelTableGetPixelCount Get the number of pixels encoded in the pixel table
- PixelTableGetRun Get the contents of a pixel table run
- PixelTableGetRuns Get the number of runs in a pixel table
- PixelTableGetSize The number of pixels in the pixel table's bounding box
- PixelTableNew Allocate a new pixel table
- PixelTableSetRun Initialises a pixel table run
- PixelTableShiftOrigin Changes the origin of the pixel table
- PixelTableToBinaryImage Convert a pixel table to a binary image

1.1.7 Data Structures

- PixelHeapFree Destroy heap structure
- PixelHeapIsEmpty Query heap
- PixelHeapNew Create a new heap structure
- PixelHeapPop Pop item onto heap
- PixelHeapPush Push item onto heap
- PixelQueueFree Destroy queue structure
- PixelQueueIsEmpty Query queue
- PixelQueueNew Create a new queue structure
- PixelQueuePop Pop item from queue
- PixelQueuePush Push item onto queue
- StablePixelHeapFree Destroy heap structure
- StablePixelHeapIsEmpty Query heap
- StablePixelHeapNew Create a new heap structure
- StablePixelHeapPop Pop item onto heap
- StablePixelHeapPush Push item onto heap

1.1.8 Numerical Algorithms

• OneDimensionalSearch - Numerical algorithm

1.1.9 Sorting

- DistributionSort Sort a block of data
- DistributionSortIndices Sort indices to block of data
- DistributionSortIndices16 Sort indices to a block of data
- GetRank Value selection function
- ImageSort Sort image data
- ImageSortIndices Sort indices to image data
- InsertionSort Sort a block of data
- InsertionSortIndices Sort indices to a block of data
- InsertionSortIndices16 Sort indices to a block of data

- QuickSort Sort a block of data
- QuickSortAnything Sort data of any type
- QuickSortIndices Sort indices to a block of data
- QuickSortIndices16 Sort indices to a block of data
- Sort Sort a block of data
- SortAnything Sort data of any type
- SortCompareFunction Typedef for comparison function (sorting)
- SortIndices Sort indices to a block of data
- SortIndices16 Sort indices to a block of data
- SortSwapFunction Typedef for swap and copy function (sorting)

1.1.10 Indexing

- CoordinateToIndex Convert coordinate to pixel index
- dip_PixelGetFloat Midlevel PixelIO function
- dip_PixelGetInteger Midlevel PixelIO function
- dip_PixelSetFloat Midlevel PixelIO function
- dip_PixelSetInteger Midlevel PixelIO function
- Get Get a pixel value
- GetComplex Get complex pixel value
- GetFloat Get float pixel value
- GetInteger Get integer pixel value
- IndexToCoordinate Convert pixel index to coordinate
- IndexToCoordinateWithSingletons Convert pixel index to coordinate
- NeighbourIndicesListMake Get indices to direct neighbours
- NeighbourListMake Get list of direct neighbours
- NeighbourListMakeChamfer Get list of neighbours based on Chamfer metric
- NeighbourListMakeImage Get list of neighbours based on metric in image
- NeighbourListToIndices Get indices to neighbours
- Set the value of a pixel
- SetComplex Set a pixel value
- SetFloat Set a pixel value
- SetInteger Set a pixel value

1.1.11 Memory Management

- MemoryCopy Copy memory blocks
- MemoryFree Free a chunk of memory
- MemoryFunctionsSet Sets memory allocation functions
- MemoryNew Allocate and track memory
- MemoryReallocate Reallocate a chunk of memory
- ResourcesFree Free resources
- ResourcesMerge Add one resource list to another
- ResourcesNew Allocate a resource tracking structure
- ResourceSubscribe Track a resource
- ResourceUnsubscribe Stop tracking a resource

1.1.12 Support Functions

- DataTypeAllowed Check whether a data type is allowed
- DataTypeGetInfo Get information about a data type
- error.h Contains error messages
- ErrorFree Free a DIPlib call tree
- Exit Clean up before exiting
- FillBoundaryArray Fill the border of array according to the boundary condition
- GetLibraryInformation Support function
- GetUniqueNumber Obtain an unique value
- GlobalBoundaryConditionGet Get global Boundary Conditions
- GlobalBoundaryConditionSet Set global boundary conditions
- GlobalFilterShapeGet Get global filter shape value
- GlobalFilterShapeSet Set the global filter shape value
- GlobalGaussianTruncationGet Get the global gaussian truncation
- GlobalGaussianTruncationSet Set the global gaussian truncation
- Initialise Initialise DIPlib
- macros.h Various macros
- ovl.h Call an overloaded function

- PhysicalDimensionsCopy Copy a Physical Dimensions
- PhysicalDimensionsFree Free a Physical Dimensions data structure
- PhysicalDimensionsIsIsotropic Checks if the Physical Dimensions are isotropic
- Physical Dimensions New Allocates a new Physical Dimensions structure
- Register Generic registry function
- RegisterClass Register a registry class
- RegistryArrayNew Allocate a registry array
- RegistryGet Get a registry item
- RegistryList Get an array of registry IDs
- RegistryValid Validate an registry item
- TimerGet Timing functions
- TimerSet Timing functions
- tpi.h Type iterator
- Unregister Remove a registry item

1.2 File I/O functions

1.2.1 File IO

- Colour2Gray Convert ND image with colour information to a (n-1)D grayvalue image (in dipIO)
- ImageFileGetInfo Get information about image in file (in dipIO)
- ImageFileInformationFree Free a Image File Information structure (in dipIO)
- ImageFileInformationNew Allocate an Image File Information structure (in dipIO)
- ImageIsGIF Confirm that a file is a GIF file (in dipIO)
- ImageIsICS Confirm that a file is an ICS file (in dipIO)
- ImageIsJPEG Confirm that a file is a JPEG file (in dipIO)
- ImageIsLSM Confirm that a file is a Zeiss LSM file (in dipIO)
- ImageIsTIFF Confirm that a file is a TIFF file (in dipIO)
- ImageRead Read grey-value image from file (in dipIO)
- ImageReadColour Read colour image from file (in dipIO)
- ImageReadCSV Read comma-separated values from file (in dipIO)
- ImageReadCSVInfo Get information about image in comma-separated values file (in dipIO)

- ImageReadGIF Read a GIF image from file (in dipIO)
- ImageReadGIFInfo Get information about image in GIF file (in dipIO)
- ImageReadICS Read ICS image from file (in dipIO)
- ImageReadICSInfo Get information about image in ICS file (in dipIO)
- ImageReadJPEG Read JPEG image from file (in dipIO)
- ImageReadJPEGInfo Get information about image in JPEG file (in dipIO)
- ImageReadLSM Read Zeiss LSM image from file (in dipIO)
- ImageReadLSMInfo Get information about image in LSM file (in dipIO)
- ImageReadPIC Read BioRad PIC image from file (in dipIO)
- ImageReadPICInfo Get information about image in BioRad PIC file (in dipIO)
- ImageReadROI Read a portion of a grey-value image from file (in dipIO)
- ImageReadTIFF Read TIFF image from file (in dipIO)
- ImageReadTIFFInfo Get information about image in TIFF file (in dipIO)
- ImageWrite Write grey-value image to file (in dipIO)
- ImageWriteColour Write colour image to file (in dipIO)
- ImageWriteCSV Write image to a comma-separated-value file (in dipIO)
- ImageWriteEPS Write image to Encapsulated PostScript file (in dipIO)
- ImageWriteFLD Write image to AVS field file (in dipIO)
- ImageWriteGIF Write image to a GIF file (in dipIO)
- ImageWriteICS Write ICS image to file (in dipIO)
- ImageWriteJPEG Write JPEG image to file (in dipIO)
- ImageWritePS Write image to PostScript file (in dipIO)
- ImageWriteTIFF Write TIFF image to file (in dipIO)
- MeasurementRead Read measurement results from a file
- MeasurementWrite Write measurement results to a file
- MeasurementWriteCSV Write measurement results to a CSV file
- MeasurementWriteHTML Write measurement results to an HTML file
- MeasurementWriteText Write measurement results as readable text

1.3 Image Processing Functions

1.3.1 Mathematics

- Abs Arithmetic function
- Acos trigonometric function
- Add arithmetic function
- AddComplex arithmetic function
- AddFloat arithmetic function
- AddInteger arithmetic function
- And logic operation
- Arith arithmetic function
- Arith_ComplexSeparated arithmetic function
- Asin trigonometric function
- Atan trigonometric function
- Atan2 arithmetic function
- BesselJO mathematical function
- BesselJ1 mathematical function
- BesselJN mathematical function
- BesselYO mathematical function
- BesselY1 mathematical function
- Besselyn mathematical function
- Ceil Arithmetic function
- Compare Compare grey values in two images
- Cos trigonometric function
- Cosh trigonometric function
- CumulativeSum statistics function
- Div arithmetic function
- DivComplex arithmetic function
- DivFloat arithmetic function
- DivInteger arithmetic function
- Equal Compare grey values in two images

- Erf mathematical function
- Erfc mathematical function
- Exp arithmetic function
- Exp10 arithmetic function
- Exp2 arithmetic function
- Floor Arithmetic function
- Fraction Arithmetic function
- GetMaximumAndMinimum statistics function
- Greater Compare grey values in two images
- IDivergence difference measure
- Imaginary Arithmetic function
- Invert logic operation
- Lesser Compare grey values in two images
- Ln arithmetic function
- LnGamma mathematical function
- LnNormError difference measure
- \bullet Log10 arithmetic function
- Log2 arithmetic function
- Max arithmetic function
- MaxFloat arithmetic function
- Maximum statistics function
- mBesselJO mathematical function
- mBesselJ1 mathematical function
- mBesselJN mathematical function
- mBesselYO mathematical function
- mBesselY1 mathematical function
- mBesselYN mathematical function
- Mean statistics function
- MeanAbsoluteError difference measure
- MeanError difference measure

- MeanModulus statistics function
- MeanSquareError difference measure
- MeanSquareModulus statistics function
- Median statistics function
- mErf mathematical function
- mErfc mathematical function
- mExp10 mathematical function
- mExp2 mathematical function
- mFraction mathematical function
- mGammaP mathematical function
- mGammaQ mathematical function
- Min arithmetic function
- MinFloat arithmetic function
- Minimum statistics function
- mLnGamma mathematical function
- mLog2 mathematical function
- mNearestInt mathematical function
- Modulo Arithmetic function
- Modulus Arithmetic function
- mReciprocal mathematical function
- mSign mathematical function
- mSinc mathematical function
- mTruncate mathematical function
- Mul arithmetic function
- MulComplex arithmetic function
- MulConjugate arithmetic function
- MulConjugateComplex arithmetic function
- MulFloat arithmetic function
- Mulinteger arithmetic function
- NearestInt Arithmetic function

- NormaliseSum Normalise the sum of the pixel values
- NotEqual Compare grey values in two images
- NotGreater Compare grey values in two images
- NotLesser Compare grey values in two images
- Or logic operation
- Percentile statistics function
- Phase Arithmetic function
- RadialMaximum statistics function
- RadialMean statistics function
- RadialMinimum statistics function
- RadialSum statistics function
- Real Arithmetic function
- Reciprocal arithmetic function
- RootMeanSquareError difference measure
- Select Configurable selection function
- Sign Arithmetic function
- Sin trigonometric function
- Sinc mathematical function
- SingularValueDecomposition Singular value decomposition
- Sinh trigonometric function
- \bullet Sqrt arithmetic function
- StandardDeviation statistics function
- Sub arithmetic function
- SubComplex arithmetic function
- SubFloat arithmetic function
- SubInteger arithmetic function
- Sum statistics function
- SumModulus statistics function
- Tan trigonometric function
- Tanh trigonometric function

- TensorImageInverse Invert tensor image
- Truncate Arithmetic function
- Variance statistics function
- WeightedAdd arithmetic function
- WeightedDiv arithmetic function
- WeightedMul arithmetic function
- WeightedSub arithmetic function
- Xor logic operation

1.3.2 Statistics

- ChordLength Compute the chord lengths of the different phases
- CumulativeSum statistics function
- GetMaximumAndMinimum statistics function
- IDivergence difference measure
- LnNormError difference measure
- Maximum statistics function
- Mean statistics function
- MeanAbsoluteError difference measure
- MeanError difference measure
- MeanModulus statistics function
- MeanSquareError difference measure
- MeanSquareModulus statistics function
- Median statistics function
- Minimum statistics function
- PairCorrelation Compute the pair correlation function
- Percentile statistics function
- ProbabilisticPairCorrelation Compute the probabilistic pair correlation function
- RadialMaximum statistics function
- RadialMean statistics function
- RadialMinimum statistics function

- RadialSum statistics function
- RootMeanSquareError difference measure
- StandardDeviation statistics function
- Sum statistics function
- SumModulus statistics function
- Variance statistics function

1.3.3 Manipulation

- Crop Remove the outer parts of an image
- dip_PixelGetFloat Midlevel PixelIO function
- dip_PixelGetInteger Midlevel PixelIO function
- dip_PixelSetFloat Midlevel PixelIO function
- dip_PixelSetInteger Midlevel PixelIO function
- ExtendRegion Image manipulation functions
- Get Get a pixel value
- GetComplex Get complex pixel value
- GetFloat Get float pixel value
- GetInteger Get integer pixel value
- GetLine Get a line from an image
- GetSlice Get a slice from an image
- Map Remaps an image
- Mirror Mirrors an image
- PutLine Put a line in an image
- PutSlice Put a slice in an image
- Resampling Interpolation function
- Rotation Interpolation function
- Rotation3d Interpolation function
- Rotation3d_Axis Interpolation function
- Set the value of a pixel
- SetComplex Set a pixel value

- SetFloat Set a pixel value
- SetInteger Set a pixel value
- Shift an image manipulation function
- Skewing Interpolation function
- Subsampling Interpolation function
- Wrap Wrap an image

1.3.4 Interpolation

- Resampling Interpolation function
- Rotation Interpolation function
- Rotation3d Interpolation function
- Rotation3d_Axis Interpolation function
- Skewing Interpolation function
- SubpixelLocation Gets coordinates of an extremum with sub-pixel precision
- SubpixelMaxima Gets coordinates of local maxima with sub-pixel precision
- SubpixelMinima Gets coordinates of local minima with sub-pixel precision
- Subsampling Interpolation function

1.3.5 Painting

- PaintBox Paint a box
- PaintDiamond Paint a diamond-shaped object
- PaintEllipsoid Paint an ellipsoid

1.3.6 Linear Filters

- Convolve1d Perform a 1D convolution
- ConvolveFT Fourier transform—based convolution filter
- Derivative Derivative filter
- FiniteDifference A linear gradient filter
- FiniteDifferenceEx A linear gradient filter
- GaborIIR Infinite impulse response filter
- Gauss Gaussian Filter

- Gaussett Gaussian Filter through the Fourier Domain
- Gaussiir Infinite impulse response filter
- GeneralConvolution Genaral convolution filter
- Laplace Second order derivative filter
- SeparableConvolution FrameWork for separable convolution filters
- Sharpen Enhance an image
- SobelGradient A linear gradient filter
- Uniform Uniform filter

1.3.7 Derivative Filters

- Derivative Derivative filter
- Dgg Second order derivative filter
- FiniteDifference A linear gradient filter
- FiniteDifferenceEx A linear gradient filter
- Gauss Gaussian Filter
- Gaussett Gaussian Filter through the Fourier Domain
- GradientDirection2D Derivative filter
- GradientMagnitude Derivative filter
- Laplace Second order derivative filter
- LaplaceMinDgg Second order derivative filter
- LaplacePlusDgg Second order derivative filter
- SobelGradient A linear gradient filter

1.3.8 Non-Linear Filters

- BiasedSigma Adaptive edge sharpening & contrast enhancing filter
- Closing Morphological closing operation
- Dilation Local maximum filter
- Erosion Local minimum filter
- GaussianSigma Adaptive Gaussian smoothing filter
- GeneralisedKuwahara Generalised Kuwahara filter
- GeneralisedKuwaharaImproved Generalised Kuwahara filter

- Kuwahara Edge perserving smoothing filter
- KuwaharaImproved Edge perserving smoothing filter
- MedianFilter Non-linear smoothing filter
- MorphologicalSmoothing Morphological smoothing filter
- Opening Morphological opening operation
- PercentileFilter Rank-order filter
- Sigma Adaptive uniform smoothing filter
- VarianceFilter Sample Variance Filter

1.3.9 Binary Filters

- BinaryClosing Binary morphological closing operation
- BinaryDilation Binary morphological dilation operation
- BinaryErosion Binary morphological erosion operation
- BinaryOpening Binary morphological opening operation
- BinaryPropagation Morphological propagation of binary objects
- EdgeObjectsRemove Remove binary edge objects
- EuclideanSkeleton binary skeleton operation
- GrowRegions Dilate the regions in a labelled image
- Label Label a binary image

1.3.10 Mathematical Morphology

- AreaOpening Morphological filter
- BinaryClosing Binary morphological closing operation
- BinaryDilation Binary morphological dilation operation
- BinaryErosion Binary morphological erosion operation
- BinaryOpening Binary morphological opening operation
- BinaryPropagation Morphological propagation of binary objects
- \bullet ${\tt Closing}$ Morphological closing operation
- Dilation Local maximum filter
- DirectedPathOpening Morphological filter
- EdgeObjectsRemove Remove binary edge objects

- Erosion Local minimum filter
- EuclideanSkeleton binary skeleton operation
- GrowRegions Dilate the regions in a labelled image
- GrowRegionsWeighted Grow labelled regions using grey-weighted distances
- Lee Morphological edge detector
- LocalMinima Marks local minima (or regional minima)
- Maxima Detects local maxima
- Minima Detects local minima
- $\bullet \ {\tt MorphologicalGradientMagnitude} \ \ {\tt Morphological\ edge\ detector}$
- MorphologicalRange Morphological edge detector
- MorphologicalReconstruction Morphological filter
- MorphologicalSmoothing Morphological smoothing filter
- MorphologicalThreshold Morphological smoothing filter
- MultiScaleMorphologicalGradient Morphological edge detector
- Opening Morphological opening operation
- PathOpening Morphological filter
- SeededWatershed Morphological segmentation
- Tophat Morphological high-pass filter
- UpperEnvelope Upper envelope transform (a flooding and an algebraic closing)
- Watershed Morphological segmentation

1.3.11 Point Operations

- Clip Point operation
- Compare Compare grey values in two images
- ContrastStretch Point operation
- Equal Compare grey values in two images
- ErfClip Point Operation
- Greater Compare grey values in two images
- HysteresisThreshold Point Operation
- IsodataThreshold Point operation

- Lesser Compare grey values in two images
- NotEqual Compare grey values in two images
- NotGreater Compare grey values in two images
- NotLesser Compare grey values in two images
- NotZero Point Operation
- RangeThreshold Point Operation
- Select Configurable selection function
- SelectValue Point Operation
- Threshold Point Operation

1.3.12 Transforms

- FourierTransform Computes the Fourier transform
- HartleyTransform Computes the Hartley transform

1.3.13 Distance Transforms

- Euclidean Distance Transform Euclidean distance transform
- GreyWeightedDistanceTransform Grey weighted distance transform
- GrowRegionsWeighted Grow labelled regions using grey-weighted distances
- VectorDistanceTransform Euclidean vector distance transform

1.4 Application Functions

1.4.1 Smoothing

- Closing Morphological closing operation
- Gauss Gaussian Filter
- GaussFT Gaussian Filter through the Fourier Domain
- Kuwahara Edge perserving smoothing filter
- KuwaharaImproved Edge perserving smoothing filter
- MedianFilter Non-linear smoothing filter
- MorphologicalSmoothing Morphological smoothing filter
- MorphologicalThreshold Morphological smoothing filter
- Opening Morphological opening operation

- PercentileFilter Rank-order filter
- Uniform Uniform filter
- UpperEnvelope Upper envelope transform (a flooding and an algebraic closing)

1.4.2 Sharpening

• Sharpen - Enhance an image

1.4.3 Line and Edge Detection

- DanielsonLineDetector Line detector
- GradientMagnitude Derivative filter
- Laplace Second order derivative filter
- Lee Morphological edge detector
- MorphologicalGradientMagnitude Morphological edge detector
- MorphologicalRange Morphological edge detector
- $\bullet \ {\tt MultiScaleMorphologicalGradient} \ \ {\tt Morphological\ edge\ detector}$

1.4.4 Extrema Detection

- LocalMinima Marks local minima (or regional minima)
- Maxima Detects local maxima
- Minima Detects local minima
- SubpixelLocation Gets coordinates of an extremum with sub-pixel precision
- SubpixelMaxima Gets coordinates of local maxima with sub-pixel precision
- SubpixelMinima Gets coordinates of local minima with sub-pixel precision

1.4.5 Object Generation

- CityBlockDistanceToPoint Distance generation function
- EllipticDistanceToPoint Distance generation function
- EuclideanDistanceToPoint Distance generation function
- FTBox Generates the Fourier transform of a box
- FTCross Generates the Fourier transform of a cross
- FTCube Generates the Fourier transform of a cube

- FTEllipsoid Generates Fourier transform of a ellipsoid
- FTGaussian Generates the Fourier transform of a Gaussian
- FTSphere Generated Fourier transform of a sphere
- IncoherentOTF Generates an incoherent OTF
- IncoherentPSF Generates an incoherent PSF
- PaintBox Paint a box
- PaintDiamond Paint a diamond-shaped object
- PaintEllipsoid Paint an ellipsoid
- TestObjectAddNoise TestObject generation function
- TestObjectBlur TestObject generation function
- TestObjectCreate TestObject generation function
- TestObjectModulate TestObject generation function

1.4.6 Noise Generation

- BinaryNoise Generates an image disturbed by binary noise
- BinaryRandomVariable Binary random variable generator
- GaussianNoise Generate an image disturbed by Gaussian noise
- GaussianRandomVariable Gaussian random variable generator
- PoissonNoise Generate an image disturbed by Poisson noise
- PoissonRandomVariable Poisson random variable generator
- RandomSeed Initialise random number generator
- RandomSeedVector Initialise random number generator
- RandomVariable Random number generator
- UniformNoise Generate an image disturbed by uniform noise
- UniformRandomVariable Uniform random variable generator

1.4.7 Image Restoration

- AttenuationCorrection Attenuation correction algorithm
- ExponentialFitCorrection Exponential fit based attenuation correction
- PseudoInverse Image restoration filter
- SimulatedAttenuation Simulation of the attenuation process

- TikhonovMiller Image restoration filter
- TikhonovRegularizationParameter Determine the value of the regularisation parameter
- Wiener Image Restoration Filter

1.4.8 Shift Estimation

- CrossCorrelationFT Normalized cross-correlation using the Fourier Transform
- FindShift Estimate the shift between images

1.4.9 Segmentation

- Canny Edge detector
- HysteresisThreshold Point Operation
- IsodataThreshold Point operation
- RangeThreshold Point Operation
- SeededWatershed Morphological segmentation
- Threshold Point Operation
- Watershed Morphological segmentation

1.4.10 Analysis

- Canny Edge detector
- ChordLength Compute the chord lengths of the different phases
- DanielsonLineDetector Line detector
- ImageChainCode Extracts all chain codes from a labeled image
- Label Label a binary image
- Measure Measure object features
- PairCorrelation Compute the pair correlation function
- ProbabilisticPairCorrelation Compute the probabilistic pair correlation function
- StructureTensor2D Two dimensional Structure Tensor
- SubpixelLocation Gets coordinates of an extremum with sub-pixel precision
- SubpixelMaxima Gets coordinates of local maxima with sub-pixel precision
- SubpixelMinima Gets coordinates of local minima with sub-pixel precision

1.4.11 Measurement

- ChainCodeArrayFree Chain code array deallocation
- ChainCodeArrayNew Chain code array allocation
- ChainCodeConvexHull Compute convex hull from chain code
- ChainCodeFree Chain code object deallocation
- ChainCodeGetChains Chain code access function
- ChainCodeGetConnectivity Chain code access function
- ChainCodeGetFeret Chain code measurement function
- ChainCodeGetLabel Chain code access function
- ChainCodeGetLength Chain code measurement function
- ChainCodeGetLongestRun Chain code measurement function
- ChainCodeGetRadius Chain code measurement function
- ChainCodeGetSize Chain code access function
- ChainCodeGetStart Chain code access function
- ChainCodeNew Chain code object allocation
- ConvexHullGetArea Convex hull measurement function
- ConvexHullGetFeret Convex hull measurement function
- ConvexHullGetPerimeter Convex hull measurement function
- FeatureAnisotropy2D Measure the anisotropy in a labeled region
- FeatureBendingEnergy Undocumented measurement function
- FeatureCenter Measure the object's center
- FeatureChainCodeBendingEnergy Undocumented measurement function
- ullet FeatureChainCodeFunction Measurement feature #measure function
- FeatureComposeFunction Measurement feature #compose function
- FeatureCompositeFunction Measurement feature #measure function
- FeatureConvertFunction Measurement feature #convert function
- FeatureConvexArea Measure the area of the object's convex hull
- FeatureConvexity Measure the object's convexity
- FeatureConvexPerimeter Measure the perimeter of the object's convex hull
- FeatureConvHullFunction Measurement feature #measure function

- FeatureCreateFunction Measurement feature #create function
- FeatureDescriptionFree Free a Feature Description
- FeatureDescriptionFunction Measurement feature #description function
- FeatureDescriptionGetDescription Get the description of the described feature
- FeatureDescriptionGetLabels Get the labels of the described feature
- FeatureDescriptionGetName Get the name of the described feature
- FeatureDescriptionGetUnits Get the Units of the described feature
- FeatureDescriptionNew Allocate a new FeatureDescription
- FeatureDescriptionSetDescription Set the description of the described feature
- FeatureDescriptionSetDimensionLabels Label set convenience function
- FeatureDescriptionSetLabel Set the name of a particular feature label
- FeatureDescriptionSetLabels Set the labels of the described feature
- FeatureDescriptionSetName Set the name of the described feature
- FeatureDescriptionSetUnit Set the units of a particular feature dimension
- FeatureDescriptionSetUnits Set the units of a described feature
- FeatureDimension Measure the object's dimensions
- FeatureExcessKurtosis Undocumented measurement function
- FeatureFeret Measure the object's Feret diameters
- FeatureGinertia Measure the object's inertia
- FeatureGmu Measure the object's inertia
- FeatureGravity Measure the object's gravity
- FeatureImageFunction Measurement feature #measure function
- FeatureInertia Measure the object's inertia
- FeatureLineFunction Measurement feature #measure function
- FeatureLongestChaincodeRun Undocumented measurement function
- FeatureMass Measure the mass of the object (sum of grey-values)
- FeatureMaximum Measure the object's maximum coordinate value
- FeatureMaxVal Measure the object's maximum intensity
- FeatureMean Measure the object's mean intensity
- FeatureMinimum Measure the object's minimum coordinate value

- FeatureMinVal Measure the object's minimum intensity
- FeatureMu Measure the object's inertia
- FeatureOrientation2D Undocumented measurement function
- FeatureP2A Measure the circularity of the object
- FeaturePerimeter Measure the object's perimeter length
- FeatureRadius Measure the object's radius statistics
- FeatureShape Measure shape parameters of the object
- FeatureSize Measure the object's area/volume
- FeatureSkewness Undocumented measurement function
- FeatureStdDev Measure the standard deviation of the object's intensity
- FeatureSum Measure the sum of the grey values of the object
- FeatureSurfaceArea Measure the area of the object's surface
- \bullet Feature ValueFunction - Measurement feature #value function
- GetObjectLabels Lists object labels in image
- ImageChainCode Extracts all chain codes from a labeled image
- Label Label a binary image
- Measure Measure object features
- MeasurementFeatureConvert Convert the data of a measurement feature
- MeasurementFeatureDescription Measurement Description access function
- MeasurementFeatureFormat Feature data format convenience function
- MeasurementFeatureRegister Register a measurement function
- MeasurementFeatureRegistryFeatureDescription Get the feature description of a registered measurement feature
- MeasurementFeatureRegistryFeatureNeedsIntensityImage Checks whether the measurement function needs an intensity image
- MeasurementFeatureRegistryGet Get the registry information of a measurement feature
- MeasurementFeatureRegistryList Obtain a list of the registered measurement features
- MeasurementFeatures Get the measurement ID array
- MeasurementFeatureSize Feature data convenience function
- MeasurementFeatureValid Verify a measurement feature ID
- MeasurementForge Allocate the data of a measurement data structure

- MeasurementFree Free a measurement data structure
- MeasurementGetName Get the name of a Measurement structure
- MeasurementGetPhysicalDimensions Get the physical dimensions info of a measurement
- MeasurementID Get the ID of a Measurement structure
- MeasurementIsValid Checks whether a measurement is valid
- MeasurementNew Create new measurement data structure
- MeasurementNumberOfFeatures Get the number of measurement feature IDs
- MeasurementNumberOfObjects Get the number of object IDs
- MeasurementObjectData Object data access function
- MeasurementObjects Get an object ID array
- MeasurementObjectValid Verify an object ID
- MeasurementObjectValue Object value access function
- MeasurementSetName Set the name of a measurement structure
- MeasurementSetPhysicalDimensions Set the physical dimensions info of the measurement
- MeasurementToHistogram Creats a histogram for a measurement
- MeasurementToImage Exports the data in a measurement structure to an image
- ObjectToMeasurement Convert object label value to measurement value
- Physical Dimensions Copy Copy a Physical Dimensions
- Physical Dimensions Free Free a Physical Dimensions data structure
- Physical Dimensions Is Isotropic Checks if the Physical Dimensions are isotropic
- Physical Dimensions New Allocates a new Physical Dimensions structure
- SmallObjectsRemove Remove small objects from an image

1.4.12 Functions for Microscopy

- AttenuationCorrection Attenuation correction algorithm
- ExponentialFitCorrection Exponential fit based attenuation correction
- IncoherentOTF Generates an incoherent OTF
- IncoherentPSF Generates an incoherent PSF
- SimulatedAttenuation Simulation of the attenuation process

Chapter 2

Function reference

Abs

Arithmetic function

SYNOPSIS

dip_Error dip_Abs (in, out)

DATA TYPES

binary, integer, integer, float, complex

FUNCTION

Computes the absolute value of the input image values.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	in	Input
$dip_{-}Image$	out	Output

SEE ALSO

 ${\tt Ceil}, \, {\tt Floor}, \, {\tt Sign}, \, {\tt Truncate}, \, {\tt Fraction}, \, {\tt NearestInt}$

Acos

trigonometric function

SYNOPSIS

dip_Error dip_Acos (in, out)

DATA TYPES

binary, integer, ${f float}$

FUNCTION

Computes the arc cosine of the input image values.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	in	Input
dip_Image	out	Output

SEE ALSO

Sin, Cos, Tan, Asin, Atan, Atan2, Sinh, Cosh, Tanh

AdaptiveBanana

Performs Gaussian filtering steered by parameter images

SYNOPSIS

#include "dip_adaptive.h"
dip_Error DIP_TPI_FUNC(dip_AdaptiveBanana)(in, out, para_images, curv_image,
filterSize, order, truncation)

DATA TYPES

sfloat

FUNCTION

This function performs Gaussian filtering steerd by the information stored in the parameter images (local orientation) and in the curvature image. The meaning of the parameter images depends on the dimensionality of the input image. Up to now only 2 and 3D images are supported for adaptive filtering. If the input image is not of type **float** it is converted to that type.

para_images: ImageArray containing orientation images.

2D: angle of the orientation.

3D: polar coordinate phi, theta for intrinsic 1D structures polar coordinates of two orientations for intrinsic 2D structures.

filterSize: Array containing the sigmas of the derivatives.

For intrinsic 1D structures, the first value is along the contour, the second perpendicular to it.

For intrinsic 2D structures, the first two are in the plane, whereas the other is perpendicular to them. If a value is zero no convolution is done is this direction.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input image
dip_Image	out	Output image
dip_ImageArray	para_images	Parameter images
$\mathtt{dip}_{-}\mathtt{Image}$	curv	Curvature image
dip_FloatArray	filterSize	Size of the filter
dip_IntegerArray	order	Order of the Gaussian derivative
dip_int	truncation	Truncation of the Gaussian

SEE ALSO

AdaptivePercentile, AdaptiveGauss, Gauss

LITERATURE

- P. Bakker, "Image structure analysis for seismic interpretation", PhD Thesis, TU Delft, The Netherlands, 2001
- L. Haglund, *Adaptive Mulitdimensional Filtering*", PhD Thesis, Link"oping University, Sweden, 1992 W.T. Freeman, "Steerable Filters and Local Analysis of Image Structure", PhD Thesis, MIT, USA, 1992

AdaptiveGauss

Performs Gaussian filtering steered by parameter images

SYNOPSIS

#include "dip_adaptive.h"
dip_Error dip_AdaptiveGauss(in,out,para_images,filterSize,order,truncation)

DATA TYPES

sfloat

FUNCTION

This function performs Gaussian filtering steerd by the information stored in the parameter images. The meaning of the parameter images depends on the dimensionality of the input image. Up to now only 2 and 3D images are supported for adaptive filtering. If the input image is not of type **float** it is converted to that type.

para_images: ImageArray containing orientation images.

2D: angle of the orientation

3D: polar coordinate phi, theta for intrinsic 1D structures polar coordinates of two orientations for intrinsic 2D structures

filterSize: Array containing the sigmas of the derivatives.

For intrinsic 1D structures, the first value is along the contour, the second perpendicular to it.

For intrinsic 2D structures, the first two are in the plane, whereas the other is perpendicular to them. If a value is zero no convolution is done is this direction.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	in	Input image
dip_Image	out	Output image
${ t dip_ImageArray}$	para_images	Parameter images
dip_FloatArray	filterSize	Size of the filter
dip_IntegerArray	order	Order of the Gaussian derivative
dip_int	truncation	Truncation of the Gaussian

SEE ALSO

AdaptivePercentile, AdaptiveBanana, Gauss

LITERATURE

- P. Bakker, "Image structure analysis for seismic interpretation", PhD Thesis, TU Delft, The Netherlands, 2001
- L. Haglund, *Adaptive Mulitdimensional Filtering*", PhD Thesis, Link" oping University, Sweden, 1992 W.T. Freeman, "Steerable Filters and Local Analysis of Image Structure", PhD Thesis, MIT, USA, 1992

AdaptivePercentile

Performs Percentile filtering steered by parameter images

SYNOPSIS

#include "dip_adaptive.h"
dip_Error DIP_TPI_FUNC(dip_AdaptivePercentile)(in, out, para_images, filterSize,
precentile)

DATA TYPES

sfloat

FUNCTION

This function performs percentile filtering steerd by the information stored in the parameter images (local orientation). The meaning of the parameter images depends on the dimensionality of the input image. Up to now only 2 and 3D images are supported for adaptive filtering. If the input image is not of type **float** it is converted to that type.

para_images: ImageArray containing orientation images.

2D: angle of the orientation.

3D: polar coordinate phi, theta for intrinsic 1D structures polar coordinates of two orientations for intrinsic 2D structures.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{\mathtt{-}}\mathtt{Image}$	in	Input image
dip_Image	out	Output image
dip_ImageArray	para_images	Parameter images
dip_FloatArray	filterSize	Size of the filter
dip_float	percentile	Percentile value

SEE ALSO

 ${\tt Adaptive Banana,\ Adaptive Gauss,\ Percentile Filter,\ Median Filter}$

LITERATURE

P. Bakker, "Image structure analysis for seismic interpretation", PhD Thesis, TU Delft, The Netherlands, 2001

L. Haglund, Adaptive Mulitdimensional Filtering", PhD Thesis, Link" oping University, Sweden, 1992 W.T. Freeman, "Steerable Filters and Local Analysis of Image Structure", PhD Thesis, MIT, USA, 1992

Add

arithmetic function

SYNOPSIS

```
dip_Error dip_Add ( in1, in2, out )
Calls Arith ( in1, in2, out, DIP_ARITHOP_ADD, DIP_DT_MINIMUM )
```

AddComplex

arithmetic function

SYNOPSIS

dip_Error dip_AddComplex (in, out, constant)

DATA TYPES

binary, integer, float, complex

FUNCTION

This function computes out = in + constant on a pixel by pixel basis. The data types of the in1 image and constant may be of different types. See Information about dyadic operations for more information about what the type of the output will be.

ARGUMENTS

Data type	Name	Description	
dip_Image	in	Input	
dip_Image	out	Output	
dip_complex	constant	Constant	

SEE ALSO

AddFloat

arithmetic function

SYNOPSIS

dip_Error dip_AddFloat (in, out, constant)

DATA TYPES

 $binary, \, \mathbf{integer}, \, \mathbf{float}, \, \mathbf{complex}$

FUNCTION

This function computes out = in + constant on a pixel by pixel basis. The data types of the in1 image and constant may be of different types. See Information about dyadic operations for more information about what the type of the output will be.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	in	Input
dip_Image	out	Output
dip_float	constant	Constant

SEE ALSO

Arith, Arith_ComplexSeparated, AddInteger, AddComplex, SubFloat, MulFloat, DivFloat

AddInteger

arithmetic function

SYNOPSIS

dip_Error AddInteger (in, out, constant)

DATA TYPES

 $binary, \, \mathbf{integer}, \, \mathbf{float}, \, \mathbf{complex}$

FUNCTION

This function computes out = in + constant on a pixel by pixel basis. The data types of the in1 image and constant may be of different types. See Information about dyadic operations for more information about what the type of the output will be.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	in	Input
dip_Image	out	Output
dip_int	constant	Constant

SEE ALSO

Arith, Arith_ComplexSeparated, AddFloat, AddComplex, SubInteger, MulInteger, DivInteger

And

logic operation

SYNOPSIS

dip_Error dip_And (in1, in2, out)

DATA TYPES

binary

FUNCTION

The function And performs the logic AND operation between the corresponding pixels in in1 and in2, and stores the result in out.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	in1	First binary input image
dip_Image	in2	Second binary input image
dip_Image	out	Output image

SEE ALSO

Arith, Xor, Or, Invert

AreaOpening

Morphological filter

SYNOPSIS

```
#include "dip_morphology.h"
dip_Error dip_AreaOpening ( grey, mask, out, filtersize, connectivity, closing )
```

DATA TYPES

integer, float

FUNCTION

The image grey will be filtered to remove local maxima (closing is DIP_FALSE) or local minima (closing is DIP_TRUE) with an area smaller than filtersize (in pixels).

Theoretically, the area opening can be written as the supremum of all the openings with each of the possible compact structuring elements of filtersize pixels. The connectivity parameter indicates which shapes are considered compact (i.e. all pixels are connected). See The connectivity parameter for more information.

ARGUMENTS

Data type	Name	Description
dip_Image	grey	Grey-value input image
$\mathtt{dip}_{-}\mathtt{Image}$	mask	Mask image for ROI processing
$\mathtt{dip}_{-}\mathtt{Image}$	out	Output image
dip_int	filtersize	Size of structuring element
dip_int	connectivity	Connectivity
dip_Boolean	closing	DIP_FALSE for area opening, DIP_TRUE for area closing

LITERATURE

L. Vincent, Grayscale area openings and closings, their efficient implementation and applications, Mathematical Morphology and Its Applications to Signal Processing, pages 22-27, 1993.

SEE ALSO

Opening, Closing, PathOpening, DirectedPathOpening, MorphologicalReconstruction

Arith

arithmetic function

SYNOPSIS

dip_Error dip_Arith (in1, in2, out, op, dt)

DATA TYPES

binary, integer, float, complex

FUNCTION

This function computes out = in1 op in2 on a pixel by pixel basis. The data types of the in1 and in2 image may be of different types. dt may be any of DIPlib's data types, or the constants DIP_DT_MINIMUM, DIP_DT_FLEX or DIP_DT_FLEXBIN, and determines what the output data type will be. DIP_DT_MINIMUM selects a data type according to the default for dyadic operations, see Information about dyadic operations for more information. DIP_DT_FLEX will choose a floating point (real or complex) type, the precision depends on the input types. DIP_DT_FLEXBIN is the same as DIP_DT_FLEX, except that two binary inputs will produce a binary output.

ARGUMENTS

Data type	Name	Description
dip_Image	in1	First input
dip_Image	in2	Second input
dip_Image	out	Output
dipf_ArithOperation	ор	Dyadic arithmetic operation
dip_DataType	dt	Data type for output

The dipf_ArithOperation flag can be one of:

Name	Description
DIP_ARITHOP_ADD	Addition (in1+in2)
DIP_ARITHOP_SUB	Subtraction (in1-in2)
DIP_ARITHOP_MUL	Multiplication (in1*in2)
DIP_ARITHOP_DIV	Division (in1/in2)
DIP_ARITHOP_MUL_CONJUGATE	Conjugate multiplication (in1*conj(in2))

For two binary inputs, and with $dt = DIP_DT_FLEXBIN$, the operations performed are equivalent to logical operations:

Name	Description
DIP_ARITHOP_ADD	Or (in1 in2)
DIP_ARITHOP_SUB	And not (in1&!in2)
DIP_ARITHOP_MUL	And (in1&in2)
DIP_ARITHOP_DIV	<pre>Xor (in1^in2)</pre>
DIP_ARITHOP_MUL_CONJUGATE	And (in1&in2)

SEE ALSO

Arith_ComplexSeparated, Add, Sub, Mul, Div, MulConjugate, AddInteger, AddFloat, AddComplex, SubInteger, SubFloat, SubComplex, MulInteger, MulFloat, MulComplex, MulConjugateComplex, DivInteger, DivFloat, DivComplex

Arith_ComplexSeparated

arithmetic function

SYNOPSIS

dip_Error dip_Arith (in1_real, in2_imag, in2_real, in2_imag, out_real, out_imag, op,
dt)

DATA TYPES

binary, integer, float

FUNCTION

This function computes out = in1 op in2 on a pixel by pixel basis. The data types of the in1 and in2 image may be of different types. The two input images and the output images have the complex portion of the data as a separate image, that is, in1 = in1_real + iin1_imag. in1_imag and in2_imag may be 0. dt may be any of DIPlib's data types, or the constants DIP_DT_MINIMUM or DIP_DT_FLEX, and determines what the output data type will be. DIP_DT_MINIMUM selects a data type according to the default for dyadic operations, see Information about dyadic operations for more information. DIP_DT_FLEX will choose a floating point (real or complex) type, the precision depends on the input types.

ARGUMENTS

Data type	Name	Description
dip_Image	in1_real	First input, real part
dip_Image	$in1_imag$	First input, imaginary part (or NULL)
dip_Image	in2_real	Second input, real part
dip_Image	$in2_imag$	Second input, imaginary part (or NULL)
dip_Image	out_real	Output, real part
dip_Image	out_imag	Output, imaginary part
dipf_ArithOperation	op	Dyadic arithmetic operation
dip_DataType	dt	Data type for output

The dipf_ArithOperation flag can be one of:

Name	Description
DIP_ARITHOP_ADD	Addition (in1+in2)
DIP_ARITHOP_SUB	Subtraction (in1-in2)
DIP_ARITHOP_MUL	Multiplication (in1*in2)
DIP_ARITHOP_DIV	Division (in1/in2)
DIP_ARITHOP_MUL_CONJUGATE	Conjugate multiplication (in1*conj(in2))

SEE ALSO

Arith, Add, Sub, Mul, MulConjugate, Div, AddInteger, AddFloat, AddComplex, SubInteger, SubFloat, SubComplex, MulInteger, MulFloat, MulComplex, MulConjugateComplex, DivInteger, DivFloat, DivComplex

ArrayFree

Array free function

SYNOPSIS

dip_Error dip_ArrayFree (array)

FUNCTION

This function frees *array, and sets array to zero.

ARGUMENTS

Data type	Name	Description
dip_Array *	array	pointer to a dip_Array

SEE ALSO

ArrayNew, ArrayFree

ArrayFree, IntegerArrayFree, FloatArrayFree, ComplexArrayFree, BoundaryArrayFree, FrameWorkProcessArrayFree, DataTypeArrayFree, ImageArrayFree, BooleanArrayFree, VoidPointerArrayFree, StringArrayFree, CoordinateArrayFree

ArrayNew

Array allocation function

SYNOPSIS

dip_Error dip_ArrayNew (array, size, elementSize, resources)

FUNCTION

This functions allocates the size elements of a dip_Array and sets the size of the array to size. The size of each element is determined by elementSize.

ARGUMENTS

Data type	Name	Description
dip_Array *	array	Array
dip_int	size	Size
dip_int	elementSize	ElementSize
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

ArrayNew, ArrayFree

ArrayNew, IntegerArrayNew, FloatArrayNew, ComplexArrayNew, BoundaryArrayNew, FrameWorkProcessArrayNew, DataTypeArrayNew, ImageArrayNew, BooleanArrayNew, VoidPointerArrayNew, StringArrayNew, CoordinateArrayNew

Asin

trigonometric function

SYNOPSIS

dip_Error dip_Asin (in, out)

DATA TYPES

binary, integer, ${f float}$

FUNCTION

Computes the arc sine of the input image values.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	in	Input
$dip_{-}Image$	out	Output

SEE ALSO

Sin, Cos, Tan, Acos, Atan, Atan2, Sinh, Cosh, Tanh

Atan

trigonometric function

SYNOPSIS

dip_Error dip_Atan (in, out)

DATA TYPES

binary, integer, ${f float}$

FUNCTION

Computes the arc tangent of the input image values.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	in	Input
dip_Image	out	Output

SEE ALSO

Sin, Cos, Tan, Asin, Acos, Atan2, Sinh, Cosh, Tanh

Atan2

arithmetic function

SYNOPSIS

dip_Error dip_Atan2 (in1, in2, out)

DATA TYPES

binary, integer, float

FUNCTION

This function computes out = atan2(in1, in2) on a pixel by pixel basis. The data types of the in1 and in2 image may be of different types. See Information about dyadic operations for more information about what the type of the output will be.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	in1	First input
dip_Image	in2	Second input
dip_Image	out	Output

SEE ALSO

Sin, Cos, Tan, Asin, Acos, Atan, Sinh, Cosh, Tanh

AttenuationCorrection

Attenuation correction algorithm

SYNOPSIS

#include "dip_microscopy.h"

dip_Error dip_AttenuationCorrection (in, out, fAttenuation, bAttenuation, background, threshold, NA, refIndex, ratio, method)

DATA TYPES

binary, integer, float

FUNCTION

This function implements an attenuation correction using three different recursive attenuation correction algorithms. The RAC-DET algorithm is the most accurate one, since it takes both forward and backward attenuation into account. It is however considerably slower that the RAC-LT2 and RAC-LT1 algorithms which take only forward attenuation into account. These last two algorithms assume a constant attenuation (background) for pixels with an intensity lower than the threshold.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input image
dip_Image	out	Output image
dip_float	fAttenuation	Forward attenuation factor
dip_float	bAttenuation	Backward attenuation factor
dip_float	background	Background attenuation factor
dip_float	threshold	Background threshold
dip_float	NA	Numerical aperture
dip_float	refIndex	Refractive index
dip_float	ratio	Z/X sampling ratio
dipf_AttenuationCorrection	method	Correction method

The dipf_AttenuationCorrection enumaration consists of the following flags:

Name	Description
DIP_ATTENUATION_RAC_LT2	Recursive Attenuation Correction algorithm using two Light Cone
	convolutions
DIP_ATTENUATION_RAC_LT1	Recursive Attenuation Correction algorithm using one Light Cone
	convolution
DIP_ATTENUATION_RAC_DET	Recursive Attenuation Correction algorithm using Directional
	Extinction Tracking

LITERATURE

K.C. Strasters, H.T.M. van der Voort, J.M. Geusebroek, and A.W.M. Smeulders, *Fast attenuation correction in fluorescence confocal imaging: a recursive approach*, BioImaging, vol. 2, no. 2, 1994, 78-92.

AUTHOR

Karel Strasters, adapted to DIPlib by Geert van Kempen.

SEE ALSO

 ${\tt SimulatedAttenuation}, {\tt ExponentialFitCorrection}$

BesselJ0

mathematical function

SYNOPSIS

dip_Error dip_BesselJO (in, out)

DATA TYPES

binary, integer, float

FUNCTION

Computes the Bessel function J0 of the input image values.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	out	Output

SEE ALSO

BesselJ1, BesselJN, BesselY0, BesselY1, BesselYN, LnGamma, Erf, Erfc, Sinc

BesselJ1

mathematical function

SYNOPSIS

dip_Error dip_BesselJ1 (in, out)

DATA TYPES

binary, integer, float

FUNCTION

Computes the Bessel function J1 of the input image values.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	in	Input
dip_Image	out	Output

SEE ALSO

BesselJO, BesselJN, BesselYO, BesselYI, BesselYN, LnGamma, Erf, Erfc, Sinc

BesselJN

mathematical function

SYNOPSIS

dip_Error dip_BesselJN (in, out, n)

DATA TYPES

binary, integer, ${f float}$

FUNCTION

Computes the Bessel function J of the order ${\tt n}$ of the input image values.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{\mathtt{I}}\mathtt{Image}$	in	Input
$\mathtt{dip}_{-}\mathtt{Image}$	out	Output
dip_int	n	Order of the Bessel function

SEE ALSO

BesselJO, BesselJO, BesselYO, BesselYO, BesselYN, LnGamma, Erf, Erfc, Sinc

BesselY0

mathematical function

SYNOPSIS

dip_Error dip_BesselYO (in, out)

DATA TYPES

binary, integer, float

FUNCTION

Computes the Bessel function Y0 of the input image values.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	in	Input
dip_Image	out	Output

SEE ALSO

BesselJO, BesselJI, BesselJN, BesselYI, BesselYN, LnGamma, Erf, Erfc, Sinc

BesselY1

mathematical function

SYNOPSIS

dip_Error dip_BesselY1 (in, out)

DATA TYPES

binary, integer, float

FUNCTION

Computes the Bessel function Y1 of the input image values.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	in	Input
dip_Image	out	Output

SEE ALSO

BesselJO, BesselJI, BesselJN, BesselYO, BesselYN, LnGamma, Erf, Erfc, Sinc

BesselYN

mathematical function

SYNOPSIS

dip_Error dip_BesselYN (in, out, n)

DATA TYPES

binary, integer, float

FUNCTION

Computes the Bessel function Y of the order ${\tt n}$ of the input image values.

ARGUMENTS

Data type	Name	Description	
$\mathtt{dip}_{\mathtt{I}}\mathtt{Image}$	in	Input	
$\mathtt{dip}_{-}\mathtt{Image}$	out	Output	
dip_int	n	Order of the Bessel function	

SEE ALSO

BesselJO, BesselJI, BesselJN, BesselYO, BesselYI, LnGamma, Erf, Erfc, Sinc

BiasedSigma

Adaptive edge sharpening & contrast enhancing filter

SYNOPSIS

```
#include "dip_filtering.h"
dip_Error dip_BiasedSigma ( in, out, se, boundary, param, shape, sigma, outputCount
)
```

DATA TYPES

integer, float

FUNCTION

The Biased Sigma filter is an adaptive edge sharpening and contrast enhancing filter. Its operation differs from the Sigma filter by separating the pixels with intensities higher than the pixel being filtered, from the pixels with lower intensities. The output for this pixel is the average closest in value to this pixel. If outputCount is DIP_TRUE, the output values represent the number of pixels over which the average has been calculated. When threshold is DIP_TRUE, the pixel intensities are thresholded at +/- 2 sigma, when it is set to DIP_FALSE, the intensities are weighted with the Gaussian difference with the intensity of the center pixel.

Only the rectangular, elliptic and diamond filter shapes are supported (DIP_FLT_SHAPE_RECTANGULAR, DIP_FLT_SHAPE_ELLIPTIC and DIP_FLT_SHAPE_DIAMOND). Other filter shapes can be implemented by setting shape to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, and passing a binary image in se. The "on" pixels define the shape of the filter window. Other values of shape are illegal.

If shape is not equal to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, se can be set to zero. When shape is set to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, param is ignored, and can be set to zero.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input image
dip_Image	out	Output image
dip_Image	se	Custom filter window (binary)
dip_BoundaryArray	boundary	Boundary conditions
dip_FloatArray	param	Filter sizes
dip_FilterShape	shape	Filter shape
dip_float	sigma	Sigma
dip_Boolean	outputCount	Output the Count

The enumerator dip_FilterShape contains the following constants:

Name	Description
DIP_FLT_SHAPE_DEFAULT	Default filter window, same as
	DIP_FLT_SHAPE_RECTANGULAR
DIP_FLT_SHAPE_RECTANGULAR	Rectangular filter window, can be even in size
DIP_FLT_SHAPE_ELLIPTIC	Elliptic filter window, always odd in size
DIP_FLT_SHAPE_DIAMOND	Diamond-shaped filter window, always odd in size
DIP_FLT_SHAPE_PARABOLIC	Parabolic filter window (morphology only)
DIP_FLT_SHAPE_DISCRETE_LINE	Rotated line structuring element (morphology only)
DIP_FLT_SHAPE_INTERPOLATED_LINE	Rotated line structuring element, through interpolation
	(morphology only)
DIP_FLT_SHAPE_PERIODIC_LINE	(not implemented)
DIP_FLT_SHAPE_STRUCTURING_ELEMENT	Use se as filter window, can be any size

Only the rectangular, elliptic and diamond filter shapes are supported (DIP_FLT_SHAPE_RECTANGULAR, DIP_FLT_SHAPE_ELLIPTIC and DIP_FLT_SHAPE_DIAMOND). Other filter shapes can be implemented by setting shape to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, and passing a binary image in se. The "on" pixels define the shape of the filter window. Other values of shape are illegal.

If shape is not equal to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, se can be set to zero. When shape is set to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, param is ignored, and can be set to zero.

LITERATURE

John-Sen Lee, Digital Image Smoothing and the Sigma Filter, Computer Vision, Graphics and Image Processing, 24, 255-269, 1983

SEE ALSO

Sigma, GaussianSigma

BinaryClosing

Binary morphological closing operation

SYNOPSIS

```
#include "dip_binary.h"
dip_Error dip_BinaryClosing ( in, out, connectivity, iterations, edge )
```

DATA TYPES

binary

FUNCTION

The connectivity parameter defines the metric, that is, the shape of the structuring element. 1 indicates city-block metric, or a diamond-shaped structuring element. 2 indicates chessboard metric, or a square structuring element. -1 and -2 indicate alternating connectivity and produce an octagonal structuring element. See The connectivity parameter for more information. The edge parameter specifies whether the border of the image should be treated as object (DIP_TRUE) or as background (DIP_FALSE). Additionally, you can set it to -1 for special handling: DIP_FALSE for the dilation, DIP_TRUE for the erosion; this avoids the border effect you can get in the corners of the image in some cases.

See section 9.6, "Morphology-based operations", in Fundamentals of Image Processing for a description of binary mathematical morphology operations.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	out	Output
dip_int	connectivity	Connectivity
dip_int	iterations	Iterations
dip_int	edge	Edge condition

KNOWN BUGS

This function is only implemented for images with a dimension up to three.

SEE ALSO

BinaryDilation, BinaryErosion, BinaryOpening, BinaryPropagation

BinaryDilation

Binary morphological dilation operation

SYNOPSIS

```
#include "dip_binary.h"
dip_Error dip_BinaryDilation ( in, out, connectivity, iterations, edge )
```

DATA TYPES

binary

FUNCTION

The connectivity parameter defines the metric, that is, the shape of the structuring element. 1 indicates city-block metric, or a diamond-shaped structuring element. 2 indicates chessboard metric, or a square structuring element. -1 and -2 indicate alternating connectivity and produce an octagonal structuring element. See The connectivity parameter for more information. The edge parameter specifies whether the border of the image should be treated as object (DIP_TRUE) or as background (DIP_FALSE).

See section 9.6, "Morphology-based operations", in Fundamentals of Image Processing for a description of binary mathematical morphology operations.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	out	Output
dip_int	connectivity	Connectivity
dip_int	iterations	Iterations
dip_Boolean	edge	Edge pixels on

KNOWN BUGS

This function is only implemented for images with a dimension up to three.

SEE ALSO

BinaryErosion, BinaryClosing, BinaryOpening, BinaryPropagation

BinaryErosion

Binary morphological erosion operation

SYNOPSIS

```
#include "dip_binary.h"
dip_Error dip_BinaryErosion ( in, out, connectivity, iterations, edge )
```

DATA TYPES

binary

FUNCTION

The connectivity parameter defines the metric, that is, the shape of the structuring element. 1 indicates city-block metric, or a diamond-shaped structuring element. 2 indicates chessboard metric, or a square structuring element. -1 and -2 indicate alternating connectivity and produce an octagonal structuring element. See The connectivity parameter for more information. The edge parameter specifies whether the border of the image should be treated as object (DIP_TRUE) or as background (DIP_FALSE).

See section 9.6, "Morphology-based operations", in Fundamentals of Image Processing for a description of binary mathematical morphology operations.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	out	Output
dip_int	connectivity	Connectivity
dip_int	iterations	Iterations
dip_Boolean	edge	Edge condition

KNOWN BUGS

This function is only implemented for images with a dimension up to three.

SEE ALSO

BinaryDilation, BinaryClosing, BinaryOpening, BinaryPropagation

BinaryImageToPixelTable

Convert a binary image to a pixel table

SYNOPSIS

```
#include "dip_pixel_table.h"
dip_Error dip_BinaryImageToPixelTable ( im, table, resources )
```

DATA TYPES

binary

FUNCTION

This functions converts a binary image to a newly allocated pixel table table.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	im	Binary image
dip_PixelTable *	table	Pixel table
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

Description of DIPlib's pixel tables

 ${\tt PixelTableCreateFilter, GreyValuesInPixelTable, PixelTableToBinaryImage}$

BinaryNoise

Generates an image disturbed by binary noise

SYNOPSIS

```
#include "dip_noise.h"
dip_Error dip_BinaryNoise ( in, out, p10, p01, random )
```

DATA TYPES

binary

FUNCTION

Generate an image disturbed by binary noise. See BinaryRandomVariable for more information on the random number generator.

ARGUMENTS

Data type	Name	Description	
dip_Image	in	Input	
dip_Image	out	Output	
dip_float	p10	Probability of a one to zero transition	
dip_float	p01	Probability of a zero to one transition	
dip_Random *	random	Pointer to a random value structure	

EXAMPLE

Get a binary noise disturbed image as follows:

```
dip_Image in, out;
dip_float p10, p01;
dip_Random random;

p10 = 0.1;
p01 = 0.2;
DIPXJ( dip_RandomSeed( &random, 0 ));
DIPXJ( dip_BinaryNoise( in, out, p10, p01, &random ));
```

SEE ALSO

 ${\tt BinaryRandomVariable, RandomVariable, RandomSeed, RandomSeedVector, UniformNoise, GaussianNoise, PoissonNoise}$

BinaryOpening

Binary morphological opening operation

SYNOPSIS

```
#include "dip_binary.h"
dip_Error dip_BinaryOpening ( in, out, connectivity, iterations, edge )
```

DATA TYPES

binary

FUNCTION

The connectivity parameter defines the metric, that is, the shape of the structuring element. 1 indicates city-block metric, or a diamond-shaped structuring element. 2 indicates chessboard metric, or a square structuring element. -1 and -2 indicate alternating connectivity and produce an octagonal structuring element. See The connectivity parameter for more information. The edge parameter specifies whether the border of the image should be treated as object (DIP_TRUE) or as background (DIP_FALSE). Additionally, you can set it to -1 for special handling: DIP_TRUE for the erosion, DIP_FALSE for the dilation; this avoids the border effect you can get in the corners of the image in some cases.

See section 9.6, "Morphology-based operations", in Fundamentals of Image Processing for a description of binary mathematical morphology operations.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	out	Output
dip_int	connectivity	Connectivity
dip_int	iterations	Iterations
dip_int	edge	Edge condition

KNOWN BUGS

This function is only implemented for images with a dimension up to three.

SEE ALSO

BinaryDilation, BinaryErosion, BinaryClosing, BinaryPropagation

BinaryPropagation

Morphological propagation of binary objects

SYNOPSIS

#include "dip_binary.h"
dip_Error dip_BinaryPropagation (seed, mask, out, connectivity, iterations, edge)

DATA TYPES

binary

FUNCTION

The connectivity parameter defines the metric, that is, the shape of the structuring element. 1 indicates city-block metric, or a diamond-shaped structuring element. 2 indicates chessboard metric, or a square structuring element. -1 and -2 indicate alternating connectivity and produce an octagonal structuring element. See The connectivity parameter for more information. The edge parameter specifies whether the border of the image should be treated as object (DIP_TRUE) or as background (DIP_FALSE).

See section 9.6, "Morphology-based operations", in Fundamentals of Image Processing for a description of binary mathematical morphology operations, and section 10.3, "Segmentation", for applications of binary propagation.

ARGUMENTS

Data type	Name	Description
dip_Image	seed	Input seed
dip_Image	mask	Input mask
dip_Image	out	Output
dip_int	connectivity	Connectivity
dip_int	iterations (0)	Iterations
dip_Boolean	edge	Edge condition

KNOWN BUGS

This function is only implemented for images with a dimension up to three.

SEE ALSO

 ${\tt BinaryDilation,\,BinaryErosion,\,BinaryClosing,\,BinaryOpening,\,EdgeObjectsRemove,\,GrowRegions}$

BinaryRandomVariable

Binary random variable generator

SYNOPSIS

```
#include "dip_noise.h"
dip_Error dip_BinaryRandomVariable ( random, input, p10, p01, output )
```

FUNCTION

The binary random variable is generated by altering the input value, if the value of a generated random variable is higher than the p10 probability, if input is DIP_TRUE, or higher than p01 otherwise.

See RandomVariable for more information on the random number generator.

ARGUMENTS

Data type	Name	Description	
dip_Random *	random	Pointer to a random value structure	
dip_Boolean	input	Input	
dip_float	p10	Probability of a one to zero transition	
dip_float	p01	Probability of a zero to one transition	

EXAMPLE

Get a binary random variable as follows:

```
dip_Random random;
dip_float p10, p01, value;

p10 = 0.1;
p01 = 0.2;
DIPXJ( dip_RandomSeed( &random, 0 ));
DIPXJ( dip_BinaryRandomVariable( &random, 1, p10, p01, &value ));
```

SEE ALSO

RandomVariable, RandomSeed, RandomSeedVector, UniformRandomVariable, GaussianRandomVariable, PoissonRandomVariable

BooleanArrayCopy

Copy an array

SYNOPSIS

dip_Error dip_BooleanArrayCopy (dest, src, resources)

FUNCTION

This function copies the boolean array **src** to **dest**. The array **dest** is created by this function as well.

ARGUMENTS

Data type	Name	Description
dip_BooleanArray *	dest	Destination array
dip_BooleanArray	src	Source array
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

BooleanArrayNew, BooleanArrayFree, BooleanArrayCopy, BooleanArrayFind
IntegerArrayCopy, FloatArrayCopy, ComplexArrayCopy, DataTypeArrayCopy, BooleanArrayCopy,
VoidPointerArrayCopy, StringArrayCopy

BooleanArrayFind

Find value in array

SYNOPSIS

dip_Error dip_BooleanArrayFind (array, value, index, found)

FUNCTION

Finds a value in an array and "returns" its index in the array. If found is zero, BooleanArrayFind will produce an error if value is not found, otherwise found obtains the search result (DIP_FALSE if value is not found).

ARGUMENTS

Data type	Name	Description	
dip_BooleanArray	array	Array to find value in	
dip_Boolean	value	Value to find	
dip_int *	index	Index of the found value	
dip_Boolean *	found	Value found or not	

SEE ALSO

BooleanArrayNew, BooleanArrayFree, BooleanArrayCopy, BooleanArrayFind
IntegerArrayFind, FloatArrayFind, ComplexArrayFind, DataTypeArrayFind, BooleanArrayFind,
VoidPointerArrayFind

BooleanArrayFree

Array free function

SYNOPSIS

dip_Error dip_BooleanArrayFree (array)

FUNCTION

This function frees *array, and sets array to zero.

ARGUMENTS

Data type	Name	Description
dip_BooleanArray *	array	Array

SEE ALSO

VoidPointerArrayNew, VoidPointerArrayFree, VoidPointerArrayCopy, VoidPointerArrayFind ArrayFree, IntegerArrayFree, FloatArrayFree, ComplexArrayFree, BoundaryArrayFree, FrameWorkProcessArrayFree, DataTypeArrayFree, ImageArrayFree, BooleanArrayFree, VoidPointerArrayFree, StringArrayFree, CoordinateArrayFree

BooleanArrayNew

Array allocation function

SYNOPSIS

dip_Error dip_BooleanArrayNew (array, size, value, resources)

FUNCTION

This function allocates the size elements of a dip_BooleanArray and sets the size of the array to size. Each array element is initialized with value.

ARGUMENTS

Data type	Name	Description
dip_BooleanArray *	array	Array
dip_int	size	Size
dip_Boolean	value	Initial value
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

BooleanArrayNew, BooleanArrayFree, BooleanArrayCopy, BooleanArrayFind ArrayNew, IntegerArrayNew, FloatArrayNew, ComplexArrayNew, BoundaryArrayNew, FrameWorkProcessArrayNew, DataTypeArrayNew, ImageArrayNew, BooleanArrayNew, VoidPointerArrayNew, StringArrayNew, CoordinateArrayNew

BoundaryArrayFree

Array free function

SYNOPSIS

dip_Error dip_BoundaryArrayFree (array)

FUNCTION

This function frees *array, and sets array to zero.

ARGUMENTS

Data type	Name	Description
dip_BoundaryArray *	array	Boundary conditions

SEE ALSO

BoundaryArrayNew, BoundaryArrayFree

ArrayFree, IntegerArrayFree, FloatArrayFree, ComplexArrayFree, BoundaryArrayFree, FrameWorkProcessArrayFree, DataTypeArrayFree, ImageArrayFree, BooleanArrayFree, VoidPointerArrayFree, StringArrayFree, CoordinateArrayFree

BoundaryArrayNew

Array allocation function

SYNOPSIS

dip_Error dip_BoundaryArrayNew (array, size, value, resources)

FUNCTION

This function allocates the size elements of a dip_BoundaryArray and sets the size of the array to size. Each array element is initialized with value.

ARGUMENTS

Data type	Name	Description
dip_BoundaryArray *	array	Boundary conditions
dip_int	size	Size
dip_Boundary	value	Initial value
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

BoundaryArrayNew, BoundaryArrayFree

ArrayNew, IntegerArrayNew, FloatArrayNew, ComplexArrayNew, BoundaryArrayNew, FrameWorkProcessArrayNew, DataTypeArrayNew, ImageArrayNew, BooleanArrayNew, VoidPointerArrayNew, StringArrayNew, CoordinateArrayNew

Canny

Edge detector

SYNOPSIS

```
#include "dip_detection.h"
dip_Error dip_Canny ( in, out, sigma, upper, lower )
```

DATA TYPES

Input is integer or float; output is **binary**.

FUNCTION

The Canny edge detector finds the ridges in the gradient magnitude, which correspond to the edges in the image. The gradient magnitude (see GradientMagnitude) is computed using Gaussian derivatives, with a sigma of sigma in both dimensions. The found ridges are pruned to remove the less salient edges. A threshold t1 is computed so that the 1-upper fraction of pixels with the highest gradient magnitude are kept. A second threshold, t2 = t1*lower, is selected that determines the minimal gradient magnitude expected of an edge. All edge pixels that exceed t2, and are in the same connected region as at least one pixel that exceeds t1, are selected as the output of this function (see HysteresisThreshold).

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	out	Output
dip_float	sigma	Sigma parameter for Gaussian derivatives
dip_float	lower	Lower threshold, as a fraction of upper threshold
dip_float	upper	Percentile used to compute upper threshold

LIMITATIONS

This function only works on 2D images.

Ceil

Arithmetic function

SYNOPSIS

dip_Error dip_Ceil (in, out)

DATA TYPES

binary, integer, ${f float}$

FUNCTION

Computes the ceil of the input image values, and outputs a signed integer typed image.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	in	Input
dip_Image	out	Output

SEE ALSO

Abs, Floor, Sign, Truncate, Fraction, NearestInt

${\tt ChainCodeArrayFree}$

Chain code array deallocation

SYNOPSIS

```
#include "dip_chaincode.h"
dip_Error dip_ChainCodeArrayFree ( array )
```

FUNCTION

This function frees *array, and sets array to zero.

ARGUMENTS

Data type	Name	Description
dip_ChainCodeArray *	array	Pointer to chain code array

SEE ALSO

 ${\tt ImageChainCode}, {\tt ChainCodeNew}, {\tt ChainCodeFree}, {\tt ChainCodeArrayNew}$

${\tt ChainCodeArrayNew}$

Chain code array allocation

SYNOPSIS

```
#include "dip_chaincode.h"
dip_Error dip_ChainCodeArrayNew ( array, size, resources )
```

FUNCTION

This function allocates the \mathtt{size} elements of a $\mathtt{dip_ChainCodeArrayNew}$ and sets the size of the array to \mathtt{size} .

ARGUMENTS

Data type	Name	Description
dip_ChainCodeArray *	array	Receives pointer to allocated structure
dip_int	size	Number of chains to allocate space for
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

 ${\tt ImageChainCode}, {\tt ChainCodeNew}, {\tt ChainCodeFree}, {\tt ChainCodeArrayFree}$

ChainCodeConvexHull

Compute convex hull from chain code

SYNOPSIS

```
#include "dip_chaincode.h"
dip_Error dip_ChainCodeConvexHull ( chaincode, polygon, resources )
```

FUNCTION

(To be documented)

We're using Melkman's algorithm to determine the convex hull of a polygonal representation of the boundary encoded by chaincode.

ARGUMENTS

Data type	Name	Description
dip_ChainCode	chaincode	Input chain code
dip_Polygon*	polygon	Output convex polygon
dip_Resources	resources	Resources tracking structure. See ResourcesNew

LITERATURE

Avraham A. Melkman, "On-line construction of the convex hull of a simple polyline," Information Processing Letters 25(1):11-12, 1987.

SEE ALSO

 ${\tt ImageChainCode,\ ConvexHullGetArea,\ ConvexHullGetPerimeter,\ ConvexHullGetFereter}$

ChainCodeFree

Chain code object deallocation

SYNOPSIS

```
#include "dip_chaincode.h"
dip_Error dip_ChainCodeFree ( chaincode )
```

FUNCTION

Deallocates the chaincode object created by ChainCodeNew, and sets the pointer to zero.

ARGUMENTS

Data type	Name	Description
dip_ChainCode *	chaincode	Pointer to chain code

SEE ALSO

 ${\tt ImageChainCode}, {\tt ChainCodeNew}, {\tt ChainCodeArrayNew}, {\tt ChainCodeArrayFree}$

ChainCodeGetChains

Chain code access function

SYNOPSIS

```
#include "dip_chaincode.h"
dip_Error dip_ChainCodeGetChains ( chaincode, chain )
```

FUNCTION

Returns a pointer to the first element of the chain. Each chain element contains a pointer to the next element. The last element has a NULL pointer. ChainCodeGetSize returns the number of elements in the chain. See ChainCodeNew for more information.

ARGUMENTS

Data type	Name	Description
dip_ChainCode	chaincode	Chain code
dip_Chain **	chain	Receives the pointer to the first element in the chain

SEE ALSO

ImageChainCode, ChainCodeNew, ChainCodeFree, ChainCodeGetSize, ChainCodeGetStart, ChainCodeGetLabel, ChainCodeGetConnectivity

ChainCodeGetConnectivity

Chain code access function

SYNOPSIS

```
#include "dip_chaincode.h"
dip_Error dip_ChainCodeGetConnectivity ( chaincode, connectivity )
```

FUNCTION

Returns the connectivity used when extracting the boundary described in chaincode. connectivity==1 indicates 4-connected neighbours, and the code uses integers 0 through 3. connectivity==2 indicates 8-connected neighbours, and the code uses values 0 through 7. See The connectivity parameter.

ARGUMENTS

Data type	Name	Description
dip_ChainCode	chaincode	Chain code
dip_int *	connectivity	Receives the connectivity value in chaincode

SEE ALSO

 ${\tt ImageChainCode}, {\tt ChainCodeNew}, {\tt ChainCodeFree}, {\tt ChainCodeGetSize}, {\tt ChainCodeGetChains}, {\tt ChainCodeGetStart}, {\tt ChainCodeGetLabel}$

ChainCodeGetFeret

Chain code measurement function

SYNOPSIS

```
#include "dip_chaincode.h"
dip_Error dip_ChainCodeGetFeret ( chaincode, stepSize, feret )
```

FUNCTION

This function measures the longest and shortest projections of the object encoded by chaincode. The chain code is rotated in stepSize degree intervals and the length of the projection on the x and y axes is computed for each orientation. The sizes of maximum and minimum projections, as well as the rotation at which they were obtained, are returned in the feret structure, which contains the following elements:

Data type	Name	Description
dip_float	maxDiameter	The widest projection of the object
dip_float	minDiameter	The narrowest projection of the object
dip_float	maxPerpendicular	The width of the projection perpendicular to minDiameter
dip_float	maxAngle	The angle of the projection for maxDiameter
dip_float	minAngle	The angle of the projection for minDiameter

ChainCodeGetFeret is the function formerly used by Measure for the FeatureFeret measurement. This measurement now uses ConvexHullGetFeret instead.

ARGUMENTS

Data type	Name	Description
dip_ChainCode	chaincode	Input chain code
dip_float	stepSize	The step size, in degrees
dip_Feret *	feret	Output measurement

CREDITS

The original code on which the current implementation is based, was donated by Gerie van der Heijden.

SEE ALSO

 ${\tt ImageChainCode}, ChainCodeNew, ChainCodeFree, ChainCodeGetLength, ChainCodeGetLongestRun, ChainCodeGetRadius}$

ChainCodeGetLabel

Chain code access function

SYNOPSIS

```
#include "dip_chaincode.h"
dip_Error dip_ChainCodeGetLabel ( chaincode, label )
```

FUNCTION

Returns the label ID of the object whose boundary is described by chaincode.

ARGUMENTS

Data type	Name	Description
dip_ChainCode	chaincode	Chain code
dip_int *	label	Receives the label ID in chaincode

SEE ALSO

${\tt ChainCodeGetLength}$

Chain code measurement function

SYNOPSIS

```
#include "dip_chaincode.h"
dip_Error dip_ChainCodeGetLength ( chaincode, length )
```

FUNCTION

Computes the length of the boundary encoded by chaincode. See FeaturePerimeter for a description of the algorithm. ChainCodeGetLength is the function used by Measure for the FeaturePerimeter measurement.

ARGUMENTS

Data type	Name	Description
dip_ChainCode	chaincode	Input chain code
dip_float *	length	Output measurement

SEE ALSO

 ${\tt ImageChainCode}, {\tt ChainCodeNew}, {\tt ChainCodeFree}, {\tt ChainCodeGetLongestRun}, {\tt ChainCodeGetFeret}, {\tt ChainCodeGetRadius}$

${\tt ChainCodeGetLongestRun}$

Chain code measurement function

SYNOPSIS

```
#include "dip_chaincode.h"
dip_Error dip_ChainCodeGetLongestRun ( chaincode, longestRun )
```

FUNCTION

Returns the number of pixels in the longest run of identical codes in chaincode. This represents the longest straight section of the boundary. ChainCodeGetLongestRun is the function used by Measure for the FeatureLongestChaincodeRun measurement.

ARGUMENTS

Data type	Name	Description
dip_ChainCode	chaincode	Input chain code
dip_int *	longestRun	Receives the pixel count for the longest run

SEE ALSO

 ${\tt ImageChainCode}, {\tt ChainCodeNew}, {\tt ChainCodeFree}, {\tt ChainCodeGetLength}, {\tt ChainCodeGetFeret}, {\tt ChainCodeGetRadius}$

ChainCodeGetRadius

Chain code measurement function

SYNOPSIS

```
#include "dip_chaincode.h"
dip_Error dip_ChainCodeGetRadius ( chaincode, radius )
```

FUNCTION

This function computes statistics on the radius of an object. The centre of gravity of the object's border pixels is used as the centre of the object. The distance from each border pixel to this centre is is computed. The maximum, minimum, mean and variance of these distances are returned in the radius structure, which contains the following elements:

Data type	Name	Description
dip_float	max	Maximum object radius
dip_float	mean	Mean object radius
dip_float	min	Minimum object radius
dip_float	var	Variance of object radius

ChainCodeGetRadius is the function used by Measure for the FeatureRadius measurement.

ARGUMENTS

Data type	Name	Description
dip_ChainCode	chaincode	Input chain code
dip_CCRadius *	radius	Output measurement

SEE ALSO

 ${\tt ImageChainCode}, {\tt ChainCodeNew}, {\tt ChainCodeFree}, {\tt ChainCodeGetLength}, {\tt ChainCodeGetLength}, {\tt ChainCodeGetLength}, {\tt ChainCodeGetLength}, {\tt ChainCodeGetFeret}, {\tt Ch$

ChainCodeGetSize

Chain code access function

SYNOPSIS

```
#include "dip_chaincode.h"
dip_Error dip_ChainCodeGetSize ( chaincode, number )
```

FUNCTION

Returns the number of elements in the chain code.

Note: this is not a correct measure for the object's perimeter, use ChainCodeGetLength instead.

ARGUMENTS

Data type	Name	Description	
dip_ChainCode	chaincode	Chain code	
dip_int *	number	Receives the number of elements in the chain.	

SEE ALSO

 ${\tt ImageChainCode}, ChainCodeNew, ChainCodeFree, ChainCodeGetChains, ChainCodeGetStart, ChainCodeGetLabel, ChainCodeGetConnectivity$

ChainCodeGetStart

Chain code access function

SYNOPSIS

```
#include "dip_chaincode.h"
dip_Error dip_ChainCodeGetStart ( chaincode, startX, startY )
```

FUNCTION

Returns the start coordinates of the chain.

ARGUMENTS

Data type	Name	Description	
dip_ChainCode	chaincode	Chain code	
dip_int *	startX	Receives the start x-coordinate in chaincode	
dip_int *	startY	Receives the start y-coordinate in chaincode	

SEE ALSO

 ${\tt ImageChainCode}, {\tt ChainCodeNew}, {\tt ChainCodeFree}, {\tt ChainCodeGetSize}, {\tt ChainCodeGetChains}, {\tt ChainCodeGetLabel}, {\tt ChainCodeGetConnectivity}$

ChainCodeNew

Chain code object allocation

SYNOPSIS

```
#include "dip_chaincode.h"
dip_Error dip_ChainCodeNew ( chaincode, resources )
```

FUNCTION

Allocates an object of type dip_ChainCode. However, since its fields are private and currently there exist only read access functions, it is of little use creating such an object.

A dip_ChainCode object stores the following data:

Data type	Name	Description
dip_int	startX	Start coordinates for chain, ChainCodeGetStart
dip_int	startY	Start coordinates for chain, ChainCodeGetStart
dip_int	label	Label ID of object, ChainCodeGetLabel
$\mathtt{dip}_{-}\mathtt{int}$	connectivity	Connectivity of chain, ChainCodeGetConnectivity
dip_int	number	Number of elements in chain, ChainCodeGetSize
dip_Chain *	chain	Pointer to first element in chain, ChainCodeGetChains

The dip_Chain structure has the following elements:

Data type	Name	Description
dip_uint8	code	Direction of step taken from previous to this pixel (Freeman code)
dip_Boolean	border	Pixel is on the border
dip_Chain *	next	Pointer to the next element in the chain

The chain parameter points to the first dip_Chain object in the chain, which points to the next through its next value. The last element in the chain has a NULL pointer.

Each chain element contains the code value (between 0 and 3 or between 0 and 7, depending on the connectivity) as well as a border value, which indicates whether the pixel is on the edge of the image or not. The border value is important because it indicates that the object is cut by the imaging window and needs to be treated differently.

The chain code for an object always has as many elements as the object has border pixels.

ARGUMENTS

Data type	Name	Description
dip_ChainCode *	chaincode	Receives pointer to allocated structure
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

ImageChainCode, ChainCodeFree, ChainCodeArrayNew, ChainCodeArrayFree, ChainCodeGetSize, ChainCodeGetChains, ChainCodeGetStart, ChainCodeGetLabel, ChainCodeGetConnectivity, ChainCodeGetLength, ChainCodeGetLongestRun, ChainCodeGetFeret

ChangeDataType

Change the data type of an image

SYNOPSIS

dip_Error dip_ChangeDataType(example, target, dataType)

FUNCTION

Inherit all properties of the input image except the data type. The data type is explicitly specified through dataType. When dataType is zero, the data type of the output image is not modified. The example image may be either "raw" or "forged".

ARGUMENTS

Data type	Name	Description
dip_Image	example	An example image
dip_Image	target	The target image
dip_DataType	dataType	The data type

SEE ALSO

DIPlib's data types

 ${\tt ImageCopyProperties,\ ImageAssimilate,\ ChangeToOd}$

ChangeDimensions

Changes the order of the dimensions in an image

SYNOPSIS

dip_Error dip_ChangeDimensions(image, neworder)

FUNCTION

Re-orders the dimensions in an image, optionally removing or adding singleton dimensions (those dimensions with size 1), without copying the data. **neworder** is a list of the dimension numbers in the new order, for example (1,0,2) will swap the first two dimensions. Setting **neworder** to 0 removes all singleton dimensions without altering the order. This is useful, for example, after calling a function such as **Maximum** to compute a maximum projection over one dimension. The output image of **Maximum** keeps the dimensionality of the input image, and thus has a singleton dimension. To add singleton dimensions, use a negative value in the **neworder** array. For example, (0,1,-1) adds a 3rd dimension of size 1 to the image.

ARGUMENTS

Data type	Name	Description
dip_Image	image	The image to modify
dip_IntegerArray	neworder	The new order of the dimensions

SEE ALSO

The image structure

ImageGetDimensions, ImageSetDimensions, ImageGetStride

ChangeTo0d

Make an image zero dimensional

SYNOPSIS

dip_Error dip_ChangeToOd(example, target, dataType)

FUNCTION

Inherit all properties of the input image except the data type and the dimensionality. The data type is explicitly specified through dataType. When dataType is zero, the data type of the output image is not modified. The dimensionality is set to zero. The example image may be either "raw" or "forged".

ARGUMENTS

Data type	Name	Description	
dip_Image	example	An example image	
dip_Image	target	The target image	
dip_DataType	dataType	The data type. See DIPlib's data types	

SEE ALSO

ImageCopyProperties, ImageAssimilate, ChangeDataType

ChordLength

Compute the chord lengths of the different phases

SYNOPSIS

```
#include "dip_analysis.h"
dip_Error dip_ChordLength ( object, mask, dist, probes, length, sampling )
```

DATA TYPES

binary, integer

FUNCTION

This function computes the chord lengths of the different phases in object. If object is a binary image, the image is a regarded as a two phase image. In case object is of the integer type, the image is regarded as a labeled image, with each integer value encoding a phase. Optionally a mask image can be provided to select which pixels in object should be used to compute the chord lengths. The probes variable specifies how many random point pairs should be drawn to compute the lengths. Length specifies the maximum correlation length. The correlation function can be computed using a random (DIP_CORRELATION_ESTIMATOR_RANDOM) or grid method (DIP_CORRELATION_ESTIMATOR_GRID), as specified by sampling.

ARGUMENTS

Data type	Name	Description
dip_Image	object	Object image
dip_Image	mask	Mask image
dip_Distribution	dist	Ouput distribution
dip_int	probes	Number of probes
dip_int	length	Maximum chord length
dipf_CorrelationEstimator	sampling	Samplings method

SEE ALSO

PairCorrelation, ProbabilisticPairCorrelation

CityBlockDistanceToPoint

Distance generation function

SYNOPSIS

```
#include "dip_generation.h"
dip_Error dip_CityBlockDistanceToPoint ( output, origin, scale )
```

DATA TYPES

Output: sfloat

FUNCTION

Computes the cityblock distance of each pixel in the output image to a point at origin. The coordinates of origin may lie outside the image. The scale parameter may be used to specify the relative distance between pixels in each dimension.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{\mathtt{-}}\mathtt{Image}$	output	Output Image
dip_FloatArray	origin	Origin
dip_FloatArray	scale	Relative scale of the pixel distances for each dimension

SEE ALSO

 ${\tt EllipticDistanceToPoint}, {\tt EuclideanDistanceToPoint}$

Clip

Point operation

SYNOPSIS

```
#include "dip_point.h"
dip_Error dip_Clip ( in, out, clipLow, clipHigh, clipFlag )
```

DATA TYPES

integer, **float**

FUNCTION

Clips in at either the minimum value clipLow of the maximum value clipHigh or both. If the flag DIP_CLIP_THRESHOLD_AND_RANGE is specified, the clip bound are defined by clipLow +/-clipHigh/2.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input image
dip_Image	out	Output image
dip_float	clipLow	Lower clip bound value
dip_float	clipHigh	Higher clip bound value
dipf_Clip	clipFlag	Clip flag

The following dipf_Clip flags are defined:

Name	Description
DIP_CLIP_BOTH	clip both the lower and upper bound
DIP_CLIP_LOW	clip lower bound only
DIP_CLIP_HIGH	clip upper bound only
DIP_CLIP_THRESHOLD_AND_RANGE	use clipLow and clipHigh as threshold and range value
DIP_CLIP_LOW_AND_HIGH_BOUNDS	same as DIP_CLIP_BOTH

SEE ALSO

 ${\tt Threshold}, {\tt RangeThreshold}, {\tt ErfClip}, {\tt ContrastStretch}$

Closing

Morphological closing operation

SYNOPSIS

```
#include "dip_morphology.h"
dip_Error dip_Closing ( in, out, se, boundary, param, shape )
```

DATA TYPES

integer, float, binary

FUNCTION

Grey-value closing with different structuring elements.

The rectangular, elliptic and diamond structuring elements are "flat", i.e. these structuring elements have a constant value. For these structuring elements, param determines the sizes of the structuring elements.

When shape is DIP_FLT_SHAPE_DISCRETE_LINE or DIP_FLT_SHAPE_INTERPOLATED_LINE, the structuring element is a line. param->array[0] determines the length, param->array[1] the angle. This is currently only supported for 2D images. Interpolated lines use interpolation to obtain a more accurate result, but loose the strict increasingness and extensivity (these properties are satisfied only by approximation).

When shape is set to DIP_FLT_SHAPE_PARABOLIC, params specifies the curvature of the parabola.

When shape is set to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, se is used as structuring element. It can be either a binary or a grey-value image. Its origin is the center, or one pixel to the left of the center if the size is even.

If shape is not equal to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, se can be set to zero. When shape is set to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, param is ignored, and can be set to zero.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	in	Input
dip_Image	out	Output
dip_Image	se	Custom structuring element
dip_BoundaryArray	boundary	Boundary conditions
dip_FloatArray	param	Filter parameters
dip_FilterShape	shape	Structuring element

The enumerator dip_FilterShape contains the following constants:

Name	Description	
DIP_FLT_SHAPE_DEFAULT	Default filter window, same as	
	DIP_FLT_SHAPE_RECTANGULAR	
DIP_FLT_SHAPE_RECTANGULAR	Rectangular filter window, can be even in size	
DIP_FLT_SHAPE_ELLIPTIC	Elliptic filter window, always odd in size	
DIP_FLT_SHAPE_DIAMOND	Diamond-shaped filter window, always odd in size	
DIP_FLT_SHAPE_PARABOLIC	Parabolic filter window (morphology only)	
DIP_FLT_SHAPE_DISCRETE_LINE	Rotated line structuring element (morphology only)	
DIP_FLT_SHAPE_INTERPOLATED_LINE	Rotated line structuring element, through interpolation	
	(morphology only)	
DIP_FLT_SHAPE_PERIODIC_LINE	(not implemented)	
DIP_FLT_SHAPE_STRUCTURING_ELEMENT	Use se as filter window, can be any size	

SEE ALSO

Opening, Dilation, Erosion

Colour2Gray

Convert ND image with colour information to a (n-1)D grayvalue image (in dipIO)

SYNOPSIS

```
#include "dipio_tools.h"
dip_Error dipio_Colour2Gray ( in, out, photometric )
```

FUNCTION

This function converts a colour image, as read by ImageReadColour, to a grayvalue intensity image. in is expected to contain the colour information along the last axis. out will be a scalar image with one less dimension than the input.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input image
dip_Image	out	Output image
dipio_PhotometricInterpretation	photometric	Photometric interpretation

The enumerator dipio_PhotometricInterpretation contains the following constants:

Name	Description
DIPIO_PHM_GREYVALUE	No colour information present; it's a grey-value image.
DIPIO_PHM_RGB	RGB image (the first three planes are red, green and blue)
DIPIO_PHM_RGB_NONLINEAR	Non-linear R'G'B' image (RGB channels to the power of 0.4)
DIPIO_PHM_CMY	CMY image (the first three planes are cyan, magenta and yellow)
DIPIO_PHM_CMYK	CMYK image (the first four planes are cyan, magenta, yellow and
	black)
DIPIO_PHM_CIELUV	CIE L*u'v' image (the first three planes are luminosity, u* and v*)
DIPIO_PHM_CIELAB	CIE L*a*b* image (the first three planes are luminosity, a* and b*)
DIPIO_PHM_CIEXYZ	CIE XYZ (the first three planes are X, Y and Z)
DIPIO_PHM_CIEYXY	CIE Yxy (the first three planes are Y, x and y)
DIPIO_PHM_HCV	HCV image (the first three planes are hue, chroma and value)
DIPIO_PHM_HSV	HSV image (the first three planes are hue, saturation and value)
DIPIO_PHM_DEFAULT	Same as DIPIO_PHM_GREYVALUE
DIPIO_PHM_GENERIC	Anything can be coded in the channels; the same as DIPIO_PHM_CMYK

Most file formats support only some of these.

KNOWN BUGS

Some colourspaces are not converted correctly. R'G'B' (DIPIO_PHM_RGB_NONLINEAR), is treated like RGB. From a CIE Lab (DIPIO_PHM_CIELAB) or Luv (DIPIO_PHM_CIELUV) the luminosity channel is

extracted, which is also a non-linear conversion away from the intensity. From HCV (DIPIO_PHM_HCV) and HSV (DIPIO_PHM_HSV) the value channel is extracted, which again is a non-linear conversion away from the intensity. CMYK (DIPIO_PHM_CMYK) and CMY (DIPIO_PHM_CMY) conversion is not implemented. Specifying these values will result in an error.

SEE ALSO

ImageRead, ImageReadColour, ImageReadROI

Compare

Compare grey values in two images

SYNOPSIS

dip_Error dip_Compare (in1, in2, out, selector)

DATA TYPES

binary, integer, float

FUNCTION

This function can perform various pixel-by-pixel comparisons (smaller, smaller- equal, equal, not equal, greater-equal, greater) between in1 ans in2. out contains the binary result. This is implemented with a call to Select whose in3 and in4 are set to binary true and false, respectively. in2 can be a 0D image for comparison of pixel values with a single scalar value. This leads to the functionality of Threshold, but with more options.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	in1	First input
dip_Image	in2	Second input
dip_Image	out	Output
dipf_Select	selector	Select flag

The $dipf_Select$ flag can be one of:

Name	Description
DIP_SELECT_LESSER	<, Lesser than
DIP_SELECT_LESSER_EQUAL	<=, Lesser or equal
DIP_SELECT_NOT_EQUAL	!=, Unequal
DIP_SELECT_EQUAL	==, Equal
DIP_SELECT_GREATER_EQUAL	>=, Greater or equal
DIP_SELECT_GREATER	>, Greater

SEE ALSO

 ${\tt Select, Threshold, Equal, Greater, Lesser, NotEqual, NotGreater, NotLesser, SelectValue, NotZero}$

ComplexArrayCopy

Copy an array

SYNOPSIS

dip_Error dip_ComplexArrayCopy (dest, src, resources)

FUNCTION

This function copies the complex array src to dest. The array dest is created by this function as well.

ARGUMENTS

Data type	Name	Description
dip_ComplexArray *	dest	Destination array
dip_ComplexArray	src	Source array
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

ComplexArrayNew, ComplexArrayFree, ComplexArrayCopy, ComplexArrayFind
IntegerArrayCopy, FloatArrayCopy, ComplexArrayCopy, DataTypeArrayCopy, BooleanArrayCopy,
VoidPointerArrayCopy, StringArrayCopy

ComplexArrayFind

Find value in array

SYNOPSIS

dip_Error dip_ComplexArrayFind (array, value, index, found)

FUNCTION

Finds a value in an array and "returns" its index in the array. If found is zero, ComplexArrayFind will produce an error if value is not found, otherwise found obtains the search result (DIP_FALSE if value is not found).

ARGUMENTS

Data type	Name	Description
dip_ComplexArray	array	Array to find value in
dip_complex	value	Value to find
dip_int *	index	Index of the found value
dip_Boolean *	found	Value found or not

SEE ALSO

ComplexArrayNew, ComplexArrayFree, ComplexArrayCopy, ComplexArrayFind
IntegerArrayFind, FloatArrayFind, ComplexArrayFind, DataTypeArrayFind, BooleanArrayFind,
VoidPointerArrayFind

ComplexArrayFree

Array free function

SYNOPSIS

dip_Error dip_ComplexArrayFree (array)

FUNCTION

This function frees *array, and sets array to zero.

ARGUMENTS

Data type	Name	Description
dip_ComplexArray *	array	Array

SEE ALSO

ComplexArrayNew, ComplexArrayFree, ComplexArrayCopy, ComplexArrayFind

ArrayFree, IntegerArrayFree, FloatArrayFree, ComplexArrayFree, BoundaryArrayFree, FrameWorkProcessArrayFree, DataTypeArrayFree, ImageArrayFree, BooleanArrayFree, VoidPointerArrayFree, StringArrayFree, CoordinateArrayFree

ComplexArrayNew

Array allocation function

SYNOPSIS

dip_Error dip_ComplexArrayNew (array, size, value, resources)

FUNCTION

This function allocates the size elements of a dip_ComplexArray and sets the size of the array to size. Each array element is initialized with value.

ARGUMENTS

Data type	Name	Description
dip_ComplexArray *	array	Array
dip_int	size	Size
dip_complex	value	Initial value
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

ComplexArrayNew, ComplexArrayFree, ComplexArrayCopy, ComplexArrayFind
ArrayNew, IntegerArrayNew, FloatArrayNew, ComplexArrayNew, BoundaryArrayNew,
FrameWorkProcessArrayNew, DataTypeArrayNew, ImageArrayNew, BooleanArrayNew,

VoidPointerArrayNew, StringArrayNew, CoordinateArrayNew

ContrastStretch

Point operation

SYNOPSIS

#include "dip_point.h"

dip_Error dip_ContrastStretch (in, out, lowerBound, upperBound, outMaximum, outMinimum, method, sigmoidSlope, sigmoidPoint, maxDecade)

DATA TYPES

integer, **float**

FUNCTION

ContrastStretch stretches the pixel values of the input image. Pixel values higher or equal to UpperBound are stretched to the OutMaximum value. A similar thing holds for LowerBound and OutMinimum. Method determines how pixel values are stretched. SigmoidSlope and SigmoidPoint are used by the DIP_CST_SIGMOID method. MaxDecade determines the maximum number of decades the method DIP_CST_DECADE will stretch (values lower than MaxDecade will be set to zero).

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	in	Input
dip_Image	out	Output
dip_float	lowerBound	LowerBound (%)
dip_float	upperBound	UpperBound (%)
dip_float	outMax	OutMaximum
dip_float	outMin	OutMinimum
dipf_ContrastStretch	method	Method
dip_float	sigmoidSlope	SigmoidSlope
dip_float	sigmoidPoint	SigmoidPoint
dip_float	maxDecade	MaxDecade

The following dipf_ContrastStretch flags are defined:

Name	Description
DIP_CST_LINEAR	linear contrast stretch
DIP_CST_SIGNED_LINEAR	linear stretch with zero at fixed value
DIP_CST_LOGARITHMIC	logarithmic contrast stretch
DIP_CST_SIGNED_LOGARITHMIC	signed logarithmic contrast stretch
DIP_CST_ERF	linear contrast stretch with erf clipping
DIP_CST_DECADE	Decade contrast stretching
DIP_CST_SIGMOID	Contrast stretched by sigmoid function
DIP_CST_CLIP	Simple clipping
DIP_CST_01	Stretching of [0,1] input values
DIP_CST_PI	Stretching of [-Pi,Pi] input values

In the explanaition of the different contrast stretch flags, the variables input, output, inMin, inMax, outMin and outMax are used. With input and output is meant the pixel being processed of respecitively the input and output image. inMin and inMax are the pixel values corresponding to the lowerBound and upperBound of the input image. outMin and outMax are parameters passed to the function dip_ContrastStretch.

The DIP_CST_LINEAR stretches the input in the following way:

```
scale = (outMax - outMin) / (inMax - inMin)
output = scale * (MIN(inMax, MAX(inMin, input )) - inMin) + outMin
```

The DIP_CST_SIGNED_LINEAR stretches the input in the following way:

```
max = MAX(inMax, ABS( inMin ));
scale = (outMax - outMin) / (2 * max)
offset = (outMax - outMin) / 2
output = scale * (MIN(inMax, MAX(inMin, input)) - offset) + outMin
```

The DIP_CST_LOGARITHMIC stretches the input in the following way:

```
scale = (outMax - outMin) / log( inMax - inMin + 1)
offset = inMin - 1
output = scale * log(MIN(inMax, MAX(inMin, input)) - offset) + outMin
```

The DIP_CST_SIGNED_LOGARITHMIC stretches the input in the following way:

```
max = MAX(inMax, ABS( inMin ))
scale = (outMax - outMin) / (2 * log( max + 1))
offset = (outMax + outMin) / 2
output = scale * log(MIN(inMax, MAX(inMin, input)) - offset) + outMin
```

The DIP_CST_ERF stretches the input in the following way:

```
scale = (outMax - outMin) / (inMax - inMin)
threshold = (inMax + inMin)/ 2
range = inMax - inMin
in = MIN(inMax, MAX(inMin, input))
out = (range / 2) * erf( SQRT_PI * (in - threshold) / range )
output = scale * (out + threshold ) + outMin
```

The DIP_CST_DECADE stretches the input in the following way:

```
inScale = inMax - inMin
outScale = outMax - outMin
in = MIN(inMax, DIP_MAX(inMin, input))
decade = log10(inScale / (in - inMin + EPSILON))
if(decade < maxDecade)
  decade -= floor(decade)
  output = outScale * (1 - decade) + outMin
else
  output = 0</pre>
```

The DIP_CST_SIGMOID stretches the input in the following way:

The DIP_CST_CLIP stretches the input in the following way:

```
output = MIN(outMax, MAX(outMin, input))
```

The DIP_CST_01 stretches the input in the following way:

```
scale = (outMax - outMin)
output = scale * input + outMin
```

The DIP_CST_01 stretches the input in the following way:

```
scale = (outMax - outMin) / 2 * Pi
output = scale * (input + Pi) + outMin
```

SEE ALSO

See section 9.1, "Histogram-based operations", in Fundamentals of Image Processing. Clip, ErfClip

ConvertArray

converts the data type of an array

SYNOPSIS

```
#include "dip_convert_array.h"
dip_Error DIP_TWO_FUNC(dip_ConvertArray)( in, inStride, inPlane, out, outStride,
outPlane, number )
```

FUNCTION

Converts the in array to the out array.

ARGUMENTS

Data type	Name	Description	
void *	in	input array	
dip_int	inStride	Stride of the input array	
dip_int	inPlane	plane number in case in is a binary array	
void *	out	output array	
$\mathtt{dip}_{-}\mathtt{int}$	outStride	Stride of the output array	
$\mathtt{dip}_{-}\mathtt{int}$	outPlane	plane number in case out is a binary array	
dip_int	number	size of the arrays	

ConvertDataType

Converts the data type of an image

SYNOPSIS

dip_Error dip_ConvertDataType (in, out, dataType)

FUNCTION

Convert the data type of the input data to dataType and stores the result in out.

Conversion from a *complex* type to another (non-complex) type, is done by taking the real part.

Conversion to a binary type from another (non-binary) type, is done as follows; any non-zero number becomes 1, zero becomes zero.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	out	Output
dip_DataType	dataType	Data type. See DIPlib's data types

ConvexHullGetArea

Convex hull measurement function

SYNOPSIS

```
#include "dip_chaincode.h"
dip_Error dip_ConvexHullGetArea ( polygon, area )
```

FUNCTION

ConvexHullGetArea is the function used by Measure for the FeatureConvexArea measurement.

ARGUMENTS

Data type	Name	Description
dip_Polygon	polygon	Input convex polygon
dip_float*	area	Output measurement

SEE ALSO

ChainCodeConvexHull, ConvexHullGetPerimeter, ConvexHullGetFeret

ConvexHullGetFeret

Convex hull measurement function

SYNOPSIS

```
#include "dip_chaincode.h"
dip_Error dip_ConvexHullGetFeret ( polygon, feret )
```

FUNCTION

This function measures the longest and shortest projections of the object encoded by polygon. A "rotating calipers" algorithm finds all antipodal edges and vertices. Then the distances between each of these pairs is computed. These distances correspond to the lengths of the projection under an orientation perpendicular to the edge used. The maximum and minimum distance, as well as the corresponding orientations, are returned in the feret structure, which contains the following elements:

Data type	Name	Description
dip_float	maxDiameter	The widest projection of the object
dip_float	minDiameter	The narrowest projection of the object
dip_float	maxPerpendicular	The width of the projection perpendicular to minDiameter
dip_float	maxAngle	The angle of the projection for maxDiameter
dip_float	minAngle	The angle of the projection for minDiameter

ConvexHullGetFeret is the function used by Measure for the FeatureFeret measurement.

ARGUMENTS

Data type	Name	Description	
dip_Polygon	polygon	Input convex polygon	
dip_Feret*	feret	Output measurements	

NOTE

This function is more accurate than ChainCodeGetFeret, given a correct polygonal representation of the convex hull of the object. Because antipodal pairs are identified, and the angle of the edges is used, this algorithm doesn't depend on an angle step size, as ChainCodeGetFeret does.

LITERATURE

Algorithm A3.7 in M.I. Shamos, "Computational geometry," Ph.D. thesis, Yale University, 1978.

SEE ALSO

 ${\tt ChainCodeConvexHull}, {\tt ConvexHullGetArea}, {\tt ConvexHullGetPerimeter}, {\tt ChainCodeGetFeret}$

ConvexHullGetPerimeter

Convex hull measurement function

SYNOPSIS

```
#include "dip_chaincode.h"
dip_Error dip_ConvexHullGetPerimeter ( polygon, length )
```

FUNCTION

 ${\tt ConvexHullGetPerimeter} \ is \ the \ function \ used \ by \ {\tt Measure} \ for \ the \ {\tt FeatureConvexPerimeter} \ measurement.$

ARGUMENTS

Data type	Name	Description		
dip_Polygon	polygon	Input convex polygon		
dip_float*	length	Output measurement		

SEE ALSO

 ${\tt ChainCodeConvexHull}, {\tt ConvexHullGetArea}, {\tt ConvexHullGetFeret}, {\tt ChainCodeGetLength}$

Convolve1d

Perform a 1D convolution

SYNOPSIS

```
#include "dip_linear.h"
dip_Error DIP_TPI_FUNC(dip_Convolve1d)( in, out, filter, size, filterSize, origin, flags, boundary )
```

DATA TYPES

integer, float

FUNCTION

This function performs a one-dimensional convolution of the input data with the given filter kernel. In general your filter will be centered around the origin. The origin is uniquely defined if the filter size is odd, but if the filter size is even you'll have to specify whether the origin of the filter lies to the left or the right. Words cannot possibly suffice here, so here is a small pictorial representation:

When the filter size is even, one of the flags <code>DIP_CNV_LEFT</code> or <code>DIP_CNV_RIGHT</code> must be specified. When the filter size is odd both flags are ignored. It is also possible to specify the origin of the filter directly by using the <code>DIP_CNV_USE_ORIGIN</code> flag in combination with the <code>origin</code> parameter. Again a small pictorial representation:

when DIP_CNV_USE_ORIGIN is NOT specified origin is computed as follows:

The input data is copied to a temporary buffer, after which the input data is extended according to the boundary condition specified. You can use the flags DIP_CNV_HAS_BORDER to indicate that the input data already has a border. In this case you must make sure that there are enough pixels on either side of the array:

```
on the left : ((filterSize - 1) - origin) pixels
on the right : (origin) pixels
```

If DIP_CNV_HAS_BORDER is specified and in != out no auxiliary storage is used.

You must also specify the symmetry of the filter as follows:

```
odd filter size
                                       DIP_CNV_EVEN
                     a b
                           c b
                        b
                          c -b -a
                                       DIP_CNV_ODD
                              d
                                       DIP_CNV_GENERAL
even filter size
                                       DIP_CNV_EVEN
                     a b
                          С
                              c b a
                                       DIP_CNV_ODD
                     a b c -c -b -a
                     a b c d e f
                                       DIP_CNV_GENERAL
```

ARGUMENTS

Data type	Name	Description	
void *	in	Pointer to the input data	
void *	out	Pointer to the output data	
void *	filter	Pointer to the filter data	
dip_int	size	Size of the input data	
dip_int	filterSize	filterSize Size of the filter	
dip_int	origin	Origin of the filter. Only valid in conjunction with	
		DIP_CNV_USE_ORIGIN	
dipf_Convolve	flags	A combination of the flags described above	
dip_Boundary	boundary	One of the standard boundary conditions. See Boundary	
		conditions	

SEE ALSO

General information about convolution

 ${\tt Separable Convolution}, {\tt Separable Frame Work}$

ConvolveFT

Fourier transform-based convolution filter

SYNOPSIS

```
#include "dip_linear.h"
dip_Error dip_ConvolveFT ( in, psf, out, inrep, psfrep, outrep )
```

DATA TYPES

binary, integer, float, complex

FUNCTION

This function convolves the in image with the point spread function psf, by multiplying their Fourier transforms. The inrep, psfrep and outrep specify whether the images are spatial images (DIP_IMAGE_REPRESENTATION_SPATIAL) or their Fourier transform. (DIP_IMAGE_REPRESENTATION_SPECTRAL).

out is cast to a real type if and only if both in and psf are real and in the spatial domain. That is, no effort is made to check for evenness of images in the Fourier domain, nor to check the values of the imaginary component of the result. To convert the output to a real-valued type, use the function ConvertDataType.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input image
dip_Image	psf	Psf image
dip_Image	out	Output image
${\tt dipf_ImageRepresentation}$	inrep	Input spatial or spectral
${ t dipf_{-}ImageRepresentation}$	psfrep	PSF spatial or spectral
dipf_ImageRepresentation	outrep	Output spatial or spectral

SEE ALSO

General information about convolution

CoordinateArrayFree

Array free function

SYNOPSIS

dip_Error dip_CoordinateArrayFree (array)

FUNCTION

This function frees *array, and sets array to zero.

ARGUMENTS

Data type	Name	Description
dip_CoordinateArray * array		Array

SEE ALSO

CoordinateArrayNew, CoordinateArrayFree

ArrayFree, IntegerArrayFree, FloatArrayFree, ComplexArrayFree, BoundaryArrayFree, FrameWorkProcessArrayFree, DataTypeArrayFree, ImageArrayFree, BooleanArrayFree, VoidPointerArrayFree, StringArrayFree, CoordinateArrayFree

CoordinateArrayNew

Array allocation function

SYNOPSIS

dip_Error dip_CoordinateArrayNew (array, ndims, size, resources)

FUNCTION

This function allocates the size elements of a dip_CoordinateArray and sets the size of the array to size. Each element has ndims values, to store coordinates of an ndims-dimensional image. Each array element is initialized to 0.

ARGUMENTS

Data type	Name	Description
dip_CoordinateArray *	array	Array
dip_int	ndims	Dimensionality
dip_int	size	Size
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

CoordinateArrayNew, CoordinateArrayFree

ArrayNew, IntegerArrayNew, FloatArrayNew, ComplexArrayNew, BoundaryArrayNew, FrameWorkProcessArrayNew, DataTypeArrayNew, ImageArrayNew, BooleanArrayNew, VoidPointerArrayNew, StringArrayNew, CoordinateArrayNew

CoordinateToIndex

Convert coordinate to pixel index

SYNOPSIS

```
#include "dip_coordsindx.h"
dip_Error dip_CoordinateToIndex ( coordinates, index, stride )
```

FUNCTION

This function converts a pixel coordinate to an pixel index which is specific for the image from which stride was obtained. coordinages and stride must have the same number of elements.

ARGUMENTS

Data type	Name	Description
dip_IntegerArray	coordinates	Coordinate array
dip_int *	index	Pointer to pixel index
dip_IntegerArray	stride	stride array

SEE ALSO

IndexToCoordinate

Cos

trigonometric function

SYNOPSIS

dip_Error dip_Cos (in, out)

DATA TYPES

binary, integer, float, complex

FUNCTION

Computes the cosine of the input image values.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	in	Input
$dip_{-}Image$	out	Output

SEE ALSO

 ${\tt Sin,\ Cos,\ Tan,\ Asin,\ Acos,\ Atan,\ Atan2,\ Sinh,\ Cosh,\ Tanh}$

Cosh

trigonometric function

SYNOPSIS

dip_Error dip_Cosh (in, out)

DATA TYPES

binary, integer, ${f float}$

FUNCTION

Computes the hyperbolic cosine of the input image values.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	out	Output

SEE ALSO

Sin, Cos, Tan, Asin, Acos, Atan, Atan2, Sinh, Tanh

Crop

Remove the outer parts of an image

SYNOPSIS

```
#include "dip_manipulation.h"
dip_Error dip_Crop ( in, out, origin, size )
```

DATA TYPES

binary, integer, float, complex

FUNCTION

Crop a part of the image. The requested part is selected by specifying its upper left corner (origin), and its size (size). If in has a different type than out, it will be converted to the type of out.

ARGUMENTS

Data type	Name	Description		
dip_Image	in	Input Image		
dip_Image	out	Output Image		
dip_IntegerArray	origin	Coordinate in in of the upper left corner of the section		
dip_IntegerArray	size Size of the new image			

SEE ALSO

GetSlice, GetLine

CrossCorrelationFT

Normalized cross-correlation using the Fourier Transform

SYNOPSIS

```
#include "dip_findshift.h"
dip_Error dip_CrossCorrelationFT ( in1, in2, out, in1rep, in2rep, outrep )
```

DATA TYPES

binary, integer, float, complex

FUNCTION

This function calculates the cross-correlation between two images of equal size. The returned image is the cross-correlation normalized in such a way that only the phase information is of importance. This results as a very sharp peak in the spatial domain. This function performs out = (Conj(in1)*in2)/((Abs(in1))^2) in the Fourier domain. It is used by FindShift. The inrep, psfrep and outrep specify whether the images are spatial images (DIP_IMAGE_REPRESENTATION_SPATIAL) or their Fourier transform. (DIP_IMAGE_REPRESENTATION_SPECTRAL).

ARGUMENTS

Data type	Name	Description
dip_Image	in1	Input image
dip_Image	in2	Input image
dip_Image	out	Output image
dipf_ImageRepresentation	in1rep	Input 1 spatial or spectral
dipf_ImageRepresentation	in2rep	Input 2 spatial or spectral
dipf_ImageRepresentation	outrep	Output spatial or spectral

SEE ALSO

FindShift

CumulativeSum

statistics function

SYNOPSIS

dip_Error dip_CumulativeSum (in, mask, out, ps)

DATA TYPES

binary, integer, float, complex

FUNCTION

Calculates the cumulative sum of the pixel values over all those dimensions which are specified by ps, i.e.:

 $\operatorname{out}(x,y)=\operatorname{sum_i=0:x}, j=0:y \text{ in}(i,j) \text{ when ps specifies both x and y out}(x,y)=\operatorname{sum_j=0:y in}(x,j) \text{ when ps specifies only y}$

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	mask (0)	Mask
$\mathtt{dip}_{-}\mathtt{Image}$	out	Output
dip_BooleanArray	ps (0)	Dimensions to project

SEE ALSO

From images to scalars

Mean, Variance, StandardDeviation, MeanModulus, SumModulus, MeanSquareModulus, Maximum, Minimum, Median, Percentile

DanielsonLineDetector

Line detector

SYNOPSIS

#include "dip_orientation.h"
dip_Error dip_DanielsonLineDetector (in, line, energy, angle, boundary, sigma,
truncation, flavour)

DATA TYPES

binary, integer, float

FUNCTION

The Danielson line dectector uses second derivatives to detect lines in 2D images and to estimate their orientation. See the literature reference for an in-depth information on this detector.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input image
dip_Image	line	Line image
dip_Image	energy	Energy image
dip_Image	angle	Angle image
dip_BoundaryArray	boundary	Boundary conditions
dip_FloatArray	sigma	Sigma of second derivatives
dip_float	truncation	Gauss Truncation, see GlobalGaussianTruncationGet
dip_DerivativeFlavour	flavour	Derivative filter flavour

LITERATURE

P.E. Danielson, Q. Lin and Q-Z Yes, i "Efficient detection of second degree variations in 2D and 3D images", Report LiTH-ISY-R-2155, Linkoping University, Linkoping, Sweden, 1999

SEE ALSO

Derivative, StructureTensor2D

DataTypeAllowed

Check whether a data type is allowed

SYNOPSIS

dip_Error dip_DataTypeAllowed(dataType, allow, allowedTypes, allowed)

FUNCTION

This function checks whether the dataType is (or is not) in the set of data types specified by allowedTypes. If allow is DIP_TRUE, the data type should be in this set. If allow is DIP_FALSE, the data type should not be in this set. If the allowed parameter is zero, the routine returns dip_errorDataTypeNotSupported if the required condition is not satisfied. If nonzero, it should point to a boolean variable. This boolean variable will be set to DIP_TRUE if the condition is satisfied, or DIP_FALSE if not.

ARGUMENTS

Data type	Name	Description
dip_DataType	dataType	The data type to check
dip_Boolean	allow	DIP_TRUE: check if the data type is included.
		DIP_FALSE: check if the data type is not included
dip_DataTypeProperties	allowedTypes	The set of data types to check against, see
		DataTypeGetInfo
dip_Boolean *	allowed	Pointer to a boolean to store the answer, or 0 to
		indicate that dip_errorDataTypeNotSupported
		should be returned if the condition is not satisfied

SEE ALSO

DIPlib's data types
DataTypeGetInfo

${\tt DataTypeArrayCopy}$

Copy an array

SYNOPSIS

dip_Error dip_DataTypeArrayCopy (dest, src, resources)

FUNCTION

This function copies the data type array **src** to **dest**. The array **dest** is created by this function as well.

ARGUMENTS

Data type	Name	Description
dip_DataTypeArray *	dest	Destination array
dip_DataTypeArray	src	Source array
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

DIPlib's data types

DataTypeArrayNew, DataTypeArrayFree, DataTypeArrayCopy, DataTypeArrayFind
IntegerArrayCopy, FloatArrayCopy, ComplexArrayCopy, DataTypeArrayCopy, BooleanArrayCopy,
VoidPointerArrayCopy, StringArrayCopy

DataTypeArrayFind

Find value in array

SYNOPSIS

dip_Error dip_DataTypeArrayFind (array, value, index, found)

FUNCTION

Finds a value in an array and "returns" its index in the array. If found is zero, DataTypeArrayFind will produce an error if value is not found, otherwise found obtains the search result (DIP_FALSE if value is not found).

ARGUMENTS

Data type	Name	Description
dip_DataTypeArray	array	Array to find value in
dip_DataType	value	Value to find
dip_int *	index	Index of the found value
dip_Boolean *	found	Value found or not

SEE ALSO

DIPlib's data types

DataTypeArrayNew, DataTypeArrayFree, DataTypeArrayCopy, DataTypeArrayFind
IntegerArrayFind, FloatArrayFind, ComplexArrayFind, DataTypeArrayFind, BooleanArrayFind,
VoidPointerArrayFind

DataTypeArrayFree

Array free function

SYNOPSIS

dip_Error dip_DataTypeArrayFree (array)

FUNCTION

This function frees *array, and sets array to zero.

ARGUMENTS

Data type	Name	Description
dip_DataTypeArray *	array	Array

SEE ALSO

DIPlib's data types

DataTypeArrayNew, DataTypeArrayFree, DataTypeArrayCopy, DataTypeArrayFind
ArrayFree, IntegerArrayFree, FloatArrayFree, ComplexArrayFree, BoundaryArrayFree,
FrameWorkProcessArrayFree, DataTypeArrayFree, ImageArrayFree, BooleanArrayFree,
VoidPointerArrayFree, StringArrayFree, CoordinateArrayFree

DataTypeArrayNew

Array allocation function

SYNOPSIS

dip_Error dip_DataTypeArrayNew (array, size, value, resources)

FUNCTION

This function allocates the size elements of a dip_DataTypeArray and sets the size of the array to size. Each array element is initialized with value.

ARGUMENTS

Data type	Name	Description
dip_DataTypeArray *	array	Array
dip_int	size	Size
dip_DataType	value	Initial value
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

DIPlib's data types

DataTypeArrayNew, DataTypeArrayFree, DataTypeArrayCopy, DataTypeArrayFind

${\tt DataTypeGetInfo}$

Get information about a data type

SYNOPSIS

dip_Error dip_DataTypeGetInfo(dataType, info, whatInfo)

ARGUMENTS

Data type	Name	Description
dip_DataType	dataType	The data type to get information about
void *	info	Pointer to a variable to put the information in
dipf_DataTypeGetInfo	whatInfo	What information should be returned

FUNCTION

Get information about a data type. Depending on the whatInfo flag this routine will return information about the data type through the info parameter. A pointer must be passed to this routine which must point to a variable of the proper type to contain the information which will be returned. This pointer is passed as a void pointer through the info parameter. Below is a table of the flags that determine what information is returned, the type of the variable that is used to store the information in and a description of the information that is returned.

dipf_DataTypeGetInfo	type	description
DIP_DT_INFO_PROPS	dip_DataTypeProperties	a set of flags as shown in the
		table below
DIP_DT_INFO_SIZEOF	dip_int	sizeof(data type)
DIP_DT_INFO_C2R	dip_DataType	for complex types returns the
		corresponding floating point
		type (i.e. dip_scomplex ->
		dip_sfloat) for other data
		types returns the data type
		itself

The following table shows which dip_DataTypeProperties flags are set for which data types:

Data type identifier group	data types
DIP_DT_IS_UINT	unsigned integer
DIP_DT_IS_UNSIGNED	unsigned integer
DIP_DT_IS_SINT	signed integer
DIP_DT_IS_INT	signed and unsigned integer
DIP_DT_IS_INTEGER	signed and unsigned intege
DIP_DT_IS_FLOAT	floating-point
DIP_DT_IS_REAL	integer and floating-point
DIP_DT_IS_COMPLEX	complex floating-point
DIP_DT_IS_SIGNED	signed integer, floating-point and complex
DIP_DT_IS_BINARY	binary
DIP_DT_IS_ANY	all

SEE ALSO

DIPlib's data types

Derivative

Derivative filter

SYNOPSIS

```
#include "dip_derivatives.h"
dip_Error dip_Derivative ( in, out, boundary, ps, sigmas, order, truncation, flavour)
```

DATA TYPES

Depends on the underlying implementation, but expect: binary, integer, **float**

FUNCTION

This function provides a common interface to different families of regularised derivative operators. Which family is used, is specified by the flavour parameter. The order of the derivative operator along each of the cartesian axes may be specified independently.

Be sure to read the documentation on the underlying implementation to learn about the properties and limitations of the various families.

For the Gaussian family of filters, sigmas must be given, but order can be 0 (only smooth, don't take the derivative).

For the finite difference filter, sigmas can be 0, in which case the non-derivative dimensions will not be processed. Any element of sigmas that is non-zero where the corresponding order is zero, indicates a dimension that will be smoothed. Note it's possible to reporduce the SobelGradient filter this way.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	out	Output
dip_BoundaryArray	bc	Boundary conditions
dip_BooleanArray	ps (0)	Dimensions to process
dip_FloatArray	sigmas	Sigma of Gaussian
dip_int	order (0)	Derivative order
dip_float	truncation	Truncation, see GlobalGaussianTruncationGet
dip_DerivativeFlavour	flavour	Derivative filter flavour

The enumerator flavour parameter is one of:

Name	Description
DIP_DF_DEFAULT	Default derivative flavour (==DIP_DF_FIRGAUSS)
DIP_DF_FIRGAUSS	Gaussian family, FIR implementation, Gauss
DIP_DF_IIRGAUSS	Gaussian family, IIR implementation, GaussIIR
DIP_DF_FTGAUSS	Gaussian family, FT implementation, GaussFT
DIP_DF_FINITEDIFF	Finite difference implementation, FiniteDifferenceEx

SEE ALSO

See section 9.5, "Derivative-based operations", in Fundamentals of Image Processing.

Gauss, GaussFT, GaussIIR, FiniteDifferenceEx, GradientMagnitude, GradientDirection2D,
Laplace, SobelGradient

Dgg

Second order derivative filter

SYNOPSIS

```
#include "dip_derivatives.h"
dip_Error dip_Dgg ( in, out, boundary, ps, sigmas, tc, flavour )
```

DATA TYPES

Depends on the underlying implementation, but expect: binary, integer, **float**

FUNCTION

Computes the second derivative in gradient direction of an image using the Derivative function.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	out	Output
dip_BoundaryArray	boundary	Boundary conditions
dip_BooleanArray	ps	Dimensions to process
dip_FloatArray	sigmas	Sigma of Gaussian
dip_float	tc	Truncation of Gaussian, see
		GlobalGaussianTruncationGet
dip_DerivativeFlavour	flavour	Derivative flavour

The enumerator flavour parameter is one of:

Name	Description
DIP_DF_DEFAULT	Default derivative flavour (==DIP_DF_FIRGAUSS)
DIP_DF_FIRGAUSS	Gaussian family, FIR implementation, Gauss
DIP_DF_IIRGAUSS	Gaussian family, IIR implementation, GaussIIR
DIP_DF_FTGAUSS	Gaussian family, FT implementation, GaussFT
DIP_DF_FINITEDIFF	Finite difference implementation, FiniteDifferenceEx

SEE ALSO

See section 9.5, "Derivative-based operations", in Fundamentals of Image Processing (Dgg is called SDGD in the text).

Derivative, GradientMagnitude, GradientDirection2D, Laplace, LaplacePlusDgg,

LaplaceMinDgg

Dilation

Local maximum filter

SYNOPSIS

```
#include "dip_morphology.h"
dip_Error dip_Dilation ( in, out, se, boundary, param, shape )
```

DATA TYPES

integer, float, binary

FUNCTION

Grey-value dilation with different structuring elements.

The rectangular, elliptic and diamond structuring elements are "flat", i.e. these structuring elements have a constant value. For these structuring elements, param determines the sizes of the structuring elements.

When shape is DIP_FLT_SHAPE_DISCRETE_LINE or DIP_FLT_SHAPE_INTERPOLATED_LINE, the structuring element is a line. param->array[0] determines the length, param->array[1] the angle. This is currently only supported for 2D images. Interpolated lines use interpolation to obtain a more accurate result, but loose the strict increasingness and extensivity (these properties are satisfied only by approximation).

When shape is set to DIP_FLT_SHAPE_PARABOLIC, params specifies the curvature of the parabola.

When shape is set to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, se is used as structuring element. It can be either a binary or a grey-value image. Its origin is the center, or one pixel to the left of the center if the size is even.

If shape is not equal to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, se can be set to zero. When shape is set to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, param is ignored, and can be set to zero.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	out	Output
dip_Image	se	Custom structuring element
dip_BoundaryArray	boundary	Boundary conditions
dip_FloatArray	param	Filter parameters
dip_FilterShape	shape	Structuring element

The enumerator dip_FilterShape contains the following constants:

Name	Description
DIP_FLT_SHAPE_DEFAULT	Default filter window, same as
	DIP_FLT_SHAPE_RECTANGULAR
DIP_FLT_SHAPE_RECTANGULAR	Rectangular filter window, can be even in size
DIP_FLT_SHAPE_ELLIPTIC	Elliptic filter window, always odd in size
DIP_FLT_SHAPE_DIAMOND	Diamond-shaped filter window, always odd in size
DIP_FLT_SHAPE_PARABOLIC	Parabolic filter window (morphology only)
DIP_FLT_SHAPE_DISCRETE_LINE	Rotated line structuring element (morphology only)
DIP_FLT_SHAPE_INTERPOLATED_LINE	Rotated line structuring element, through interpolation
	(morphology only)
DIP_FLT_SHAPE_PERIODIC_LINE	(not implemented)
DIP_FLT_SHAPE_STRUCTURING_ELEMENT	Use se as filter window, can be any size

SEE ALSO

Closing, Opening, Erosion

${\tt dip_PixelGetFloat}$

Midlevel PixelIO function

SYNOPSIS

dip_Error dip_PixelGetFloat (vptr, type, position, stride, plane, val)

FUNCTION

The dip_PixelGet/SetInteger and dip_PixelGet/SetFloat functions provide midlevel access to image pixel values. These functions are faster than the highlevel Get and Set functions, but are easier to use than the lowlevel DIP_PIXEL_GET and DIP_PIXEL_SET macros as defined in dip_macros.h.

ARGUMENTS

Data type	Name	Description
void *	vptr	Void pointer to the image data
dip_DataType	type	Image data type. See DIPlib's data types
dip_IntegerArray	position	Position of the pixel in the image
dip_IntegerArray	stride	Image data stride
dip_int	plane	Plane of the pixel (binary images)
dip_float *	val	Pointer to the variable receiving the obtained pixel value

SEE ALSO

dip__PixelGetInteger, dip__PixelSetInteger, dip__PixelSetFloat, Get, Set, GetInteger,
SetInteger, GetFloat

dip__PixelGetInteger

Midlevel PixelIO function

SYNOPSIS

dip_Error dip_PixelGetInteger (vptr, type, position, stride, plane, val)

FUNCTION

The dip_PixelGet/SetInteger and dip_PixelGet/SetFloat functions provide midlevel access to image pixel values. These functions are faster than the highlevel Get and Set functions, but are easier to use than the lowlevel DIP_PIXEL_GET and DIP_PIXEL_SET macros as defined in dip_macros.h.

ARGUMENTS

Data type	Name	Description
void *	vptr	Void pointer to the image data
dip_DataType	type	Image data type. See DIPlib's data types
dip_IntegerArray	position	Position of the pixel in the image
dip_IntegerArray	stride	Image data stride
dip_int	plane	Plane of the pixel (binary images)
dip_int *	val	Pointer to the variable receiving the obtained pixel value

SEE ALSO

dip__PixelGetFloat, dip__PixelSetInteger, dip__PixelSetFloat, Get, Set, GetInteger,
SetInteger, GetFloat

$dip_PixelSetFloat$

Midlevel PixelIO function

SYNOPSIS

dip_Error dip_PixelSetFloat (val, vptr, type, position, stride, plane)

FUNCTION

The dip_PixelGet/SetInteger and dip_PixelGet/SetFloat functions provide midlevel access to image pixel values. These functions are faster than the highlevel Get and Set functions, but are easier to use than the lowlevel DIP_PIXEL_GET and DIP_PIXEL_SET macros as defined in dip_macros.h.

ARGUMENTS

Data type	Name	Description
dip_float	val	Value to write to the pixel
void *	vptr	Void pointer to the image data
dip_DataType	type	Image data type. See DIPlib's data types
dip_IntegerArray	position	Position of the pixel in the image
dip_IntegerArray	stride	Image data stride
dip_int	plane	Plane of the pixel (binary images)

SEE ALSO

$$\label{linear_pixel} \begin{split} & \texttt{dip_PixelGetInteger}, \ \texttt{dip_PixelSetInteger}, \ \texttt{Get}, \ \texttt{Set}, \ \texttt{GetInteger}, \\ & \texttt{SetInteger}, \ \texttt{GetFloat}, \ \texttt{SetFloat} \end{split}$$

dip__PixelSetInteger

Midlevel PixelIO function

SYNOPSIS

dip_Error dip_PixelSetInteger (val, vptr, type, position, stride, plane)

FUNCTION

The dip_PixelGet/SetInteger and dip_PixelGet/SetFloat functions provide midlevel access to image pixel values. These functions are faster than the highlevel Get and Set functions, but are easier to use than the lowlevel DIP_PIXEL_GET and DIP_PIXEL_SET macros as defined in dip_macros.h.

ARGUMENTS

Data type	Name	Description
dip_int	val	Value to write to the pixel
void *	vptr	Void pointer to the image data
dip_DataType	type	Image data type. See DIPlib's data types
dip_IntegerArray	position	Position of the pixel in the image
dip_IntegerArray	stride	Image data stride
dip_int	plane	Plane of the pixel (binary images)

SEE ALSO

dip__PixelGetInteger, dip__PixelGetFloat, dip__PixelSetFloat, Get, Set, GetInteger,
SetInteger, GetFloat

DirectedPathOpening

Morphological filter

SYNOPSIS

```
#include "dip_morphology.h"
dip_Error DirectedPathOpening ( grey, mask, out, param, closing, constrained )
```

DATA TYPES

binary, integer, float

FUNCTION

Theoretically, the path opening can be written as the supremum of all the openings with each of the possible linear structuring elements of composed of a set number of pixels in the general orientation of param. The param parameter is interpreted as follows: length is set to max(param). direction is set to round(param/length). direction now contains only values 0, -1 or 1. A 90 degree cone is defined around the given direction, and this cone gives the neighbourhood connectivity. The structuring element is formed by length pixels connected according to this neighbourhood connectivity. For example, in 2D, if param is [10,0], the structuring element will be formed by 10 pixels connected either diagonally or horizontally. It will extend across exactly 10 horizontal pixels, but can vary in shape to adapt to local image content.

If closing is DIP_TRUE, the path closing will be computed instead of the opening.

If constrained is DIP_TRUE, the algorithm is modified as follows: Only one consecutive step is allowed in a direction other than the exact direction specified. For example, following the [10,0] example above, a diagonal step must be followed by at least one horizontal step. This avoids zig-zag lines, especially if the main direction is diagonal. It also reduces the maximal angle that a straight line can deviate from the chosen direction. The unconstrained algorithm will keep lines rotated by up to 45 degrees; the constrained algorithm limits this to 22.5 degrees.

The algorithm uses a boundary condition such that any line connected to the border is considered infinite in length. To constrain the lines to the image domain, set a **two** pixel border around the image to the minimum value, e.g. 0 (for the opening), or the maximum value, e.g. 255 (for the closing).

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	grey	Grey-value input image
$\mathtt{dip}_{\mathtt{-}}\mathtt{Image}$	mask	Mask image for ROI processing
dip_Image	out	Output image
dip_FloatArray	param	Size of structuring element
dip_Boolean	closing	DIP_FALSE for path opening, DIP_TRUE for path closing
dip_Boolean	constrained	DIP_TRUE for constrained paths, DIP_FALSE for the original
		path opening algorithm

LITERATURE

- H. Talbot and B. Appleton, Efficient complete and incomplete path openings and closings, Image and Vision Computing 25:416-425, 2007.
- C.L. Luengo Hendriks, Constrained and dimensionality-independent path openings, IEEE Transactions on Image Processing 19(6):1587-1595, 2010.

SEE ALSO

Opening, Closing, PathOpening, AreaOpening

DistributionSort

Sort a block of data

SYNOPSIS

```
#include "dip_sort.h"
dip_Error dip_DistributionSort ( data, size, dataType )
```

FUNCTION

Sorts a block of data (of size size and data type dataType) using the distribution sort algorithm.

ARGUMENTS

Data type	Name	Description
void *	data	Data
dip_int	size	Size
dip_DataType	dataType	Data type. See DIPlib's data types

SEE ALSO

General information about sorting

 $\label{lem:continuity} DistributionSortIndices, DistributionSortIndices 16, Sort, ImageSort, SortIndices, SortIndices 16, ImageSortIndices$

DistributionSortIndices

Sort indices to block of data

SYNOPSIS

```
#include "dip_sort.h"
dip_Error dip_DistributionSortIndices ( data, indices, size, dataType )
```

FUNCTION

Sorts a list of indices rather than the data itself using the distribution sort algorithm.

ARGUMENTS

Data type	Name	Description
void *	data	Data
dip_sint32 *	indices	Indices
dip_int	size	Size
dip_DataType	dataType	Data type. See DIPlib's data types

SEE ALSO

General information about sorting

DistributionSort, DistributionSortIndices16, Sort, ImageSort, SortIndices, SortIndices16, ImageSortIndices

DistributionSortIndices16

Sort indices to a block of data

SYNOPSIS

```
#include "dip_sort.h"
dip_Error dip_DistributionSortIndices16 ( data, indices, size, dataType )
```

FUNCTION

Sorts a list of (16 bit) indices rather than the data itself using the distribution sort algorithm.

ARGUMENTS

Data type	Name	Description
void *	data	Data
dip_sint16 *	indices	Indices
dip_int	size	Size
dip_DataType	dataType	Data type. See DIPlib's data types

SEE ALSO

General information about sorting

 ${\tt DistributionSort, DistributionSortIndices, Sort, ImageSort, SortIndices, SortIndices16, ImageSortIndices}$

Div

arithmetic function

SYNOPSIS

```
dip_Error dip_Div ( in1, in2, out )
Calls Arith ( in1, in2, out, DIP_ARITHOP_DIV, DIP_DT_MINIMUM )
```

DivComplex

arithmetic function

SYNOPSIS

dip_Error dip_DivComplex (in, out, constant)

DATA TYPES

binary, integer, float, complex

FUNCTION

This function computes out = in / constant on a pixel by pixel basis. If constant is zero, out will be set to zero. The data types of the in1 image and constant may be of different types. See Information about dyadic operations for more information about what the type of the output will be.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	in	Input
dip_Image	out	Output
dip_complex	constant	Constant

SEE ALSO

Arith, Arith_ComplexSeparated, DivInteger, DivFloat, AddComplex, SubComplex, MulComplex, MulConjugateComplex

DivFloat

arithmetic function

SYNOPSIS

dip_Error dip_DivFloat (in, out, constant)

DATA TYPES

binary, integer, float, complex

FUNCTION

This function computes out = in / constant on a pixel by pixel basis. If constant is zero, out will be set to zero. The data types of the in1 image and constant may be of different types. See Information about dyadic operations for more information about what the type of the output will be.

ARGUMENTS

Data type	Name	Description	
$dip_{-}Image$	in	Input	
dip_Image	out	Output	
dip_float	constant	Constant	

SEE ALSO

Arith, Arith_ComplexSeparated, DivInteger, DivComplex, AddFloat, SubFloat, MulFloat

DivInteger

arithmetic function

SYNOPSIS

dip_Error dip_DivInteger (in, out, constant)

DATA TYPES

 $binary, \, \mathbf{integer}, \, \mathbf{float}, \, \mathbf{complex}$

FUNCTION

This function computes out = in / constant on a pixel by pixel basis. If constant is zero, out will be set to zero. The data types of the in1 image and constant may be of different types. See Information about dyadic operations for more information about what the type of the output will be.

ARGUMENTS

Data type	Name	Description	
dip_Image	in	Input	
dip_Image	out	Output	
dip_int	constant	Constant	

SEE ALSO

Arith, Arith_ComplexSeparated, DivFloat, DivComplex, AddInteger, SubInteger, MulInteger

EdgeObjectsRemove

Remove binary edge objects

SYNOPSIS

```
#include "dip_binary.h"
dip_Error dip_EdgeObjectsRemove ( in, out, connectivity )
```

DATA TYPES

binary

FUNCTION

The function EdgeObjectsRemove removes those binary objects from in which are connected to the edges of the image. The connectivity of the objects is determined by connectivity. This function is a front-end to BinaryPropagation. It calls BinaryPropagation with no seed image and the edge pixels turned on. The result of the propagation is xor-ed with the input image. The connectivity parameter defines the metric, that is, the shape of the structuring element. 1 indicates city-block metric, or a diamond-shaped structuring element. 2 indicates chessboard metric, or a square structuring element. -1 and -2 indicate alternating connectivity and produce an octagonal structuring element. See The connectivity parameter for more information. The edge parameter specifies whether the border of the image should be treated as object (DIP_TRUE) or as background (DIP_FALSE).

See section 10.3, "Segmentation", in Fundamentals of Image Processing for a description of the edge object removal operation.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Binary input image
dip_Image	out	Output
dip_int	connectivity	Pixel connectivity

KNOWN BUGS

This function is only implemented for images with a dimension up to three.

SEE ALSO

BinaryPropagation, Xor

${\tt EllipticDistanceToPoint}$

Distance generation function

SYNOPSIS

```
#include "dip_generation.h"
dip_Error dip_EllipticDistanceToPoint ( output, origin, scale )
```

DATA TYPES

Output: sfloat

FUNCTION

Computes the elliptic distance of each pixel in the output image to a point at origin. The coordinates of origin may lie outside the image. The scale parameter may be used to specify the relative distance between pixels in each dimension.

ARGUMENTS

Data type	Name	Description	
$\mathtt{dip}_{\mathtt{-}}\mathtt{Image}$	output	Output Image	
dip_FloatArray	origin	Coordinates of the Origin	
dip_FloatArray	scale	Relative scale of the pixel distances for each dimension	

SEE ALSO

 ${\tt Euclidean Distance To Point, City Block Distance To Point}$

Equal

Compare grey values in two images

SYNOPSIS

dip_Error dip_Equal (in1, in2, out)

DATA TYPES

binary, integer, float

FUNCTION

This function sets each pixel in out to "true" when corresponding pixels in in1 and in2 are equal. This is the same as Compare with the DIP_SELECT_EQUAL selector flag.

in 2 can be a 0D image for comparison of pixel values with a single scalar value. This leads to a functionality similar to SelectValue.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{\mathtt{-}}\mathtt{Image}$	in1	First input
$\mathtt{dip}_{-}\mathtt{Image}$	in2	Second input
dip_Image	out	Output

SEE ALSO

Compare, Threshold, Greater, Lesser, NotEqual, NotGreater, NotLesser, SelectValue, NotZero

Erf

mathematical function

SYNOPSIS

dip_Error dip_Erf (in, out)

DATA TYPES

binary, integer, ${f float}$

FUNCTION

Computes the error function of the input image values.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	in	Input
$dip_{-}Image$	out	Output

SEE ALSO

BesselJO, BesselJI, BesselJN, BesselYO, BesselYI, BesselYN, LnGamma, Erfc, Sinc

Erfc

mathematical function

SYNOPSIS

dip_Error dip_Erfc (in, out)

DATA TYPES

binary, integer, ${f float}$

FUNCTION

Computes the complementary error function of the input image values.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{\mathtt{-}}\mathtt{Image}$	in	Input
$\mathtt{dip}_{-}\mathtt{Image}$	out	Output

SEE ALSO

BesselJO, BesselJO, BesselJN, BesselYO, BesselYO, BesselYN, LnGamma, Erf, Sinc

ErfClip

Point Operation

SYNOPSIS

```
#include "dip_point.h"
dip_Error dip_ErfClip ( in, out, threshold, range, clipFlag )
```

DATA TYPES

integer, **float**

FUNCTION

Clips in using the erf function at either or both the values threshold +/- range/2. If the flag DIP_CLIP_LOW_AND_HIGH_BOUNDS is specified, threshold and range are used as lower and upper bounds respectively.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	in	Input image
dip_Image	out	Output image
dip_float	threshold	Threshold value
dip_float	range	Range value
dipf_Clip	clipFlag	clipFlag

The following dipf_Clip flags are defined:

Name	Description
DIP_CLIP_BOTH	clip both the lower and upper bound
DIP_CLIP_LOW	clip lower bound only
DIP_CLIP_HIGH	clip upper bound only
DIP_CLIP_THRESHOLD_AND_RANGE	same as DIP_CLIP_BOTH
DIP_CLIP_LOW_AND_HIGH_BOUNDS	use threshold and range as lower and upper bounds

LITERATURE

L.J. van Vliet, *Grey-Scale Measurements in Multi-Dimensional Digitized Images*, Ph.D. thesis Delft University of Technology, Delft University Press, Delft, 1993

SEE ALSO

Clip, ContrastStretch

Erosion

Local minimum filter

SYNOPSIS

```
#include "dip_morphology.h"
dip_Error dip_Erosion ( in, out, se, boundary, param, shape )
```

DATA TYPES

integer, float, binary

FUNCTION

Grey-value erosion with different structuring elements.

The rectangular, elliptic and diamond structuring elements are "flat", i.e. these structuring elements have a constant value. For these structuring elements, param determines the sizes of the structuring elements.

When shape is DIP_FLT_SHAPE_DISCRETE_LINE or DIP_FLT_SHAPE_INTERPOLATED_LINE, the structuring element is a line. param->array[0] determines the length, param->array[1] the angle. This is currently only supported for 2D images. Interpolated lines use interpolation to obtain a more accurate result, but loose the strict increasingness and extensivity (these properties are satisfied only by approximation).

When shape is set to DIP_FLT_SHAPE_PARABOLIC, params specifies the curvature of the parabola.

When shape is set to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, se is used as structuring element. It can be either a binary or a grey-value image. Its origin is the center, or one pixel to the left of the center if the size is even.

If shape is not equal to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, se can be set to zero. When shape is set to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, param is ignored, and can be set to zero.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	in	Input
dip_Image	out	Output
dip_Image	se	Custom structuring element
dip_BoundaryArray	boundary	Boundary conditions
dip_FloatArray	param	Filter parameters
dip_FilterShape	shape	Structuring element

The enumerator dip_FilterShape contains the following constants:

Name	Description
DIP_FLT_SHAPE_DEFAULT	Default filter window, same as
	DIP_FLT_SHAPE_RECTANGULAR
DIP_FLT_SHAPE_RECTANGULAR	Rectangular filter window, can be even in size
DIP_FLT_SHAPE_ELLIPTIC	Elliptic filter window, always odd in size
DIP_FLT_SHAPE_DIAMOND	Diamond-shaped filter window, always odd in size
DIP_FLT_SHAPE_PARABOLIC	Parabolic filter window (morphology only)
DIP_FLT_SHAPE_DISCRETE_LINE	Rotated line structuring element (morphology only)
DIP_FLT_SHAPE_INTERPOLATED_LINE	Rotated line structuring element, through interpolation
	(morphology only)
DIP_FLT_SHAPE_PERIODIC_LINE	(not implemented)
DIP_FLT_SHAPE_STRUCTURING_ELEMENT	Use se as filter window, can be any size

SEE ALSO

Closing, Opening, Dilation

error.h

Contains error messages

SYNOPSIS

#include "dip_error.h"

FUNCTION

Contains a lot of definitions to do with DIPlib's error mechanism. In particular, this include file contains definitions for a number of error messages. These are all of the type extern const char *. A list of the error sorted by category follows below:

Memory allocation

Name	Description
dip_errorCouldNotAllocateMemory	No memory could be allocated

Image creation errors

Name	Description
dip_errorImageIsLocked	Image is locked
dip_errorImageNotRaw	Image is not in the RAW state
dip_errorImageNotValid	Image is not in the VALID state
dip_errorImagesNotUnique	Image is used as an output image more than once
${ t dip_errorImageLockInvalidKey}$	Cannot unlock. Wrong key

Image type errors

Name	Description
dip_errorIllegalImageType	Illegal image type
dip_errorImageTypeDoesNotExist	Image type does not exist
dip_errorImageTypeAlreadyExists	Adding image type failed. Type already exists
dip_errorImageTypeNotSupported	Image type not supported
dip_errorImageTypeHandlerMissing	No type handler for image type

Image data type errors

Name	Description
dip_errorDataTypeNotSupported	Data type not supported
dip_errorIllegalDataType	Illegal data type

Image dimension(ality) errors

Name	Description
$ ext{dip_errorIllegalDimensionality}$	Illegal dimensionality
dip_errorDimensionalityNotSupported	Dimensionality not supported
dip_errorIllegalDimension	Illegal dimension

ErrorFree

Free a DIPlib call tree

SYNOPSIS

void dip_ErrorFree(error)

FUNCTION

Free a DIPlib call tree.

ARGUMENTS

Data type	Name	Description
dip_Error	error	DIPlib call tree

RETURNS

Nothing

EuclideanDistanceToPoint

Distance generation function

SYNOPSIS

```
#include "dip_generation.h"
dip_Error dip_EuclideanDistanceToPoint ( output, origin )
```

DATA TYPES

Output: sfloat

FUNCTION

Computes the Euclidean distance of each pixel in the output image to a point at origin. The coordinates of origin may lie outside the image.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{\mathtt{-}}\mathtt{Image}$	output	Output Image
dip_FloatArray	origin	Coordinates of the Origin

SEE ALSO

 ${\tt EllipticDistanceToPoint}, {\tt CityBlockDistanceToPoint}$

EuclideanDistanceTransform

Euclidean distance transform

SYNOPSIS

#include "dip_distance.h"
dip_Error dip_EuclideanDistanceTransform (in, out, distance, border, method)

DATA TYPES

binary

FUNCTION

This function computes the Euclidean distance transform of an input binary image using the vector-based method as opposed to the chamfer method. This method computes distances from the objects (binary 1's) to the nearest background (binary 0's) of in and stored the result in out. The out image is a sfloat type image.

The distance parameter can be used to specify anisotropic sampling densities. If it is set to zero, the sampling density is assumed to be 1.0 along all axes.

The border parameter specifies whether the edge of the image should be treated as objects (border = DIP_TRUE) or as background (border = DIP_FALSE).

Individual vector components of the Euclidean distance transform can be obtained with the VectorDistanceTransform.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input image
dip_Image	out	Output image
dip_FloatArray	distance	Sampling distances
dip_Boolean	border	Image border type
dipf_DistanceTransform	method	Transform method

dipf_DistanceTransform defines the following distance transform types:

Name	Description
DIP_EDT_FAST	fastest, but most errors
DIP_EDT_TIES	slower, but fewer errors
DIP_EDT_TRUE	slow, uses lots of memory, but is "error free"
DIP_EDT_BRUTE_FORCE	gives a result from which errors are calculated for the other methods. This
	method is extremly slow and should only be used for testing purposes.

LITERATURE

Danielsson, P.E. (1980). "Euclidean distance mapping." Computer Graphics and Image Processing 14: 227-248.

Mullikin, J.C. (1992). "The vector distance transform in two and three dimensions." CVGIP: Graphical Models and Image Processing 54(6): 526-535.

Ragnemalm, I. (1990). Generation of Euclidean Distance Maps, Thesis No. 206. Licentiate thesis. Linkoing University, Sweden.

Ye, Q.Z. (1988). "The signed Euclidean distance transform and its applications." in Proceedings, 9th International Conference on Pattern Recognition, Rome, 495-499.

KNOWN BUGS

The EDT_TRUE transform type is prone to produce an internal buffer overflow when applied to larger (almost) spherical objects. It this cases use EDT_TIES or EDT_BRUTE_FORCE instead.

The option border = DIP_FALSE is not supported for EDT_BRUTE_FORCE.

This function supports 2 and 3-dimensional images.

AUTHOR

James C. Mullikin, adapted to DIPlib by Geert M.P. van Kempen

SEE ALSO

VectorDistanceTransform, GreyWeightedDistanceTransform

EuclideanSkeleton

binary skeleton operation

SYNOPSIS

#include "dip_binary.h"
dip_Error dip_EuclideanSkeleton (in, out, endpixelCondition, edgeCondition)

DATA TYPES

binary

FUNCTION

This function calculates an accurate (euclidean)skeleton. It tests Hilditch conditions to preserve topology. The algorithms uses the following distance metrics:

2D

5	4-connected neighbor
7	8-connected neighbor
11	neighbors reachable with a knight's move
$\overline{3D}$	
4	6-connected neighbors

$\mid 4$	6-connected neighbors
6	18-connected neighbors
7	26-connected neighbors
9	neighbors reachable with knight's move
10	(2,1,1) neighbors
12	(2,2,1) neighbors

The edge parameter specifies whether the border of the image should be treated as object (DIP_TRUE) or as background (DIP_FALSE). See section 9.6, "Morphology-based operations", in Fundamentals of Image Processing for a description of the skeleton operation.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Binary input image
dip_Image	out	Output image
dip_EndpixelCondition	endpixelCondition	Endpixel condition
dip_Boolean	edgeCondition	Edge condition

The dip_EndpixelCondition enumeration consists of the following flags:

Name	Description
DIP_ENDPIXEL_CONDITION_LOOSE_ENDS_AWAY	Loose ends are eaten away
DIP_ENDPIXEL_CONDITION_NATURAL	"natural" endpixel condition of this
	algorithm
DIP_ENDPIXEL_CONDITION_KEEP_WITH_ONE_NEIGHBOR	Keep endpoint if it has a neighbor
DIP_ENDPIXEL_CONDITION_KEEP_WITH_TWO_NEIGHBORS	Keep endpoint if it has two neighbors
DIP_ENDPIXEL_CONDITION_KEEP_WITH_THREE_NEIGHBORS	Keep endpoint if it has three neighbors

KNOWN BUGS

EuclideanSkeleton is only implemented for 2 and 3 D images.

EuclideanSkeleton does not process pixels in a 2-pixel border around the edge. If this is an issue, consider adding 2 pixels on each side of your image.

The function is buggy for 3D images. DIP_ENDPIXEL_CONDITION_LOOSE_ENDS_AWAY and DIP_ENDPIXEL_CONDITION_KEEP_WITH_ONE_NEIGHBOR produce the same result as DIP_ENDPIXEL_CONDITION_KEEP_WITH_THREE_NEIGHBORS. Both DIP_ENDPIXEL_CONDITION_NATURAL and DIP_ENDPIXEL_CONDITION_KEEP_WITH_TWO_NEIGHBORS produce resonable results under most circumstances, but don't count on it!

LITERATURE

"Improved metrics in image processing applied to the Hilditch skeleton", B.J.H. Verwer, 9th ICPR, Rome, November 14-17, 1988.

AUTHOR

Ben Verwer, adapted to DIPlib by Geert van Kempen.

SEE ALSO

BinaryPropagation

Exit

Clean up before exiting

SYNOPSIS

```
dip_Error dip_Exit( void )
dip_Error dipio_Exit( void )
```

FUNCTION

Free all memory used internally by DIPlib. Call this function when you stop using DIPlib (before exiting your program).

SEE ALSO

Initialise

Exp

arithmetic function

SYNOPSIS

dip_Error dip_Exp (in, out)

DATA TYPES

binary, integer, **float**

FUNCTION

Computes the natural exponent of the input image values.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	in	Input
$dip_{-}Image$	out	Output

SEE ALSO

Sqrt, Exp2, Exp10, Ln, Log2, Log10

Exp10

arithmetic function

SYNOPSIS

dip_Error dip_Exp10 (in, out)

DATA TYPES

binary, integer, **float**

FUNCTION

Computes the base ten exponent of the input image values.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	in	Input
$dip_{-}Image$	out	Output

SEE ALSO

Sqrt, Exp, Exp2, Ln, Log2, Log10

Exp2

arithmetic function

SYNOPSIS

dip_Error dip_Exp2 (in, out)

DATA TYPES

binary, integer, **float**

FUNCTION

Computes the base two exponent of the input image values.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	in	Input
dip_Image	out	Output

SEE ALSO

Sqrt, Exp, Exp10, Ln, Log2, Log10

ExponentialFitCorrection

Exponential fit based attenuation correction

SYNOPSIS

#include "dip_microscopy.h"
dip_Error dip_ExponentialFitCorrection (in, out, method, percentile, fromWhere,
hysteresis, varWeighted)

DATA TYPES

binary, integer, float

FUNCTION

This routine implements a simple absorption, reflection and bleaching correction based upon the assumption that the sum of these effects result in a exponential extinction of the signal as a function of depth. Only pixels that are non-zero are taken into account. Depending upon the chosen method, the mean or a percentile of all the non-zero pixels are calculated as a function of the slice number (depth). Then an exponential function is fitted through these slice-representing values. The starting point of the fit is determined by fromWhere. The first maximum is found with point[z+1] > hysteresis * point[z]. If the mean variant is chosen one can chose to apply a variance weighting to the fit.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input image
dip_Image	out	Output image
dipf_ExpFitData	method	Data statistic to fit on
dip_float	percentile	Percentile
dipf_ExpFitStart	fromWhere	From where to start the fit
dip_float	hysteresis	First maximum hysteresis
dip_Boolean	varWeighted	Fit with variance weights

The dipf_ExpFitData enumaration consists of the following flags:

Name	Description
DIP_ATTENUATION_EXP_FIT_DATA_MEAN	Fit on the mean values
DIP_ATTENUATION_EXP_FIT_DATA_PERCENTILE	Fit on the specified percentile of the data

The dipf_ExpFitStart enumaration consists of the following flags:

Name	Description
DIP_ATTENUATION_EXP_FIT_START_FIRST_PIXEL	Start fit on first pixel
DIP_ATTENUATION_EXP_FIT_START_GLOBAL_MAXIMUM	Start fit on global maximum
DIP_ATTENUATION_EXP_FIT_START_FIRST_MAXIMUM	Start fit on first maximum

LITERATURE

K.C. Strasters, H.T.M. van der Voort, J.M. Geusebroek, and A.W.M. Smeulders, "Fast attenuation correction in fluorescence confocal imaging: a recursive approach", BioImaging, vol. 2, no. 2, 1994, 78-92.

AUTHOR

Karel Strasters, adapted to DIPlib by Geert van Kempen.

SEE ALSO

 ${\tt Attenuation Correction}, {\tt Simulated Attenuation}$

ExtendRegion

Image manipulation functions

SYNOPSIS

```
#include "dip_manipulation.h"
dip_Error dip_ExtendRegion ( image, origin, regDims, bc, ordering, imValues )
```

FUNCTION

This functions extends a region in an image, defined by origin and regDims, with a specified boundary condition bc. The pixels outside the region are modified according to bc. ordering changes the order in which the dimensions are processed, set to 0 to use default process order.

ARGUMENTS

Data type	Name	Description	
$\mathtt{dip}_{-}\mathtt{Image}$	image	Image, will be modified	
dip_IntegerArray	origin	Origin of region	
dip_IntegerArray	regDims	Size of region	
dip_BoundaryArray	bc	Boundary conditions	
dip_IntegerArray	ordering	Ordering of dimensions	
dip_Image *	imValues	Unused, set to 0.	

NOTE

Boundary conditions DIP_BC_ZERO_ORDER_EXTRAPOLATE, DIP_BC_FIRST_ORDER_EXTRAPOLATE and DIP_BC_SECOND_ORDER_EXTRAPOLATE are not supported.

FeatureAnisotropy2D

Measure the anisotropy in a labeled region

SYNOPSIS

```
#include "dip_measurement.h"
dip_int dip_FeatureAnisotropy2DID ( void )
```

OUTPUT DATA TYPE

dip_float

FUNCTION

dip_FeatureAnisotropy2DID returns the ID value of this measurement function, that is registered by Initialise.

The grey value input image should contain an orientation field. For each labeled region, a tensor is constructed at each of the region's pixels. This tensor is as follow:

```
cos^2(phi) cos(phi)sin(phi)
[
cos(phi)sin(phi) sin^2(phi)
```

The next step is to compute a new tensor, each element computed by averaging the corresponding elements of all the individual tensors. This average tensor represents the orientation information of the region as a whole. Eigenvalue analysis of this tensor yields two eigenvalues, the largest 10, the smallest 11. The anisotropy measure is:

```
(10 - 11) / (10 + 11)
```

which is zero for a fully isotropic regions (i.e. one where there is no preferred orientation), and one for a fully anisotropic region (i.e. when there is a single orientation).

NOTE

This function ignores any physical dimensions passed through the Measure function.

SEE ALSO

Measure

 $\label{lem:convex} Feature Anisotropy 2D, Feature Bending Energy, Feature Center, \\ Feature Chain Code Bending Energy, Feature Convex Area, Feature Convex Perimeter, \\ \\ \end{tabular}$

FeatureConvexity, FeatureDimension, FeatureExcessKurtosis, FeatureFeret,
FeatureGinertia, FeatureGmu, FeatureGravity, FeatureInertia, FeatureLongestChaincodeRun,
FeatureMass, FeatureMaxVal, FeatureMaximum, FeatureMean, FeatureMinVal, FeatureMinimum,
FeatureMu, FeatureOrientation2D, FeatureP2A, FeaturePerimeter, FeatureRadius,
FeatureShape, FeatureSize, FeatureSkewness, FeatureStdDev, FeatureSum,
FeatureSurfaceArea

FeatureBendingEnergy

Undocumented measurement function

FUNCTION

This measurement function is undocumented and not meant for public use. Measure

FeatureCenter

Measure the object's center

SYNOPSIS

```
#include "dip_measurement.h"
dip_int dip_FeatureCenterID ( void )
```

OUTPUT DATA TYPE

dip_FloatArray

FUNCTION

dip_FeatureCenterID returns the ID value of this measurement function, that is registered by Initialise.

This functions measures the centre of an object by calculating the first moments of the object using the object labels as binary mask. The intensity information is not taken into account.

SEE ALSO

Measure

FeatureChainCodeBendingEnergy

Undocumented measurement function

FUNCTION

This measurement function is undocumented and not meant for public use.

NOTE: this function uses chain codes. It expects each measured object to be compact, that is, to have only one chain code. Additional chain codes are ignored, meaning that non-compact objects are not measured properly. Take care in providing the correct connectivity value: if you object is compact only with 2-connectivity, this measure will fail if you call Measure with a value of 1 for the connectivity.

Measure

FeatureChainCodeFunction

Measurement feature #measure function

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error (*dip_FeatureChainCodeFunction) ( measurement, featureID, objectID, chaincode, iterations )
```

FUNCTION

The chaincode measure function is meant for 2-D measurement functions that only require information on the shape of the object's contour, such as FeaturePerimeter. The chaincode function is called for each object separately, with the contour of that object stored in chaincode.

ARGUMENTS

Data type	Name	Description	
dip_Measurement	measurement	Measurement data structure	
dip_int	featureID	Measurement function ID	
dip_int	objectID	ID of the object to be measured	
dip_ChainCode	chaincode	Chaincode data structure encoding the object's contour	
dip_int	iterations	Number of iterations the measure function needs to scan the	
		data	

SEE ALSO

 $\label{lem:measurement} {\tt MeasurementFeatureRegister}, {\tt FeatureLineFunction}, {\tt FeatureConvHullFunction}, {\tt FeatureCompositeFunction}, {\tt FeatureCreateFunction}$

${\tt Feature Compose Function}$

Measurement feature #compose function

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error (*dip_FeatureComposeFunction) ( measurement, featureID, label, intensity,
compositeFeatureID, resources )
```

FUNCTION

The compose function is called to obtain a list of measurement features. These features are measured before the measure function of a composite feature is called (FeatureCompositeFunction). This parameter is ignored for other measurement types. The compose function is called after the create function (FeatureCreateFunction).

ARGUMENTS

Data type	Name	Description
dip_Measurement	measurement	Measurement data structure
dip_int	featureID	Measurement feature function ID
$\mathtt{dip}_{-}\mathtt{Image}$	label	Image with pixel intensities represending object
		IDs
$\mathtt{dip}_{-}\mathtt{Image}$	intensity	Image containing corresponding intensity values
dip_IntegerArray *	compositeFeatureID	Pointer to an integer array containing the the IDs
		of the measurement features this function requires
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

 ${\tt MeasurementFeatureRegister}, {\tt FeatureCompositeFunction}$

FeatureCompositeFunction

Measurement feature #measure function

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error (*dip_FeatureCompositeFunction) ( measurement, featureID, objectID,
composite, iterations )
```

FUNCTION

The composite measure function is meant for features that derive their measurements from the results of other measurement functions. The measurement IDs this function is based on is obtained by calling the compose function (FeatureComposeFunction). The composite measure function obtains the results of these measurements through its composite function parameter. Use the regular measurement structure access method to read the values in this parameter (i.e. MeasurementObjectValue).

ARGUMENTS

Data type	Name	Description	
dip_Measurement	measurement	Measurement data structure	
dip_int	featureID	Measurement function ID	
dip_int	objectID	ID of the object to be measured	
dip_Measurement	composite	Measurement structure containing the measurement data this	
		function is based on	
dip_int	iterations	Number of iterations the measure function needs to scan the	
		data	

SEE ALSO

 $\label{lem:measurement} {\tt MeasurementFeatureRegister}, \ {\tt FeatureLineFunction}, \ {\tt FeatureChainCodeFunction}, \ {\tt FeatureConvHullFunction}, \ {\tt FeatureConvection}, \ {\tt FeatureComposeFunction}$

FeatureConvertFunction

Measurement feature #convert function

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error (*dip_FeatureConvertFunction) ( in, featureID, inID, out, outID, resources
)
```

FUNCTION

The convert function should convert the measurement data of the feature featureID for the object inID in the measurement in to object outID of measurement out. This function is called by MeasurementFeatureConvert.

ARGUMENTS

Data type	Name	Description	
dip_Measurement	in	Input measurement data structure	
dip_int	featureID	Measurement function ID	
dip_int	inID	ID of the object in in	
dip_Measurement	out	Output measurement data structure	
dip_int	outID	ID of the object in out	
dip_Resources	resources	Resources tracking structure. See ResourcesNew	

SEE ALSO

 ${\tt MeasurementFeatureRegister}, {\tt MeasurementFeatureConvert}$

FeatureConvexArea

Measure the area of the object's convex hull

SYNOPSIS

#include "dip_measurement.h"
dip_int dip_FeatureConvexAreaID (void)

OUTPUT DATA TYPE

dip_float

FUNCTION

dip_FeatureConvexAreaID returns the ID value of this measurement function, that is registered by Initialise.

This function measures the area of the convex hull of the object. The convex hull is a polygon derived from the border pixels, and thus its area is not necessarily an integer. This function supports 2D images only.

NOTE

This function uses chain codes. It expects each measured object to be compact, that is, to have only one chain code. Additional chain codes are ignored, meaning that non-compact objects are not measured properly. Take care in providing the correct connectivity value: if you object is compact only with 2-connectivity, this measure will fail if you call Measure with a value of 1 for the connectivity.

SEE ALSO

 ${\tt Measure, Image Chain Code, Chain Code Convex Hull, Convex Hull Get Area}$

FeatureConvexity

Measure the object's convexity

SYNOPSIS

#include "dip_measurement.h"
dip_int dip_FeatureConvexityID (void)

OUTPUT DATA TYPE

dip_float

FUNCTION

dip_FeatureConvexityID returns the ID value of this measurement function, that is registered by Initialise

This function is a composite measurement function, that returns the ratio between FeatureSize and FeatureConvexArea. A convex object will have a convexity of 1 (or slightly smaller due to discretization issues). Convexity values smaller than 1 indicate that the object boundary has concavities or that the object has holes. This function supports 2D images only.

NOTE

This function uses chain codes. It expects each measured object to be compact, that is, to have only one chain code. Additional chain codes are ignored, meaning that non-compact objects are not measured properly. Take care in providing the correct connectivity value: if you object is compact only with 2-connectivity, this measure will fail if you call Measure with a value of 1 for the connectivity.

SEE ALSO

Measure, ImageChainCode, ChainCodeConvexHull

FeatureConvexPerimeter

Measure the perimeter of the object's convex hull

SYNOPSIS

#include "dip_measurement.h"
dip_int dip_FeatureConvexPerimeterID (void)

OUTPUT DATA TYPE

dip_float

FUNCTION

dip_FeatureConvexPerimeterID returns the ID value of this measurement function, that is registered by Initialise.

This function measures the perimeter of the convex hull of the object. This function supports 2D images only.

NOTE

If any physical dimensions are passed to this function through Measure, only the sample distance along the first dimension is used. All other dimensions are assumed to be sampled the same way. This produces incorrect results for anisotropically sampled images.

NOTE

This function uses chain codes. It expects each measured object to be compact, that is, to have only one chain code. Additional chain codes are ignored, meaning that non-compact objects are not measured properly. Take care in providing the correct connectivity value: if you object is compact only with 2-connectivity, this measure will fail if you call Measure with a value of 1 for the connectivity.

SEE ALSO

 ${\tt Measure,\ Image Chain Code,\ Chain Code Convex Hull,\ Convex Hull Get Perimeter}$

FeatureAnisotropy2D, FeatureBendingEnergy, FeatureCenter,
FeatureChainCodeBendingEnergy, FeatureConvexArea, FeatureConvexPerimeter,
FeatureConvexity, FeatureDimension, FeatureExcessKurtosis, FeatureFeret,
FeatureGinertia, FeatureGmu, FeatureGravity, FeatureInertia, FeatureLongestChaincodeRun,
FeatureMass, FeatureMaxVal, FeatureMaximum, FeatureMean, FeatureMinimum,

 $\label{lem:peatureMu} FeatureOrientation2D, FeatureP2A, FeaturePerimeter, FeatureRadius, FeatureShape, FeatureSize, FeatureSkewness, FeatureStdDev, FeatureSum, FeatureSurfaceArea$

FeatureConvHullFunction

Measurement feature #measure function

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error (*dip_FeatureConvHullFunction) ( measurement, featureID, objectID,
convhull, iterations )
```

FUNCTION

The convhull measure function is meant for 2-D measurement functions that only require information on the convex hull of the object, such as FeatureFeret. The convhull function is called for each object seperately, with the convex hull of that object stored in convhull.

ARGUMENTS

Data type	Name	Description	
dip_Measurement	measurement	Measurement data structure	
dip_int	featureID	Measurement function ID	
dip_int	objectID	ID of the object to be measured	
dip_Polygon	convhull	Polygon data structure representing the object's convex hull	
dip_int	iterations	Number of iterations the measure function needs to scan the	
		data	

SEE ALSO

 $\label{lem:measurement} {\tt MeasurementFeatureRegister}, \ {\tt FeatureLineFunction}, \ {\tt FeatureChainCodeFunction}, \ {\tt FeatureCompositeFunction}, \ {\tt FeatureCreateFunction}$

FeatureCreateFunction

Measurement feature #create function

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error (*dip_FeatureCreateFunction) ( measurement, featureID, label, intensity,
physDims, params, data, resources )
```

FUNCTION

The create function is called to initialise the measurement function. It should allocate and initialise a memory block for internal use, assign this block to the pointer *data, and register it in resources.

ARGUMENTS

Data type	Name	Description
dip_Measurement	measurement	Measurement data structure
dip_int	featureID	Measurement feature function ID
dip_Image	label	Image with pixel intensities represending object IDs
dip_Image	intensity	Image containing corresponding intensity values
dip_PhysicalDimensions	physDims	Physical dimensions data structure
void *	params	For future expansion, is currently always NULL
void **	data	Pointer to a data block that can later be accessed
		using MeasurementObjectData
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

MeasurementFeatureRegister, FeatureLineFunction, FeatureImageFunction, FeatureChainCodeFunction, FeatureConvHullFunction, FeatureCompositeFunction

FeatureDescriptionFree

Free a Feature Description

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error dip_FeatureDescriptionFree ( description )
```

FUNCTION

This function frees a Feature Description data structure. This is not the preferred way of freeing a Feature Description. Use the resources mechanism instead (Resources tracking structure. See ResourcesNew).

ARGUMENTS

Data type	Name	Description
dip_FeatureDescription *	description	Feature Description to be freed

SEE ALSO

FeatureDescriptionFunction

Measurement feature #description function

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error (*dip_FeatureDescriptionFunction) ( measurement, featureID, physDims,
decription, resources )
```

FUNCTION

The description function should return a dip_FeatureDescription structure containing information on the measurement function, such as its name, a short description, labels for each value measured, and units of its measurement. This function is called by MeasurementFeatureDescription.

The description structure should be allocated by this function using FeatureDescriptionNew, and registered in resources. The functions FeatureDescriptionSetName, FeatureDescriptionSetDescription, FeatureDescriptionSetDimensionLabels and

FeatureDescriptionSetUnits should be used to populate the structure.

ARGUMENTS

Data type	Name	Description
dip_Measurement	measurement	Measurement data structure
dip_int	featureID	Measurement function ID
dip_PhysicalDimensions	physDims	Physical dimensions data structure
dip_FeatureDescription *	description	Pointer to a structure containing descriptive
		information of the measurement feature function
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

 $\label{lem:measurement} \begin{tabular}{ll} Measurement Feature Description, Feature Description Set Name, Feature Description Set Description, Feature Description Set Dimension Labels, Feature Description Set Units \\ \begin{tabular}{ll} Feature Description Set Units \\ \begin$

FeatureDescriptionGetDescription

Get the description of the described feature

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error dip_FeatureDescriptionGetDescription ( description, text, resources )
```

FUNCTION

Gets the description of the feature described by description.

ARGUMENTS

Data type	Name	Description
dip_FeatureDescription	description	Feature description data structure
dip_String *	text	Description text
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

FeatureDescriptionGetLabels

Get the labels of the described feature

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error dip_FeatureDescriptionGetLabels ( description, labels, resources )
```

FUNCTION

Gets the labels of the data of the feature descripted by description.

ARGUMENTS

Data type	Name	Description
dip_FeatureDescription	description	Feature description data structure
dip_StringArray *	labels	Feature Labels
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

FeatureDescriptionGetName

Get the name of the described feature

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error dip_FeatureDescriptionGetName ( description, name, resources )
```

FUNCTION

Gets the name of the feature described by description.

ARGUMENTS

Data type	Name	Description
dip_FeatureDescription	description	Feature description data structure
dip_String *	name	Name of the measurement feature
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

FeatureDescriptionGetUnits

Get the Units of the described feature

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error dip_FeatureDescriptionGetUnits ( description, units, resources )
```

FUNCTION

Gets the units of the data of the feature descripted by description.

ARGUMENTS

Data type	Name	Description
dip_FeatureDescription	description	Feature description data structure
dip_StringArray *	units	Array of Unit texts
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

FeatureDescriptionNew

Allocate a new FeatureDescription

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error dip_FeatureDescriptionNew ( description, resources )
```

FUNCTION

This function allocates a new dip_FeatureDescription data structure. A feature description contains the name, a short description of a measurement feature, as well as the labels and units of the data measured by the feature.

ARGUMENTS

Data type	Name	Description
dip_FeatureDescription	description	Feature description data structure
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

${\tt Feature Description Set Description}$

Set the description of the described feature

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error dip_FeatureDescriptionSetDescription ( description, text )
```

FUNCTION

Sets the description of the feature descripted by description.

ARGUMENTS

Data type	Name	Description
dip_FeatureDescription	description	Feature description data structure
char *	text	Description text

SEE ALSO

FeatureDescriptionSetDimensionLabels

Label set convenience function

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error dip_FeatureDescriptionSetDimensionLabels ( description, measurement,
featureID, baseLabel )
```

FUNCTION

This function set the labels of the feature, described by description, by adding for each label a dimension indicator to baseLabel. For dimensions 0 to 3, X, Y or Z is added. For dimensions higher, the numerical value of the dimension is added.

ARGUMENTS

Data type	Name	Description
dip_FeatureDescription	description	Feature description data structure
dip_Measurement	measurement	Measurement data structure
dip_int	featureID	ID of the measurement feature
char *	baseLabel	Base label

SEE ALSO

FeatureDescriptionSetLabel

Set the name of a particular feature label

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error dip_FeatureDescriptionSetLabel ( description, number, label )
```

FUNCTION

This function sets the name of a particular label of the described feature.

ARGUMENTS

Data type	Name	Description
dip_FeatureDescription	description	Feature description data structure
dip_int	number	Index of the label
char *	label	Label text

SEE ALSO

FeatureDescriptionSetLabels

Set the labels of the described feature

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error dip_FeatureDescriptionSetLabels ( description, measurement, featureID, labels, label )
```

FUNCTION

Sets the labels of the data of the feature descripted by description.

ARGUMENTS

Data type	Name	Description
dip_FeatureDescription	description	Feature description data structure
dip_Measurement	measurement	Measurement data structure
dip_int	featureID	ID of the measurement feature
dip_StringArray	labels	Array of label describing strings, one for each label
char *	label	Single description of all feature labels

SEE ALSO

FeatureDescriptionSetName

Set the name of the described feature

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error dip_FeatureDescriptionSetName ( description, name )
```

FUNCTION

Sets the name of the feature descripted by description.

ARGUMENTS

Data type	Name	Description
dip_FeatureDescription	description	Feature description data structure
char *	name	Name of the measurement feature

SEE ALSO

FeatureDescriptionSetUnit

Set the units of a particular feature dimension

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error dip_FeatureDescriptionSetUnit ( description, number, unit )
```

FUNCTION

This function sets the name of the unit along a particular dimension of the described feature.

ARGUMENTS

Data type	Name	Description
dip_FeatureDescription	description	Feature description data structure
dip_int	number	Index of the dimension
char *	unit	Units text

SEE ALSO

FeatureDescriptionSetUnits

Set the units of a described feature

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error dip_FeatureDescriptionSetUnits ( description, measurement, featureID,
units, unit )
```

FUNCTION

Sets the units of the data of the feature descripted by description.

ARGUMENTS

Data type	Name	Description
dip_FeatureDescription	description	Feature description data structure
dip_Measurement	measurement	Measurement data structure
dip_int	featureID	ID of the measurement feature
dip_StringArray	units	Array of Unit texts, one for each unit
char *	unit	Single text for all units

SEE ALSO

FeatureDimension

Measure the object's dimensions

SYNOPSIS

```
#include "dip_measurement.h"
dip_int dip_FeatureDimensionID ( void )
```

OUTPUT DATA TYPE

dip_FloatArray

FUNCTION

dip_FeatureDimensionID returns the ID value of this measurement function, that is registered by Initialise

This function measures the length of an object along the principal axes of the label image (e.g. the length object along the X, Y & Z axes).

SEE ALSO

Measure

FeatureExcessKurtosis

Undocumented measurement function

FUNCTION

This measurement function is undocumented and not meant for public use. Measure

FeatureFeret

Measure the object's Feret diameters

SYNOPSIS

#include "dip_measurement.h"
dip_int dip_FeatureFeretID (void)

OUTPUT DATA TYPE

dip_FloatArray

FUNCTION

dip_FeatureFeretID returns the ID value of this measurement function, that is registered by Initialise.

This function measures the Feret maximum and minimum diameters of an object. The Feret diameters are found by a "rotating caliper" algorithm on the convex hull polygon, see ConvexHullGetFeret. This function supports 2D images only.

The values returned are:

FeretMax	The widest projection of the object
FeretMin	The narrowest projection of the object
FeretPerpMin	The width of the projection perpendicular to
	"FeretMin"
FeretMaxAng	The angle of the projection for "FeretMax"
FeretMinAng	The angle of the projection for "FeretMin"

NOTE

If any physical dimensions are passed to this function through Measure, only the sample distance along the first dimension is used. All other dimensions are assumed to be sampled the same way. This produces incorrect results for anisotropically sampled images.

NOTE

This function uses chain codes. It expects each measured object to be compact, that is, to have only one chain code. Additional chain codes are ignored, meaning that non-compact objects are not measured properly. Take care in providing the correct connectivity value: if you object is compact only with 2-connectivity, this measure will fail if you call Measure with a value of 1 for the connectivity.

NOTE

In DIPlib version 2.3 and earlier, this measure was computed from the chain code directly, using the function ChainCodeGetFeret.

SEE ALSO

Measure, ImageChainCode, ChainCodeConvexHull, ConvexHullGetFeret

FeatureGinertia

Measure the object's inertia

SYNOPSIS

#include "dip_measurement.h"
dip_int dip_FeatureGinertiaID (void)

OUTPUT DATA TYPE

dip_FloatArray

FUNCTION

dip_FeatureGinertiaID returns the ID value of this measurement function, that is registered by Initialise.

This function calculates the inertia (weighted by its grey values) of an object by calculating the eigenvalues of the object's second order moments tensor. This measure only supports 2D and 3D objects. FeatureGinteria supports a measurement parameter (see Measure). If a pointer to a non-zero Boolean is supplied, this function will not only measure the eigenvalues of the second order moments tensor, but also the angles of its eigenvectors.

ARGUMENTS

Data type	Name	Description
dip_Boolean *	angles	Pointer to a Boolean specifying that "eigenangles" should be
		measured (not yet implemented)

NOTE

If any physical dimensions are passed to this function through Measure, only the sample distance along the first dimension is used. All other dimensions are assumed to be sampled the same way. This produces incorrect results for anisotropically sampled images.

LITERATURE

"Practical Handbook on Image Processing for Scientific Applications, chapter 16", Bernd Jahne, CRC Press, 1999.

SEE ALSO

Measure

FeatureGmu

Measure the object's inertia

SYNOPSIS

```
#include "dip_measurement.h"
dip_int dip_FeatureGmuID ( void )
```

OUTPUT DATA TYPE

dip_FloatArray

FUNCTION

dip_FeatureGmuID returns the ID value of this measurement function, that is registered by Initialise.

This function calculates the inertia (weighted by its grey values) of an object by calculating the object's second order moments tensor. This measure only supports 2D and 3D objects. The output tensor is ordered as follows:

```
2D: xx, xy, yy3D: xx, xy, xz, yy, yz, zz
```

NOTE

If any physical dimensions are passed to this function through Measure, only the sample distance along the first dimension is used. All other dimensions are assumed to be sampled the same way. This produces incorrect results for anisotropically sampled images.

LITERATURE

"Mechanics", Florian Scheck, Springer, 1999.

SEE ALSO

Measure

 ${\tt FeatureShape,\,FeatureSize,\,FeatureSkewness,\,FeatureStdDev,\,FeatureSum,}$ FeatureSurfaceArea

FeatureGravity

Measure the object's gravity

SYNOPSIS

```
#include "dip_measurement.h"
dip_int dip_FeatureGravityID ( void )
```

OUTPUT DATA TYPE

dip_FloatArray

FUNCTION

dip_FeatureGravityID returns the ID value of this measurement function, that is registered by Initialise

This function measures the point of gravity of the object, by calculating the object's first moment weighted by the intensity of each object pixel.

SEE ALSO

Measure

FeatureImageFunction

Measurement feature #measure function

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error (*dip_FeatureImageFunction) ( measurement, featureID, label, intensity,
objectID, iterations )
```

FUNCTION

The image measurement function is meant for measurement operation that need neighborhood or global object shape information for its operation (e.g. the FeatureSurfaceArea function needs to evaluate the 6 connected neighborhood of each boundary voxel). The object ID image label can contain values that are not present in objectID. These labels should be ignored.

ARGUMENTS

Data type	Name	Description
dip_Measurement	measurement	Measurement data structure
dip_int	featureID	Measurement function ID
dip_Image	label	Image with pixel intensities represending object IDs
$\mathtt{dip}_{-}\mathtt{Image}$	intensity	Image containing corresponding intensity values
dip_IntegerArray	objectID	Array of objectIDs to be measured
dip_int	iterations	Number of iterations the measure function needs to scan the
		image

SEE ALSO

MeasurementFeatureRegister, FeatureLineFunction, FeatureChainCodeFunction, FeatureConvHullFunction, FeatureCompositeFunction, FeatureCreateFunction

FeatureInertia

Measure the object's inertia

SYNOPSIS

```
#include "dip_measurement.h"
dip_int dip_FeatureInertiaID ( void )
```

OUTPUT DATA TYPE

dip_FloatArray

FUNCTION

dip_FeatureInertiaID returns the ID value of this measurement function, that is registered by Initialise.

This function calculates the inertia of an object by calculating the eigenvalues of the object's second order moments tensor. This measure only supports 2D and 3D objects. FeatureInteria supports a measurement parameter (see Measure). If a pointer to a non-zero Boolean is supplied, this function will not only measure the eigenvalues of the second order moments tensor, but also the angles of its eigenvectors.

ARGUMENTS

Data type	Name	Description
dip_Boolean *	angles	Pointer to a Boolean specifying that "eigenangles" should be
		measured (not yet implemented)

NOTE

If any physical dimensions are passed to this function through Measure, only the sample distance along the first dimension is used. All other dimensions are assumed to be sampled the same way. This produces incorrect results for anisotropically sampled images.

LITERATURE

"Practical Handbook on Image Processing for Scientific Applications, chapter 16", Bernd Jahne, CRC Press, 1999.

SEE ALSO

Measure

FeatureLineFunction

Measurement feature #measure function

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error (*dip_FeatureLineFunction) ( measurement, featureID, label,intensity,
size, objectID, dim, iterations )
```

FUNCTION

The line measure function obtains two arrays (label and intensity) with label and intensity information of the objects to be measured. The line measurement function is called for every line in the image (the scan dimension is determined at run time to be optimal). Since label can contain more than one different label, line itself is responsible for storing the measurement results for the appropriate object (using, for example, MeasurementObjectData). The object ID array label can contain values that are not present in objectID. These labels should be ignored.

ARGUMENTS

Data type	Name	Description
dip_Measurement measurement		Measurement data structure
dip_int	featureID	Measurement function ID
dip_sint32 *	label	Pointer to a list (image line) of object IDs
dip_float *	intensity	Pointer to a list of corresponding intensity values
dip_int	size	Size of the label and intensity list
dip_IntegerArray	objectID	Array of objectIDs to be measured
$\mathtt{dip}_{-}\mathtt{int}$	dim	Dimension of the line, see ScanFrameWork
dip_int	iterations	Number of iterations the measure function needs to scan the
		line

SEE ALSO

 ${\tt MeasurementFeatureRegister, FeatureImageFunction, FeatureChainCodeFunction, FeatureConvHullFunction, FeatureCompositeFunction, FeatureCreateFunction}$

FeatureLongestChaincodeRun

Undocumented measurement function

FUNCTION

This measurement function is undocumented and not meant for public use.

NOTE: this function uses chain codes. It expects each measured object to be compact, that is, to have only one chain code. Additional chain codes are ignored, meaning that non-compact objects are not measured properly. Take care in providing the correct connectivity value: if you object is compact only with 2-connectivity, this measure will fail if you call Measure with a value of 1 for the connectivity.

Measure

FeatureMass

Measure the mass of the object (sum of grey-values)

SYNOPSIS

```
#include "dip_measurement.h"
dip_int dip_FeatureMassID ( void )
```

OUTPUT DATA TYPE

FUNCTION

dip_FeatureMassID returns the ID value of this measurement function, that is registered by Initialise. This function is just an alias for dip_FeatureSumID.

This function measures the sum of the grey-value in the intensity image (see Measure) of pixels inside the object, and is equivalent to FeatureSum

SEE ALSO

Measure

FeatureMaximum

Measure the object's maximum coordinate value

SYNOPSIS

```
#include "dip_measurement.h"
dip_int dip_FeatureMaximumID ( void )
```

OUTPUT DATA TYPE

dip_IntegerArray, dip_FloatArray

FUNCTION

dip_FeatureMaximumID returns the ID value of this measurement function, that is registered by Initialise.

This function measures the maximum coordinate value of each dimension of the object.

If a dip_PhysicalDimensions parameter is given to Measure, the maximum coordinate of the object is given in physical units, and is a dip_FloatArray rather than a dip_IntegerArray.

SEE ALSO

Measure

FeatureMaxVal

Measure the object's maximum intensity

SYNOPSIS

```
#include "dip_measurement.h"
dip_int dip_FeatureMaxValID ( void )
```

OUTPUT DATA TYPE

dip_float

FUNCTION

dip_FeatureMaxValID returns the ID value of this measurement function, that is registered by Initialise.

This function measures the maximum intensity in the intensity image (see Measure) of pixels inside the object.

SEE ALSO

Measure

FeatureMean

Measure the object's mean intensity

SYNOPSIS

```
#include "dip_measurement.h"
dip_int dip_FeatureMeanID ( void )
```

OUTPUT DATA TYPE

dip_float

FUNCTION

dip_FeatureMeanID returns the ID value of this measurement function, that is registered by Initialise

This function measures the mean intensity in the intensity image (see Measure) of pixels inside the object.

SEE ALSO

Measure

FeatureMinimum

Measure the object's minimum coordinate value

SYNOPSIS

```
#include "dip_measurement.h"
dip_int dip_FeatureMinimumID ( void )
```

OUTPUT DATA TYPE

dip_IntegerArray, dip_FloatArray

FUNCTION

dip_FeatureMinimumID returns the ID value of this measurement function, that is registered by Initialise.

This function measures the minimum coordinate value of each dimension of the object.

If a dip_PhysicalDimensions parameter is given to Measure, the minimum coordinate of the object is given in physical units, and is a dip_FloatArray rather than a dip_IntegerArray.

SEE ALSO

Measure

FeatureMinVal

Measure the object's minimum intensity

SYNOPSIS

```
#include "dip_measurement.h"
dip_int dip_FeatureMinValID ( void )
```

OUTPUT DATA TYPE

dip_float

FUNCTION

dip_FeatureMinValID returns the ID value of this measurement function, that is registered by Initialise

This function measures the minimum intensity in the intensity image (see Measure) of pixels inside the object.

SEE ALSO

Measure

FeatureMu

Measure the object's inertia

SYNOPSIS

```
#include "dip_measurement.h"
dip_int dip_FeatureMuID ( void )
```

OUTPUT DATA TYPE

dip_FloatArray

FUNCTION

dip_FeatureMuID returns the ID value of this measurement function, that is registered by Initialise

This function calculates the inertia of an object by calculating the object's second order moments tensor. This measure only supports 2D and 3D objects. The output tensor is ordered as follows:

```
2D: xx, xy, yy
3D: xx, xy, xz, yy, yz, zz
```

NOTE

If any physical dimensions are passed to this function through Measure, only the sample distance along the first dimension is used. All other dimensions are assumed to be sampled the same way. This produces incorrect results for anisotropically sampled images.

LITERATURE

"Mechanics", Florian Scheck, Springer, 1999.

SEE ALSO

Measure

FeatureSurfaceArea

FeatureOrientation2D

Undocumented measurement function

FUNCTION

This measurement function is undocumented and not meant for public use. Measure

FeatureP2A

Measure the circularity of the object

SYNOPSIS

```
#include "dip_measurement.h"
dip_int dip_FeatureP2AID ( void )
```

OUTPUT DATA TYPE

dip_float

FUNCTION

dip_FeatureP2AID returns the ID value of this measurement function, that is registered by Initialise.

This function is a composite measurement function, that uses FeatureSize, FeaturePerimeter, and FeatureSurfaceArea to determine the circularity of an object by calculating: 2D: P2A = perimeter^2 / (4Pi * size) 3D: P2A = surface-area^1.5 / (6 Sqrt(Pi) * size)

NOTE

This function ignores any physical dimensions passed through the Measure function. The units are always pixels.

NOTE

This function uses chain codes. It expects each measured object to be compact, that is, to have only one chain code. Additional chain codes are ignored, meaning that non-compact objects are not measured properly. Take care in providing the correct connectivity value: if you object is compact only with 2-connectivity, this measure will fail if you call Measure with a value of 1 for the connectivity.

SEE ALSO

Measure

FeatureAnisotropy2D, FeatureBendingEnergy, FeatureCenter,
FeatureChainCodeBendingEnergy, FeatureConvexArea, FeatureConvexPerimeter,
FeatureConvexity, FeatureDimension, FeatureExcessKurtosis, FeatureFeret,
FeatureGinertia, FeatureGmu, FeatureGravity, FeatureInertia, FeatureLongestChaincodeRun,
FeatureMass, FeatureMaxVal, FeatureMaximum, FeatureMean, FeatureMinimum,

FeatureMu, FeatureOrientation2D, FeatureP2A, FeaturePerimeter, FeatureRadius, FeatureShape, FeatureSize, FeatureSkewness, FeatureStdDev, FeatureSum, FeatureSurfaceArea

FeaturePerimeter

Measure the object's perimeter length

SYNOPSIS

```
#include "dip_measurement.h"
dip_int dip_FeaturePerimeterID ( void )
```

OUTPUT DATA TYPE

dip_float

FUNCTION

dip_FeaturePerimeterID returns the ID value of this measurement function, that is registered by Initialise.

This measures the perimeter of 2D objects by calculating the length of the chain code of its enclosing border. This function assumes that each object has a single connected border. The used method for measuring the length of the chain code is optimal for circles, and for a collection of objects that are randomly oriented, see the referenced literature for details. This function supports 2D images only.

NOTE

If any physical dimensions are passed to this function through Measure, only the sample distance along the first dimension is used. All other dimensions are assumed to be sampled the same way. This produces incorrect results for anisotropically sampled images.

NOTE

This function uses chain codes. It expects each measured object to be compact, that is, to have only one chain code. Additional chain codes are ignored, meaning that non-compact objects are not measured properly. Take care in providing the correct connectivity value: if you object is compact only with 2-connectivity, this measure will fail if you call Measure with a value of 1 for the connectivity.

LITERATURE

A.M. Vossepoel and A.W.M. Smeulders (1982), "Vector Code Probability and Metrication Error in the Representation of Straight Lines of Finite Length", Computer Graphics and Image Processing 20: 347-364

SEE ALSO

 ${\tt Measure, ImageChainCode, ChainCodeGetLength}$

FeatureRadius

Measure the object's radius statistics

SYNOPSIS

#include "dip_measurement.h"
dip_int dip_FeatureRadiusID (void)

OUTPUT DATA TYPE

dip_FloatArray

FUNCTION

dip_FeatureRadiusID returns the ID value of this measurement function, that is registered by Initialise.

This function obtains various statistics from the distance of each boundary pixel to the object's centre. The centre of the object is obtained from the centre of gravity of the border pixels only. See ChainCodeGetRadius for more information. This function supports 2D images only.

The values returned are:

RadiusMax	The maximum radius
RadiusMean	The average radius
RadiusMin	The minimum radius
RadiusStD	The standard deviation of the radii

NOTE

If any physical dimensions are passed to this function through Measure, only the sample distance along the first dimension is used. All other dimensions are assumed to be sampled the same way. This produces incorrect results for anisotropically sampled images.

NOTE

This function uses chain codes. It expects each measured object to be compact, that is, to have only one chain code. Additional chain codes are ignored, meaning that non-compact objects are not measured properly. Take care in providing the correct connectivity value: if you object is compact only with 2-connectivity, this measure will fail if you call Measure with a value of 1 for the connectivity.

SEE ALSO

 ${\tt Measure,\ ImageChainCode,\ ChainCodeGetRadius}$

FeatureShape

Measure shape parameters of the object

SYNOPSIS

```
#include "dip_measurement.h"
dip_int dip_FeatureShapeID ( void )
```

OUTPUT DATA TYPE

dip_FloatArray

FUNCTION

dip_FeatureShapeID returns the ID value of this measurement function, that is registered by Initialise

This function is a composite measurement function, that uses FeatureSize, FeaturePerimeter, and FeatureFeret to measure the following shape characteristics of 2D objects:

Squarity	area / (s * sp)
Circularity	area / (Pi/4 * sp^2)
Triangularity	area / (1/2 * s * sp)
Elliticity	area / (Pi/4 * s * sp)
Elongation	p / 1

with area the size, s the shortest Feret diameter, 1 the longest Feret diameter, sp the Feret diameter perpendicular to s, and p the perimeter of the object. The values in the output array are given in this order.

When the measured object is either a perfect square, circle, triangle or ellipse, the values obtained by FeatureShape will be 1.0.

NOTE

This function assumes isotropic sampling, even if the physical dimensions given through the Measure function say otherwise.

NOTE

This function uses chain codes. It expects each measured object to be compact, that is, to have only one chain code. Additional chain codes are ignored, meaning that non-compact objects are not measured properly. Take care in providing the correct connectivity value: if you object is compact only with 2-connectivity, this measure will fail if you call Measure with a value of 1 for the connectivity.

SEE ALSO

Measure

FeatureSize

Measure the object's area/volume

SYNOPSIS

```
#include "dip_measurement.h"
dip_int dip_FeatureSizeID ( void )
```

OUTPUT DATA TYPE

dip_int, dip_float

FUNCTION

dip_FeatureSizeID returns the ID value of this measurement function, that is registered by Initialise.

This function measures the object's size by counting the number of pixels having the same object ID. This measure is the optimal procedure for estimating the area (2D) or volume (3D) of an object with an arbitrary size. The measurement value's unit are in pixels (pixels^2 in 2D, pixels^3 in 3D).

If a dip_PhysicalDimensions parameter is given to Measure, the size of the object is given in physical units, and is a dip_float rather than a dip_int.

SEE ALSO

Measure

FeatureSkewness

Undocumented measurement function

FUNCTION

This measurement function is undocumented and not meant for public use. Measure

FeatureStdDev

Measure the standard deviation of the object's intensity

SYNOPSIS

```
#include "dip_measurement.h"
dip_int dip_FeatureStdDevID ( void )
```

OUTPUT DATA TYPE

dip_float

FUNCTION

dip_FeatureStdDevID returns the ID value of this measurement function, that is registered by Initialise

This function measures the standard deviation of the intensity in the intensity image (see Measure) of pixels inside the object.

SEE ALSO

Measure

FeatureSum

Measure the sum of the grey values of the object

SYNOPSIS

```
#include "dip_measurement.h"
dip_int dip_FeatureSumID ( void )
```

OUTPUT DATA TYPE

dip_float

FUNCTION

dip_FeatureSumID returns the ID value of this measurement function, that is registered by Initialise.

This function measures the sum of the grey-value in the intensity image (see Measure) of pixels inside the object.

SEE ALSO

Measure

FeatureSurfaceArea

Measure the area of the object's surface

SYNOPSIS

```
#include "dip_measurement.h"
dip_int dip_FeatureSurfaceAreaID ( void )
```

OUTPUT DATA TYPE

dip_float

FUNCTION

dip_FeatureSurfaceAreaID returns the ID value of this measurement function, that is registered by Initialise.

This function measures the area of a 3D object's surface using six-connected boundary voxels.

NOTE

If any physical dimensions are passed to this function through Measure, only the sample distance along the first dimension is used. All other dimensions are assumed to be sampled the same way. This produces incorrect results for anisotropically sampled images.

LITERATURE

J.C. Mullikin and P.W. Verbeek (1993), "Surface area estimation of digitized planes.", bioimaging 1(1): 6-16.

SEE ALSO

Measure

FeatureValueFunction

Measurement feature #value function

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error (*dip_FeatureValueFunction) ( measurement, featureID, objectID, physDims,
data, format, resources )
```

FUNCTION

The value function should return the measurement values produced by the measurement function, for one specific object. This function is called by MeasurementObjectValue.

ARGUMENTS

Data type	Name	Description
dip_Measurement	measurement	Measurement data structure
dip_int	featureID	Measurement function ID
dip_int	objectID	ID of the object to be measured
dip_PhysicalDimensions	physDims	Physical dimensions data structure
void **	data	Pointer to a measurement-specific internal
		data block
dipf_MeasurementValueFormat *	format	Pointer to a data format label, See
		MeasurementObjectValue
dip_Resources	resources	Resources tracking structure. See
		ResourcesNew

SEE ALSO

 ${\tt MeasurementFeatureRegister}, {\tt MeasurementObjectValue}$

FillBoundaryArray

Fill the border of array according to the boundary condition

SYNOPSIS

dip_Error DIP_TPI_FUNC(dip_FillBoundaryArray)(in, out, size, border, boundary)

FUNCTION


Set the values of the border pixels of an array. The pixels of out outside the range of the array in are set to a value determined by the boundary condition and the pixel values of in.

ARGUMENTS

Data type	Name	Description
void *	in	input array
void *	out	output array
dip_int	size	size of input array
dip_int	border	size of the extended borders
dip_Boundary	boundary	Boundary conditions

NOTE

The out array has to be allocated before this function is called, and should at least has the size of (size + 2 * border). Thus, border specifies the length of the border on both sides of the in array. Furthermore, the out pointer should point to that element in the out array that corresponds to the first element in the in array:



The enumerator dip_boundary contains the following constants:

Name	Description
DIP_BC_SYM_MIRROR	Symmetric mirroring
DIP_BC_ASYM_MIRROR	Asymmetric mirroring
DIP_BC_PERIODIC	Periodic copying
DIP_BC_ASYM_PERIODIC	Asymmetric periodic copying
DIP_BC_ADD_ZEROS	Extending the image with zeros
DIP_BC_ADD_MAX_VALUE	Extending the image with +infinity
DIP_BC_ADD_MIN_VALUE	Extending the image with -infinity

SEE ALSO

SeparableFrameWork

FindShift

Estimate the shift between images

SYNOPSIS

```
#include "dip_findshift.h"
dip_Error dip_FindShift ( in1, in2, out, method, parameter )
```

DATA TYPES

binary, integer, float, complex

FUNCTION

This function estimates the (sub-pixel) global shift between in1 and in2. The numbers found represent the shift of in1 with respect to in2, or the position of the first pixel of in2 in the coordinate system of in1. There are two methods that can be used: CPF and MTS. Both methods require that the shift be small. Therefore, first the integer pixel is calculated, and both images are cropped to the common part.

If method is 0, DIP_FSM_MTS is used.

If method is DIP_FSM_INTEGER_ONLY, integer shifts are calculated using cross correlation. This works for images of any dimensionality.

CPF

The CPF method (marked as FFTS in the literature below) uses the phase of the cross-correlation (as calculated by CrossCorrelationFT) to estimate the shift. parameter sets the amount of frequencies used in this estimation. The maximum value that makes sense is sqrt(1/2). Any larger value will give the same result. Choose smaller values to ignore the higher frequencies, which have a smaller SNR and are more affected by aliasing. If parameter is set to 0, the optimal found for images sub-sampled by a factor four will be used (parameter = 0.2).

This method only supports 2-D images.

MTS

The MTS method (marked as GRS in the literature below) uses a first order Taylor approximation of the equation in1(t) = in2(t-s) at scale parameter. Setting parameter to zero, a scale of 1 will be used. This means that the images will be smoothed with a Gaussian kernel of 1. This method is more accurate than CPF.

This method supports images with a dimensionality between 1 and 3.

ITER

The ITER method is an iterative version of the MTS method. It is known that a single estimation with MTS has a bias due to truncation of the Taylor expansion series (Pham et al., 2005) The bias can be expressed as a polynomial of the subpixel displacements. As a result, if the MTS method is applied iteratively, and the shift is refined after each iteration, the bias eventually becomes negligible. By using just 3 iterations, and noticing that $\log(\text{bias_increment})$ is a linear sequence, it is possible to correct for the bias up to O(1e-6).

Set parameter to 0 for normal behaviour. Other bahaviour is not supported, but: a parameter in the range (0,0.1] speficies the desired accuracy; parameter<0 causes round(-parameter) iterations to be run.

This method supports images with a dimensionality between 1 and 3.

PROJ

The PROJ method computes the shift in each dimension separately, applying the ITER method on the various projections of the image onto a single axis. It is fast and fairly accurate for high SNR. Should not be used for low SNR.

This method supports images with any number of dimensions.

ARGUMENTS

Data type	Name	Description
dip_Image	in1	Input image
dip_Image	in2	Input image
dip_FloatArray	out	Estimated shift
dipf_FindShiftMethod	method	Estimation method
dip_float	parameter	Parameter

The dipf_FindShiftMethod enumeration consists of the following flags:

Name	Description
DIP_FSM_DEFAULT	Default method (MTS)
DIP_FSM_INTEGER_ONLY	Find only integer shift
DIP_FSM_CPF	Cross-correlation method
DIP_FSM_FFTS	Same
DIP_FSM_MTS	Taylor series method
DIP_FSM_GRS	Same
DIP_FSM_ITER	Iterative version of MTS
DIP_FSM_PROJ	One-dimensional projectoin method
DIP_FSM_NCC	Undocumented

LITERATURE

C.L. Luengo Hendriks, *Improved Resolution in Infrared Imaging Using Randomly Shifted Images*, M.Sc. Thesis, Delft University of Technology, 1998

T.Q. Pham, M. Bezuijen, L.J. van Vliet, K. Schutte, C.L. Luengo Hendriks, Performance of Optimal

Registration Estimators, In Proc. of SPIE 5817 - Visual Information Processing XIV, Defense and Security Symposium, Orlando, 2005

SEE ALSO

 ${\tt CrossCorrelationFT}$

FiniteDifference

A linear gradient filter

SYNOPSIS

```
#include "dip_linear.h"
dip_Error dip_FiniteDifference ( in, out, boundary, processDim, filter )
```

DATA TYPES

binary, integer, float, complex

FUNCTION

The FiniteDifference filter implements several basic one dimensional FIR convolution filters. The dimension in which the operation is to be performed is specified by processDim. The operation itself is selected with filter. The $(1\ 0\ -1)/2$, $(1\ -1\ 0)$ & $(0\ 1\ -1)$ are difference filters that approximate a first order derivative, the $(1\ -2\ 1)$ filter approximates a second order derivative operation. The triangular $(1\ 2\ 1)/4$ filter is a local smoothing filter.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input image
dip_Image	out	Output image
dip_BoundaryArray	boundary	Boundary conditions
dip_int	processDim	ProcessDim
dipf_FiniteDifference	filter	Filter selection

The dipf_FiniteDifference enumeration consists of the following flags:

Name	Description
DIP_FINITE_DIFFERENCE_M101	out[ii] = (in[ii+1] - in[ii-1])/2
DIP_FINITE_DIFFERENCE_OM11	out[ii] = in[ii+1] - in[ii]
DIP_FINITE_DIFFERENCE_M110	$\operatorname{out}[ii] = \operatorname{in}[ii] - \operatorname{in}[ii-1]$
DIP_FINITE_DIFFERENCE_1M21	out[ii] = in[ii-1] - 2*in[ii] + in[ii+1]
DIP_FINITE_DIFFERENCE_121	out[ii] = (in[ii-1] + 2*in[ii] + in[ii+1])/4

SEE ALSO

General information about convolution

FiniteDifferenceEx, SobelGradient, Uniform, Gauss, SeparableConvolution, Convolve1d, Derivative

FiniteDifferenceEx

A linear gradient filter

SYNOPSIS

```
#include "dip_linear.h"
dip_Error dip_FiniteDifferenceEx ( in, out, boundary, process, parOrder, smoothflag
)
```

DATA TYPES

binary, integer, float, complex

FUNCTION

The FiniteDifferenceEx filter implements several basic one dimensional FIR convolution filters. The difference between this function and FiniteDifference is that this one has an interface more similar to Gauss and Derivative: it can process different derivatives along different dimensions at the same time. The first derivative is a convolution with $(1\ 0\ -1)/2$, and the second derivative is a convolution with $(1\ -2\ 1)$. When parOrder is 0 for a dimension, either the triangular smoothing filter $(1\ 2\ 1)/4$ is applied (smoothflag set to DIP_TRUE), or the dimension is not processed at all (smoothflag set to DIP_FALSE).

Setting all process to DIP_TRUE, all parOrder to 0 except one dimension to 1, and smoothflag to DIP_TRUE yields the SobelGradient.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input image
dip_Image	out	Output image
dip_BoundaryArray	boundary	Boundary conditions
dip_BooleanArray	process	Dimensions to process
dip_IntegerArray	parOrder	Order of Derivative along each dimension
dip_Boolean	smoothflag	Whether or not to smooth in the non-derivative directions

SEE ALSO

General information about convolution

 ${\tt FiniteDifference,\,SobelGradient,\,Uniform,\,Gauss,\,Derivative}$

FloatArrayCopy

Copy an array

SYNOPSIS

dip_Error dip_FloatArrayCopy (dest, src, resources)

FUNCTION

This function copies the float array src to dest. The array dest is created by this function as well.

ARGUMENTS

Data type	Name	Description	
dip_FloatArray *	dest	Destination array	
dip_FloatArray src		Source array	
dip_Resources resources		Resources tracking structure. See ResourcesNew	

SEE ALSO

 ${\tt FloatArrayNew,\,FloatArrayFree,\,FloatArrayCopy,\,FloatArrayFind}$

FloatArrayFind

Find value in array

SYNOPSIS

dip_Error dip_FloatArrayFind (array, value, index, found)

FUNCTION

Finds a value in an array and "returns" its index in the array. If found is zero, FloatArrayFind will produce an error if value is not found, otherwise found obtains the search result (DIP_FALSE if value is not found).

ARGUMENTS

Data type	Name	Description	
dip_FloatArray	array	Array to find value in	
dip_float	value	Value to find	
dip_int *	index	Index of the found value	
dip_Boolean *	found	Value found or not	

SEE ALSO

FloatArrayNew, FloatArrayFree, FloatArrayCopy, FloatArrayFind

FloatArrayFree

Array free function

SYNOPSIS

dip_Error dip_FloatArrayFree (array)

FUNCTION

This function frees *array, and sets array to zero.

ARGUMENTS

Data type	Name	Description
dip_FloatArray *	array	Array

SEE ALSO

FloatArrayNew, FloatArrayFree, FloatArrayCopy, FloatArrayFind

ArrayFree, IntegerArrayFree, FloatArrayFree, ComplexArrayFree, BoundaryArrayFree, FrameWorkProcessArrayFree, DataTypeArrayFree, ImageArrayFree, BooleanArrayFree, VoidPointerArrayFree, StringArrayFree, CoordinateArrayFree

FloatArrayNew

Array allocation function

SYNOPSIS

dip_Error dip_FloatArrayNew (array, size, value, resources)

FUNCTION

This function allocates the size elements of a dip_FloatArray and sets the size of the array to size. Each array element is initialized with value.

ARGUMENTS

Data type	Name	Description	
dip_FloatArray * array		Array	
dip_int	size	Size	
dip_float	value	Initial value	
dip_Resources	resources	Resources tracking structure. See ResourcesNew	

SEE ALSO

FloatArrayNew, FloatArrayFree, FloatArrayCopy, FloatArrayFind

ArrayNew, IntegerArrayNew, FloatArrayNew, ComplexArrayNew, BoundaryArrayNew, FrameWorkProcessArrayNew, DataTypeArrayNew, ImageArrayNew, BooleanArrayNew, VoidPointerArrayNew, StringArrayNew, CoordinateArrayNew

Floor

Arithmetic function

SYNOPSIS

dip_Error dip_Floor (in, out)

DATA TYPES

binary, integer, ${f float}$

FUNCTION

Computes the floor of the input image values, and outputs a signed integer typed image.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	in	Input
$dip_{-}Image$	out	Output

SEE ALSO

Abs, Ceil, Sign, Truncate, Fraction, NearestInt

FourierTransform

Computes the Fourier transform

SYNOPSIS

```
#include "dip_transform.h"
dip_Error dip_FourierTransform ( in, out, trFlags, process, theFuture )
```

DATA TYPES

binary, integer, float, complex

FUNCTION

Performs a Fourier transform on in and places the result in out.

Normalisation: 1/sqrt(dimension) for each dimension.

Defaults: process may be zero, indicating that all dimensions should be processed.

Sampling in Fourier Domain (FD): Let one pixel in the spatial domain (SD) be Delta_SD [m], then one pixel in the FD is Delta_FD = $1/(Delta_SD * N) [m^-1]$, where N is the width of the image in pixels. As a consequence the maximal frequency in the FD image is $N/2 * 1/(Delta_SD * N) = 1/(2 * Delta_SD) [m^-1]$ and is thus independent of the image width N and only related to the Nyquist frequency. The frequency of one FD pixel is therefore related to the image width N.

Note: In consequence of the above the FD resolution will not be isotropic if the image size are not square.

Note: Spatial zero-padding of the image increases the FD resolution only apparently (empty magnification).

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
$\mathtt{dip}_{-}Image$	out	Output
dipf_FourierTransform	trFlags	Transform flags
dip_BooleanArray	process (0)	Dimensions to process
void *	theFuture	For future use, should be set to zero

The dipf_FourierTransform enumeration consists of the following flags:

Name	Description	
DIP_TR_FORWARD	Forward transformation	
DIP_TR_INVERSE	Inverse transformation	

SEE ALSO

 ${\tt HartleyTransform}$

Fraction

Arithmetic function

SYNOPSIS

dip_Error dip_Fraction (in, out)

DATA TYPES

binary, integer, float

FUNCTION

Computes the fraction of the input image values, and outputs a float typed image.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	out	Output

SEE ALSO

Abs, Ceil, Floor, Sign, Truncate, NearestInt

FrameWorkProcessArrayFree

Array free function

SYNOPSIS

dip_Error dip_FrameWorkProcessArrayFree (array)

FUNCTION

This function frees *array, and sets array to zero.

ARGUMENTS

Data type	Name	Description
dip_FrameWorkProcessArray *	array	Array

SEE ALSO

 ${\tt FrameWorkProcessArrayNew,\,FrameWorkProcessArrayFree}$

ArrayFree, IntegerArrayFree, FloatArrayFree, ComplexArrayFree, BoundaryArrayFree, FrameWorkProcessArrayFree, DataTypeArrayFree, ImageArrayFree, BooleanArrayFree, VoidPointerArrayFree, StringArrayFree, CoordinateArrayFree

FrameWorkProcessArrayNew

Array allocation function

SYNOPSIS

dip_Error dip_FrameWorkProcessArrayNew (array, size, value, resources)

FUNCTION

This function allocates the size elements of a dip_FrameWorkProcessArray and sets the size of the array to size. Each array element is initialized with value.

ARGUMENTS

Data type	Name	Description
dip_FrameWorkProcessArray *	array	Array
dip_int	size	Size
dip_FrameWorkProcess	value	Initial value
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

FrameWorkProcessArrayNew, FrameWorkProcessArrayFree

ArrayNew, IntegerArrayNew, FloatArrayNew, ComplexArrayNew, BoundaryArrayNew, FrameWorkProcessArrayNew, DataTypeArrayNew, ImageArrayNew, BooleanArrayNew, VoidPointerArrayNew, StringArrayNew, CoordinateArrayNew

FTBox

Generates the Fourier transform of a box

SYNOPSIS

```
#include "dip_generation.h"
dip_Error dip_FTBox ( image, length, scale, amplitude )
```

DATA TYPES

Output: sfloat, scomplex

FUNCTION

Generates the Fourier transform of a box with the half length of its sides specified by length and scale.

ARGUMENTS

Data type	Name	Description
dip_Image	image	Output Image
dip_float	length	Length
dip_FloatArray	scale	Scale
dip_float	amplitude	Amplitude

SEE ALSO

FTEllipsoid, FTSphere, FTCube, FTCross, FTGaussian

FTCross

Generates the Fourier transform of a cross

SYNOPSIS

```
#include "dip_generation.h"
dip_Error dip_FTCross ( image, length, scale, amplitude )
```

DATA TYPES

Output: sfloat, scomplex

FUNCTION

Generates the Fourier transform of a cross with the length of its sides specified by length and radius.

ARGUMENTS

Data type	Name	Description
dip_Image	image	Output Image
dip_float	length	Length of the cross' axes
dip_FloatArray	scale	Scale
dip_float	amplitude	Amplitude

SEE ALSO

FTEllipsoid, FTSphere, FTBox, FTCube, FTGaussian

FTCube

Generates the Fourier transform of a cube

SYNOPSIS

```
#include "dip_generation.h"
dip_Error dip_FTCube ( image, length, amplitude )
```

DATA TYPES

Output: sfloat, scomplex

FUNCTION

Generates the Fourier transform of a cube with the length of its sides equal to two times length.

ARGUMENTS

Data type	Name	Description	
$dip_{-}Image$	image	Output Image	
dip_float	length	Length	
dip_float	amplitude	Amplitude	

SEE ALSO

FTEllipsoid, FTSphere, FTBox, FTCross, FTGaussian

FTEllipsoid

Generates Fourier transform of a ellipsoid

SYNOPSIS

```
#include "dip_generation.h"
dip_Error dip_FTEllipsoid ( image, radius, scale, amplitude )
```

DATA TYPES

Output: sfloat, scomplex

FUNCTION

Generates the Fourier transform of an ellipsoid with the length of its axes specified by radius and scale.

ARGUMENTS

Data type	Name	Description
dip_Image	image	Output Image
dip_float	radius	Radius
dip_FloatArray	scale	Scale
dip_float	amplitude	Amplitude

LITERATURE

L.J. van Vliet, *Grey-Scale Measurements in Multi-Dimensional Digitized Images*, Ph.D. thesis Delft University of Technology, Delft University Press, Delft, 1993

KNOWN BUGS

This function is only implemented for images with a dimensionality up to three.

SEE ALSO

FTSphere, FTBox, FTCube, FTCross, FTGaussian

FTGaussian

Generates the Fourier transform of a Gaussian

SYNOPSIS

```
#include "dip_generation.h"
dip_Error dip_FTGaussian ( output, sigma, volume, cutoff )
```

DATA TYPES

Output: sfloat, scomplex

FUNCTION

Generates the Fourier transform of a Gaussian with sigma's sigma. (The Fourier transform of a Gaussian, is a Gaussian.) volume is the integral of the Gaussian in the spatial domain. The cutoff variable can be used to avoid the calculation of the exponent of large negative values, which is can be very time consuming. Values of the exponent that are below cutoff yield a 0 value for the exponent. When cutoff is set to 0 or a positive value, DIP_GENERATION_EXP_CUTOFF is used (it is defined as -50).

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	output	Output Image
dip_FloatArray	sigma	Sigma of the Gaussian
dip_float	volume	Total intensity of the Gaussian
dip_float	cutoff	Cutoff value for the exponent

SEE ALSO

FTEllipsoid, FTSphere, FTBox, FTCube, FTCross

${\tt FTSphere}$

Generated Fourier transform of a sphere

SYNOPSIS

```
#include "dip_generation.h"
dip_Error dip_FTSphere ( image, radius, amplitude )
```

DATA TYPES

Output: sfloat, scomplex

FUNCTION

Generates the Fourier transform of a sphere with radius radius and an amplitude of amplitude.

ARGUMENTS

Data type	Name	Description	
dip_Image	image	Output Image	
dip_float	radius	Radius	
dip_float	amplitude	Amplitude	

KNOWN BUGS

This function is only implemented for images with a dimensionality up to three.

SEE ALSO

FTEllipsoid, FTBox, FTCube, FTCross, FTGaussian

GaborIIR

Infinite impulse response filter

SYNOPSIS

```
#include "dip_iir.h"
dip_Error dip_GaborIIR ( in, out, boundary, ps, sigmas, frequencies, order,
truncation )
```

DATA TYPES

binary, integer, float

FUNCTION

Recursive infinite impulse response implementation of the Gabor filter.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
$\mathtt{dip}_{\mathtt{-}}\mathtt{Image}$	out	Output
dip_BoundaryArray	boundary	Boundary conditions
dip_BooleanArray	ps	Dimensions to process
dip_FloatArray	sigmas	Sigma of Gaussian
dip_FloatArray	frequencies	frequencies
dip_IntegerArray	order	order
dip_float	truncation	Truncation, see GlobalGaussianTruncationGet

SEE ALSO

GaussIIR

Gauss

Gaussian Filter

SYNOPSIS

```
#include "dip_linear.h"
dip_Error dip_Gauss ( in, out, boundary, process, sigmas, order, truncation )
```

DATA TYPES

binary, integer, float

FUNCTION

Finite impulse response implementation of a Gaussian convolution filter and Gaussian derivative convolution filters.

The Gaussian kernel is cut off at truncation times the sigma of the filter (in each dimension). The sum of the Gaussian's coefficients is normalised to one. A truncation of zero or less indicates that the global preferred truncation ought to be used, see GlobalGaussianTruncationGet. For the derivatives, the truncation value is increased slightly: the actual value for truncation used is truncation + 0.5*order. The minimum filter size is 3 pixels, or 5 pixels for the 3rd order derivative.

Both the process and the order parameter may be zero. If process is zero all dimensions are processed. If order is zero no derivatives are taken.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	out	Output
dip_BoundaryArray	boundary	Boundary conditions
dip_BooleanArray	process (0)	Dimensions to process
dip_FloatArray	sigmas	Sigma of Gaussian
dip_IntegerArray	order (0)	Order of Derivative along each dimension
dip_float	truncation	Truncation of Gaussian

LIMITATIONS

The order of the derivative is limited to the interval 0-3. Sigmas considerably smaller than 1.0 will yield nonsensical results.

SEE ALSO

See sections 9.4, "Smoothing operations", and 9.5, "Derivative-based operations", in Fundamentals of Image Processing.

General information about convolution

GaussFT, GaussIIR, Derivative, GlobalGaussianTruncationGet

GaussFT

Gaussian Filter through the Fourier Domain

SYNOPSIS

```
#include "dip_linear.h"
dip_Error dip_GaussFT ( in, out, sigmas, order, truncation )
```

DATA TYPES

binary, integer, float, complex

FUNCTION

Fourier Domain implementation of a Gaussian convolution filter and Gaussian derivative convolution filters. The Gaussian kernel in the Fourier Domain is cut off at the equivalent of truncation times sigmas. If truncation is smaller or equal to 0, it is cut off where the argument to exp is smaller than -50, as in FTGaussian.

The order parameter may be zero, in which case no derivatives are taken.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	out	Output
dip_FloatArray	sigmas	Sigma of Gaussian
dip_IntegerArray	order (0)	Order of Derivative along each dimension
dip_float	truncation	Truncation of Gaussian kernel, see
		GlobalGaussianTruncationGet

SEE ALSO

See sections 9.4, "Smoothing operations", and 9.5, "Derivative-based operations", in Fundamentals of Image Processing.

General information about convolution

Gauss, GaussIIR, Derivative

GaussianNoise

Generate an image disturbed by Gaussian noise

SYNOPSIS

```
#include "dip_noise.h"
dip_Error dip_GaussianNoise ( in, out, variance, random )
```

DATA TYPES

integer, **float**

FUNCTION

Generate an image disturbed by additive Gaussian noise. See GaussianRandomVariable for more information on the random number generator.

ARGUMENTS

Data type	Name	Description	
dip_Image	in	Input	
dip_Image	out	Output	
dip_float	variance	Variance of the Gaussian distribution the noise is drawn from	
dip_Random *	random	Pointer to a random value structure	

EXAMPLE

Get a image with additive Gaussian noise as follows:

```
dip_Image in, out;
dip_float variance;
dip_Random random;

variance = 1.0;
DIPXJ( dip_RandomSeed( &random, 0 ));
DIPXJ( dip_GaussianNoise( in, out, variance, &random ));
```

SEE ALSO

 ${\tt Gaussian Random Variable, Random Seed, Random Seed Vector, Uniform Noise, Poisson Noise, Binary Noise}$

GaussianRandomVariable

Gaussian random variable generator

SYNOPSIS

```
#include "dip_noise.h"
dip_Error dip_GaussianRandomVariable ( random, mean, variance, output1, output2 )
```

FUNCTION

GaussianRandomVariable uses the algorithm described by D.E. Knuth as the Polar Method to generate two Gaussian distributed random variables. See RandomVariable for more information on the random number generator.

ARGUMENTS

Data type	Name	Description	
dip_Random *	random	Pointer to a random value structure	
dip_float	mean	Mean of the distribution, the samples are drawn from	
dip_float	variance	Variance of the distribution, the samples are drawn from	
dip_float *	output1	First output value	
dip_float *	output2	Second output value	

EXAMPLE

Get two Gaussian random variable as follows:

```
dip_Random random;
dip_float mean, variance, value1, value2;

mean = 0.0;
variance = 1.0;
DIPXJ( dip_RandomSeed( &random, 0 ));
DIPXJ( dip_GaussianRandomVariable( &random, mean, variance, &value1, &value2 ));
```

LITERATURE

Knuth, D.E., Seminumerical algorithms, The art of computer programming, vol. 2, second edition Addison-Wesley, Menlo Park, California, 1981.

SEE ALSO

RandomVariable, RandomSeed, RandomSeedVector, UniformRandomVariable, PoissonRandomVariable, BinaryRandomVariable

GaussianSigma

Adaptive Gaussian smoothing filter

SYNOPSIS

```
#include "dip_filtering.h"
dip_Error dip_GaussianSigma ( in, out, boundary, sigma, gaussSigma, outputCount,
truncation )
```

DATA TYPES

integer, float

FUNCTION

The GaussianSigma filter is an adaptive Gauss-ian smoothing filter. The value of the pixel under investigation is replaced by the Gaussian-weighted average of the pixelvalues in the filter region which lie in the interval +/- 2 sigma from the value of the pixel that is filtered. The filter region is specified by gaussSigma and truncation. If outputCount is DIP_TRUE, the output values represent the number of pixel over which the average has been calculated. When threshold is DIP_TRUE, the pixel intensities are thresholded at +/- 2 sigma, when it is set to DIP_FALSE, the intensities are weighted with the Gaussian difference with the intensity of the center pixel.

With threshold set to DIP_FALSE, this filter is also known as the bilateral filter.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	in	Input image
dip_Image	out	Output image
dip_BoundaryArray	boundary	Boundary conditions
dip_float	sigma	Sigma
dip_FloatArray	gaussSigma	Sigma of Gaussian
dip_Boolean	outputCount	Output the Count
dip_float	truncation	Truncation of Gaussian, see
		GlobalGaussianTruncationGet

LITERATURE

John-Sen Lee, *Digital Image Smoothing and the Sigma Filter*, Computer Vision, Graphics and Image Processing, 24, 255-269, 1983

SEE ALSO

Sigma, BiasedSigma, Gauss

GaussIIR

Infinite impulse response filter

SYNOPSIS

```
#include "dip_iir.h"
dip_Error dip_GaussIIR ( in, out, boundary, process, sigmas, order, truncation )
```

DATA TYPES

binary, integer, float

FUNCTION

Recursive infinite impulse response implementation of the Gauss filter.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	out	Output
dip_BoundaryArray	boundary	Boundary conditions
dip_BooleanArray	process	Dimensions to process
dip_FloatArray	sigmas	Sigma of Gaussian
dip_IntegerArray	order	Order of Derivative
dip_IntegerArray	order	Order of the IIR Filter
dip_int	designMethod	Method of IIR design
dip_float	truncation	Truncation of Gaussian, see
		GlobalGaussianTruncationGet

SEE ALSO

Gauss, Derivative

GeneralConvolution

Genaral convolution filter

SYNOPSIS

```
#include "dip_linear.h"
dip_Error dip_GeneralConvolution ( in, psf, out, boundary )
```

DATA TYPES

integer, float, complex

FUNCTION

This function convolves the in image with the point spread function psf, directly in the spatial domain. If the kernel psf is separable, use the function SeparableConvolution instead. If psf is large (and not separable), use the function ConvolveFT instead.

If the image psf is even in size, the origin is taken as the pixel to the right of the middle.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input image
dip_Image	psf	Psf image
dip_Image	out	Output image
dip_BoundaryArray	boundary	Boundary conditions

SEE ALSO

General information about convolution

 ${\tt SeparableConvolution}, {\tt ConvolveFT}, {\tt Uniform}$

GeneralisedKuwahara

Generalised Kuwahara filter

SYNOPSIS

```
#include "dip_filtering.h"
dip_Error dip_GeneralisedKuwahara ( in, selection, out, se, boundary, param, shape,
minimum )
```

DATA TYPES

binary, integer, float

FUNCTION

This function is a generalisation of the Kuwahara filter in the sense that is does not use the variance criterion to select the smoothed value, but instead accepts an image with the selection values. The algorithm finds, for every pixel, the minimum or maximum (as specified with minimum) value of selection within the filter window (its size specified by param), and outputs the corresponding value in in. When in is the output of Uniform, and selection is the output of VarianceFilter, this function produces the same result as Kuwahara.

Only the rectangular, elliptic and diamond filter shapes are supported (DIP_FLT_SHAPE_RECTANGULAR, DIP_FLT_SHAPE_ELLIPTIC and DIP_FLT_SHAPE_DIAMOND). Other filter shapes can be implemented by setting shape to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, and passing a binary image in se. The "on" pixels define the shape of the filter window. Other values of shape are illegal.

If shape is not equal to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, se can be set to zero. When shape is set to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, param is ignored, and can be set to zero.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	in	Input
dip_Image	selection	Selection
dip_Image	out	Output
dip_Image	se	Custom filter window (binary)
dip_BoundaryArray	boundary	Boundary conditions
dip_FloatArray	param	Filter sizes
dip_FilterShape	shape	Filter shape
dip_Boolean	minimum	Select minimum or maximum?

The enumerator dip_FilterShape contains the following constants:

Name	Description
DIP_FLT_SHAPE_DEFAULT	Default filter window, same as
	DIP_FLT_SHAPE_RECTANGULAR
DIP_FLT_SHAPE_RECTANGULAR	Rectangular filter window, can be even in size
DIP_FLT_SHAPE_ELLIPTIC	Elliptic filter window, always odd in size
DIP_FLT_SHAPE_DIAMOND	Diamond-shaped filter window, always odd in size
DIP_FLT_SHAPE_PARABOLIC	Parabolic filter window (morphology only)
DIP_FLT_SHAPE_DISCRETE_LINE	Rotated line structuring element (morphology only)
DIP_FLT_SHAPE_INTERPOLATED_LINE	Rotated line structuring element, through interpolation
	(morphology only)
DIP_FLT_SHAPE_PERIODIC_LINE	(not implemented)
DIP_FLT_SHAPE_STRUCTURING_ELEMENT	Use se as filter window, can be any size

SEE ALSO

 ${\tt KuwaharaImproved,\,GeneralisedKuwaharaImproved,\,VarianceFilter,\,Uniform}$

GeneralisedKuwaharaImproved

Generalised Kuwahara filter

SYNOPSIS

#include "dip_filtering.h"

dip_Error dip_GeneralisedKuwaharaImproved (in, selection, out, se, boundary, param,
shape, threshold, minimum)

DATA TYPES

binary, integer, float

FUNCTION

This function implements an improved version of <code>GeneralisedKuwahara</code>, see that function's description for more information. This function adds a <code>threshold</code> parameter that avoids false edges in uniform regions. If the difference between maximal and minimal values within the filter window is smaller or equal to <code>threshold</code>, the centre pixel is taken, instead of the minimum (or maximum). Setting <code>threshold</code> to zero yields the same result as <code>GeneralisedKuwahara</code>.

Only the rectangular, elliptic and diamond filter shapes are supported (DIP_FLT_SHAPE_RECTANGULAR, DIP_FLT_SHAPE_ELLIPTIC and DIP_FLT_SHAPE_DIAMOND). Other filter shapes can be implemented by setting shape to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, and passing a binary image in se. The "on" pixels define the shape of the filter window. Other values of shape are illegal.

If shape is not equal to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, se can be set to zero. When shape is set to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, param is ignored, and can be set to zero.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	selection	Selection
dip_Image	out	Output
dip_Image	se	Custom filter window (binary)
dip_BoundaryArray	boundary	Boundary conditions
dip_FloatArray	param	Filter sizes
dip_FilterShape	shape	Filter shape
dip_float	threshold	Minimal value difference within window
dip_Boolean	minimum	Select minimum or maximum?

The enumerator dip_FilterShape contains the following constants:

Name	Description
DIP_FLT_SHAPE_DEFAULT	Default filter window, same as
	DIP_FLT_SHAPE_RECTANGULAR
DIP_FLT_SHAPE_RECTANGULAR	Rectangular filter window, can be even in size
DIP_FLT_SHAPE_ELLIPTIC	Elliptic filter window, always odd in size
DIP_FLT_SHAPE_DIAMOND	Diamond-shaped filter window, always odd in size
DIP_FLT_SHAPE_PARABOLIC	Parabolic filter window (morphology only)
DIP_FLT_SHAPE_DISCRETE_LINE	Rotated line structuring element (morphology only)
DIP_FLT_SHAPE_INTERPOLATED_LINE	Rotated line structuring element, through interpolation
	(morphology only)
DIP_FLT_SHAPE_PERIODIC_LINE	(not implemented)
DIP_FLT_SHAPE_STRUCTURING_ELEMENT	Use se as filter window, can be any size

SEE ALSO

 ${\tt Kuwahara}, {\tt GeneralisedKuwahara}, {\tt KuwaharaImproved}, {\tt VarianceFilter}, {\tt Uniform}$

Get

Get a pixel value

SYNOPSIS

```
#include "dip_manipulation.h"
dip_Error dip_Get ( in, const, cor, adjust )
```

FUNCTION

This functions get the value of a pixel in image in at the coordinate cor. If cor is zero, the first pixel value is retrieved.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	in	Input
dip_Image	const	0-D output image
dip_IntegerArray	cor	Pixel coordinate
dip_Boolean	adjust	Adjust data type of output image

SEE ALSO

GetInteger, GetFloat, GetComplex, dip__PixelGetInteger, dip__PixelGetFloat, Set

GetComplex

Get complex pixel value

SYNOPSIS

```
#include "dip_manipulation.h"
dip_Error dip_GetComplex ( in, value, cor )
```

FUNCTION

This functions get the value of a pixel in image in at the coordinate cor. If cor is zero, the first pixel value is retrieved.

ARGUMENTS

Data type	Name	Description	
$\mathtt{dip}_{-}\mathtt{Image}$	in	Input	
dip_complex *	value	Value	
dip_IntegerArray	cor	Pixel coordinate	

SEE ALSO

Get, GetInteger, GetFloat, dip__PixelGetInteger, dip__PixelGetFloat, Set

GetFloat

Get float pixel value

SYNOPSIS

```
#include "dip_manipulation.h"
dip_Error dip_GetFloat ( in, value, cor )
```

FUNCTION

This functions get the value of a pixel in image in at the coordinate cor. If cor is zero, the first pixel value is retrieved.

ARGUMENTS

Data type	Name	Description	
$\mathtt{dip}_{-}\mathtt{Image}$	in	Input	
dip_float *	value	Value	
dip_IntegerArray	cor	Pixel coordinate	

SEE ALSO

Get, GetInteger, GetComplex, dip_PixelGetInteger, dip_PixelGetFloat, Set

GetInteger

Get integer pixel value

SYNOPSIS

```
#include "dip_manipulation.h"
dip_Error dip_GetInteger ( in, value, cor )
```

FUNCTION

This functions get the value of a pixel in image in at the coordinate cor. If cor is zero, the first pixel value is retrieved.

ARGUMENTS

Data type	Name	Description	
$\mathtt{dip}_{-}\mathtt{Image}$	in	Input	
dip_int *	value	Value	
dip_IntegerArray	cor	Pixel coordinate	

SEE ALSO

Get, GetFloat, GetComplex, dip__PixelGetInteger, dip__PixelGetFloat, Set

${\tt GetLibraryInformation}$

Support function

SYNOPSIS

```
#include "dip_information.h"
dip_Error dip_GetLibraryInformation ( info )
#include "dipio_image.h"
dip_Error dipio_GetLibraryInformation ( info )
```

FUNCTION

This function fills the given dip_LibraryInformation structure with information about the release version and date, copyright information and author information of the DIPlib library.

ARGUMENTS

Data type	Name	Description
dip_LibraryInformation*	info	DIPlib library information

GetLine

Get a line from an image

SYNOPSIS

```
#include "dip_manipulation.h"
dip_Error dip_GetLine ( in, out, cor, dimension )
```

DATA TYPES

binary, integer, float, complex

FUNCTION

Get a orthogonal line form an image. The position of the line in the image is specified by the coordinates at which its left most pixel (cor) should be placed and on which dimension of the image, the dimension of the line maps (dimension). If in has If in has a different type than out, it will be converted to the type of out.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input Image
$\mathtt{dip}_{-}\mathtt{Image}$	out	Output Line Image
dip_IntegerArray	cor	Coordinate in the image of the left most pixel of the line
dip_int	dimension	Dimension of the image on which the line's dimension maps

SEE ALSO

GetSlice, PutSlice, PutLine

GetMaximumAndMinimum

statistics function

SYNOPSIS

dip_Error dip_GetMaximumAndMinimum (in, mask, max, min)

DATA TYPES

 $integer,\,float\\$

FUNCTION

This function gets both the maximum and minimum of all the pixel values in the in image. Optionally, a mask image can be specified to exclude pixels from this search.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	in	Input image
$dip_{-}Image$	mask (0)	Mask image
dip_float	*max	Pointer to maximum variable
dip_float	*min	Pointer to minimum variable

SEE ALSO

Maximum, Minimum

GetObjectLabels

Lists object labels in image

SYNOPSIS

dip_Error dip_GetObjectLabels (in, mask, labels, nullIsObject, resources)

DATA TYPES

binary, integer

FUNCTION

This function produces an array of object labels present in the image in. Optionally, mask can mask the regions in in where to search for labels. The boolean nullIsObject specifies whether or not to treat the value zero as an object label.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	in	Input label image
$\mathtt{dip}_{-}\mathtt{Image}$	mask	Mask image
dip_IntegerArray *	labels	Array of labels
dip_Boolean	nullIsObject	treat the value zero ad an object label
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

Label, IntegerArrayFind

GetRank

Value selection function

SYNOPSIS

```
#include "dip_sort.h"
dip_Error dip_GetRank ( array, datatype, min, max, rank, value )
```

FUNCTION

GetRank gets the value at rank rank in the array array. min should be set to the first index of array, max to the last. dip_GetRank will use array for temporary storage, so the values in the array will be changed are this function is ready.

ARGUMENTS

Data type	Name	Description
dip_float *	array	Array to searched in
dip_DataType	datatype	
dip_int	min	minimal array index
dip_int	max	maximal array index
dip_int	rank	Rank
dip_float *	value	Value of the rank element

EXAMPLE

This example finds the median value for the array.

```
dip_float array[ SIZE ], median;
dip_int rank;

/* fill the array with values */

rank = SIZE/2;
DIPXX( dip_GetRank( array, DIP_DT_FLOAT, 0, (SIZE - 1), rank, &median ));
```

SEE ALSO

General information about sorting

DistributionSort, DistributionSortIndices, DistributionSortIndices16, InsertionSort, InsertionSortIndices, InsertionSortIndices16, QuickSort, QuickSortIndices, QuickSortIndices, SortIndices, SortIndices SortIndices SortIndices

GetSlice

Get a slice from an image

SYNOPSIS

```
#include "dip_manipulation.h"
dip_Error dip_GetSlice ( in, out, cor, dim1, dim2 )
```

DATA TYPES

binary, integer, float, complex

FUNCTION

Get a orthogonal slice from a image. The requested slice is selected by specifying its upper left corner (cor) and on which dimensions of the image, the dimensions of the slice map (dim1, dim2). If in has a different type than out, it will be converted to the type of out.

ARGUMENTS

Data type	Name	Description	
dip_Image	in	3D Input Image	
$\mathtt{dip}_{-}\mathtt{Image}$	out	2D Output Image	
dip_IntegerArray	cor	Coordinate in in of the upper left corner of the slice	
dip_int	dim1	Dimension of in on which the slice's first dimension maps	
dip_int	dim2	Dimension of in on which the slice's second dimensionmaps	

SEE ALSO

 ${\tt PutSlice},\,{\tt GetLine},\,{\tt PutLine}$

${\tt GetUniqueNumber}$

Obtain an unique value

SYNOPSIS

```
#include "dip_globals.h"
dip_Error dip_GetUniqueNumber ( number )
```

FUNCTION

This function gives an unique integer value. The value is unique is the sense that its value has not yet been returned by this function nor will it be returned by subsequent calls.

ARGUMENTS

Data type	Name	Description
dip_int *	number	Pointer to an integer in which the number is stored

${\tt Global Boundary Condition Get}$

Get global Boundary Conditions

SYNOPSIS

```
#include "dip_globals.h"
dip_Error dip_GlobalBoundaryConditionGet ( boundary, size, resources )
```

FUNCTION

This function allocates the boundary array array of size size with the global default boundary conditions for each dimension of the image. The initial values of this global array is DIP_BC_SYMMETRIC_MIRROR.

ARGUMENTS

Data type	Name	Description
dip_BoundaryArray *	boundary	Pointer to Boundary conditions
dip_int	size	Size of the new array
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

Boundary conditions

 ${\tt Global Boundary Condition Set, Global Gaussian Truncation Get, Global Gaussian Truncation Set, Global Filter Shape Get, Global Filter Shape Set}$

${\tt Global Boundary Condition Set}$

Set global boundary conditions

SYNOPSIS

```
#include "dip_globals.h"
dip_Error dip_GlobalBoundaryConditionSet ( boundary )
```

FUNCTION

This function sets the global boundary conditions equal to boundary.

ARGUMENTS

Data type	Name	Description
dip_BoundaryArray	boundary	Boundary conditions

SEE ALSO

 ${\tt Global Boundary Condition Get, Global Gaussian Truncation Get, Global Gaussian Truncation Set, Global Filter Shape Get, Global Filter Shape Set}$

GlobalFilterShapeGet

Get global filter shape value

SYNOPSIS

```
#include "dip_globals.h"
dip_Error dip_GlobalFilterShapeGet ( shape )
```

FUNCTION

This function gets the global default of the filter shape used by DIPlib's linear and morphology filters. The initial value of this global is DIP_FLT_SHAPE_RECTANGULAR.

This setting currently has no effect on any of the filters in DIPlib.

ARGUMENTS

Data type	Name	Description
dip_FilterShape *	shape	Filter shape

SEE ALSO

 ${\tt Global Boundary Condition Set, Global Boundary Condition Set, Global Gaussian Truncation Get, Global Gaussian Truncation Set, Global Filter Shape Set}$

GlobalFilterShapeSet

Set the global filter shape value

SYNOPSIS

```
#include "dip_globals.h"
dip_Error dip_GlobalFilterShapeSet ( shape )
```

FUNCTION

This function sets the global default of the filter shape used by DIPlib's linear and morphology filters. The initial value of this global is DIP_FLT_SHAPE_RECTANGULAR.

This setting currently has no effect on any of the filters in DIPlib.

ARGUMENTS

Data type	Name	Description	
dip_FilterShape	shape	Filter shape	

SEE ALSO

 ${\tt Global Boundary Condition Get, Global Boundary Condition Set, Global Gaussian Truncation Get, Global Gaussian Truncation Set, Global Filter Shape Get}$

GlobalGaussianTruncationGet

Get the global gaussian truncation

SYNOPSIS

```
#include "dip_globals.h"
dip_Error dip_GlobalGaussianTruncationGet ( truncation )
```

FUNCTION

This function gets the global default of the truncation used by the finite impluse response implementation of the Gauss (derivative) filter. The initial value of this global is 3.0.

ARGUMENTS

Data type	Name	Description
dip_float *	truncation	Gaussian truncation

SEE ALSO

 ${\tt Global Boundary Condition Get, Global Boundary Condition Set, Global Gaussian Truncation Set, Global Filter Shape Get, Global Filter Shape Set}$

GlobalGaussianTruncationSet

Set the global gaussian truncation

SYNOPSIS

```
#include "dip_globals.h"
dip_Error dip_GlobalGaussianTruncationSet ( truncation )
```

FUNCTION

This function sets the global default of the truncation used by the finite impluse response implementation of the Gauss (derivative) filter. The initial value of this global is 3.0.

ARGUMENTS

Data type	Name	Description
dip_float	truncation	Truncation

SEE ALSO

 ${\tt Global Boundary Condition Set, Global Boundary Condition Set, Global Gaussian Truncation Get, Global Filter Shape Get, Global Filter Shape Set}$

GradientDirection2D

Derivative filter

SYNOPSIS

```
#include "dip_derivatives.h"
dip_Error dip_GradientDirection2D ( in, out, boundary, ps, sigmas, tc, atanFlavour,
flavour )
```

DATA TYPES

Depends on the underlying implementation, but expect: binary, integer, **float**

FUNCTION

Computes the gradient direction of an image using the Derivative function. This functions supports only two dimensional images.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	out	Output
dip_BoundaryArray	boundary	Boundary conditions
dip_BooleanArray	ps	Dimensions to process
dip_FloatArray	sigmas	Sigma of Gaussian
dip_float	tc	Truncation of Gaussian, see
		GlobalGaussianTruncationGet
dip_GradientDirectionAtanFlavour	atanFlavour	Atan flavour
dip_DerivativeFlavour	flavour	Derivative flavour

The enumerator flavour parameter is one of:

Name	Description
DIP_DF_DEFAULT	Default derivative flavour (==DIP_DF_FIRGAUSS)
DIP_DF_FIRGAUSS	Gaussian family, FIR implementation, Gauss
DIP_DF_IIRGAUSS	Gaussian family, IIR implementation, GaussIIR
DIP_DF_FTGAUSS	Gaussian family, FT implementation, GaussFT
DIP_DF_FINITEDIFF	Finite difference implementation, FiniteDifferenceEx

SEE ALSO

See section 9.5, "Derivative-based operations", in Fundamentals of Image Processing.

DIP*lib* function reference 307

Derivative, GradientMagnitude, Laplace, Dgg, LaplacePlusDgg, LaplaceMinDgg

GradientMagnitude

Derivative filter

SYNOPSIS

```
#include "dip_derivatives.h"
dip_Error dip_GradientMagnitude ( in, out, boundary, ps, sigmas, tc, flavour )
```

DATA TYPES

Depends on the underlying implementation, but expect: binary, integer, **float**

FUNCTION

Computes the gradient magnitude of an image using the Derivative function.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	out	Output
dip_BoundaryArray	boundary	Boundary conditions
dip_BooleanArray	ps (0)	Dimensions to process
dip_FloatArray	sigmas	Sigma of the Gaussian
dip_float	tc	Gaussian truncation, see GlobalGaussianTruncationGet
dip_DerivativeFlavour	flavour	Derivative flavour

The enumerator flavour parameter is one of:

Name	Description	
DIP_DF_DEFAULT	Default derivative flavour (==DIP_DF_FIRGAUSS)	
DIP_DF_FIRGAUSS	Gaussian family, FIR implementation, Gauss	
DIP_DF_IIRGAUSS	Gaussian family, IIR implementation, GaussIIR	
DIP_DF_FTGAUSS	Gaussian family, FT implementation, GaussFT	
DIP_DF_FINITEDIFF	Finite difference implementation, FiniteDifferenceEx	

SEE ALSO

See section 9.5, "Derivative-based operations", in Fundamentals of Image Processing. Derivative, GradientDirection2D, Laplace, Dgg, LaplacePlusDgg, LaplaceMinDgg

Greater

Compare grey values in two images

SYNOPSIS

dip_Error dip_Greater (in1, in2, out)

DATA TYPES

binary, integer, float

FUNCTION

This function sets each pixel in out to "true" when for corresponding pixels in1 > in2. This is the same as Compare with the DIP_SELECT_GREATER selector flag.

in 2 can be a 0D image for comparison of pixel values with a single scalar value. This leads to a functionality similar to that of Threshold.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{\mathtt{-}}\mathtt{Image}$	in1	First input
$\mathtt{dip}_{-}\mathtt{Image}$	in2	Second input
dip_Image	out	Output

SEE ALSO

Compare, Threshold, Equal, Lesser, NotEqual, NotGreater, NotLesser, SelectValue, NotZero

GreyValuesInPixelTable

Copy greyvalues from image in pixel table

SYNOPSIS

```
#include "dip_pixel_table.h"
dip_Error dip_GreyValuesInPixelTable ( table, image, ptgreyvalues, resources )
```

DATA TYPES

integer, float

FUNCTION

This functions converts a grey-value image to a newly allocated floating-point array, in which each element is the grey value associated to a pixel in the pixel table. The image must have the same size and dimensionality as the pixel table's bounding box. For example:

```
dip_Image kernel, binkernel;
dip_PixelTable table;
dip_FloatArray values;
...
dip_NotZero( kernel, binkernel );
dip_BinaryImageToPixelTable( binkernel, &table, resources );
dip_GreyValuesInPixelTable( table, kernel, &values, resources );
...
process->filter->array[0].parameters = values;
dip_PixelTableFrameWork( in, out, boundary, process, table );
```

ARGUMENTS

Data type	Name	Description
dip_PixelTable	table	Pixel table
dip_Image	image	Grey-value image
dip_FloatArray *	ptgreyvalues	Array to which to write pixel grey values
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

Description of DIPlib's pixel tables

BinaryImageToPixelTable, PixelTableCreateFilter

GreyWeightedDistanceTransform

Grey weighted distance transform

SYNOPSIS

```
#include "dip_distance.h"
dip_Error dip_GreyWeightedDistanceTransform ( in, seed, out, distance, chamfer, neighborhood, metric )
```

DATA TYPES

in: integer, float
seed: binary

FUNCTION

GreyWeightedDistanceTransform determines the grey weighted distance transform of the object elements in the in image and returns the result in the out image. The implemented algorithm uses a heap sort for sorting the pixels to be processed.

The images in and seed must have the same dimensions. The out image will be converted to a sfloat typed image. The seed image defines the elements that are part of the object for which the GDT is determined. It can be any type of image where all image elements not equal to 0 are considered to be part of the object(s). Those elements that are neighboring an object element in the output image are considered seeds. Before any seeds are detected the borders of the out image are set to 0. The size of the border is determined by the chamfer metric size (see below). In case of a 3 by 3 chamfer metric the image border is one element, in case of a 5 by 5 chamfer it is 2 elements. Elements in the border are not considered seeds. If no valid seeds are found the routine will terminate with an Illegal value error code.

The chamfer metric is defined by two parameters: neighborhood and metric. neigborhood should supply the different relative addresses of the neighboring elements according to the chamfer metric. The first element neighborhood[0] contains the number of elements in the chamfer neighborhood. The next three elements contain the maximum number of elements a chamfer metric exceeds the central element. The rest of the elements (starting from the fifth element) contain addresses of the different chamfer elements relative to the central element. The metric array contains the corresponding chamfer metric value. An example of a 3x3 neighborhood array with the corresponding metric is:

where imagewidth represents the width of the image in image pixels. If both neighborhood and metric pointers are NULL, the chamfer variable can be set to either 1 (indicating a 3x3 or 3x3x3 chamfer using only 4 or 6 direct neighbors), 3 (indicating a 3x3 or 3x3x3 chamfer, using all neighbors) or 5 (indicating a 5x5 or 5x5x5 chamfer). In these cases a preset neighborhood and metric arrays will be used.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	in	Input image
dip_Image	seed	Seed image
dip_Image	out	Integrated grey-value over least-resistance path (output
		image)
$\mathtt{dip}_{-}\mathtt{Image}$	distance	Metric distance over least-resistance path (output image)
dip_int	chamfer	Chamfer distance metric
dip_IntegerArray	neighborhood	Neighborhood
dip_FloatArray	metric	Metric

LITERATURE

"An efficient uniform cost algorithm applied to distance transforms", B.J.H. Verwer, P.W. Verbeek, and S.T. Dekker, IEEE Transactions on Pattern Analysis and Machine Intelligence, vol. 11, no. 4, 1989, 425-429.

"Shading from shape, the eikonal equation solved by grey-weighted distance transform", P.W. Verbeek and B.J.H. Verwer, Pattern Recognition Letters, vol. 11, no. 10, 1990, 681-690.

"Local distances for distance transformations in two and three dimensions", B.J.H. Verwer, Pattern Recognition Letters, vol. 12, no. 11, 1991, 671-682.

"Distance Transforms, Metrics, Algorithms, and Applications", B.J.H. Verwer, Ph.D. thesis Delft University of Technology, Delft University Press, Delft, 1991.

"3-D Texture characterized by Accessibility measurements, based on the grey weighted distance transform", K.C. Strasters, A.W.M. Smeulders, and H.T.M. van der Voort, BioImaging, vol 2, no. 1, 1994, p. 1-21.

"Quantitative Analysis in Confocal Image Cytometry", Karel C. Strasters, Delft University Press, Delft, 1994. ISBN 90-407-1038-4, NUGI 841

KNOWN BUGS

GreyWeightedDistanceTransform works only on 2 or 3-dimensional images. It will not work if any of the images has different strides.

GreyWeightedDistanceTransform produces incomplete results in a 2-pixel border around the edge (4 for chamfer = 5). If this is an issue, consider adding 2 pixels on each side of your image. Make sure that in has high grey values in the border to avoid unexpected output.

The function <code>GrowRegionsWeighted</code> produces a grey-weighted distance transform without these limitations and with some other possibilities.

AUTHOR

Karel C. Strasters, adapted to DIPlib by Geert M.P. van Kempen

SEE ALSO

 ${\tt GrowRegionsWeighted}, {\tt EuclideanDistanceTransform}, {\tt VectorDistanceTransform}$

GrowRegions

Dilate the regions in a labelled image

SYNOPSIS

#include "dip_regions.h"
dip_Error dip_GrowRegions (in, grey, mask, out, connectivity, iterations, order)

DATA TYPES

in: binary, integer

grey: interger, float (converted to dip_sfloat)

mask: dip_uint8

FUNCTION

The regions in the input image in are grown with several options:

If grey is NULL, the regions are dilated iterations steps, according to connectivity (see The connectivity parameter), and optionally constrained by mask. This is the labelled equivalent to BinaryPropagation. If iterations is 0, the objects are dilated until no further change is possible. order is ignored.

If an image grey is given, the labels are grown in order of the grey-values in grey. order indicates whether pixels with high grey-values are added first or last. iterations is ignored, and mask is an optional constraint. This is a watershed algorithm with initial labels. The function Watershed does not accept an initial segmentation, so these two functions complement each other. Note that GrowRegions does not leave any watershed pixels in between the regions.

ARGUMENTS

Data type	Name	Description
$ ext{dip}_{-} ext{Image}$	in	Input binary or labelled image
dip_Image	grey	Input grey-value image
dip_Image	mask	Mask image
dip_Image	out	Output binary or labelled image
dip_int	connectivity	Connectivity
dip_int	iterations	Number of iterations
dipf_GreyValueSortOrder	order	Whether to grow from low to high or high to low

The dipf_GreyValueSortOrder enumeration consists of the following values:

Name	Description
DIP_GVSO_HIGH_FIRST	Process the pixels from high grey-value to low grey-value.
DIP_GVSO_LOW_FIRST	Process the pixels from low grey-value to high grey-value.

SEE	ALSO
-----	------

 ${\tt GrowRegionsWeighted}, \, {\tt Watershed}, \, {\tt BinaryPropagation}, \, {\tt Label}$

GrowRegionsWeighted

Grow labelled regions using grey-weighted distances

SYNOPSIS

#include "dip_regions.h"
dip_Error GrowRegionsWeighted (in, grey, mask, out, distance, pixelsize, chamfer,
metric)

DATA TYPES

in: binary, integer

grey: interger, float (converted to dip_sfloat)

 $mask: dip_uint8$

FUNCTION

The regions in the input image in are grown according to a grey-weighted distance metric; the weights are given by grey. The optional mask image mask limits the growing. out contains the grown regions, and distance, if not 0, contains the grey-weighted distance of each pixel in mask to the nearest pixel in in. Non-isotropic sampling is supported through pixelsize, which can be set to 0 to assume isotropic sampling. chamfer selects the size of the chamfer metric: 3 or 5. Set chamfer to 0 to use a custom metric given by the image metric. This image should be odd in size, and each pixel gives the distance to the center pixel. The pixels set to 0 will not be considered as neighbors.

The chamfer metric used is the following for chamfer==3 (with ps0=pixelsize->array[0] and ps1=pixelsize->array[1]):

sqrt(ps0*ps0+ps1*ps1)	ps1	sqrt(ps0*ps0+ps1*ps1)
ps0	0	ps0
sqrt(ps0*ps0+ps1*ps1)	ps1	sqrt(ps0*ps0+ps1*ps1)

and the following for chamfer==5:

0		0		0	
	sqrt(ps0*ps0+4*p	s1*ps1)	sqrt(ps0*ps0+4*p	s1*ps1)	
		ps1			
sqrt(4*ps0*ps0+	ps 1qp5(p)s0*ps0+ps1	*ps1)	sqrt(ps0*ps0+ps1	* pq11) (4*ps0*ps0+p	s1*ps1)
0	ps0	0	ps0	0	
		ps1			
sqrt(4*ps0*ps0+	ps 1qp5(p)s0*ps0+ps1	*ps1)	sqrt(ps0*ps0+ps1	* pq11) (4*ps0*ps0+p	s1*ps1)
0		0		0	
	sqrt(ps0*ps0+4*p	s1*ps1)	sqrt(ps0*ps0+4*p	s1*ps1)	

Setting chamfer to 0 and metric to an image with these values produces the same results as setting chamfer to 3 or 5.

The output image distance is comparable to the out image of GreyWeightedDistanceTransform, except that that function uses optimal chamfer distances whereas this one uses the (sub-optimal) true distance. In return, this function works on images of any dimensionality, allows for non-isotropic sampling, does not skip pixels close to the edge of the image, and can be used with a mask image to constrain the propagation. Note that the seed image in GreyWeightedDistanceTransform corresponds to the zero pixels of in for this function.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	in	Input binary or labelled image
dip_Image	grey	Input grey-value image
dip_Image	mask	Mask image
dip_Image	out	Output binary or labelled image
dip_Image	distance	Output distance image
dip_FloatArray	pixelsize	Pixel size
dip_int	chamfer	Chamfer distance
dip_Image	metric	Custom metric

LITERATURE

"3-D Texture characterized by Accessibility measurements, based on the grey weighted distance transform", K.C. Strasters, A.W.M. Smeulders, and H.T.M. van der Voort, BioImaging, vol 2, no. 1, 1994, p. 1-21.

"An efficient uniform cost algorithm applied to distance transforms", B.J.H. Verwer, P.W. Verbeek, and S.T. Dekker, IEEE Transactions on Pattern Analysis and Machine Intelligence, vol. 11, no. 4, 1989, 425-429.

SEE ALSO

GrowRegions, GreyWeightedDistanceTransform, Label

HartleyTransform

Computes the Hartley transform

SYNOPSIS

```
#include "dip_transform.h"
dip_Error dip_HartleyTransform ( in, out, trFlags, process )
```

DATA TYPES

binary, integer, float

FUNCTION

This function computes a Hartley transform on in and places the result in out.

Normalisation: 1/sqrt(dimension) for each dimension.

The main advantage of the Hartley transform over the Fourier transform is that is requires half the storage for real valued images. Note, that is also possible to directly reduce the storage requirements of the Fourier transform by just storing the right half plane, since for real valued images the left half plane can be derived from the right half using the symmetry properties of the Fourier transform.

Unfortunately there seem to be two definitions of the multi-dimensional Hartley transform (they are identical in the 1-D case). DIPlib implements the Bracewell (see below) variant, since this one is easy to implement and inherits the storage advantage from the 1-D case. The following are references which each use a different variant (all scaling factors have been dropped):

Bracewell, "Discrete Hartley Transform", J. Opt. Soc. Am, vol. 73, no. 12, December 1983:

Kenneth R. Castleman, "Digital image processing", Prentice Hall, 1996:

$$DHT(u,v) = Sum Sum I(x,y) cas(ux + vy)$$

$$y x$$

Using cas(a) = cos(a) + sin(a):

```
cas(ux)cas(vy) = cos(ux)cos(vy)+cos(ux)sin(vy)+sin(ux)cos(vy)+sin(ux)sin(vy)
cas(ux+vy) = cos(ux)cos(vy)+cos(ux)sin(vy)+sin(ux)cos(vy)-sin(ux)sin(vy)
```

A subtle difference. The two definitions have very similar properties, for example the convolution property.

In implementation terms, Bracewell is equivalent to perform the one-dimensional Hartley transform along each dimension. The Castleman variant is equivalent to the definition: DHT = re(DFT) - im(DFT). On a final note, I've not noticed mention of the difference between the two variants, so the indications Bracewell's and Castleman's variant are not and should not be accepted "labels" to refer to the variants (For both variants I have selected the first reference I came across, not chronologically the first reference to use the variant).

Defaults: process may be zero, indicating that all dimensions should be processed.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	out	Output
dipf_FourierTransform	trFlags	Transformation flags
dip_BooleanArray	process (0)	Dimensions to process

The dipf_FourierTransform enumeration consists of the following flags:

Name	Description	
DIP_TR_FORWARD	Forward transformation	
DIP_TR_INVERSE	Inverse transformation	

SEE ALSO

FourierTransform

HasContiguousData

Determines whether an image has all data contiguous in memory

SYNOPSIS

dip_Error dip_HasContiguousData(image, &answer)

FUNCTION

Determines whether an image has all data contiguous in memory. This can potentially not be the case if the image is an ROI, for example, or if it was allocated with strides that cause unused gaps in the image's memory block. If answer is not zero, the verdict is passed in this variable. Otherwise, HasContiguousData returns an error in case image does not have contiguous data.

ARGUMENTS

Data type	Name	Description
dip_Image	image	The image under investigation
dip_Boolean *	answer	The verdict

SEE ALSO

The image structure

HasNormalStride, ImageGetStride, IsScalar

HasNormalStride

Determines whether an image has a normal stride

SYNOPSIS

dip_Error dip_HasNormalStride(image, &answer)

FUNCTION

Determines whether an image has a normal stride. Normal stride is defined as a stride of 1 in the first dimension, a stride of image width in the second dimension, etc. If answer is not zero, the verdict is passed in this variable. Otherwise, HasNormalStride returns an error in case image does not have a normal stride.

ARGUMENTS

Data type	Name	Description
dip_Image	image	The image under investigation
dip_Boolean *	answer	The verdict

SEE ALSO

The image structure

HasContiguousData, ImageGetStride, IsScalar

${\tt HysteresisThreshold}$

Point Operation

SYNOPSIS

```
#include "dip_point.h"
dip_Error dip_HysteresisThreshold ( in, out, low, high )
```

DATA TYPES

integer, float

FUNCTION

Performs hysteresis thresholding. From the binary image (in>low) only those regions are selected for which at least one location also has (in>high). The output image will be a binary image with foreground pixel 1 and background pixel 0;

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	in	Input image
$\mathtt{dip}_{-}\mathtt{Image}$	out	Output image
dip_float	low	Lower threshold
dip_float	high	Higher threshold

SEE ALSO

Threshold, RangeThreshold, IsodataThreshold

IDivergence

difference measure

SYNOPSIS

dip_Error dip_IDivergence (in1, in2, mask, out)

DATA TYPES

binary, integer, float

FUNCTION

Calculates the I-divergence between each pixel value of in1 and in2. Optionally the mask image can be used to exclude pixels from the calculation by setting the value of these pixels in mask to zero.

The I-Divergence is defined as: $I(x,y) = x \ln(x/y) - (x - y)$ and is divided by the number of pixels. It is the -log of a possion distribution $p(x,y)=e^{-(-y)/x!-y^2}x$ with the stirling approximation for $\ln x!$. For x=0, the stirling approximation would fail, y is returned.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	in1	First input, Data:x
$\mathtt{dip}_{-}\mathtt{Image}$	in2	Second input, Model:y
dip_Image	mask	Mask
dip_Image	out	Output

LITERATURE

Why Least Squares and Maximum Entropy? An axiomatic approach to inference for linear inverse problems, I. Csiszar, The Annals of Statistics, 19, 2032-2066, 1991.

SEE ALSO

 ${\tt MeanError}, \, {\tt MeanSquareError}, \, {\tt RootMeanSquareError}, \, {\tt MeanAbsoluteError}, \, {\tt LnNormError}$

ImageArrayFree

Array free function

SYNOPSIS

dip_Error dip_ImageArrayFree (array)

FUNCTION

This function frees *array, and sets array to zero.

ARGUMENTS

Data type	Name	Description
dip_ImageArray *	array	Array

SEE ALSO

ImageArrayNew, ImageArrayFree

ArrayFree, IntegerArrayFree, FloatArrayFree, ComplexArrayFree, BoundaryArrayFree, FrameWorkProcessArrayFree, DataTypeArrayFree, ImageArrayFree, BooleanArrayFree, VoidPointerArrayFree, StringArrayFree, CoordinateArrayFree

ImageArrayNew

Array allocation function

SYNOPSIS

dip_Error dip_ImageArrayNew (array, size, resources)

FUNCTION

This function allocates the size elements of a dip_ImageArray and sets the size of the array to size.

ARGUMENTS

Data type	Name	Description
dip_ImageArray *	array	Array
dip_int	size	Size
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

ImageArrayNew, ImageArrayFree

ArrayNew, IntegerArrayNew, FloatArrayNew, ComplexArrayNew, BoundaryArrayNew, FrameWorkProcessArrayNew, DataTypeArrayNew, ImageArrayNew, BooleanArrayNew, VoidPointerArrayNew, StringArrayNew, CoordinateArrayNew

ImageAssimilate

Inherit properties of another image

SYNOPSIS

dip_Error dip_ImageAssimilate(example, target)

FUNCTION

Give the target image the same properties (type, data type, etc...) as the example image. The example image may be either "raw" or "forged". The target image is forged.

If the target was forged before calling this function, and it exactly matches the example, nothing happens. If it doesn't match the example, it is stripped before the properties are copied.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	example	An example image
dip_Image	target	The target image

SEE ALSO

 ${\tt ImageCopyProperties, ChangeDataType, ChangeToOd}$

ImageChainCode

Extracts all chain codes from a labeled image

SYNOPSIS

```
#include "dip_chaincode.h"
dip_Error dip_ImageChainCode ( objectIm, connectivity, objectID, chaincodearray,
resources )
```

DATA TYPES

integer

FUNCTION

Extracts the chain codes for the objects in objectIm (only 2D images supported) that are listed in objectID, assuming that each object is compact (i.e. it returns the chain code for only one border for each label ID in objectID). Chain codes are constructed according to connectivity, which can only be 1 or 2 (see The connectivity parameter). The output structure chaincodearray is allocated by this function and registered in resources.

The dip_ChainCodeArray structure, like all arrays in DIPlib, contains a size and an array element. Each element is of type dip_ChainCode, and accessed by chaincodearray->array[ii], where ii is between 0 and chaincodearray->size-1. Data in the dip_ChainCode structures can only be accessed through the corresponding access functions, see ChainCodeNew.

ARGUMENTS

Data type	Name	Description
dip_Image	objectIm	Labeled input image
dip_int	connectivity	Pixel connectivity of the objects
dip_IntegerArray	objectID	Array containing object label values
dip_ChainCodeArray *	chaincodearray	Output chain codes
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

ChainCodeNew, ChainCodeFree, ChainCodeArrayNew, ChainCodeArrayFree, ChainCodeGetSize, ChainCodeGetChains, ChainCodeGetStart, ChainCodeGetLabel, ChainCodeGetConnectivity, ChainCodeGetLength, ChainCodeGetLongestRun, ChainCodeGetFeret

ImageCheckBooleanArray

Check a boolean array

SYNOPSIS

dip_Error dip_ImageCheckBooleanArray (im, array, answer)

FUNCTION

This functions check whether the size of array is equal to the dimensionality of im. If answer is not zero, it will contain the result of the test, otherwise the DIP_E_ARRAY_ILLEGAL_SIZE will be set when the test has failed.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	im	Image
dip_BooleanArray	array	Array
dip_Boolean *	answer	Answer

SEE ALSO

 $\label{lem:lemma$

${\tt ImageCheckBoundaryArray}$

Check a boundary array

SYNOPSIS

dip_Error dip_ImageCheckBoundaryArray (im, array, answer)

FUNCTION

This functions check whether the size of array is equal to the dimensionality of im. If answer is not zero, it will contain the result of the test, otherwise the DIP_E_ARRAY_ILLEGAL_SIZE will be set when the test has failed.

ARGUMENTS

Data type	Name	Description
dip_Image	im	Image
dip_BoundaryArray	array	Boundary conditions
dip_Boolean *	answer	Answer

SEE ALSO

 $\label{lem:lemma$

${\tt ImageCheckComplexArray}$

Check a complex array

SYNOPSIS

dip_Error dip_ImageCheckComplexArray (im, array, answer)

FUNCTION

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	im	Image
dip_ComplexArray	array	Array
dip_Boolean *	answer	Answer

SEE ALSO

 $\label{lem:lemmageCheckIntegerArray} I mageCheckFloatArray, I mageCheckComplexArray, I mageCheckBoundaryArray$

${\tt ImageCheckFloatArray}$

Check a float array

SYNOPSIS

dip_Error dip_ImageCheckFloatArray (im, array, answer)

FUNCTION

This functions check whether the size of array is equal to the dimensionality of im. If answer is not zero, it will contain the result of the test, otherwise the DIP_E_ARRAY_ILLEGAL_SIZE will be set when the test has failed.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{\mathtt{-}}\mathtt{Image}$	im	Image
dip_FloatArray	array	Array
dip_Boolean *	answer	Answer

SEE ALSO

 $\label{lem:lemma$

ImageCheckIntegerArray

Check an integer array

SYNOPSIS

dip_Error dip_ImageCheckIntegerArray (im, array, answer)

FUNCTION

This functions check whether the size of array is equal to the dimensionality of im. If answer is not zero, it will contain the result of the test, otherwise the DIP_E_ARRAY_ILLEGAL_SIZE will be set when the test has failed.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	im	Image
dip_IntegerArray	array	Integer rray
dip_Boolean *	answer	Answer

SEE ALSO

 ${\tt ImageCheckIntegerArray, ImageCheckFloatArray, ImageCheckComplexArray, ImageCheckBoundaryArray}$

ImageCopyProperties

Copy the properties of an image

SYNOPSIS

dip_Error dip_ImageCopyProperties(example, target)

FUNCTION

Give the target image the same properties (type, data type, etc...) as the example image. The example image may be either "raw" or "forged", whereas the target image must be "raw". See <code>ImageAssimilate</code>.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	example	An example image
dip_Image	target	The target image

SEE ALSO

The image structure

ImageFileGetInfo

Get information about image in file (in dipIO)

SYNOPSIS

dip_Error dipio_ImageFileGetInfo (imInfo, filename, format, addExtensions, recognised, resources)

FUNCTION

This function opens an image file and fills a dipio_ImageFileInformation structure with the information from that file. imInfo is allocated by this function. Use ImageFileInformationFree to free this structure, or set the resources parameter for automatic deallocation. If format is 0, all different ImageRead functions are called in sequence until the correct format has been found. If you know the format, get the correct format ID through the registry functions. See File formats recognized by dipIO for a list of currently supported formats.

The boolean addExtensions specifies whether ImageFileGetInfo should try to add file format extensions to filename, if the registered file format reader fails to recognise filename straight away. The extensions are provided by the registered file readers.

If recognised is not zero, ImageFileGetInfo will set it to DIP_TRUE when it has been able to read filename, and it will set it to DIP_FALSE when it is not able to read the file. No error will be generated in this case.

ARGUMENTS

Data type	Name	Description
dipio_ImageFileInformation *	imInfo	Output image file information. See
		ImageFileInformationNew
dip_String	filename	File name
dip_int	format	ID of file format
dip_Boolean	addExtensions	Add file format extensions to filename
dip_Boolean *	recognised	Pointer to boolean containing the file read
		status
dip_Resources	resources	Resources tracking structure. See
		ResourcesNew

SEE ALSO

ImageReadCSVInfo, ImageReadGIFInfo, ImageReadICSInfo, ImageReadJPEGInfo, ImageReadLSMInfo, ImageReadPICInfo, ImageReadPNGInfo, ImageReadTIFFInfo, ImageRead, ImageReadColour, ImageReadR0I

${\tt ImageFileInformationFree}$

Free a Image File Information structure (in dipIO)

SYNOPSIS

dip_Error dipio_ImageFileInformationFree (imInfo)

FUNCTION

Frees a dipio_ImageFileInformation structure allocated through ImageFileInformationNew or by ImageFileGetInfo.

ARGUMENTS

Data type	Name	Description
$ ext{dipio}_{-} ext{ImageFileInformation} *$	imInfo	Structure to free

SEE ALSO

ImageFileInformationNew, ImageFileGetInfo

ImageFileInformationNew

Allocate an Image File Information structure (in dipIO)

SYNOPSIS

dip_Error dipio_ImageFileInformationNew (newImInfo, name, filetype, datatype, dims,
resources)

FUNCTION

Allocates a dipio_ImageFileInformation structure. It must be freed through ImageFileInformationFree, unless a resources parameter is given, in which case it will be freed automatically when freeing the resources. This structure is usually allocated by ImageFileGetInfo.

This function will fill out some of the values in the structure with the values given on the command line. All of these can be 0.

ARGUMENTS

Data type	Name	Description
dipio_ImageFileInformation *	newImInfo	Output structure
dip_String	name	Initial value for name
dip_String	filetype	Initial value for filetype
dip_DataType	datatype	Initial value for datatype
dip_IntegerArray	dims	Initial value for dimensions
dip_Resources	resources	Resources tracking structure. See ResourcesNew

The structure dipio_ImageFileInformation contains the following elements:

Data type	Name	Description
dip_String	name	File name
dip_String	filetype	File format string
dip_DataType	datatype	Data type of image
dip_int	sigbits	Significant bits
dip_IntegerArray	dimensions	Dimensions of image
dipio_PhotometricInterpretation	photometric	Color space
dip_PhysicalDimensions	physDims	Physical dimensions structure. See
		PhysicalDimensionsNew
dip_int	numberOfImages	Number of images in a TIFF file. If
		filetype is not "TIFF", this number
		is not set
dip_StringArray	history	History tags
dip_Resources	resources	Resource tracking; all elements within
		this structure are tracked here

The enumerator dipio_PhotometricInterpretation contains the following constants:

Name	Description
DIPIO_PHM_GREYVALUE	No colour information present; it's a grey-value image.
DIPIO_PHM_RGB	RGB image (the first three planes are red, green and blue)
DIPIO_PHM_RGB_NONLINEAR	Non-linear R'G'B' image (RGB channels to the power of 0.4)
DIPIO_PHM_CMY	CMY image (the first three planes are cyan, magenta and yellow)
DIPIO_PHM_CMYK	CMYK image (the first four planes are cyan, magenta, yellow and
	black)
DIPIO_PHM_CIELUV	CIE L*u'v' image (the first three planes are luminosity, u* and v*)
DIPIO_PHM_CIELAB	CIE L*a*b* image (the first three planes are luminosity, a* and b*)
DIPIO_PHM_CIEXYZ	CIE XYZ (the first three planes are X, Y and Z)
DIPIO_PHM_CIEYXY	CIE Yxy (the first three planes are Y, x and y)
DIPIO_PHM_HCV	HCV image (the first three planes are hue, chroma and value)
DIPIO_PHM_HSV	HSV image (the first three planes are hue, saturation and value)
DIPIO_PHM_DEFAULT	Same as DIPIO_PHM_GREYVALUE
DIPIO_PHM_GENERIC	Anything can be coded in the channels; the same as DIPIO_PHM_CMYK

Most file formats support only some of these.

SEE ALSO

 ${\tt ImageFileInformationFree,\ ImageFileGetInfo,\ Physical DimensionsNew}$

ImageForge

Allocate pixel data for an image

SYNOPSIS

dip_Error dip_ImageForge(image)

FUNCTION

Allocates a block of memory to store pixel data for an image. The image must be "raw", and will be "forged" afterwards. The routine will fail if the image fields do not contain a valid combination of values for the image type.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	image	The image for which the pixel data must be allocated

SEE ALSO

The image structure

 ${\tt ImageNew,\ ImageFree,\ ImageStrip,\ ImageCopyProperties}$

ImageFree

Free an image

SYNOPSIS

dip_Error dip_ImageFree(image)

FUNCTION

Free any pixel data associated with the image and return all fields to their initial ("raw") state by calling ImageStrip. Then the image structure itself is freed. Notice that you must pass a pointer to the image instead of the image itself. This allows ImageFree to set your image variable to zero, preventing further use of the now freed image.

Because ImageNew accepts a resources structure to keep track of allocated images, direct calls to ImageFree should be unnecessary.

ARGUMENTS

Data type	Name	Description
dip_Image *	image	A pointer to the image to be freed

SEE ALSO

The image structure

ImageNew, ImageForge, ImageStrip, ImageCopyProperties

ImageGetData

Get the data pointers of a set of images

SYNOPSIS

dip_Error dip_ImageGetData(in, idp, iflags, out, odp, oflags, flags, resources)

FUNCTION

Get the data pointers of a set of images. This function should not be called before the clean up of the previous invocation (by ResourcesFree) has been performed. Currently no clean up is required by ImageGetData, but any data pointers obtained by a previous call to this function should be considered invalid when calling this function. The iflags, oflags, and flags parameters are not used in the current version. These fields should be set to zero. The resources parameter is mandatory. Any of the image arrays' elements may be set to zero, indicating that it is to be ignored.

No functions that will possibly modify an image should be called after the call to ImageGetData and before its clean up. The proper time to call ImageGetPlane and ImageGetStride is right after the call to ImageGetData.

ARGUMENTS

Data type	Name	Description
dip_ImageArray	in	Array of input images
<pre>dip_VoidPointerArray *</pre>	idp	Returns input data pointers
dipf_ImageGetDataArray	iflags	Flags for input images
dip_ImageArray	out	Array of output images
dip_VoidPointerArray *	odp	Returns output data pointers
dipf_ImageGetDataArray	oflags	Flags for output images
${ t dipf_ImageGetData}$	flags	Flags for all images
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

The image structure

ImageGetPlane, ImageGetStride

${\tt ImageGetDataType}$

Read the data type field

SYNOPSIS

dip_Error dip_ImageGetDataType(image, dataType)

FUNCTION

Read the dip_Image data type field.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	image	An image
<pre>dip_DataType *</pre>	dataType	Returns the data type field

SEE ALSO

The image structure DIPlib's data types ImageSetDataType

${\tt ImageGetDimensionality}$

Read the dimensionality field

SYNOPSIS

dip_Error dip_ImageGetDimensionality(image, dimensionality)

FUNCTION

Read the dip_Image dimensionality field.

ARGUMENTS

Data type	Name	Description
dip_Image	image	An image
dip_int *	dimensionality	Returns the dimensionality field

SEE ALSO

The image structure

 ${\tt ImageGetDimensions}$

${\tt ImageGetDimensions}$

Read the dimensions array

SYNOPSIS

dip_Error dip_ImageGetDimensions(image, dimensions, resources)

FUNCTION

Read the dip_Image dimensions Array. The array that is used to return the dimensions in, is allocated by this routine using IntegerArrayNew.

ARGUMENTS

Data type	Name	Description
dip_Image	image	An image
dip_IntegerArray *	dimensions	Returns the dimensions Array
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

The image structure

 ${\tt ImageGetDimensionality}$

ImageGetPlane

Read the plane number

SYNOPSIS

dip_Error dip_ImageGetPlane(image, plane)

FUNCTION

Read the dip_Image plane number. For binary images this is the number of the bit in which the data is stored. For other data types it is meaningless. The proper time to call this function is right after ImageGetData.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	image	An image
dip_Int *	plane	Returns the plane number

SEE ALSO

The image structure

 ${\tt ImageGetData, ImageGetStride}$

ImageGetStride

Read the stride array

SYNOPSIS

dip_Error dip_ImageGetStride(image, &stride, resources)

FUNCTION

Read the dip_Image stride array. The array that is used to return the dimensions in, is allocated by this routine using IntegerArrayNew. The proper time to call this function is right after ImageGetData.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	image	An image
dip_IntegerArray *	stride	Returns the stride array
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

The image structure

ImageGetData, ImageGetPlane

${\tt ImageGetType}$

Read the type field

SYNOPSIS

dip_Error dip_ImageGetType(image, type)

FUNCTION

Read the dip_Image type field.

ARGUMENTS

Data type	Name	Description
dip_Image	image	An image
dip_ImageType *	type	Returns the type field

SEE ALSO

The image structure

ImageSetType

ImageIsGIF

Confirm that a file is a GIF file (in dipIO)

SYNOPSIS

```
#include "dipio_gif.h"
dip_Error dipio_ImageIsGIF ( filename, veredict )
```

FUNCTION

This function verifies that the file is an GIF file. veredict is set to DIP_TRUE if it is, and to DIP_FALSE if it isn't.

ARGUMENTS

Data type	Name	Description
dip_String	filename	File name
dip_Boolean *	veredict	Set to DIP_TRUE or DIP_FALSE

SOFTWARE

This function uses GifLib (version 4.1.0 or later), which supports GIF 87a & 98a. Copyright (c)1997 Eric S. Raymond

SEE ALSO

ImageWriteGIF, ImageReadGIF

ImageIsICS

Confirm that a file is an ICS file (in dipIO)

SYNOPSIS

```
#include "dipio_ics.h"
dip_Error dipio_ImageIsICS ( filename, veredict )
```

FUNCTION

This function verifies that the file is an ICS file. veredict is set to DIP_TRUE if it is, and to DIP_FALSE if it isn't.

ARGUMENTS

Data type	Name	Description	
dip_String	filename	File name	
dip_Boolean *	veredict	Set to DIP_TRUE or DIP_FALSE	

SOFTWARE

This function uses libics (version 1.3), which supports the ICS specification revision 2.0. Copyright (c)2000-2002 Cris L. Luengo Hendriks, Dr. Hans T.M. van der Voort and many others.

SEE ALSO

ImageWriteICS, ImageReadICS

ImageIsJPEG

Confirm that a file is a JPEG file (in dipIO)

SYNOPSIS

```
#include "dipio_jpeg.h"
dip_Error dipio_ImageIsJPEG ( filename, veredict )
```

FUNCTION

This function verifies that the file is a JPEG file. veredict is set to DIP_TRUE if it is, and to DIP_FALSE if it isn't.

ARGUMENTS

Data type	Name	Description	
dip_String	filename	File name	
dip_Boolean *	veredict	Set to DIP_TRUE or DIP_FALSE	

SOFTWARE

This function uses libjpeg (version 6b or later). Copyright (c)1994-1998, Thomas G. Lane.

SEE ALSO

 ${\tt ImageWriteJPEG,\ ImageReadJPEG,\ ImageReadJPEGInfo}$

ImageIsLSM

Confirm that a file is a Zeiss LSM file (in dipIO)

SYNOPSIS

```
#include "dipio_ics.h"
dip_Error dipio_ImageIsLSM ( filename, veredict )
```

FUNCTION

This function verifies that the file is a Zeiss LSM file. veredict is set to DIP_TRUE if it is, and to DIP_FALSE if it isn't.

ARGUMENTS

Data type	Name	Description	
dip_String	filename	File name	
dip_Boolean *	veredict	Set to DIP_TRUE or DIP_FALSE	

SOFTWARE

This function uses libtiff (version 3.6.1 or later), which supports the TIFF specification revision 6.0. Copyright (c)1988-1997 Sam Leffler and Copyright (c)1991-1997 Silicon Graphics, Inc.

SEE ALSO

 ${\tt ImageReadLSM}$

ImageIsTIFF

Confirm that a file is a TIFF file (in dipIO)

SYNOPSIS

```
#include "dipio_tiff.h"
dip_Error dipio_ImageIsTIFF ( filename, veredict )
```

FUNCTION

This function verifies that the file is a TIFF file. veredict is set to DIP_TRUE if it is, and to DIP_FALSE if it isn't.

ARGUMENTS

Data type	Name	Description	
dip_String	filename	File name	
dip_Boolean *	veredict	Set to DIP_TRUE or DIP_FALSE	

SOFTWARE

This function uses libtiff (version 3.6.1 or later), which supports the TIFF specification revision 6.0. Copyright (c)1988-1997 Sam Leffler and Copyright (c)1991-1997 Silicon Graphics, Inc.

SEE ALSO

 ${\tt ImageWriteTIFF}, \, {\tt ImageReadTIFF}$

ImageNew

Allocate a structure

SYNOPSIS

dip_Error dip_ImageNew(image, resources)

FUNCTION

Allocates a dip_Image structure and initializes all fields to their default values. The resulting image is in the "raw" state, see The image structure. By using ImageCopyProperties and the "ImageSet" access functions, the image fields can be set to their desired values. Pixel data for the image can be allocated using the ImageForge function, which will will put the image in the "forged" state.

ARGUMENTS

Data type	Name	Description
dip_Image *	image	Used to return the newly allocated image
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

The image structure

ImageFree, ImageForge, ImageStrip, ImageCopyProperties

ImageRead

Read grey-value image from file (in dipIO)

SYNOPSIS

dip_Error dipio_ImageRead (image, filename, format, addExtensions, recognised)

FUNCTION

This function reads an image from a file and puts it in image. image must be allocated before calling this function. If format is 0, all different ImageRead functions are called in sequence until the correct format has been found. If you know the format, get the correct format ID through the registry functions. See File formats recognized by dipIO for a list of currently supported formats.

The boolean addExtensions specifies whether ImageRead should try to add file format extensions to filename, if the registered file format reader fails to recognise filename straight away. The extensions are provided by the registered file readers.

If recognised is not zero, ImageRead will set it to DIP_TRUE when it has been able to read filename, and it will set it to DIP_FALSE when it is not able to read the file. No error will be generated in this case.

If the file contains a colour image, Colour2Gray is called. That is, this function always returns a grey-value image.

ARGUMENTS

Data type	Name	Description
dip_Image	image	Output image
dip_String	filename	File name
dip_int	format	ID of file format
dip_Boolean	addExtensions	Add file format extensions to filename
dip_Boolean *	recognised	Pointer to boolean containing the file read status

SEE ALSO

ImageReadColour, ImageReadROI, ImageFileGetInfo, ImageReadCSV, ImageReadGIF, ImageReadICS, ImageReadJPEG, ImageReadLSM, ImageReadPIC, ImageReadPNG, ImageReadTIFF, ImageWrite, Colour2Gray

ImageReadColour

Read colour image from file (in dipIO)

SYNOPSIS

dip_Error dipio_ImageReadColour (image, filename, photometric, format, addExtensions, recognised)

FUNCTION

This function reads an image from a file and puts it in image. image must be allocated before calling this function. It works the same as ImageRead, except that, if the file contains a colour image, Colour2Gray is not called. The returned image has an extra dimension with colours (always the last dimension), and photometric is set to the colour space.

ARGUMENTS

Data type	Name	Description
dip_Image	image	Output image
dip_String	filename	File name
dipio_PhotometricInterpretation*	photometric	Photometric interpretation (==colour
		space)
dip_int	format	ID of file format
dip_Boolean	addExtensions	Add file format extensions to filename
dip_Boolean *	recognised	Pointer to boolean containing the file
		read status

The enumerator dipio_PhotometricInterpretation contains the following constants:

Name	Description
DIPIO_PHM_GREYVALUE	No colour information present; it's a grey-value image.
DIPIO_PHM_RGB	RGB image (the first three planes are red, green and blue)
DIPIO_PHM_RGB_NONLINEAR	Non-linear R'G'B' image (RGB channels to the power of 0.4)
DIPIO_PHM_CMY	CMY image (the first three planes are cyan, magenta and yellow)
DIPIO_PHM_CMYK	CMYK image (the first four planes are cyan, magenta, yellow and
	black)
DIPIO_PHM_CIELUV	CIE L*u'v' image (the first three planes are luminosity, u* and v*)
DIPIO_PHM_CIELAB	CIE L*a*b* image (the first three planes are luminosity, a* and b*)
DIPIO_PHM_CIEXYZ	CIE XYZ (the first three planes are X, Y and Z)
DIPIO_PHM_CIEYXY	CIE Yxy (the first three planes are Y, x and y)
DIPIO_PHM_HCV	HCV image (the first three planes are hue, chroma and value)
DIPIO_PHM_HSV	HSV image (the first three planes are hue, saturation and value)
DIPIO_PHM_DEFAULT	Same as DIPIO_PHM_GREYVALUE
DIPIO_PHM_GENERIC	Anything can be coded in the channels; the same as DIPIO_PHM_CMYK

Most file formats support only some of these.

SEE ALSO

ImageRead, ImageReadROI, ImageFileGetInfo, ImageReadCSV, ImageReadGIF, ImageReadICS, ImageReadJPEG, ImageReadLSM, ImageReadPIC, ImageReadPNG, ImageReadTIFF, ImageWrite, Colour2Gray

${\tt ImageReadCSV}$

Read comma-separated values from file (in dipIO)

SYNOPSIS

```
#include "dipio_csv.h"
dip_Error dipio_ImageReadCSV ( image, filename, separator )
```

FUNCTION

This function reads the comma-separated values from a file and puts it in image. image must be allocated before calling this function.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	image	Output image
dip_String	filename	File name
char	separator	Separator character

SEE ALSO

ImageRead, ImageWriteCSV

${\tt ImageReadCSVInfo}$

Get information about image in comma-separated values file (in dipIO)

SYNOPSIS

```
#include "dipio_csv.h"
dip_Error dipio_ImageReadCSVInfo ( imInfo, filename )
```

FUNCTION

Opens a comma-separated values (CSV) file and fills a dipio_ImageFileInformation structure with the information from that file. imInfo must be allocated before calling this function.

ARGUMENTS

Data type	Name	Description
dipio_ImageFileInformation	imInfo	Output image file information. See
		ImageFileInformationNew
dip_String	filename	File name

SEE ALSO

 ${\tt ImageFileGetInfo,\,ImageReadCSV,\,ImageWriteCSV,\,ImageFileInformationNew}$

ImageReadGIF

Read a GIF image from file (in dipIO)

SYNOPSIS

```
#include "dipio_gif.h"
dip_Error dipio_ImageReadGIF ( image, filename, photometric )
```

FUNCTION

This function reads an image from a GIF file and puts it in image. image must be allocated before calling this function.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	image	Output image
dip_String	filename	File name
dipio_PhotometricInterpretation *	photometric	Photometric interpretation

The enumerator ${\tt dipio_PhotometricInterpretation}$ contains the following constants:

Name	Description
DIPIO_PHM_GREYVALUE	No colour information present; it's a grey-value image.
DIPIO_PHM_RGB	RGB image (the first three planes are red, green and blue)
DIPIO_PHM_RGB_NONLINEAR	Non-linear R'G'B' image (RGB channels to the power of 0.4)
DIPIO_PHM_CMY	CMY image (the first three planes are cyan, magenta and yellow)
DIPIO_PHM_CMYK	CMYK image (the first four planes are cyan, magenta, yellow and
	black)
DIPIO_PHM_CIELUV	CIE L*u'v' image (the first three planes are luminosity, u* and v*)
DIPIO_PHM_CIELAB	CIE L*a*b* image (the first three planes are luminosity, a* and b*)
DIPIO_PHM_CIEXYZ	CIE XYZ (the first three planes are X, Y and Z)
DIPIO_PHM_CIEYXY	CIE Yxy (the first three planes are Y, x and y)
DIPIO_PHM_HCV	HCV image (the first three planes are hue, chroma and value)
DIPIO_PHM_HSV	HSV image (the first three planes are hue, saturation and value)
DIPIO_PHM_DEFAULT	Same as DIPIO_PHM_GREYVALUE
DIPIO_PHM_GENERIC	Anything can be coded in the channels; the same as DIPIO_PHM_CMYK

Most file formats support only some of these.

SOFTWARE

This function uses GifLib (version 4.1.0 or later), which supports GIF 87a & 98a. Copyright (c)1997 Eric S. Raymond

SEE A	LSO
-------	-----

 ${\tt ImageReadColour, ImageWriteGIF, ImageIsGIF}$

ImageReadGIFInfo

Get information about image in GIF file (in dipIO)

SYNOPSIS

```
#include "dipio_gif.h"
dip_Error dipio_ImageReadGIFInfo ( imInfo, filename )
```

FUNCTION

Opens a GIF file and fills a dipio_ImageFileInformation structure with the information from that file. imInfo must be allocated before calling this function.

ARGUMENTS

Data type	Name	Description
dipio_ImageFileInformation	imInfo	Output image file information. See
		ImageFileInformationNew
dip_String	filename	File name

SOFTWARE

This function uses GifLib (version 4.1.0 or later), which supports GIF 87a & 98a. Copyright (c)1997 Eric S. Raymond

SEE ALSO

 ${\tt ImageFileGetInfo, ImageIsGIF, ImageReadGIF, ImageWriteGIF, ImageFileInformationNew}$

ImageReadICS

Read ICS image from file (in dipIO)

SYNOPSIS

```
#include "dipio_ics.h"
dip_Error dipio_ImageReadICS ( image, filename, photometric, offset, roisize,
sampling )
```

FUNCTION

This function reads the image in the ICS file and puts it in image. image must be allocated before calling this function. photometric is set to match the photometric interpretation of the data in the file, if it is recognised. The colour dimension is always the last dimension of the image (no matter how it was saved in the ICS file). offset, roisize and sampling define a ROI to read in. See the comments in ImageReadROI for more information on this.

ARGUMENTS

Data type	Name	Description
dip_Image	image	Output image
dip_String	filename	File name
dipio_PhotometricInterpretation*	photometric	Photometric interpretation
dip_IntegerArray	offset	ROI offset
dip_IntegerArray	roisize	ROI size
dip_IntegerArray	sampling	ROI sampling rate

The enumerator dipio_PhotometricInterpretation contains the following constants:

Name	Description
DIPIO_PHM_GREYVALUE	No colour information present; it's a grey-value image.
DIPIO_PHM_RGB	RGB image (the first three planes are red, green and blue)
DIPIO_PHM_RGB_NONLINEAR	Non-linear R'G'B' image (RGB channels to the power of 0.4)
DIPIO_PHM_CMY	CMY image (the first three planes are cyan, magenta and yellow)
DIPIO_PHM_CMYK	CMYK image (the first four planes are cyan, magenta, yellow and
	black)
DIPIO_PHM_CIELUV	CIE L*u'v' image (the first three planes are luminosity, u* and v*)
DIPIO_PHM_CIELAB	CIE L*a*b* image (the first three planes are luminosity, a* and b*)
DIPIO_PHM_CIEXYZ	CIE XYZ (the first three planes are X, Y and Z)
DIPIO_PHM_CIEYXY	CIE Yxy (the first three planes are Y, x and y)
DIPIO_PHM_HCV	HCV image (the first three planes are hue, chroma and value)
DIPIO_PHM_HSV	HSV image (the first three planes are hue, saturation and value)
DIPIO_PHM_DEFAULT	Same as DIPIO_PHM_GREYVALUE
DIPIO_PHM_GENERIC	Anything can be coded in the channels; the same as DIPIO_PHM_CMYK

Most file formats support only some of these.

SOFTWARE

This function uses libics (version 1.3 or later), which supports the ICS specification revision 2.0. Copyright (c)2000-2002 Cris L. Luengo Hendriks, Dr. Hans T.M. van der Voort and many others. This function uses zlib (version 1.1.4 or later). Copyright (c)1995-2002 Jean-loup Gailly and Mark Adler

SEE ALSO

ImageRead, ImageReadColour, ImageReadROI, ImageWriteICS, ImageIsICS

ImageReadICSInfo

Get information about image in ICS file (in dipIO)

SYNOPSIS

```
#include "dipio_ics.h"
dip_Error dipio_ImageReadICSInfo ( imInfo, filename )
```

FUNCTION

Opens a ICS file and fills a dipio_ImageFileInformation structure with the information from that file. imInfo must be allocated before calling this function.

ARGUMENTS

Data type	Name	Description
dipio_ImageFileInformation	imInfo	Output image file information. See
		ImageFileInformationNew
dip_String	filename	File name

SOFTWARE

This function uses libics (version 1.3), which supports the ICS specification revision 2.0. Copyright (c)2000-2002 Cris L. Luengo Hendriks, Dr. Hans T.M. van der Voort and many others.

SEE ALSO

ImageFileGetInfo, ImageIsICS, ImageReadICS, ImageWriteICS, ImageFileInformationNew

ImageReadJPEG

Read JPEG image from file (in dipIO)

SYNOPSIS

```
#include "dipio_jpeg.h"
dip_Error dipio_ImageReadJPEG ( image, filename, imageNumber, photometric )
```

FUNCTION

This function reads an image from the JPEG file and puts it in image. image must be allocated before calling this function. photometric is set to either DIPIO_PHM_RGB or DIPIO_PHM_GREYVALUE. If photometric is 0, the image will be read in as grey-value, even if color information is present in the file. Color images are allocated as 3D images, with the different samples along the 3rd. dimension.

ARGUMENTS

Data type	Name	Description
dip_Image	image	Output image
dip_String	filename	File name
dipio_PhotometricInterpretation *	photometric	Photometric interpretation

The enumerator dipio_PhotometricInterpretation contains the following constants:

Name	Description
DIPIO_PHM_GREYVALUE	No colour information present; it's a grey-value image.
DIPIO_PHM_RGB	RGB image (the first three planes are red, green and blue)
DIPIO_PHM_RGB_NONLINEAR	Non-linear R'G'B' image (RGB channels to the power of 0.4)
DIPIO_PHM_CMY	CMY image (the first three planes are cyan, magenta and yellow)
DIPIO_PHM_CMYK	CMYK image (the first four planes are cyan, magenta, yellow and
	black)
DIPIO_PHM_CIELUV	CIE L*u'v' image (the first three planes are luminosity, u* and v*)
DIPIO_PHM_CIELAB	CIE L*a*b* image (the first three planes are luminosity, a* and b*)
DIPIO_PHM_CIEXYZ	CIE XYZ (the first three planes are X, Y and Z)
DIPIO_PHM_CIEYXY	CIE Yxy (the first three planes are Y, x and y)
DIPIO_PHM_HCV	HCV image (the first three planes are hue, chroma and value)
DIPIO_PHM_HSV	HSV image (the first three planes are hue, saturation and value)
DIPIO_PHM_DEFAULT	Same as DIPIO_PHM_GREYVALUE
DIPIO_PHM_GENERIC	Anything can be coded in the channels; the same as DIPIO_PHM_CMYK

Most file formats support only some of these.

SOFTWARE

This function uses libjpeg (version 6b or later). Copyright (c)1994-1998, Thomas G. Lane.

SEE AL

 ${\tt ImageReadColour, ImageWriteJPEG, ImageIsJPEG, ImageReadJPEGInfo, Colour2Gray}$

ImageReadJPEGInfo

Get information about image in JPEG file (in dipIO)

SYNOPSIS

```
#include "dipio_jpeg.h"
dip_Error dipio_ImageReadJPEGInfo ( imInfo, filename, imageNumber )
```

FUNCTION

Opens a JPEG file and fills a dipio_ImageFileInformation structure with the information from that file. imInfo must be allocated before calling this function.

ARGUMENTS

Data type	Name	Description
dipio_ImageFileInformation	imInfo	Output image file information. See
		ImageFileInformationNew
dip_String	filename	File name

SOFTWARE

This function uses libjpeg (version 6b or later). Copyright (c)1994-1998, Thomas G. Lane.

SEE ALSO

 ${\tt ImageFileGetInfo,\ ImageIsJPEG,\ ImageReadJPEG,\ ImageWriteJPEG,\ ImageFileInformationNeward and a support of the property of the propert$

ImageReadLSM

Read Zeiss LSM image from file (in dipIO)

SYNOPSIS

```
#include "dipio_lsm.h"
dip_Error dipio_ImageReadLSM ( image, filename, offset, roisize, sampling, imInfo,
resources )
```

FUNCTION

This function reads the image in the Zeiss LSM file and puts it in image. image must be allocated before calling this function. Depending on the recording mode and the number of channels recorded, an image with 2 to 5 dimensions is returned. If multiple channels were recorded, they will be put along the last dimension (which can be either the third, fourth or fifth). The "stack", "time series plane" ans "time series z-scan" recording modes return a 3D image, the "time series stack" returns a 4D image, all other modes return a 2D image (including the "line" mode).

imInfo->physDims contains information on the distance between pixels. resources is only used to allocate the imInfo structure, so if imInfo is 0, resources can be 0 too.

offset, roisize and sampling define a region of interest to read in. See the comments in ImageReadROI for more information on this. Note that the channel dimension is part of this ROI.

ARGUMENTS

Data type	Name	Description
dip_Image	image	Output image
dip_String	filename	File name
dip_IntegerArray	offset	ROI offset
dip_IntegerArray	roisize	ROI size
dip_IntegerArray	sampling	ROI sampling rate
dipio_ImageFileInformation*	imInfo	Image file information structure
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SOFTWARE

This function uses libtiff (version 3.6.1 or later), which supports the TIFF specification revision 6.0. Copyright (c)1988-1997 Sam Leffler and Copyright (c)1991-1997 Silicon Graphics, Inc.

This function uses zlib (version 1.1.4 or later). Copyright (c)1995-2002 Jean-loup Gailly and Mark Adler

SEE ALSO

 ${\tt ImageReadROI, ImageIsLSM}$

ImageReadLSMInfo

Get information about image in LSM file (in dipIO)

SYNOPSIS

```
#include "dipio_lsm.h"
dip_Error dipio_ImageReadLSMInfo ( imInfo, filename )
```

FUNCTION

Opens a LSM file and fills a dipio_ImageFileInformation structure with the information from that file. imInfo must be allocated before calling this function.

ARGUMENTS

Data type	Name	Description
dipio_ImageFileInformation	imInfo	Output image file information. See
		ImageFileInformationNew
dip_String	filename	File name

SOFTWARE

This function uses libtiff (version 3.6.1 or later), which supports the TIFF specification revision 6.0. Copyright (c)1988-1997 Sam Leffler and Copyright (c)1991-1997 Silicon Graphics, Inc.

SEE ALSO

 ${\tt ImageFileGetInfo, ImageIsLSM, ImageReadLSM, ImageFileInformationNew}$

ImageReadPIC

Read BioRad PIC image from file (in dipIO)

SYNOPSIS

```
#include "dipio_pic.h"
dip_Error dipio_ImageReadPIC ( image, filename, offset, roisize, sampling, info,
resources )
```

FUNCTION

This function reads the image in the BioRAD PIC file and puts it in image. image must be allocated before calling this function. The information stored in the file is put in info.

offset and roisize define a region of interest to be read in. The ROI is clipped to the actual image data, so it is safe to specify a ROI that is too large. sampling can be used to read in a subset of the pixels of the chosen ROI. Any or all of these three parameters can be NULL.

ARGUMENTS

Data type	Name	Description
dip_Image	image	Output image
dip_String	filename	File name
dip_IntegerArray	offset	ROI offset
dip_IntegerArray	roisize	ROI size
dip_IntegerArray	sampling	ROI sampling rate
dipio_ImageFileInformation *	info	File information
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

ImageRead, ImageReadROI

${\tt ImageReadPICInfo}$

Get information about image in BioRad PIC file (in dipIO)

SYNOPSIS

```
#include "dipio_pic.h"
dip_Error dipio_ImageReadPICInfo ( imInfo, filename )
```

FUNCTION

Opens a BioRAD PIC file and fills a dipio_ImageFileInformation structure with the information from that file. imInfo must be allocated before calling this function.

ARGUMENTS

Data type	Name	Description
dipio_ImageFileInformation	imInfo	Output image file information. See
		ImageFileInformationNew
dip_String	filename	File name

SEE ALSO

 ${\tt ImageFileGetInfo,\,ImageReadPIC,\,ImageFileInformationNew}$

ImageReadR0I

Read a portion of a grey-value image from file (in dipIO)

SYNOPSIS

dip_Error dipio_ImageReadROI (image, filename, offset, roisize, sampling, format, addExtensions, recognised)

FUNCTION

This function reads an image from a file and puts it in image. image must be allocated before calling this function. It works the same as ImageRead, except that the user is allowed to specify a region of the image to read. This is done through the offset and roisize parameters. The ROI is clipped to the image size, so it is safe to specify a ROI that is too large. sampling can be used to read in a subset of the pixels of the chosen ROI. Any or all of these three parameters can be NULL.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	image	Output image
dip_String	filename	File name
dip_IntegerArray	offset	ROI offset
dip_IntegerArray	roisize	ROI size
dip_IntegerArray	sampling	ROI sampling rate
dip_int	format	ID of file format
dip_Boolean	addExtensions	Add file format extensions to filename
dip_Boolean *	recognised	Pointer to boolean containing the file read status

SEE ALSO

ImageRead, ImageReadColour, ImageFileGetInfo, ImageReadCSV, ImageReadGIF, ImageReadICS, ImageReadJPEG, ImageReadLSM, ImageReadPIC, ImageReadPNG, ImageReadTIFF, ImageWrite, Colour2Gray

ImageReadTIFF

Read TIFF image from file (in dipIO)

SYNOPSIS

```
#include "dipio_tiff.h"
dip_Error dipio_ImageReadTIFF ( image, filename, imageNumber, photometric )
```

FUNCTION

This function reads an image from the TIFF file and puts it in image. image must be allocated before calling this function. imageNumber indicates which image from the multi-page TIFF file to read. 0 is the first image. photometric is set to match the photometric interpretation of the TIFF file. Colour images and multi-sample images are allocated as 3D images, with the different samples along the 3rd dimension.

Multi-page TIFF files in which all pages contain an image of the same size and type, can be read as a 3D or 4D (Colour along the 4th dimension) image by setting imageNumber to -1. If the images are not of the same size and type, an error will be generated.

ARGUMENTS

Data type	Name	Description
dip_Image	image	Output image
dip_String	filename	File name
dip_int	imageNumber	Image number to read
dipio_PhotometricInterpretation *	photometric	Photometric interpretation

The enumerator dipio_PhotometricInterpretation contains the following constants:

Name	Description
DIPIO_PHM_GREYVALUE	No colour information present; it's a grey-value image.
DIPIO_PHM_RGB	RGB image (the first three planes are red, green and blue)
DIPIO_PHM_RGB_NONLINEAR	Non-linear R'G'B' image (RGB channels to the power of 0.4)
DIPIO_PHM_CMY	CMY image (the first three planes are cyan, magenta and yellow)
DIPIO_PHM_CMYK	CMYK image (the first four planes are cyan, magenta, yellow and
	black)
DIPIO_PHM_CIELUV	CIE L*u'v' image (the first three planes are luminosity, u* and v*)
DIPIO_PHM_CIELAB	CIE L*a*b* image (the first three planes are luminosity, a* and b*)
DIPIO_PHM_CIEXYZ	CIE XYZ (the first three planes are X, Y and Z)
DIPIO_PHM_CIEYXY	CIE Yxy (the first three planes are Y, x and y)
DIPIO_PHM_HCV	HCV image (the first three planes are hue, chroma and value)
DIPIO_PHM_HSV	HSV image (the first three planes are hue, saturation and value)
DIPIO_PHM_DEFAULT	Same as DIPIO_PHM_GREYVALUE
DIPIO_PHM_GENERIC	Anything can be coded in the channels; the same as DIPIO_PHM_CMYK

Most file formats support only some of these.

SOFTWARE

This function uses libtiff (version 3.6.1 or later), which supports the TIFF specification revision 6.0. Copyright (c)1988-1997 Sam Leffler and Copyright (c)1991-1997 Silicon Graphics, Inc.

This function uses zlib (version 1.1.4 or later). Copyright (c)1995-2002 Jean-loup Gailly and Mark Adler

KNOWN BUGS

TIFF is a very flexible file format. We have to limit the types of images that can be read to the more common ones, and to the ones dipIO writes. These are the most obvious limitations:

Tiled images are not supported.

Only 1, 4, 8, 16 and 32 bits per pixel integer grayvalues are read, as well as 32-bit and 64-bit floating point.

Only 4 and 8 bits per pixel colourmapped images are read. Colourmapped images contain 16-bit gray-values: stretching of the display will be necessary.

Class Y images (YCbCr) and Log-compressed images (LogLuv or LogL) are not supported.

SEE ALSO

 ${\tt ImageRead}, {\tt ImageReadColour}, {\tt ImageWriteTIFF}, {\tt ImageIsTIFF}, {\tt Colour2Gray}$

ImageReadTIFFInfo

Get information about image in TIFF file (in dipIO)

SYNOPSIS

```
#include "dipio_tiff.h"
dip_Error dipio_ImageReadTIFFInfo ( imInfo, filename, imageNumber )
```

FUNCTION

Opens a TIFF file and fills a dipio_ImageFileInformation structure with the information from that file. imInfo must be allocated before calling this function. imageNumber indicates which image from the multi-page TIFF file to get info on. 0 is the first image. imInfo->numberOfImages gives the number of pages in the file.

ARGUMENTS

Data type	Name	Description
$ ext{dipio}_{-} ext{ImageFileInformation}$	imInfo	Output image file information. See
		ImageFileInformationNew
dip_String	filename	File name
dip_int	imageNumber	Image number to query

SOFTWARE

This function uses libtiff (version 3.6.1 or later), which supports the TIFF specification revision 6.0. Copyright (c)1988-1997 Sam Leffler and Copyright (c)1991-1997 Silicon Graphics, Inc.

SEE ALSO

ImageFileGetInfo, ImageIsTIFF, ImageReadTIFF, ImageWriteTIFF, ImageFileInformationNew

ImagesCheck

Check properties of several images

SYNOPSIS

dip_Error dip_ImagesCheck(images, imageType, dataType, compareFlag, checkFlag)

FUNCTION

This function checks whether the image type and the data type of all the images in the array match with the imageType and dataType variables, and compares selected properties of the first image with those of the other images in the array. This comparison is done by calling ImagesCompareTwo. The checkFlag can be used to compare properties not supported by ImagesCompare. An error is returned by ImagesCheck if a check or comparison fails.

ARGUMENTS

Data type	Name	Description
dip_ImageArray *	images	Array of Images
$ ext{dip_ImageType}$	imageType	Image type of the first Image
dip_DataTypeProperties	dataType	Data type of the first Image. See DataTypeGetInfo
${ t dipf_ImagesCompare}$	compareFlag	Properties to compare. See ImagesCompare
dipf_ImagesCheck	checkFlag	Extra properties to be compared

dipf_ImagesCheck

Name	Description
DIP_CKIM_MAX_PRECISION_MATCH	Check whether data types match or match to the
	DIP_GTP_MAX_PRECISION DataType
DIP_CKIM_CASTING_TYPE_MATCH	Check whether data types match or match to the
	DIP_GTP_CAST_R2C or DIP_GTP_CAST_C2R types of the first
	image in image
DIP_CKIM_IGNORE_NULL_DIM_IMAGES	Ignore images with a zero dimensionality, this flag is usefull
	when 0d images are used as generic data containers of
	constants

SEE ALSO

ImagesCompareTwo, ImagesCompare, ImagesCheckTwo

ImagesCheckTwo

Check properties of two images

SYNOPSIS

 $\label{lem:compare} {\tt dip_ImagesCheckTwo(image1, image2, imageType, dataType, compareFlag, checkFlag)}$

FUNCTION

This function checks whether the image type and the data type of the two images match with the imageType and dataType variables, and compares selected properties of the two images. This comparison is done by calling ImagesCompareTwo. The checkFlag can be used to compare properties not supported by ImagesCompare. ImagesCheckTwo returns an error code if a check or comparison fails.

ARGUMENTS

Data type	Name	Description
dip_Image	image1	First Image
dip_Image	image2	Second Image
dip_ImageType	imageType	Image type of the first Image
dip_DataTypeProperties	dataType	Data type of the first Image. See DataTypeGetInfo
dipf_ImagesCompare	compareFlag	Properties to compare. See ImagesCompare
dipf_ImagesCheck	checkFlag	Extra properties to be compared

$dipf_ImagesCheck$

Name	Description
DIP_CKIM_MAX_PRECISION_MATCH	Check whether data types match or match to the
	DIP_GTP_MAX_PRECISION DataType
DIP_CKIM_CASTING_TYPE_MATCH	Check whether data types match or match to the
	DIP_GTP_CAST_R2C or DIP_GTP_CAST_C2R types
DIP_CKIM_IGNORE_NULL_DIM_IMAGES	Ignore images with a zero dimensionality, this flag is usefull
	when 0d images are used as generic data containers of
	constants

SEE ALSO

 ${\tt ImagesCompareTwo,\,ImagesCompare,\,ImagesCheck}$

ImagesCompare

Compare properties of several images

SYNOPSIS

dip_Error dip_ImagesCompare(images, condition, result)

FUNCTION

This function compares some standard fields of a number of Images or performs a full comparison. Only if the comparison result is true between each of the Images, will the final comparison result be true. The condition parameter specifies which properties should be tested. If 0, a full comparison of the Images is performed. Otherwise it should be a logical OR of the dipf_ImagesCompare flags. DIP_CPIM_MATCH_ALL_STANDARD is equivalent to all the flags OR'ed together. The difference between DIP_CPIM_MATCH_ALL_STANDARD and the full comparison specified by 0, is that the first will compare all the standard fields (type, data type, dimensions), whereas the other compares all fields relevant to a particular DIPlib Image type. This may exclude some of the standard fields and include some fields particular to the type of DIPlib Image in question. There are two modes of operation. If the result parameter is set, it is used to store the result of the comparison, a set of OR'ed dipf_ImagesCompare flags. If the result parameter is 0, an error is returned if the condition parameter and the resulting set of flags are not the same.

ARGUMENTS

Data type	Name	Description
$ exttt{dip_ImageArray}$	images	Array of Images
${ t dipf_ImagesCompare}$	condition	Properties to compare. 0 indicates full comparison
dipf_ImagesCompare *	result	Result: flags to indicate if the properties were the same. 0
		indicates that an error should be returned if the requested
		properties do not match

$dipf_ImagesCompare$

Name	Description
DIP_CPIM_DIMENSIONALITIES_MATCH	Dimensionalities match
DIP_CPIM_DIMENSIONS_MATCH	Dimensions match. The comparison is done up to the lower
	of the of the two dimensionalities
DIP_CPIM_SIZE_MATCH	Combination of DIP_CPIM_DIMENSIONALITIES_MATCH and
	DIP_CPIM_DIMENSIONS_MATCH
DIP_CPIM_TYPES_MATCH	Types match
DIP_CPIM_DATA_TYPES_MATCH	Data types match
DIP_CPIM_MATCH_ALL_STANDARD	All flags above OR'ed together
DIP_CPIM_STRIDES_MATCH	Strides match
DIP_CPIM_FULL_MATCH	Full match. Returned in result. To test for a full match
	use 0. Note: This is NOT equivalent to the other flags
	OR'ed together, and it cannot be used as condition

SEE ALSO

 ${\tt ImagesCompareTwo,\ ImagesCheckTwo,\ ImagesCheck}$

ImagesCompareTwo

Compare properties of two images

SYNOPSIS

dip_Error dip_ImagesCompareTwo(image1, image2, condition, result)

FUNCTION

This function compares some standard fields of two Images or performs a full comparison. The condition parameter specifies which properties should be tested. See ImagesCompare for more information. There are two modes of operation. If the result parameter is set, it is used to store the result of the comparison, a set of OR'ed dipf_ImagesCompare flags. If the result parameter is 0, an error is returned if the condition parameter and the resulting set of flags are not the same.

ARGUMENTS

Data type	Name	Description
dip_Image	image1	First Image
dip_Image	image2	Second Image
${ t dipf_ImagesCompare}$	condition	Properties to compare. See ImagesCompare
dipf_ImagesCompare*	result	Result: flags to indicate if the properties were the same. 0
		indicates that an error should be returned if the requested
		properties do not match

SEE ALSO

ImagesCompare, ImagesCheckTwo, ImagesCheck

${\tt ImageSetDataType}$

Set the data type field

SYNOPSIS

dip_Error dip_ImageSetDataType(image, dataType)

FUNCTION

Set the dip_Image data type field. The image must be "raw".

ARGUMENTS

Data type	Name	Description
dip_Image	image	An image
dip_DataType	type	The image data type

SEE ALSO

DIPlib's data types

The image structure

 ${\tt ImageGetDataType}$

${\tt ImageSetDimensions}$

Set the dimensions array

SYNOPSIS

dip_Error dip_ImageSetDimensions(image, dimensions)

FUNCTION

Set the dip_Image dimensions array. The image must be "raw".

ARGUMENTS

Data type	Name	Description
dip_Image	image	An image
dip_IntegerArray	dimensions	The image dimensions

SEE ALSO

The image structure

ImageGetDimensions, ChangeDimensions

ImageSetType

Set the image type field

SYNOPSIS

dip_Error dip_ImageSetType(image, type)

FUNCTION

Set the dip_Image type field. The image must be "raw".

ARGUMENTS

Data type	Name	Description
dip_Image	image	An image
dip_ImageType	type	The image type

SEE ALSO

The image structure ${\tt ImageGetType}$

${\tt ImageSort}$

Sort image data

SYNOPSIS

```
#include "dip_sort.h"
dip_Error dip_ImageSort ( in, out, algorithm )
```

FUNCTION

Produces an image (out) with the sorted pixel values of in.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	in	Input
dip_Image	out	Output
dip_Sort	algorithm	Sort algorithm

The sortType parameter is one of:

Name	Description
DIP_SORT_DEFAULT	Default sort algorithm
DIP_SORT_QUICK_SORT	Quick sort
DIP_SORT_DISTRIBUTION_SORT	Distribution sort
DIP_SORT_INSERTION_SORT	Insertion sort

SEE ALSO

General information about sorting

 ${\tt DistributionSort, InsertionSort, QuickSort, Sort, SortIndices, SortIndices16, ImageSortIndices}$

${\tt ImageSortIndices}$

Sort indices to image data

SYNOPSIS

```
#include "dip_sort.h"
dip_Error dip_ImageSortIndices ( in, indices, algorithm, flags )
```

FUNCTION

Sorts a list of indices rather than the data itself using the algorithm specified by algorithm. Unless the DIP_ISI_USE_INDICES, the indices image will be initialised with one index for each pixel in the image.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	indices	Indices
dip_Sort	algorithm	Sort algorithm
${ t dipf_{-}ImageSortIndices}$	flags	Flags

The sortType parameter is one of:

Name	Description
DIP_SORT_DEFAULT	Default sort algorithm
DIP_SORT_QUICK_SORT	Quick sort
DIP_SORT_DISTRIBUTION_SORT	Distribution sort
DIP_SORT_INSERTION_SORT	Insertion sort

The dipf_ImageSortIndices enumeration consists of the following flags:

Name	Description
DIP_ISI_USE_INDICES	Use the indices as given in the indices image

SEE ALSO

General information about sorting

 ${\tt DistributionSort, InsertionSort, QuickSort, Sort, ImageSort, SortIndices, SortIndices16}$

ImagesSeparate

Take care of in-place operations

SYNOPSIS

dip_Error dip_ImagesSeparate(in, out, newOut, saved, resources)

FUNCTION

First the list of output images is checked to see if any output image is used more than once. If this is the case an error is returned. Then the input and output images are examined. If any of the output images is also used as an input image, the function allocates a new image. This image is returned through the newOut array. For each output image a corresponding image is returned in this array. Either the original output image itself, or either a new image as discussed above. After the call to dip_ImagesSeparate, the images in the newOut array should be used instead of the original output images. After you are done processing the images, a call to ResourcesFree will perform the necessary post-processing. The post-processing consists of copying the data from the temporary output images to the original output images and freeing the temporary images. Because the post-processing is called through ResourcesFree, the resources parameter is mandatory. Any of the image arrays' elements may be set to zero, indicating that it is to be ignored.

The boolean saved array can be used to indicate that an input image has been stored in a safe place. In this case dip_ImagesSeparate will not have to allocate a temporary image if the input image is also used as an output image. The saved parameter may either be zero, which indicates that none of the input images has been saved, or it must be an array containing booleans corresponding each of the input images. DIP_TRUE indicates that the image has been saved.

ARGUMENTS

Data type	Name	Description
dip_ImageArray	in	An array of input images
dip_ImageArray	out	An array of output images
dip_ImageArray *	newOut	Returns an array containing the replacement output images
dip_BooleanArray	saved	An array of booleans indicating which input images are safely
		stored
dip_Resources	resources	Resources tracking structure. See ResourcesNew. May not be
		zero

SEE ALSO

ImageGetData

${\tt ImageStrip}$

Restore an image to its initial ("raw") state

SYNOPSIS

dip_Error dip_ImageStrip(image)

FUNCTION

Free any pixel data associated with the image and return all fields to their initial ("raw") state. Essentially the image is returned to the state it was in right after it was allocated with ImageNew.

ARGUMENTS

Data type	Name	Description	
dip_Image	image	The image to be stripped	

SEE ALSO

The image structure

ImageNew, ImageForge, ImageFree, ImageCopyProperties

ImageWrite

Write grey-value image to file (in dipIO)

SYNOPSIS

dip_Error dipio_ImageWrite (image, filename, physDims, format, compression)

FUNCTION

This function writes a grey-vlaue image to a file, overwriting any other file with the same name. physDims gives physical dimensions of the image, and can be set to 0 for default values. Not all file formats are able to store physical dimensions. Get the format ID through the registry functions. See File formats recognized by dipIO for a list of currently supported formats. If format is 0, ICSv2 is used.

ARGUMENTS

Data type	Name	Description
dip_Image	image	Output image
dip_String	filename	File name
dip_PhysicalDimensions	physDims	Physical dimensions structure. See
		PhysicalDimensionsNew
dip_int	format	ID of file format
dipio_Compression	compression	Compression method and level. See Compression
		methods for image files

SEE ALSO

ImageWriteColour, ImageWriteCSV, ImageWriteEPS, ImageWriteFLD, ImageWriteGIF, ImageWriteICS, ImageWriteJPEG, ImageWritePNG, ImageWritePS, ImageWriteTIFF, ImageRead

ImageWriteColour

Write colour image to file (in dipIO)

SYNOPSIS

 $\label{lem:compression} \verb|dip_Error dipio_ImageWriteColour (image, filename, photometric, physDims, format, compression)| \\$

FUNCTION

This function writes a colour image to a file, overwriting any other file with the same name. photometric must be set to the correct value. Not all file formats support all photometric values, and some don't support colour at all. physDims gives physical dimensions of the image, and can be set to 0 for default values. Not all file formats are able to store physical dimensions. Get the format ID through the registry functions. See File formats recognized by dipIO for a list of currently supported formats. If format is 0, ICSv2 is used.

ARGUMENTS

Data type	Name	Description
dip_Image	image	Output image
dip_String	filename	File name
dipio_PhotometricInterpretation	photometric	Photometric interpretation (==colour
		space)
dip_PhysicalDimensions	physDims	Physical dimensions structure. See
		PhysicalDimensionsNew
dip_int	format	ID of file format
dipio_Compression	compression	Compression method and level. See
		Compression methods for image files

The enumerator dipio_PhotometricInterpretation contains the following constants:

Name	Description
DIPIO_PHM_GREYVALUE	No colour information present; it's a grey-value image.
DIPIO_PHM_RGB	RGB image (the first three planes are red, green and blue)
DIPIO_PHM_RGB_NONLINEAR	Non-linear R'G'B' image (RGB channels to the power of 0.4)
DIPIO_PHM_CMY	CMY image (the first three planes are cyan, magenta and yellow)
DIPIO_PHM_CMYK	CMYK image (the first four planes are cyan, magenta, yellow and
	black)
DIPIO_PHM_CIELUV	CIE L*u'v' image (the first three planes are luminosity, u* and v*)
DIPIO_PHM_CIELAB	CIE L*a*b* image (the first three planes are luminosity, a* and b*)
DIPIO_PHM_CIEXYZ	CIE XYZ (the first three planes are X, Y and Z)
DIPIO_PHM_CIEYXY	CIE Yxy (the first three planes are Y, x and y)
DIPIO_PHM_HCV	HCV image (the first three planes are hue, chroma and value)
DIPIO_PHM_HSV	HSV image (the first three planes are hue, saturation and value)
DIPIO_PHM_DEFAULT	Same as DIPIO_PHM_GREYVALUE
DIPIO_PHM_GENERIC	Anything can be coded in the channels; the same as DIPIO_PHM_CMYK

Most file formats support only some of these.

SEE ALSO

ImageWrite, ImageWriteCSV, ImageWriteEPS, ImageWriteFLD, ImageWriteGIF, ImageWriteICS, ImageWriteJPEG, ImageWritePNG, ImageWritePS, ImageWriteTIFF, ImageRead, Colour2Gray

ImageWriteCSV

Write image to a comma-separated-value file (in dipIO)

SYNOPSIS

```
#include "dipio_csv.h"
dip_Error dipio_ImageWriteCSV ( image, filename, separator )
dip_Error dipio_ImageWriteCSV ( dip_Image, dip_String, char );
```

FUNCTION

This function writes the image to a comma-separated-values file, overwriting any other file with the same name. Optionally, an other separator than the comma can be specified using **separator**. Sometimes a space, a tab or a colon are used instead. Each line of image data is ended by a newline.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	image	Output image
dip_String	filename	File name
char	separator	Optional alternative separator character

SEE ALSO

ImageWrite, ImageReadCSV

ImageWriteEPS

Write image to Encapsulated PostScript file (in dipIO)

SYNOPSIS

```
#include "dipio_ps.h"
dip_Error dipio_ImageWriteEPS ( image, filename, photometric, xcm, ycm, border )
```

FUNCTION

This function writes the image to an Encapsulated PostScript file, overwriting any other file with the same name. Set the image size in xcm and ycm. border sets the size of the border around the image. If border is 0, no border is drawn. For colour images, set photometric (supported are RGB and CMYK) and write the colour channels along the third image dimension.

ARGUMENTS

Data type	Name	Description
dip_Image	image	Output image
dip_String	filename	File name
dipio_PhotometricInterpretation	photometric	Photometric interpretation
dip_float	xcm	X-size of image in cm.
dip_float	ycm	Y-size of image in cm.
dip_int	border	Thickness of border, zero is no border

The enumerator dipio_PhotometricInterpretation contains the following constants:

Name	Description
DIPIO_PHM_GREYVALUE	No colour information present; it's a grey-value image.
DIPIO_PHM_RGB	RGB image (the first three planes are red, green and blue)
DIPIO_PHM_RGB_NONLINEAR	Non-linear R'G'B' image (RGB channels to the power of 0.4)
DIPIO_PHM_CMY	CMY image (the first three planes are cyan, magenta and yellow)
DIPIO_PHM_CMYK	CMYK image (the first four planes are cyan, magenta, yellow and
	black)
DIPIO_PHM_CIELUV	CIE L*u'v' image (the first three planes are luminosity, u* and v*)
DIPIO_PHM_CIELAB	CIE L*a*b* image (the first three planes are luminosity, a* and b*)
DIPIO_PHM_CIEXYZ	CIE XYZ (the first three planes are X, Y and Z)
DIPIO_PHM_CIEYXY	CIE Yxy (the first three planes are Y, x and y)
DIPIO_PHM_HCV	HCV image (the first three planes are hue, chroma and value)
DIPIO_PHM_HSV	HSV image (the first three planes are hue, saturation and value)
DIPIO_PHM_DEFAULT	Same as DIPIO_PHM_GREYVALUE
DIPIO_PHM_GENERIC	Anything can be coded in the channels; the same as DIPIO_PHM_CMYK

Most file formats support only some of these.

SEE	ΑL	_SO
~		

 ${\tt ImageWrite}, {\tt ImageWriteColour}, {\tt ImageWritePS}$

${\tt ImageWriteFLD}$

Write image to AVS field file (in dipIO)

SYNOPSIS

```
#include "dipio_fld.h"
dip_Error dipio_ImageWriteFLD ( image, filename )
```

FUNCTION

This function writes the image to an AVS Field file.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	image	Output image
dip_String	filename	File name

SEE ALSO

 ${\tt ImageWrite}$

ImageWriteGIF

Write image to a GIF file (in dipIO)

SYNOPSIS

```
#include "dipio_gif.h"
dip_Error dipio_ImageWriteGIF ( image, filename, labelImage )
```

FUNCTION

This function writes the gray-value image to a GIF file, overwriting any other file with the same name. Optionally, an integer-typed image can be identified as a labeled image using labelImage. In that case a colour GIF image will be saved.

ARGUMENTS

Data type	Name	Description
dip_Image	image	Output image
dip_String	filename	File name
dip_Boolean	labelImage	Regard an integer image as a labeled image

SOFTWARE

This function uses GifLib (version 4.1.0 or later), which supports GIF 87a & 98a. Copyright (c)1997 Eric S. Raymond

SEE ALSO

ImageWrite, ImageReadGIF, ImageIsGIF

ImageWriteICS

Write ICS image to file (in dipIO)

SYNOPSIS

```
#include "dipio_ics.h"
dip_Error dipio_ImageWriteICS ( image, filename, photometric, physDims, history,
sigbits, version, compression )
```

FUNCTION

This function writes the image to an ICS file, overwriting any other file with the same name. version can set to 1 to use the ICS v.1.0 file format (the 2-file version), instead of ICS v.2.0. For colour images, set photometric and write the colour channels along the last image dimension. Set sigbits only if the number of significant bits is different from the full range of the data type of image (use 0 otherwise). physDims can be set to 0 to fill out default values. history can be 0 if you do not want to bother.

ARGUMENTS

Data type	Name	Description
dip_Image	image	Output image
dip_String	filename	File name
dipio_PhotometricInterpretation	photometric	Photometric interpretation
dip_PhysicalDimensions	physDims	Physical dimensions structure. See
		PhysicalDimensionsNew
dip_StringArray	history	Tags that are written to the history in the
		ICS header
dip_int	sigbits	Number of significant bits.
dip_int	version	ICS version
dipio_Compression	compression	Compression method and level. See
		Compression methods for image files

The enumerator dipio_PhotometricInterpretation contains the following constants:

Name	Description
DIPIO_PHM_GREYVALUE	No colour information present; it's a grey-value image.
DIPIO_PHM_RGB	RGB image (the first three planes are red, green and blue)
DIPIO_PHM_RGB_NONLINEAR	Non-linear R'G'B' image (RGB channels to the power of 0.4)
DIPIO_PHM_CMY	CMY image (the first three planes are cyan, magenta and yellow)
DIPIO_PHM_CMYK	CMYK image (the first four planes are cyan, magenta, yellow and
	black)
DIPIO_PHM_CIELUV	CIE L*u'v' image (the first three planes are luminosity, u* and v*)
DIPIO_PHM_CIELAB	CIE L*a*b* image (the first three planes are luminosity, a* and b*)
DIPIO_PHM_CIEXYZ	CIE XYZ (the first three planes are X, Y and Z)
DIPIO_PHM_CIEYXY	CIE Yxy (the first three planes are Y, x and y)
DIPIO_PHM_HCV	HCV image (the first three planes are hue, chroma and value)
DIPIO_PHM_HSV	HSV image (the first three planes are hue, saturation and value)
DIPIO_PHM_DEFAULT	Same as DIPIO_PHM_GREYVALUE
DIPIO_PHM_GENERIC	Anything can be coded in the channels; the same as DIPIO_PHM_CMYK

Most file formats support only some of these.

SOFTWARE

This function uses libics (version 1.3 or later), which supports the ICS specification revision 2.0. Copyright (c)2000-2002 Cris L. Luengo Hendriks, Dr. Hans T.M. van der Voort and many others.

This function uses ${\tt zlib}$ (version 1.1.4 or later). Copyright (c)1995-2002 Jean-loup Gailly and Mark Adler

SEE ALSO

ImageWrite, ImageWriteColour, ImageReadICS, ImageIsICS

ImageWriteJPEG

Write JPEG image to file (in dipIO)

SYNOPSIS

```
#include "dipio_jpeg.h"
dip_Error ImageWriteJPEG ( image, filename, photometric, physDims, complevel )
```

FUNCTION

This function writes the image to a JPEG file, overwriting any other file with the same name. photometric can set to let the function know how to write the JPEG image (supported colour space is RGB).

If photometric is not DIPIO_PHM_GRAYVALUE, a 3D image is expected, in which the different planes are stored along the 3rd dimension.

physDims gives physical dimensions of the image, which will be used to set the dots per inch property of the JPEG file. It can be set to 0 for default values (300 dpi). If the physDims->dimensionUnits is not given, meters are assumed.

complevel is a number between 1 (worst quality, smallest files) and 100 (best quality, largest files). Setting complevel to 0 uses the default compression level, which is 90.

ARGUMENTS

Data type	Name	Description
dip_Image	image	Output image
dip_String	filename	File name
dipio_PhotometricInterpretation	photometric	Photometric interpretation
dip_PhysicalDimensions	physDims	Physical dimensions structure. See
		PhysicalDimensionsNew
dipio_uint	complevel	Compression level

The enumerator dipio_PhotometricInterpretation contains the following constants:

Name	Description
DIPIO_PHM_GREYVALUE	No colour information present; it's a grey-value image.
DIPIO_PHM_RGB	RGB image (the first three planes are red, green and blue)
DIPIO_PHM_RGB_NONLINEAR	Non-linear R'G'B' image (RGB channels to the power of 0.4)
DIPIO_PHM_CMY	CMY image (the first three planes are cyan, magenta and yellow)
DIPIO_PHM_CMYK	CMYK image (the first four planes are cyan, magenta, yellow and
	black)
DIPIO_PHM_CIELUV	CIE L*u'v' image (the first three planes are luminosity, u* and v*)
DIPIO_PHM_CIELAB	CIE L*a*b* image (the first three planes are luminosity, a* and b*)
DIPIO_PHM_CIEXYZ	CIE XYZ (the first three planes are X, Y and Z)
DIPIO_PHM_CIEYXY	CIE Yxy (the first three planes are Y, x and y)
DIPIO_PHM_HCV	HCV image (the first three planes are hue, chroma and value)
DIPIO_PHM_HSV	HSV image (the first three planes are hue, saturation and value)
DIPIO_PHM_DEFAULT	Same as DIPIO_PHM_GREYVALUE
DIPIO_PHM_GENERIC	Anything can be coded in the channels; the same as DIPIO_PHM_CMYK

Most file formats support only some of these.

SOFTWARE

This function uses libjpeg (version 6b or later). Copyright (c)1994-1998, Thomas G. Lane.

SEE ALSO

ImageWrite, ImageWriteColour, ImageReadJPEG, ImageIsJPEG, ImageReadJPEGInfo

ImageWritePS

Write image to PostScript file (in dipIO)

SYNOPSIS

```
#include "dipio_ps.h"
dip_Error dipio_ImageWritePS ( image, filename, photometric, caption, xcm, ycm, border )
```

FUNCTION

This function writes the image to a PostScript file, overwriting any other file with the same name. Set the image size in xcm and ycm. border sets the size of the border around the image. If border is 0, no border is drawn. You can give the page a title through caption. For colour images, set photometric (supported are RGB and CMYK) and write the colour channels along the third image dimension.

ARGUMENTS

Data type	Name	Description
dip_Image	image	Output image
dip_String	filename	File name
dipio_PhotometricInterpretation	photometric	Photometric interpretation
dip_String	caption	Title for page
dip_float	xcm	X-size of image on page, in cm.
dip_float	ycm	Y-size of image on page, in cm.
dip_int	border	Thickness of border, zero is no border

The enumerator ${\tt dipio_PhotometricInterpretation}$ contains the following constants:

Name	Description
DIPIO_PHM_GREYVALUE	No colour information present; it's a grey-value image.
DIPIO_PHM_RGB	RGB image (the first three planes are red, green and blue)
DIPIO_PHM_RGB_NONLINEAR	Non-linear R'G'B' image (RGB channels to the power of 0.4)
DIPIO_PHM_CMY	CMY image (the first three planes are cyan, magenta and yellow)
DIPIO_PHM_CMYK	CMYK image (the first four planes are cyan, magenta, yellow and
	black)
DIPIO_PHM_CIELUV	CIE L*u'v' image (the first three planes are luminosity, u* and v*)
DIPIO_PHM_CIELAB	CIE L*a*b* image (the first three planes are luminosity, a* and b*)
DIPIO_PHM_CIEXYZ	CIE XYZ (the first three planes are X, Y and Z)
DIPIO_PHM_CIEYXY	CIE Yxy (the first three planes are Y, x and y)
DIPIO_PHM_HCV	HCV image (the first three planes are hue, chroma and value)
DIPIO_PHM_HSV	HSV image (the first three planes are hue, saturation and value)
DIPIO_PHM_DEFAULT	Same as DIPIO_PHM_GREYVALUE
DIPIO_PHM_GENERIC	Anything can be coded in the channels; the same as DIPIO_PHM_CMYK

Most file formats support only some of these.

SEE ALSO

ImageWrite, ImageWriteColour, ImageWriteEPS

ImageWriteTIFF

Write TIFF image to file (in dipIO)

SYNOPSIS

```
#include "dipio_tiff.h"
dip_Error ImageWriteTIFF ( image, filename, photometric, physDims, compression )
```

FUNCTION

This function writes the image to a TIFF file, overwriting any other file with the same name. photometric can set to let the function know how to write the TIFF image (supported colour spaces are RGB, CIE Lab and CMYK).

If photometric is not DIPIO_PHM_GRAYVALUE, a 3D image is expected, in which the different planes are stored along the 3rd dimension.

physDims gives physical dimensions of the image, which will be used to set the dots per inch property of the TIFF file. It can be set to 0 for default values (300 dpi). If the physDims->dimensionUnits is not given, meters are assumed.

ARGUMENTS

Data type	Name	Description
dip_Image	image	Output image
dip_String	filename	File name
dipio_PhotometricInterpretation	photometric	Photometric interpretation
dip_PhysicalDimensions	physDims	Physical dimensions structure. See
		PhysicalDimensionsNew
dipio_Compression	compression	Compression method and level. See
		Compression methods for image files

The enumerator dipio_PhotometricInterpretation contains the following constants:

Name	Description
DIPIO_PHM_GREYVALUE	No colour information present; it's a grey-value image.
DIPIO_PHM_RGB	RGB image (the first three planes are red, green and blue)
DIPIO_PHM_RGB_NONLINEAR	Non-linear R'G'B' image (RGB channels to the power of 0.4)
DIPIO_PHM_CMY	CMY image (the first three planes are cyan, magenta and yellow)
DIPIO_PHM_CMYK	CMYK image (the first four planes are cyan, magenta, yellow and
	black)
DIPIO_PHM_CIELUV	CIE L*u'v' image (the first three planes are luminosity, u* and v*)
DIPIO_PHM_CIELAB	CIE L*a*b* image (the first three planes are luminosity, a* and b*)
DIPIO_PHM_CIEXYZ	CIE XYZ (the first three planes are X, Y and Z)
DIPIO_PHM_CIEYXY	CIE Yxy (the first three planes are Y, x and y)
DIPIO_PHM_HCV	HCV image (the first three planes are hue, chroma and value)
DIPIO_PHM_HSV	HSV image (the first three planes are hue, saturation and value)
DIPIO_PHM_DEFAULT	Same as DIPIO_PHM_GREYVALUE
DIPIO_PHM_GENERIC	Anything can be coded in the channels; the same as DIPIO_PHM_CMYK

Most file formats support only some of these.

SOFTWARE

This function uses libtiff (version 3.6.1 or later), which supports the TIFF specification revision 6.0. Copyright (c)1988-1997 Sam Leffler and Copyright (c)1991-1997 Silicon Graphics, Inc.

This function uses ${\tt zlib}$ (version 1.1.4 or later). Copyright (c)1995-2002 Jean-loup Gailly and Mark Adler

SEE ALSO

ImageWrite, ImageWriteColour, ImageReadTIFF, ImageIsTIFF

Imaginary

Arithmetic function

SYNOPSIS

dip_Error dip_Imaginary (in, out)

DATA TYPES

binary, integer, float, $\mathbf{complex}$

FUNCTION

Computes the imaginary part of the input image values, and outputs a float typed image.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	out	Output

SEE ALSO

Modulus, Phase, Real

IncoherentOTF

Generates an incoherent OTF

SYNOPSIS

```
#include "dip_microscopy.h"
dip_Error dip_IncoherentOTF ( out, defocus, xNyquist, amplitude, otf )
```

DATA TYPES

Output: sfloat

FUNCTION

This function implements the formulae for a (defocused) incoherent OTF as described by Castleman. When defocus is unequal to zero, either the Stokseth approximation or the Hopkins approximation is used. The defocus is defined a the maximum defocus path length error divided by the wave length (See Castleman for details). The summation over the Bessel functions in the Hopkins formluation, is stopped when the change is smaller than DIP_MICROSCOPY_HOPKINS_OTF_CUTOFF.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	out	Output
dip_float	defocus	Defocus
dip_float	xNyquist	Oversampling
dip_float	amplitude	Amplitude
dipf_IncoherentOTF	otf	Otf approximation

The dipf_IncoherentOTF enumeration supports the following flags:

Name	Description
DIP_MICROSCOPY_OTF_STOKSETH	Stokseth OTF approximation
DIP_MICROSCOPY_OTF_HOPKINS	Hopkins OTF approximation

LITERATURE

K.R. Castleman, "Digital image processing, second edition", Prentice Hall, Englewood Cliffs, 1996.

SEE ALSO

IncoherentPSF

IncoherentPSF

Generates an incoherent PSF

SYNOPSIS

```
#include "dip_microscopy.h"
dip_Error dip_IncoherentPSF ( output, xNyquist, amplitude )
```

DATA TYPES

Output: sfloat

FUNCTION

This function generates an incoherent in-focus point spread function of a diffraction limited objective.

ARGUMENTS

Data type	Name	Description	
$\mathtt{dip}_{\mathtt{-}}\mathtt{Image}$	output	Output Image	
dip_float	xNyquist	Oversampling Factor	
dip_float	amplitude	Amplitude	

LITERATURE

K.R. Castleman, "Digital image processing, second edition", Prentice Hall, Englewood Cliffs, 1996.

SEE ALSO

IncoherentOTF

IndexToCoordinate

Convert pixel index to coordinate

SYNOPSIS

```
#include "dip_coordsindx.h"
dip_Error dip_IndexToCoordinate ( index, coordinate, stride )
```

FUNCTION

This function is identical to IndexToCoordinateWithSingletons, but does not handle images with singleton dimensions (dimensions where the size is 1). Please use the other function instead, this one is provided for backwards compatability only.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{int}$	index	lineair index
dip_IntegerArray	coordinate	output coordinates
dip_IntegerArray	stride	stride array

SEE ALSO

 ${\tt IndexToCoordinateWithSingletons, CoordinateToIndex}$

${\tt IndexToCoordinateWithSingletons}$

Convert pixel index to coordinate

SYNOPSIS

```
#include "dip_coordsindx.h"
dip_Error dip_IndexToCoordinateWithSingletons ( index, coordinate, size, stride )
```

FUNCTION

This function converts an pixel index of an image to a coordinate array. The conversion is done by calculating the modulus of the index with the stride and size arrays obtained from the image. coordinate has to be an allocated integer array with its size equal to the size of stride and size.

A set of macros can be used instead of this function to avoid some overhead when repeatedly converting linear indices to coordinates for the same image:

DIP_INDEX_TO_COORDINATE(index, coordinates, stride, ix);

```
DIP_FNR_DECLARE; /* Declares dip_Resources rg */

dip_Image image;
dip_int index;
dip_IntegerArray coordinates;

dip_IntegerArray size;
dip_IntegerArray stride;

DIP_INDEX_TO_COORDINATE_DECL( ix ); /* This macro declares variable "ix", name it whatever

DIP_FNR_INITIALISE;

/* ... */

DIPXJ( dip_ImageGetDimensions( image, &size, rg ));
DIPXJ( dip_ImageGetStride( image, &stride, rg ));

DIP_INDEX_TO_COORDINATE_INIT( size, stride, ix, rg ); /* This macro initialises variable "ix"

DIPXJ( dip_ImageGetStride( image, &stride, ix, rg ); /* This macro initialises variable "ix"

DIPXJ( dip_IntegerArrayNew( &coordinates, stride->size, 0, rg ));

/* Now, every time you need to obtain the coordinates for an index, do: */
```

ARGUMENTS

Data type	Name	Description
dip_int	index	lineair index
dip_IntegerArray	coordinate	output coordinates
dip_IntegerArray	size	image size array
dip_IntegerArray	stride	stride array

SEE ALSO

 ${\tt IndexToCoordinate}, {\tt CoordinateToIndex}$

Initialise

Initialise DIPlib

SYNOPSIS

```
dip_Error dip_Initialise( void )
dip_Error dipio_Initialise( void )
```

FUNCTION

Initialise the DIPlib library. Must be called before using any of the other DIPlib functions. This function can be invoked more than once; all but the first invocation are ignored.

SEE ALSO

Exit

InsertionSort

Sort a block of data

SYNOPSIS

```
#include "dip_sort.h"
dip_Error dip_InsertionSort ( data, size, dataType )
```

FUNCTION

Sorts a block of data (of size size and data type dataType) using the insertion sort algorithm.

ARGUMENTS

Data type	Name	Description
void *	data	Data
dip_int	size	Size
dip_DataType	dataType	Data type. See DIPlib's data types

SEE ALSO

General information about sorting

InsertionSortIndices, InsertionSortIndices16, Sort, ImageSort, SortIndices,
SortIndices16, ImageSortIndices

InsertionSortIndices

Sort indices to a block of data

SYNOPSIS

```
#include "dip_sort.h"
dip_Error dip_InsertionSortIndices ( data, indices, size, dataType )
```

FUNCTION

Sorts a list of indices rather than the data itself using the insertion sort algorithm.

ARGUMENTS

Data type	Name	Description
void *	data	Data
dip_sint32 *	indices	Indices
dip_int	size	Size
dip_DataType	dataType	Data type, See DIPlib's data types

SEE ALSO

General information about sorting

InsertionSort, InsertionSortIndices16, Sort, ImageSort, SortIndices, SortIndices16,
ImageSortIndices

InsertionSortIndices16

Sort indices to a block of data

SYNOPSIS

```
#include "dip_sort.h"
dip_Error dip_InsertionSortIndices16 ( data, indices, size, dataType )
```

FUNCTION

Sorts a list of (16 bit) indices rather than the data itself using the insertion sort algorithm.

ARGUMENTS

Data type	Name	Description
void *	data	Data
dip_sint16 *	indices	Indices
dip_int	size	Size
dip_DataType	dataType	Data type. See DIPlib's data types

SEE ALSO

General information about sorting

InsertionSort, InsertionSortIndices, Sort, ImageSort, SortIndices, SortIndices16,
ImageSortIndices

IntegerArrayCopy

Copy an array

SYNOPSIS

dip_Error dip_IntegerArrayCopy (dest, src, resources)

FUNCTION

This function copies the integer array src to dest. The array dest is created by this function as well.

ARGUMENTS

Data type	Name	Description
dip_IntegerArray *	dest	Destination array
dip_IntegerArray	src	Source array
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

IntegerArrayNew, IntegerArrayFree, IntegerArrayCopy, IntegerArrayFind
IntegerArrayCopy, FloatArrayCopy, ComplexArrayCopy, DataTypeArrayCopy, BooleanArrayCopy,
VoidPointerArrayCopy, StringArrayCopy

IntegerArrayFind

Find value in array

SYNOPSIS

dip_Error dip_IntegerArrayFind (array, value, index, found)

FUNCTION

Finds a value in an array and "returns" its index in the array. If found is zero, IntegerArrayFind will produce an error if value is not found, otherwise found obtains the search result (DIP_FALSE if value is not found).

ARGUMENTS

Data type	Name	Description	
dip_IntegerArray	array	Array to find value in	
dip_int	value	Value to find	
dip_int *	index	Index of the found value	
dip_Boolean *	found	Value found or not	

SEE ALSO

IntegerArrayNew, IntegerArrayFree, IntegerArrayCopy, IntegerArrayFind
IntegerArrayFind, FloatArrayFind, ComplexArrayFind, DataTypeArrayFind, BooleanArrayFind,
VoidPointerArrayFind

IntegerArrayFree

Array free function

SYNOPSIS

dip_Error dip_IntegerArrayFree (array)

FUNCTION

This function frees *array, and sets array to zero.

ARGUMENTS

Data type	Name	Description
dip_IntegerArray *	array	Array

SEE ALSO

 ${\tt IntegerArrayNew,\ IntegerArrayFree,\ IntegerArrayCopy,\ IntegerArrayFind}$

ArrayFree, IntegerArrayFree, FloatArrayFree, ComplexArrayFree, BoundaryArrayFree, FrameWorkProcessArrayFree, DataTypeArrayFree, ImageArrayFree, BooleanArrayFree, VoidPointerArrayFree, StringArrayFree, CoordinateArrayFree

IntegerArrayNew

Array allocation function

SYNOPSIS

dip_Error dip_IntegerArrayNew (array, size, value, resources)

FUNCTION

This function allocates the size elements of a dip_IntegerArray and sets the size of the array to size. Each array element is initialized with value.

ARGUMENTS

Data type	Name	Description
dip_IntegerArray *	array	Array
dip_int	size	Size
dip_int	value	Initial value
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

IntegerArrayNew, IntegerArrayFree, IntegerArrayCopy, IntegerArrayFind
ArrayNew, IntegerArrayNew, FloatArrayNew, ComplexArrayNew, BoundaryArrayNew,
FrameWorkProcessArrayNew, DataTypeArrayNew, ImageArrayNew, BooleanArrayNew,
VoidPointerArrayNew, StringArrayNew, CoordinateArrayNew

Invert

logic operation

SYNOPSIS

dip_Error dip_Invert (in, out)

DATA TYPES

binary, integer

FUNCTION

The function Invert inverts the pixel value in in1 and stores the result in out.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	in	Binary input image
$dip_{-}Image$	out	Output image

SEE ALSO

And, Xor, Or

IsodataThreshold

Point operation

SYNOPSIS

#include "dip_analysis.h"
dip_Error dip_IsodataThreshold (in, out, mask, numbthresholds, values)

DATA TYPES

integer, float

FUNCTION

Thresholds in with the isodata method. Several threholds can be supplied, their value is returned in values. The different regions are label in out with different grey-values. A make image mask can be given to compute the isodata only there.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input image
$\mathtt{dip}_{-}\mathtt{Image}$	out	Output image
$\mathtt{dip}_{-}\mathtt{Image}$	mask	Mask image
dip_int	numbthresholds	Number of Thresholds
dip_FloatArray	values	Values

SEE ALSO

Threshold, RangeThreshold, HysteresisThreshold

IsScalar

Determines whether an image is a scalar

SYNOPSIS

dip_Error dip_IsScalar(image, answer)

FUNCTION

Determines whether an image is of the DIP_IMTP_SCALAR type. If answer is not zero, the verdict is passed in this variable. Otherwise, dip_IsScalar returns an error in case image fails to be a scalar.

ARGUMENTS

Data type	Name	Description	
$\mathtt{dip}_{\mathtt{I}}\mathtt{Image}$	image	The image under investigation	
dip_Boolean *	answer	The verdict	

Kuwahara

Edge perserving smoothing filter

SYNOPSIS

```
#include "dip_filtering.h"
dip_Error dip_Kuwahara ( in, out, se, boundary, param, shape )
```

DATA TYPES

binary, integer, float

FUNCTION

This function implements the kuwahara edge-preserving smoothing function. See section 9.4, "Smoothing operations", in Fundamentals of Image Processing for a description of the algorithm. However, this function does not implement the classical kuwahara filter, which only compares the variance of four regions in the filter window. Instead, it compares the variance of every region specified by the filter shape and size centered within the filter window.

Only the rectangular, elliptic and diamond filter shapes are supported (DIP_FLT_SHAPE_RECTANGULAR, DIP_FLT_SHAPE_ELLIPTIC and DIP_FLT_SHAPE_DIAMOND). Other filter shapes can be implemented by setting shape to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, and passing a binary image in se. The "on" pixels define the shape of the filter window. Other values of shape are illegal.

If shape is not equal to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, se can be set to zero. When shape is set to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, param is ignored, and can be set to zero.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	in	Input image
$\mathtt{dip}_{-}\mathtt{Image}$	out	Output image
dip_Image	se	Custom filter window (binary)
dip_BoundaryArray	boundary	Boundary conditions
dip_FloatArray	param	Filter sizes
dip_FilterShape	shape	Filter shape

The enumerator dip_FilterShape contains the following constants:

Name	Description
DIP_FLT_SHAPE_DEFAULT	Default filter window, same as
	DIP_FLT_SHAPE_RECTANGULAR
DIP_FLT_SHAPE_RECTANGULAR	Rectangular filter window, can be even in size
DIP_FLT_SHAPE_ELLIPTIC	Elliptic filter window, always odd in size
DIP_FLT_SHAPE_DIAMOND	Diamond-shaped filter window, always odd in size
DIP_FLT_SHAPE_PARABOLIC	Parabolic filter window (morphology only)
DIP_FLT_SHAPE_DISCRETE_LINE	Rotated line structuring element (morphology only)
DIP_FLT_SHAPE_INTERPOLATED_LINE	Rotated line structuring element, through interpolation
	(morphology only)
DIP_FLT_SHAPE_PERIODIC_LINE	(not implemented)
DIP_FLT_SHAPE_STRUCTURING_ELEMENT	Use se as filter window, can be any size

SEE ALSO

 ${\tt Generalised Kuwahara Improved, Generalised Kuwahara Improved, Variance Filter, Uniform}$

KuwaharaImproved

Edge perserving smoothing filter

SYNOPSIS

#include "dip_filtering.h"
dip_Error dip_KuwaharaImproved (in, out, se, boundary, param, shape, threshold)

DATA TYPES

binary, integer, float

FUNCTION

This function implements an improved version of Kuwahara, see that function's description for more information. This function adds a threshold parameter that avoids false edges in uniform regions. If the difference between maximal and minimal variance within the filter window is smaller or equal to threshold, the centre pixel is taken, instead of the minimum. Setting threshold to zero yields the same result as Kuwahara.

Only the rectangular, elliptic and diamond filter shapes are supported (DIP_FLT_SHAPE_RECTANGULAR, DIP_FLT_SHAPE_ELLIPTIC and DIP_FLT_SHAPE_DIAMOND). Other filter shapes can be implemented by setting shape to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, and passing a binary image in se. The "on" pixels define the shape of the filter window. Other values of shape are illegal.

If shape is not equal to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, se can be set to zero. When shape is set to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, param is ignored, and can be set to zero.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	in	Input image
$\mathtt{dip}_{-}\mathtt{Image}$	out	Output image
dip_Image	se	Custom filter window (binary)
dip_BoundaryArray	boundary	Boundary conditions
dip_FloatArray	param	Filter sizes
dip_FilterShape	shape	Filter shape
dip_float	threshold	Minimal variance difference within window

The enumerator dip_FilterShape contains the following constants:

Name	Description
DIP_FLT_SHAPE_DEFAULT	Default filter window, same as
	DIP_FLT_SHAPE_RECTANGULAR
DIP_FLT_SHAPE_RECTANGULAR	Rectangular filter window, can be even in size
DIP_FLT_SHAPE_ELLIPTIC	Elliptic filter window, always odd in size
DIP_FLT_SHAPE_DIAMOND	Diamond-shaped filter window, always odd in size
DIP_FLT_SHAPE_PARABOLIC	Parabolic filter window (morphology only)
DIP_FLT_SHAPE_DISCRETE_LINE	Rotated line structuring element (morphology only)
DIP_FLT_SHAPE_INTERPOLATED_LINE	Rotated line structuring element, through interpolation
	(morphology only)
DIP_FLT_SHAPE_PERIODIC_LINE	(not implemented)
DIP_FLT_SHAPE_STRUCTURING_ELEMENT	Use se as filter window, can be any size

SEE ALSO

 ${\tt Kuwahara}, {\tt GeneralisedKuwahara}, {\tt GeneralisedKuwaharaImproved}, {\tt VarianceFilter}, {\tt Uniform}$

Label

Label a binary image

SYNOPSIS

```
#include "dip_regions.h"
dip_Error dip_Label ( in, out, connectivity, flags, minsize, maxsize, nol, boundary
)
```

DATA TYPES

binary

FUNCTION

The output is an integer image. Each object (respecting the connectivity, see The connectivity parameter) in the input image receives a unique number. This number ranges from 1 to the number of objects in the image. The pixels in the output image corresponding to a given object are set to this number (label). The remaining pixels in the output image are set to 0. The minsize and maxsize set limits on the size of the objects, if the flag DIP_LB_THRESHOLD_ON_SIZE is set: Objects smaller than minsize or larger than maxsize do not receive a label and the corresponding pixels in the output image are set to zero. Setting minsize to zero implies that there is no check with respect to the minimum size of the object, and the same holds for maxsize and the maximum size of the object. If the flag DIP_LB_LABEL_IS_SIZE is set, the objects' labels are set to the objects' sizes. The boundary conditions are generally ignored (labeling stops at the boundary). The exception is DIP_BC_PERIODIC, which is the only one that makes sense for this algorithm.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input binary image
dip_Image	out	Output label image
dip_int	connectivity	Connectivity
dip_int	flags	0, or a logical OR of the flags described above
dip_int	minsize	Minimum size of the objects (0=do not check)
dip_int	maxsize	Maximum size of the objects (0=do not check)
dip_int *	nol	Pointer to dip_int. Used for returning the number of
		objects. May be set to 0.
dip_BoundaryArray	boundary	Boundary conditions

Laplace

Second order derivative filter

SYNOPSIS

```
#include "dip_derivatives.h"
dip_Error dip_Laplace ( in, out, boundary, ps, sigmas, tc, flavour )
```

DATA TYPES

Depends on the underlying implementation, but expect: binary, integer, **float**

FUNCTION

Computes the Laplace of an image using the Derivative function.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	out	Output
dip_BoundaryArray	boundary	Boundary conditions
dip_BooleanArray	ps	Dimensions to process
dip_FloatArray	sigmas	Sigma of Gaussian
dip_float	tc	Truncation of Gaussian, see
		GlobalGaussianTruncationGet
dip_DerivativeFlavour	flavour	Derivative flavour

The enumerator flavour parameter is one of:

Name	Description
DIP_DF_DEFAULT	Default derivative flavour (==DIP_DF_FIRGAUSS)
DIP_DF_FIRGAUSS	Gaussian family, FIR implementation, Gauss
DIP_DF_IIRGAUSS	Gaussian family, IIR implementation, GaussIIR
DIP_DF_FTGAUSS	Gaussian family, FT implementation, GaussFT
DIP_DF_FINITEDIFF	Finite difference implementation, FiniteDifferenceEx

SEE ALSO

See section 9.5, "Derivative-based operations", in Fundamentals of Image Processing.

Derivative, GradientMagnitude, GradientDirection2D, Dgg, LaplacePlusDgg, LaplaceMinDgg

LaplaceMinDgg

Second order derivative filter

SYNOPSIS

```
#include "dip_derivatives.h"
dip_Error dip_LaplaceMinDgg ( in, out, boundary, ps, sigmas, tc, flavour )
```

DATA TYPES

Depends on the underlying implementation, but expect: binary, integer, **float**

FUNCTION

Computes Laplace - Dgg. For two-dimensional images this is equivalent to the second order derivative in the direction perpendicular to the gradient direction.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	out	Output
dip_BoundaryArray	boundary	Boundary conditions
dip_BooleanArray	ps	Dimensions to process
dip_FloatArray	sigmas	Sigma of Gaussian
dip_float	tc	Truncation of Gaussian, see
		GlobalGaussianTruncationGet
dip_DerivativeFlavour	flavour	Derivative flavour

The enumerator flavour parameter is one of:

Name	Description
DIP_DF_DEFAULT	Default derivative flavour (==DIP_DF_FIRGAUSS)
DIP_DF_FIRGAUSS	Gaussian family, FIR implementation, Gauss
DIP_DF_IIRGAUSS	Gaussian family, IIR implementation, GaussIIR
DIP_DF_FTGAUSS	Gaussian family, FT implementation, GaussFT
DIP_DF_FINITEDIFF	Finite difference implementation, FiniteDifferenceEx

SEE ALSO

Derivative, GradientMagnitude, GradientDirection2D, Laplace, Dgg, LaplacePlusDgg

LaplacePlusDgg

Second order derivative filter

SYNOPSIS

```
#include "dip_derivatives.h"
dip_Error dip_LaplacePlusDgg ( in, out, boundary, ps, sigmas, tc, flavour )
```

DATA TYPES

Depends on the underlying implementation, but expect: binary, integer, **float**

FUNCTION

Computes the laplace and the second derivative in gradient direction of an image using the Derivative function and adds the results. The zero-crossings of the result correspond to the edges in the image, just as for the individual Laplace and Dgg operators. The localization is improved by an order of magnitude with respect to the individual operators.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	out	Output
dip_BoundaryArray	boundary	Boundary conditions
dip_BooleanArray	ps	Dimensions to process
dip_FloatArray	sigmas	Sigma of Gaussian
dip_float	tc	Truncation of Gaussian, see
		GlobalGaussianTruncationGet
dip_DerivativeFlavour	flavour	Derivative flavour

The enumerator flavour parameter is one of:

Name	Description
DIP_DF_DEFAULT	Default derivative flavour (==DIP_DF_FIRGAUSS)
DIP_DF_FIRGAUSS	Gaussian family, FIR implementation, Gauss
DIP_DF_IIRGAUSS	Gaussian family, IIR implementation, GaussIIR
DIP_DF_FTGAUSS	Gaussian family, FT implementation, GaussFT
DIP_DF_FINITEDIFF	Finite difference implementation, FiniteDifferenceEx

LITERATURE

Lucas J. van Vliet, "Grey-Scale Measurements in Multi-Dimensional Digitized Images", Delft University of Technology, 1993

P.W. Verbeek and L.J. van Vliet, "On the location error of curved edges in low-pass filtered 2-D and 3-D images", IEEE Transactions on Pattern Analysis and Machine Intelligence, vol. 16, no. 7, 1994, 726-733.

SEE ALSO

Derivative, GradientMagnitude, GradientDirection2D, Laplace, Dgg, LaplaceMinDgg

Lee

Morphological edge detector

SYNOPSIS

```
#include "dip_morphology.h"
dip_Error dip_Lee ( in, out, se, boundary, param, shape, edgeType, flags )
```

DATA TYPES

integer, float

FUNCTION

Implements a morphological edge detector based on the minimum of two complementary morphological operations. These can be chosen through the edgeType parameter.

The rectangular, elliptic and diamond structuring elements are "flat", i.e. these structuring elements have a constant value. For these structuring elements, param determines the sizes of the structuring elements.

When shape is DIP_FLT_SHAPE_DISCRETE_LINE or DIP_FLT_SHAPE_INTERPOLATED_LINE, the structuring element is a line. param->array[0] determines the length, param->array[1] the angle. This is currently only supported for 2D images. Interpolated lines use interpolation to obtain a more accurate result, but loose the strict increasingness and extensivity (these properties are satisfied only by approximation).

When shape is set to DIP_FLT_SHAPE_PARABOLIC, params specifies the curvature of the parabola.

When shape is set to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, se is used as structuring element. It can be either a binary or a grey-value image. Its origin is the center, or one pixel to the left of the center if the size is even.

If shape is not equal to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, se can be set to zero. When shape is set to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, param is ignored, and can be set to zero.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	out	Output
dip_Image	se	Custom structuring element
dip_BoundaryArray	boundary	Boundary conditions
dip_FloatArray	param	Filter parameters
dip_FilterShape	shape	Structuring element
dip_MphEdgeType	edgeType	Edge type
dipf_LeeSign	flags	Lee sign flag

The enumerator dip_FilterShape contains the following constants:

Name	Description
DIP_FLT_SHAPE_DEFAULT	Default filter window, same as
	DIP_FLT_SHAPE_RECTANGULAR
DIP_FLT_SHAPE_RECTANGULAR	Rectangular filter window, can be even in size
DIP_FLT_SHAPE_ELLIPTIC	Elliptic filter window, always odd in size
DIP_FLT_SHAPE_DIAMOND	Diamond-shaped filter window, always odd in size
DIP_FLT_SHAPE_PARABOLIC	Parabolic filter window (morphology only)
DIP_FLT_SHAPE_DISCRETE_LINE	Rotated line structuring element (morphology only)
DIP_FLT_SHAPE_INTERPOLATED_LINE	Rotated line structuring element, through interpolation
	(morphology only)
DIP_FLT_SHAPE_PERIODIC_LINE	(not implemented)
DIP_FLT_SHAPE_STRUCTURING_ELEMENT	Use se as filter window, can be any size

The enumerator dip_MphEdgeType contains the following constants:

Name	Description
DIP_MPH_TEXTURE	Response is limited to edges in texture
DIP_MPH_OBJECT	Response is limited to object edges
DIP_MPH_BOTH	All edges produce equal response

The enumerator dipf_LeeSign contains the following constants:

Name	Description
DIP_LEE_UNSIGNED	Absolute edge strength
DIP_LEE_SIGNED	Signed edge strength

SEE ALSO

 ${\tt Morphological Gradient Magnitude,\ Morphological Range,\ MultiScale Morphological Gradient,\ Tophat}$

Lesser

Compare grey values in two images

SYNOPSIS

dip_Error dip_Lesser (in1, in2, out)

DATA TYPES

binary, integer, float

FUNCTION

This function sets each pixel in out to "true" when for corresponding pixels in1 < in2. This is the same as Compare with the DIP_SELECT_LESSER selector flag.

in 2 can be a 0D image for comparison of pixel values with a single scalar value. This leads to a functionality similar to that of Threshold.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{\mathtt{-}}\mathtt{Image}$	in1	First input
$\mathtt{dip}_{-}\mathtt{Image}$	in2	Second input
dip_Image	out	Output

SEE ALSO

Compare, Threshold, Equal, Greater, NotEqual, NotGreater, NotLesser, SelectValue, NotZero

Ln

arithmetic function

SYNOPSIS

dip_Error dip_Ln (in, out)

DATA TYPES

binary, integer, ${f float}$

FUNCTION

Computes the natural logarithm of the input image values.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	out	Output

SEE ALSO

Sqrt, Exp, Exp2, Exp10, Log2, Log10

LnGamma

mathematical function

SYNOPSIS

dip_Error dip_LnGamma (in, out)

DATA TYPES

binary, integer, float

FUNCTION

Computes the natural logarithm of the gamma function of the input image values.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	out	Output

SEE ALSO

BesselJO, BesselJI, BesselJN, BesselYO, BesselYI, BesselYN, Erf, Erfc, Sinc

LnNormError

difference measure

SYNOPSIS

dip_Error dip_LnNormError (in1, in2, mask, out, order)

DATA TYPES

binary, integer, float, complex

FUNCTION

Calculates the order norm difference between each pixel value of in1 and in2. Optionally the mask image can be used to exclude pixels from the calculation by setting the value of these pixels in mask to zero.

ARGUMENTS

Data type	Name	Description	
dip_Image	in1	First input	
dip_Image	in2	Second input	
dip_Image	mask	Mask	
dip_Image	out	Output	
dip_float	order	Order	

SEE ALSO

 ${\tt MeanError}, \, {\tt MeanSquareError}, \, {\tt RootMeanSquareError}, \, {\tt MeanAbsoluteError}, \, {\tt IDivergence}$

LocalMinima

Marks local minima (or regional minima)

SYNOPSIS

```
#include "dip_morphology.h"
dip_Error dip_LocalMinima ( in, mask, out, connectivity, max_depth, max_size,
binaryOutput )
```

DATA TYPES

integer, float

FUNCTION

The binary output image is true on all pixels belonging to the minima of a region (as defined by the watershed). To find local maxima, use the inverse of the image as input to this function (see Invert). If binaryOutput is DIP_FALSE, the output is a labelled image instead of a binary one. In this case, pixels belonging to the same local minimum are assigned the same value.

The algorithm is based on the watershed transform, see Watershed for information on the parameters.

Minima is a different algorithm to obtain local minima; Maxima yields the local maxima.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	in	Input
$dip_{-}Image$	mask	Mask
dip_Image	out	Output (binary)
dip_int	connectivity	Connectivity
dip_float	max_depth	Maximum depth of a region that can be merged
dip_int	max_size	Maximum size of a region that can be merged
dip_Boolean	binaryOutput	DIP_FALSE if the output should be a labelled image

SEE ALSO

Watershed, SeededWatershed, UpperEnvelope, Minima, Maxima

Log10

arithmetic function

SYNOPSIS

dip_Error dip_Log10 (in, out)

DATA TYPES

binary, integer, **float**

FUNCTION

Computes the base ten logarithm of the input image values.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{\mathtt{-}}\mathtt{Image}$	in	Input
$\mathtt{dip}_{-}\mathtt{Image}$	out	Output

SEE ALSO

Sqrt, Exp, Exp2, Exp10, Ln, Log2

Log2

arithmetic function

SYNOPSIS

dip_Error dip_Log2 (in, out)

DATA TYPES

binary, integer, ${f float}$

FUNCTION

Computes the base two logarithm of the input image values.

ARGUMENTS

Data type	Name	Description	
dip_Image	in	Input	
dip_Image	out	Output	

SEE ALSO

Sqrt, Exp, Exp2, Exp10, Ln, Log10

macros.h

Various macros

DESCRIPTION

The include files dip_macros.h contains a number of useful macros.

Math macros

DIP_ABS(x)	Absolute value of x
DIP_MAX(x, y)	Maximum of x and y
DIP_MIN(x, y)	Minimum of x and y
DIP_FUNC(funcName, suffix)	Attaches the suffix to the function name, and
	puts and underscore in between.
DIP_SWAP(x, y, z)	Swaps variables x and y, using temporary
	variable z. Must be followed by a trailing ";"

Macros for handling complex numbers:

DIP_REAL(x)	Real part of complex number x
DIP_IMAGINARY(x)	Imaginary part of complex number x
DIP_SQUARE_MODULUS(x)	Square modulus of complex number x
DIP_MODULUS(x)	Modulus of complex number x
DIP_PHASE(x)	Phase of complex number x

Binary I/O macros

DIP_BINARY_MASK(mask, plane)	Computes a binary mask from the plane value
DIP_BINARY_READ(in, mask)	Returns the binary value from in
DIP_BINARY_WRITE(out, val, mask)	Writes the value of val to out

Random access I/O macros:

```
DIP_PIXEL_GET( ip, pos, stride, value )
DIP_PIXEL_SET( ip, pos, stride, value )
```

get/set the value of the pixel at position pos from data pointer ip with strides stride. Both pos and stride are dip_IntegerArrays.

```
DIP_PIXEL_ADD( ip, pos, stride, value )
DIP_PIXEL_SUB( ip, pos, stride, value )
DIP_PIXEL_MUL( ip, pos, stride, value )
DIP_PIXEL_DIV( ip, pos, stride, value )
```

add/subtract/multlipy/divide the value with the pixel-value at position pos from data pointer ip with strides stride. Both pos and stride are dip_IntegerArrays.

Map

Remaps an image

SYNOPSIS

```
#include "dip_manipulation.h"
dip_Error dip_Map ( in, out, map, mirror )
```

DATA TYPES

binary, integer, float, complex

FUNCTION

This function maps the dimensions of the output image to (different) dimensions of the input image. The array index of map specifies the dimension of the output image, the value of the array element of map specifies to which dimension in the input image it corresponds. Optionally, the dimensions can be mirrored, when the value of the corresponding array element in mirror is set to DIP_TRUE. The mirror operation is performed after the mapping operation.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	in	Input image
dip_Image	out	Output image
dip_IntegerArray	map	Map array
dip_BooleanArray	mirror	Mirror array

SEE ALSO

Mirror

Max

arithmetic function

SYNOPSIS

dip_Error dip_Max (in1, in2, out)

DATA TYPES

binary, integer, float

FUNCTION

This function computes out = max (in1, in2) on a pixel by pixel basis. The data types of the in1 and in2 image may be of different types. See Information about dyadic operations for more information about what the type of the output will be.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	in1	First input
dip_Image	in2	Second input
dip_Image	out	Output

SEE ALSO

MaxFloat, Min, MinFloat

MaxFloat

arithmetic function

SYNOPSIS

dip_Error dip_MaxFloat (in, out, constant)

DATA TYPES

binary, integer, float

FUNCTION

This function computes $\mathtt{out} = \max(\mathtt{in}\ , \mathtt{constant})$ on a pixel by pixel basis. The data types of the $\mathtt{in1}$ image and $\mathtt{constant}$ may be of different types. See Information about dyadic operations for more information about what the type of the output will be.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	in	Input
dip_Image	out	Output
dip_float	constant	Constant

SEE ALSO

Max, Min, MinFloat

Maxima

Detects local maxima

SYNOPSIS

```
#include "dip_filtering.h"
dip_Error dip_Maxima ( in, mask, out, connectivity, booleanOutput )
```

DATA TYPES

integer, float

FUNCTION

This function detects local maxima.

The algorithm finds a connected set of pixels with identical value, an no neighbours with higher value. This set is a local maximum and its pixels are set to 1 in the output image. If booleanOutput is false, the output image is a labelled image.

For images that have large plateaus (regions of constant value) that are not local maxima, this function can be quite slow. For example, an image that is zero everywhere except for a small peak somewhere. For such an image it is recommended to use the mask input, for example with the output of a threshold operation.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input image
$dip_{-}Image$	mask	Mask image
$\mathtt{dip}_{-}\mathtt{Image}$	out	Binary output image
dip_int	connectivity	Connectivity
dip_Boolean	booleanOutput	Give a binary output image?

NOTE

If you are looking for the old version of Maxima, it is still available through the following combination of commands:

```
dip_Dilation( in, out, se, boundary, param, shape );
dip_Equal( in, out, out );
```

SEE ALSO

 ${\tt Minima, Subpixel Maxima, Local Minima, Seeded Watershed, GrowRegions}$

Maximum

statistics function

SYNOPSIS

dip_Error dip_Maximum (in, mask, out, ps)

DATA TYPES

binary, integer, float

FUNCTION

Calculates the maximum of the pixel values over all those dimensions which are specified by ps.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	mask (0)	Mask
$\mathtt{dip}_{-}\mathtt{Image}$	out	Output
dip_BooleanArray	ps (0)	Dimensions to project

SEE ALSO

From images to scalars

 ${\tt Sum,\ Mean,\ Variance,\ Standard Deviation,\ Mean Modulus,\ Sum Modulus,\ Mean Square Modulus,\ Minimum,\ Median,\ Percentile}$

mBesselJ0

mathematical function

SYNOPSIS

dip_float dipm_BesselJO (x)

FUNCTION

Computes the Bessel function J0 of the input value.

ARGUMENTS

Data type	Name	Description	
dip_float	х	Input value	

SEE ALSO

mTruncate, mFraction, mNearestInt, mSign, mExp2, mExp10, mLog2, mSinc, mReciprocal, mBesselJ1, mBesselJN, mBesselY0, mBesselY1, mBesselYN, mLnGamma, mErf, mErfc, mGammaP, mGammaQ

mBesselJ1

mathematical function

SYNOPSIS

dip_float dipm_BesselJ1 (x)

FUNCTION

Computes the Bessel function J1 of the input value.

ARGUMENTS

Data type	Name	Description	
dip_float	х	Input value	

SEE ALSO

mTruncate, mFraction, mNearestInt, mSign, mExp2, mExp10, mLog2, mSinc, mReciprocal, mBesselJ0, mBesselJN, mBesselY0, mBesselY1, mBesselYN, mLnGamma, mErf, mErfc, mGammaP, mGammaQ

mBesselJN

mathematical function

SYNOPSIS

dip_float dipm_BesselJN (x, n)

FUNCTION

Computes the Bessel function J of the order ${\tt n}$ of the input value.

ARGUMENTS

Data type	Name	Description	
dip_float	х	Input value	
dip_int	n	Order of Bessel function	

SEE ALSO

mTruncate, mFraction, mNearestInt, mSign, mExp2, mExp10, mLog2, mSinc, mReciprocal, mBesselJ0, mBesselJ1, mBesselY0, mBesselY1, mBesselYN, mLnGamma, mErf, mErfc, mGammaP, mGammaQ

mBesselY0

mathematical function

SYNOPSIS

dip_float dipm_BesselYO (x)

FUNCTION

Computes the Bessel function Y0 of the input value.

ARGUMENTS

Data type	Name	Description	
dip_float	х	Input value	

SEE ALSO

mTruncate, mFraction, mNearestInt, mSign, mExp2, mExp10, mLog2, mSinc, mReciprocal, mBesselJ0, mBesselJ1, mBesselJN, mBesselY1, mBesselYN, mLnGamma, mErf, mErfc, mGammaP, mGammaQ

mBesselY1

mathematical function

SYNOPSIS

dip_float dipm_BesselY1 (x)

FUNCTION

Computes the Bessel function Y1 of the input value.

ARGUMENTS

Data type	Name	Description
dip_float	х	Input value

SEE ALSO

mTruncate, mFraction, mNearestInt, mSign, mExp2, mExp10, mLog2, mSinc, mReciprocal, mBesselJ0, mBesselJ1, mBesselJN, mBesselY0, mBesselYN, mLnGamma, mErf, mErfc, mGammaP, mGammaQ

mBesselYN

mathematical function

SYNOPSIS

dip_float dipm_BesselYN (x, n)

FUNCTION

Computes the Bessel function Y of the order ${\tt n}$ of the input value.

ARGUMENTS

Data type	Name	Description	
dip_float	х	Input value	
dip_int	n	Order of Bessel function	

SEE ALSO

mTruncate, mFraction, mNearestInt, mSign, mExp2, mExp10, mLog2, mSinc, mReciprocal, mBesselJ0, mBesselJ1, mBesselJN, mBesselY0, mBesselY1, mLnGamma, mErf, mErfc, mGammaP, mGammaQ

Mean

statistics function

SYNOPSIS

dip_Error dip_Mean (in, mask, out, ps)

DATA TYPES

binary, integer, float, complex

FUNCTION

Calculates the mean of the pixel values over all those dimensions which are specified by ps.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	mask (0)	Mask
$\mathtt{dip}_{-}\mathtt{Image}$	out	Output
dip_BooleanArray	ps (0)	Dimensions to project

SEE ALSO

From images to scalars

 ${\tt Sum, Variance, Standard Deviation, Mean Modulus, Sum Modulus, Mean Square Modulus, Maximum, Minimum, Median, Percentile}$

MeanAbsoluteError

difference measure

SYNOPSIS

dip_Error dip_MeanAbsoluteError (in1, in2, mask, out)

DATA TYPES

binary, integer, float, complex

FUNCTION

Calculates the mean absolute error difference between each pixel value of in1 and in2. Optionally the mask image can be used to exclude pixels from the calculation by setting the value of these pixels in mask to zero.

ARGUMENTS

Data type	Name	Description
dip_Image	in1	First input
dip_Image	in2	Second input
dip_Image	mask	Mask
dip_Image	out	Output

SEE ALSO

MeanError, MeanSquareError, RootMeanSquareError, LnNormError, IDivergence

MeanError

difference measure

SYNOPSIS

dip_Error dip_MeanError (in1, in2, mask, out)

DATA TYPES

binary, integer, float, complex

FUNCTION

Calculates the mean error difference between all pixel values of in1 and in2. Optionally the mask image can be used to exclude pixels from the calculation by setting the value of these pixels in mask to zero.

ARGUMENTS

Data type	Name	Description
dip_Image	in1	First input
dip_Image	in2	Second input
dip_Image	mask	Mask
dip_Image	out	Output

SEE ALSO

 ${\tt MeanSquareError,\,RootMeanSquareError,\,MeanAbsoluteError,\,LnNormError,\,IDivergence}$

MeanModulus

statistics function

SYNOPSIS

dip_Error dip_MeanModulus (in, mask, out, ps)

DATA TYPES

binary, integer, float, complex

FUNCTION

Calculates the mean modulus of the pixel values over all those dimensions which are specified by ps.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	in	Input
dip_Image	mask (0)	Mask
dip_Image	out	Output
dip_BooleanArray	ps (0)	Dimensions to project

SEE ALSO

From images to scalars

Sum, Mean, Variance, Standard Deviation, Sum Modulus, Mean Square Modulus, Maximum, Minimum, Median, Percentile

${\tt MeanSquareError}$

difference measure

SYNOPSIS

dip_Error dip_MeanSquareError (in1, in2, mask, out)

DATA TYPES

binary, integer, float, complex

FUNCTION

Calculates the mean square error difference between all pixel values of in1 and in2. Optionally the mask image can be used to exclude pixels from the calculation by setting the value of these pixels in mask to zero.

ARGUMENTS

Data type	Name	Description
dip_Image	in1	First input
dip_Image	in2	Second input
dip_Image	mask	Mask
dip_Image	out	Output

SEE ALSO

 ${\tt MeanError}, {\tt RootMeanSquareError}, {\tt MeanAbsoluteError}, {\tt LnNormError}, {\tt IDivergence}$

MeanSquareModulus

statistics function

SYNOPSIS

dip_Error dip_MeanSquareModulus (in, mask, out, ps)

DATA TYPES

binary, integer, float, complex

FUNCTION

Calculates the mean square modulus of the pixel values over all those dimensions which are specified by ps.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	in	Input
$\mathtt{dip}_{-}\mathtt{Image}$	mask (0)	Mask
dip_Image	out	Output
dip_BooleanArray	ps (0)	Dimensions to project

SEE ALSO

From images to scalars

Sum, Mean, Variance, Standard Deviation, Mean Modulus, Sum Modulus, Maximum, Minimum, Median, Percentile

Measure

Measure object features

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error dip_Measure ( measurement, featureID, featureParams, objectID, objectIm,
intensityIm, connectivity, physDims )
```

DATA TYPES

objectIm: integer

intensityIm: integer, float

FUNCTION

The Measure function is the top-level function of DIPlib's measurement library. This function performs measurements of the objects in the specified objectIm image. The measurements to be performed are specified by the featureID array of measurement function IDs. If featureParams is non-zero, its size should equal that of featureID. Although the current implementation of Measure does not make use of this argument, future versions will pass the data pointers of the featureParams to the corresponding measurement functions. featureParams should be set to zero for now.

The list of object IDs on which the measurements have to be performed is specified by objectID. If it is zero, Measure will call GetObjectLabels to obtain a list of all non-zero values in objectIm. The objectID values should be unequal to zero.

The state of measurement should be raw (see MeasurementNew), since Measure will forge the measurement data structure by calling MeasurementForge.

The intensityIm image defines the pixel intensity of the objects, whose shape is defined by objectIm. If none of the measurements specified in featureID require the grey-value image, it can be set to NULL.

The physDims parameter defines the physical dimensions of the pixel sizes and pixel intensity. See PhysicalDimensionsNew for more information.

ARGUMENTS

Data type	Name	Description
dip_Measurement	measurement	Measurement data structure
dip_IntegerArray	featureID	Array of measurement function IDs
dip_VoidPointerArray	featureParams (0)	Set to zero
dip_IntegerArray	objectID (0)	Array of Object IDs
$dip_{-}Image$	objectIm	Image containing object IDs, i.e. object labels
dip_Image	intensityIm	Intensity image
dip_int	connectivity	Connectivity of object's contour pixels, see
		The connectivity parameter
dip_PhysicalDimensions	physDims	Structure specifying the physical dimensions
		of the image pixels

SEE ALSO

Label, ObjectToMeasurement, MeasurementToImage, MeasurementToHistogram, MeasurementWrite, MeasurementNew, MeasurementFree, MeasurementIsValid

 ${\tt Feature Anisotropy 2D, Feature Bending Energy, Feature Center},$

FeatureChainCodeBendingEnergy, FeatureConvexArea, FeatureConvexPerimeter,

FeatureConvexity, FeatureDimension, FeatureExcessKurtosis, FeatureFeret,

 ${\tt Feature Ginertia}, \ {\tt Feature Gmu}, \ {\tt Feature Gravity}, \ {\tt Feature Inertia}, \ {\tt Feature Longest Chain code Run}, \\$

FeatureMass, FeatureMaxVal, FeatureMaximum, FeatureMean, FeatureMinVal, FeatureMinimum,

FeatureMu, FeatureOrientation2D, FeatureP2A, FeaturePerimeter, FeatureRadius,

FeatureShape, FeatureSize, FeatureSkewness, FeatureStdDev, FeatureSum,

FeatureSurfaceArea

MeasurementFeatureConvert

Convert the data of a measurement feature

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error dip_MeasurementFeatureConvert ( in, featureID, inID, out, outID, resources
)
```

FUNCTION

This function convert the data of object inID in measurement in measured by feature featureID to object outID in out.

ARGUMENTS

Data type	Name	Description	
dip_Measurement	in	Input measurement data structure	
dip_int	featureID	ID of the measurement feature	
dip_int	inID	ID of the object in in	
dip_Measurement	out	Output measurement data structure	
dip_int	outID	ID of the object in out	
dip_Resources	resources	Resources tracking structure. See ResourcesNew	

SEE ALSO

Measure, MeasurementNew

MeasurementFeatureDescription

Measurement Description access function

SYNOPSIS

#include "dip_measurement.h"
dip_Error dip_MeasurementFeatureDescription (measurement, featureID, description,
resources)

FUNCTION

The MeasurementObjectData, MeasurementObjectValue and MeasurementFeatureDescription functions provide access to the functions that are registered by each measurement function. See also MeasurementFeatureRegister.

This function gives access to a structure containing the name, a short description of the measurement feature, as well as the labels and units of the data measured by the feature specified with featureID. Use the functions FeatureDescriptionGetName,

FeatureDescriptionGetDescription, FeatureDescriptionGetLabels and FeatureDescriptionGetUnits to access the values in the description structure.

ARGUMENTS

Data type	Name	Description
dip_Measurement	measurement	Measurement data structure
dip_int	featureID	Measurement feature ID
dip_FeatureDescription *	description	Pointer to a dip_FeatureDescription structure
		containing a description of the specified feature
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

MeasurementNew, MeasurementFeatures, MeasurementFeatureValid,
MeasurementObjectData, MeasurementObjectValue, MeasurementFeatureRegister,
MeasurementFeatureRegistryList, MeasurementFeatureRegistryGet,
MeasurementFeatureRegistryFeatureDescription,
MeasurementFeatureRegistryFeatureNeedsIntensityImage FeatureDescriptionNew,
FeatureDescriptionFree, FeatureDescriptionSetName, FeatureDescriptionGetName,
FeatureDescriptionSetDescription, FeatureDescriptionGetDescription,
FeatureDescriptionSetLabels, FeatureDescriptionGetLabels, FeatureDescriptionSetLabel,
FeatureDescriptionSetDimensionLabels, FeatureDescriptionSetUnits,
FeatureDescriptionGetUnits, FeatureDescriptionSetUnit

MeasurementFeatureFormat

Feature data format convenience function

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error dip_MeasurementFeatureFormat ( measurement, featureID, format )
```

FUNCTION

This function is a convenience function on top of MeasurementObjectValue, providing an easy access to the data format of the measurement values of the featureID measurement function.

ARGUMENTS

Data type	Name	Description
dip_Measurement	measurement	Measurement data structure
dip_int	featureID	Measurement function ID
dipf_MeasurementValueFormat *	format	Pointer to measurement value data format

SEE ALSO

MeasurementObjectValue, MeasurementNew

MeasurementFeatureRegister

Register a measurement function

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error dip_MeasurementFeatureRegister ( registry )
```

FUNCTION

This function registers a measurement function, specified by registry. Once a function is registered, it can be used through the Measure function by specifying registry.id.rtid as the measurement ID. registry contains pointers to a series of functions related to making the measurement, and contains information on how these functions should be called. See below for more information.

ARGUMENTS

Data type	Name	Description
dip_MeasurementFeatureRegistry	registry	Registry

THE REGISTRY STRUCTURE

The dip_MeasurementFeatureRegistry structure contains the following fields:

Data type	Name	Description
dip_Identifier	id	Unique identifier
dipf_FeatureMeasureFunction	type	Type of function measure points to
dip_FeatureCreateFunction	create	Function pointer, see
		FeatureCreateFunction
dip_FeatureComposeFunction	compose	Function pointer, see
		FeatureComposeFunction
dip_FeatureMeasureFunction	measure	Union of function pointers
dip_FeatureValueFunction	value	Function pointer, see
		FeatureValueFunction
dip_FeatureDescriptionFunction	description	Function pointer, see
		FeatureDescriptionFunction
dip_FeatureConvertFunction	convert	Function pointer, see
		FeatureConvertFunction
dip_int	iterations	Currently ignored (set to 1)
dip_Boolean	needIntensityIm	Whether or not a grey-value image is
		needed

dip_Identifier is a struct with two values: uuid and rtid. rtid is of type dip_int, and needs to be set to a unique number (use GetUniqueNumber for that). uuid is curently ignored, but should be set to a universally unique number by using the time, date and processor ID at the time of writing

the code. The UNIX command uuidgen should be used for this.

measure points to the main measuring function, and can be of four different types, based on how it does the measuring. measure is a union with the following fields:

Data type	Name	Description
dip_FeatureLineFunction	line	Takes one image line at the time, see
		FeatureLineFunction
dip_FeatureImageFunction	image	Takes the whole image at once, see
		FeatureImageFunction
dip_FeatureChainCodeFunction	chaincode	Takes one chain code at the time, see
		FeatureChainCodeFunction
dip_FeatureConvHullFunction	convhull	Takes one convex polygon at the time, see
		FeatureConvHullFunction
dip_FeatureCompositeFunction	composite	Combines the results of various other
		measurements, see FeatureCompositeFunction

The type flag should match the function type pointed to, and can be one of the following:

Name	Description
DIP_MSR_FUNCTION_LINE_BASED	measure.line is set
DIP_MSR_FUNCTION_IMAGE_BASED	measure.image is set
DIP_MSR_FUNCTION_CHAINCODE_BASED	measure.chaincode is set
DIP_MSR_FUNCTION_CONVHULL_BASED	measure.convhull is set
DIP_MSR_FUNCTION_COMPOSITE	measure.composite is set

create points to a function that allocates and initialises any data before the measurement can start. value points to a function that returns the measurement result (called by MeasurementObjectValue). description points to a function that returns information on the measurement performed (called by MeasurementFeatureDescription). convert points to a function that copies the collected measurement data to a second measurement object (called MeasurementFeatureConvert). Finally, the compose element points to a function that returns the list of measurement IDs that the DIP_MSR_FUNCTION_COMPOSITE function depends on. This value is ignored for other types of measurement functions.

needIntensityIm should be set if the measurement function expects a grey-value input as well as the labeled image.

SEE ALSO

Measure, MeasurementFeatureRegistryList, MeasurementFeatureRegistryGet,
MeasurementFeatureRegistryFeatureDescription,
MeasurementFeatureRegistryFeatureNeedsIntensityImage, MeasurementFeatureValid,
MeasurementFeatureDescription, MeasurementObjectData, MeasurementObjectValue,
MeasurementFeatureConvert, FeatureLineFunction, FeatureImageFunction,
FeatureChainCodeFunction, FeatureConvHullFunction, FeatureCompositeFunction,
FeatureCreateFunction, FeatureComposeFunction, FeatureValueFunction,
FeatureConvertFunction, FeatureDescriptionFunction

${\tt MeasurementFeatureRegistryFeatureDescription}$

Get the feature description of a registered measurement feature

SYNOPSIS

#include "dip_measurement.h"
dip_Error dip_MeasurementFeatureRegistryFeatureDescription (featureID, description,
resources)

FUNCTION

This function obtains the feature description information of the measurement feature specified by feature ID.

ARGUMENTS

Data type	Name	Description
dip_int	featureID	Measurement feature ID
dip_FeatureDescription *	description	pointer to a dip_FeatureDescription structure
		containing descriptive information of the
		measurement feature. This data can be accessed
		with MeasurementFeatureDescription
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

Measure, MeasurementFeatureRegister, MeasurementFeatureRegistryList,
MeasurementFeatureRegistryGet, MeasurementFeatureRegistryFeatureNeedsIntensityImage,
MeasurementFeatureValid, MeasurementFeatureDescription, MeasurementObjectData,
MeasurementObjectValue

${\tt MeasurementFeatureRegistryFeatureNeedsIntensityImage}$

Checks whether the measurement function needs an intensity image

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error dip_MeasurementFeatureRegistryFeatureNeedsIntensityImage ( featureID,
    veredict )
```

FUNCTION

This function sets veredict to DIP_TRUE if the measurement feature specified by featureID requires a grey-value image, or DIP_FALSE otherwise.

ARGUMENTS

Data type	Name	Description
dip_int	featureID	Measurement feature ID
dip_Boolean *	veredict	Return value

SEE ALSO

 $\label{thm:measurement} Measurement Feature Register, Measurement Feature Registry List, \\ Measurement Feature Registry Get, Measurement Feature Registry Feature Description, \\ Measurement Feature Valid, Measurement Feature Description, Measurement Object Data, \\ Measurement Object Value$

MeasurementFeatureRegistryGet

Get the registry information of a measurement feature

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error dip_MeasurementFeatureRegistryGet ( featureID, registry )
```

FUNCTION

This function obtains (a copy of) the registry structure of the measurement feature function specified by featureID.

ARGUMENTS

Data type	Name	Description
dip_int	featureID	Measurement function ID
dip_MeasurementFeatureRegistry *	registry	Pointer to a measurement feature registry
		structure

SEE ALSO

Measure, MeasurementFeatureRegister, MeasurementFeatureRegistryList,
MeasurementFeatureRegistryFeatureDescription,
MeasurementFeatureRegistryFeatureNeedsIntensityImage, MeasurementFeatureValid,
MeasurementFeatureDescription, MeasurementObjectData, MeasurementObjectValue

MeasurementFeatureRegistryList

Obtain a list of the registered measurement features

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error dip_MeasurementFeatureRegistryList ( featureID, resources )
```

FUNCTION

This functions obtains an array of registered measurement feature IDs.

ARGUMENTS

Data type	Name	Description	
dip_IntegerArray *	featureID	Pointer to an array of measurement feature IDs	
dip_Resources	resources	Resources tracking structure. See ResourcesNew	

SEE ALSO

Measure, MeasurementFeatureRegister, MeasurementFeatureRegistryGet,
MeasurementFeatureRegistryFeatureDescription,
MeasurementFeatureRegistryFeatureNeedsIntensityImage, MeasurementFeatureValid,
MeasurementFeatureDescription, MeasurementObjectData, MeasurementObjectValue

MeasurementFeatures

Get the measurement ID array

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error dip_MeasurementFeatures ( measurement, featureID, resources )
```

FUNCTION

This function obtains an array of measurement function IDs in the measurement structure. See MeasurementForge for a (brief) explination of the measurement data structure.

ARGUMENTS

Data type	Name	Description
dip_Measurement	measurement	Measurement data structure
dip_IntegerArray *	featureID	pointer to an array of measurement function IDs
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

 ${\tt Measure}, {\tt MeasurementNew}, {\tt MeasurementNumberOfFeatures}, {\tt MeasurementFeatureValid}, {\tt MeasurementFeatureDescription}$

MeasurementFeatureSize

Feature data convenience function

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error dip_MeasurementFeatureSize ( measurement, featureID, size )
```

FUNCTION

This function is a convenience function on top of MeasurementObjectValue, providing an easy access to the number of the measurement values of the featureID measurement function.

ARGUMENTS

Data type	Name	Description
dip_Measurement	measurement	Measurement data structure
dip_int	featureID	ID of the measurement feature
dip_int *	size	Number of measurement values

SEE ALSO

MeasurementObjectValue, MeasurementNew

MeasurementFeatureValid

Verify a measurement feature ID

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error dip_MeasurementFeatureValid ( measurement, featureID, verdict )
```

FUNCTION

This function determines whether featureID is a valid measurement feature, by verifying whether featureID equals the ID of one of the registered measurement features. If verdict is not zero, the result (DIP_TRUE or DIP_FALSE) is stored in verdict, otherwise an error is returned in case the verification fails.

ARGUMENTS

Data type	Name	Description
dip_Measurement	measurement	Measurement data structure
dip_int	featureID	Measurement feature ID to validated
dip_Boolean *	verdict	Pointer to the boolean verdict

SEE ALSO

Measure, MeasurementNew, MeasurementFeatures, MeasurementFeatureDescription, MeasurementFeatureRegister, MeasurementFeatureRegistryList, MeasurementFeatureRegistryGet, MeasurementFeatureRegistryFeatureDescription, MeasurementFeatureRegistryFeatureNeedsIntensityImage

MeasurementForge

Allocate the data of a measurement data structure

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error dip_MeasurementForge ( measurement, featureID, objectID )
```

FUNCTION

This function forges a measurement data structure, that has been created with MeasurementNew. The featureID array should contain the IDs of the features to be performed. The vadility of these IDs is checked by comparing them with the IDs of registered measurement functions (see MeasurementFeatureRegister). The objectID array contains the IDs (i.e. labels) of the objects on which the features are to be performed. (For example, the Measure function accepts as one of its arguments a label image, of which the intensity of each individual pixel represents the ID of the object to which that pixel belongs. These label values should in that case correspond to the values of objectID.)

The measurement structure can be regarded as a matrix spanned by the number of features along one axis, and the number of objects along the other. MeasurementForge allocates and initialises the internal structures to contain this matrix and the data required for each conbimation of measurement and object ID.

ARGUMENTS

Data type	Name	Description
dip_Measurement	measurement	Measurement structure
dip_IntegerArray	featureID	Array of measurement function IDs
dip_IntegerArray	objectID	Array of Object IDs

SEE ALSO

Measure, MeasurementNew, MeasurementFree, MeasurementIsValid,
MeasurementFeatureRegister, MeasurementFeatureRegistryList,
MeasurementFeatureRegistryGet, MeasurementFeatureRegistryFeatureDescription,
MeasurementFeatureRegistryFeatureNeedsIntensityImage

MeasurementFree

Free a measurement data structure

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error dip_MeasurementFree ( measurement )
```

FUNCTION

This function frees a Measurement data structure. After the Measurement has been freed, the pointer measurement is set to zero.

ARGUMENTS

Data type	Name	Description
dip_Measurement *	measurement	pointer to the measurement structure to be freed

SEE ALSO

Measure, MeasurementNew, MeasurementForge, MeasurementIsValid

MeasurementGetName

Get the name of a Measurement structure

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error dip_MeasurementGetName ( measurement, name, resources )
```

FUNCTION

This function gets the name of a measurement structure

ARGUMENTS

Data type	Name	Description
dip_Measurement	measurement	Measurement
dip_String *	name	Name
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

 ${\tt Measure,\,MeasurementNew,\,MeasurementID,\,MeasurementSetName,}\\ {\tt MeasurementGetPhysicalDimensions,\,MeasurementSetPhysicalDimensions}\\$

${\tt MeasurementGetPhysicalDimensions}$

Get the physical dimensions info of a measurement

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error dip_MeasurementGetPhysicalDimensions ( measurement, physDims, resources )
```

FUNCTION

This function obtains a copy of the physical dimensions information associated with the measurement data structure. The physical dimensions data structure informs measurement features about the physical sizes and position of the pixels of the measured image.

ARGUMENTS

Data type	Name	Description
dip_Measurement	measurement	Measurement data structure
dip_PhysicalDimensions *	physDims	Pointer to a Physical Dimensions data structure
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

 ${\tt MeasurementNew,\,MeasurementID,\,MeasurementSetName,\,MeasurementGetName,\,MeasurementSetPhysicalDimensions}$

MeasurementID

Get the ID of a Measurement structure

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error dip_MeasurementID ( measurement, id )
```

FUNCTION

This function obtains the ID of a the measurement structure. The ID is a DIPlib wide unique number (see GetUniqueNumber).

ARGUMENTS

Data type	Name	Description
dip_Measurement	measurement	Measurement structure
dip_int *	id	Pointer to the id

SEE ALSO

 ${\tt Measure,\,MeasurementNew,\,MeasurementSetName,\,MeasurementGetPhysicalDimensions,\,MeasurementSetPhysicalDimensions}$

MeasurementIsValid

Checks whether a measurement is valid

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error dip_MeasurementIsValid ( measurement, verdict )
```

FUNCTION

This function determines whether measurement is forged. If verdict is not zero, the result (DIP_TRUE or DIP_FALSE) is stored in verdict, otherwise an error is returned in case the verification fails.

ARGUMENTS

Data type	Name	Description
dip_Measurement	measurement	Measurement data structure
dip_Boolean *	verdict	The validation verdict

SEE ALSO

Measure, MeasurementNew, MeasurementFree, MeasurementForge

MeasurementNew

Create new measurement data structure

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error dip_MeasurementNew ( measurement, resources )
```

FUNCTION

This function creates, by allocating and initialising it, a new Measurement data structure. After this function has been used to create a new measurement structure, the state of it is raw. It needs to be passed through MeasurementForge before it can be used to store measurement results.

ARGUMENTS

Data type	Name	Description
dip_Measurement *	measurement	pointer to the measurement structure to be created
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

MeasurementFree, MeasurementForge, MeasurementIsValid, MeasurementID, MeasurementSetName, MeasurementGetName, MeasurementGetPhysicalDimensions, MeasurementSetPhysicalDimensions, MeasurementNumberOfFeatures, MeasurementFeatureValid, MeasurementFeatureDescription, MeasurementNumberOfObjects, MeasurementObjects, MeasurementObjectValid, MeasurementObjectData, MeasurementObjectValue, MeasurementFeatureConvert, MeasurementWrite

MeasurementNumberOfFeatures

Get the number of measurement feature IDs

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error dip_MeasurementNumberOfFeatures ( measurement, features )
```

FUNCTION

This function obtains the number of measurement feature IDs in the measurement structure. See MeasurementForge for a (brief) explination of the measurement data structure.

ARGUMENTS

Data type	Name	Description
dip_Measurement	measurement	Measurement data structure
dip_int *	features	pointer to the number of measurement feature IDs

SEE ALSO

Measure, MeasurementNew, MeasurementFeatures

MeasurementNumberOfObjects

Get the number of object IDs

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error dip_MeasurementNumberOfObjects ( measurement, objects )
```

FUNCTION

This function obtains the number of object IDs belonging to the featureID measurement function ID in the measurement structure. See MeasurementForge for a (brief) explination of the measurement data structure.

ARGUMENTS

Data type	Name	Description
dip_Measurement	measurement	Measurement data structure
dip_int *	objects	Pointer to an integer containing the number of object IDs

SEE ALSO

Measure, MeasurementNew, MeasurementObjects

MeasurementObjectData

Object data access function

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error dip_MeasurementObjectData ( measurement, featureID, objectID, data,
verdict )
```

FUNCTION

The MeasurementObjectData, MeasurementObjectValue and MeasurementFeatureDescription functions provide access to the functions that are registered by each measurement function. See also MeasurementFeatureRegister.

The Object data is the data allocated by a measurement function for internal purposes, for example to store intermediate results. Its format is free. Therefore, the use of this function is only meaningful for a particular measurement function itself. To access the measurement values of a measurement function, use MeasurementObjectValue.

The verdict parameter provides a means to test whether featureID or objectID are valid within the context of measurement. If one of them is invalid, and verdict is not zero, *verdict is set to DIP_FALSE, otherwise its value is DIP_TRUE. If verdict is zero, MeasurementObjectData produces an error when either featureID or objectID is invalid.

ARGUMENTS

Data type	Name	Description
dip_Measurement	measurement	Measurement data structure
dip_int	featureID	Measurement function ID
dip_int	objectID	Object ID
void **	data	Pointer to the internal measurement data pointer
dip_Boolean *	verdict	Pointer to a boolean containing validation information

SEE ALSO

Measure, MeasurementNew, MeasurementObjects, MeasurementObjectData,
MeasurementObjectValue, MeasurementFeatureDescription, MeasurementFeatureRegister,
MeasurementFeatureRegistryList, MeasurementFeatureRegistryGet,
MeasurementFeatureRegistryFeatureDescription,
MeasurementFeatureRegistryFeatureNeedsIntensityImage

MeasurementObjects

Get an object ID array

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error dip_MeasurementObjects ( measurement, featureID, objectID, resources )
```

FUNCTION

This function obtains an array of object IDs belonging to the featureID measurement function in the measurement structure. See MeasurementForge for a (brief) explination of the measurement data structure.

ARGUMENTS

Data type	Name	Description
dip_Measurement	measurement	Measurement data structure
dip_int	featureID	Measurement function ID
dip_IntegerArray *	objectID	Pointer to an object ID array
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

Measure, MeasurementNew, MeasurementNumberOfObjects, MeasurementObjectValid, MeasurementObjectData, MeasurementObjectValue

MeasurementObjectValid

Verify an object ID

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error dip_MeasurementObjectValid ( measurement, featureID, objectID, verdict )
```

FUNCTION

This function determines whether the object ID objectID, belonging to the measurement function ID featureID, is a valid object ID, by comparing objectID to the object IDs belonging to the featureID in measurement. If verdict is not zero, the result (DIP_TRUE or DIP_FALSE) is stored in verdict, otherwise an error is returned in case the verification fails.

ARGUMENTS

Data type	Name	Description
dip_Measurement	measurement	Measurement data structure
dip_int	featureID	Measurement function ID
dip_int	objectID	Object ID
dip_Boolean *	verdict	Pointer to a boolean containing the validation verdict

SEE ALSO

Measure, MeasurementNew, MeasurementObjects, MeasurementObjectData, MeasurementObjectValue

MeasurementObjectValue

Object value access function

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error dip_MeasurementObjectValue ( measurement, featureID, objectID, data,
format, resources )
```

FUNCTION

The MeasurementObjectData, MeasurementObjectValue and MeasurementFeatureDescription functions provide access to the functions that are registered by each measurement function. See also MeasurementFeatureRegister.

The MeasurementObjectValue function provides access to the measurement values produced by the featureID measurement function measured on the objectID labeled object. The format of data is specified by format.

ARGUMENTS

Data type	Name	Description
dip_Measurement	measurement	Measurement data structure
dip_int	featureID	Measurement function ID
dip_int	objectID	Object ID
void **	data	Pointer to data pointer
dipf_MeasurementValueFormat *	format	Pointer to the data format label
dip_Resources	resources	Resources tracking structure. See
		ResourcesNew

Measurement data formats

Name	Description
DIP_MSR_VALUE_FORMAT_INTEGER	Integer scalar data format
DIP_MSR_VALUE_FORMAT_FLOAT	Float scalar data format
DIP_MSR_VALUE_FORMAT_INTEGER_ARRAY	Integer array data format
DIP_MSR_VALUE_FORMAT_FLOAT_ARRAY	Float array data format
DIP_MSR_VALUE_FORMAT_IMAGE	Data is formatted as an dip_Image

SEE ALSO

Measure, MeasurementNew, MeasurementObjects, MeasurementObjectData,
MeasurementObjectValue, MeasurementFeatureDescription, MeasurementFeatureFormat,
MeasurementFeatureSize, MeasurementFeatureRegister, MeasurementFeatureRegistryList,
MeasurementFeatureRegistryGet, MeasurementFeatureRegistryFeatureDescription,

 ${\tt MeasurementFeatureRegistryFeatureNeedsIntensityImage}$

485

DIP*lib* function reference

MeasurementRead

Read measurement results from a file

SYNOPSIS

#include "dipio_measurement.h"
dip_Error dipio_MeasurementRead (measurement, filename, format, addExtensions, recognised)

FUNCTION

This function reads measurement data from a file and puts it in measurement. measurement must be allocated before calling this function. If format is 0, all different MeasurementRead functions are called in sequence until the correct format has been found. If you know the format, get the correct format ID through the registry functions.

The boolean addExtensions specifies whether MeasurementRead should try to add file format extensions to filename, if the registered file format reader fails to recognise filename straight away. The extensions are provided by the registered file readers.

If recognised is not zero, MeasurementRead will set it to DIP_TRUE when it has been able to read filename, and it will set it to DIP_FALSE when it is not able to read the file. No error will be generated in this case.

NOTE

There are currently no measurement reading functions, so this function will always fail.

ARGUMENTS

Data type	\mathbf{Name}	Description
dip_Measurement	measurement	Measurement data structure
dip_String	filename	File name to read from
dip_int	format	ID of file format
dip_Boolean	addExtensions	Add extensions when looking for the file
dip_Boolean *	recognised	Set to DIP_TRUE if the file was found

SEE ALSO

Measure, MeasurementWrite

MeasurementSetName

Set the name of a measurement structure

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error dip_MeasurementSetName ( measurement, name )
```

FUNCTION

This function sets the name of measurement to name.

ARGUMENTS

Data type	Name	Description
dip_Measurement	measurement	Measurement
dip_String	name	Name

SEE ALSO

 ${\tt Measure,\,MeasurementNew,\,MeasurementID,\,MeasurementGetName,}\\ {\tt MeasurementGetPhysicalDimensions,\,MeasurementSetPhysicalDimensions}\\$

${\tt MeasurementSetPhysicalDimensions}$

Set the physical dimensions info of the measurement

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error dip_MeasurementSetPhysicalDimensions ( measurement, physDims )
```

FUNCTION

This function sets the physical dimensions information for the measurement data structure. The physical dimensions data structure informs measurement features about the physical sizes and position of the pixels of the measured image.

ARGUMENTS

Data type	Name	Description
dip_Measurement	measurement	Measurement data structure
dip_PhysicalDimensions	physDims	Physical Dimensions data structure

SEE ALSO

 ${\tt MeasurementNew,\,MeasurementID,\,MeasurementSetName,\,MeasurementGetName,\,MeasurementGetPhysicalDimensions}$

MeasurementToHistogram

Creats a histogram for a measurement

SYNOPSIS

#include "dip_measurement.h" #include "dip_distribution.h"
dip_Error dip_MeasurementToHistogram (histogram, measurement, featureID, binSize,
maximum, minimum, percentage, addMeasurement)

DATA TYPES

integer

FUNCTION

This function creates a (possibly multi-dimensional) histogram with the measurement results of one feature. If addMeasurement is DIP_TRUE, new data points are added to the existing histogram, and binSize, maximum, minimum and percentage input arguments are ignored. Otherwise, histogram is destroyed and recreated according to the chosen values for binSize, maximum, minimum and percentage. If percentage is DIP_TRUE, maximum and minimum represent a percentage of the data range, otherwise they represent absolute values. If maximum or minimum are NULL, the maximum or minimum of the data is used.

ARGUMENTS

Data type	Name	Description
dip_Distribution	histogram	Output histogram
dip_Measurement	measurement	Measurement data
dip_IntegerArray	featureID	List of feature IDs to use
dip_FloatArray	binSize	Size of the histogram bins
dip_FloatArray	maximum	Maximum value represented in the histogram
dip_FloatArray	minimum	Minimum value represented in the histogram
dip_Boolean	percentage	Whether maximum and minimum are percentages
dip_Boolean	addMeasurement	Whether to add data to histogram or create a new one

SEE ALSO

Measure, MeasurementToImage, ObjectToMeasurement

MeasurementToImage

Exports the data in a measurement structure to an image

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error dip_MeasurementToImage ( measurement, out, measurementIDs, objects )
```

DATA TYPES

float

FUNCTION

This function creates an image and writes the measurement data to it as if it were a table, measurements along the first dimension, objects along the second dimension.

ARGUMENTS

Data type	Name	Description
dip_Measurement	measurement	Input measurement
dip_Image	out	Output image
dip_IntegerArray	measurementIDs	List of measurement IDs to export
dip_IntegerArray	objects	List of object IDs to export

SEE ALSO

Measure, ObjectToMeasurement, MeasurementToHistogram

MeasurementWrite

Write measurement results to a file

SYNOPSIS

```
#include "dipio_measurement.h"
dip_Error dipio_MeasurementWrite ( measurement, filename, format, labels )
```

FUNCTION

This function writes measurement data to a file, overwriting any other file with the same name. Get the format ID through the registry functions. If format is 0, CSV is used.

ARGUMENTS

Data type	Name	Description
dip_Measurement	measurement	Measurement data structure
dip_String	filename	File name to write to
dip_int	format	ID of file format
dip_Boolean	labels	DIP_TRUE to write labels to file

SEE ALSO

 ${\tt Measure}, {\tt MeasurementRead}, {\tt MeasurementWriteCSV}, {\tt MeasurementWriteHTML}, {\tt MeasurementWriteText}$

MeasurementWriteCSV

Write measurement results to a CSV file

SYNOPSIS

#include "dipio_measurement.h"
dip_Error dipio_MeasurementWriteCSV (measurement, filename, separator, labels)

FUNCTION

This function writes the measurement results to a comma separated values (CSV) file, overwriting any other file with the same name.

This function calls MeasurementWriteText with the proper settings.

ARGUMENTS

Data type	Name	Description
dip_Measurement	measurement	Measurement data structure
dip_String	filename	File name to write to
char *	separator	Characters to separate values
dip_Boolean	labels	DIP_TRUE to write labels to file

SEE ALSO

Measure, MeasurementWrite, MeasurementWriteText

MeasurementWriteHTML

Write measurement results to an HTML file

SYNOPSIS

#include "dipio_measurement.h"
dip_Error dipio_MeasurementWriteHTML (measurement, filename, separator, labels)

FUNCTION

This function writes the measurement results to a formatted HTML file, overwriting any other file with the same name.

ARGUMENTS

Data type	Name	Description
dip_Measurement	measurement	Measurement data structure
dip_String	filename	File name to write to
char *	separator	Characters to separate values
dip_Boolean	labels	DIP_TRUE to write labels to file

SEE ALSO

Measure, MeasurementWrite

MeasurementWriteText

Write measurement results as readable text

SYNOPSIS

```
#include "dipio_measurement.h"
dip_Error dipio_MeasurementWriteText ( measurement, fp, options )
```

FUNCTION

This function saves/prints the results of a measurement stored in the measurement data structure. Since it will save the results to the fp FILE pointer (which has to be opened before this function is called, and closed afterwards), the results can be printed to a screen (specify stdout as fp) or to a file.

The results are saved in a matrix, with a column for each measurement, and a row for each object. The first column contains the object ID. The options structure provides a means to adjust the formatting of the measurement data. Its separator variable specifies the column separator character, the rows are separated by a newline. If the labelAlign variable is DIP_TRUE, the separator is repeated such that the columns are aligned. If the labels variable is DIP_TRUE, the first row contains measurement labels, and info specifies whether or not the short description of each measurement function should be printed before the result matrix. If results is DIP_FALSE, the measurement values are not printed.

ARGUMENTS

Data type	Name	Description
dip_Measurement	measurement	Measurement data structure
FILE *	fp	FILE pointer to which the results are saved
dipio_WriteTextFormat	options	Text formatting options

The structure dipio_WriteTextFormat contains the following elements:

Data type	Name	Description
char *	separator	Column separator character
dip_Boolean	info	Write descriptio
dip_Boolean	labels	Write labels
dip_Boolean	results	Write values
dip_Boolean	labelAlign	Align columns

SEE ALSO

Measure, MeasurementWrite

Median

statistics function

SYNOPSIS

dip_Error dip_Median (in, mask, out, ps)

DATA TYPES

binary, integer, float, complex

FUNCTION

Calculates the median of the pixel values over all those dimensions which are specified by ps.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	mask (0)	Mask
dip_Image	out	Output
dip_BooleanArray	ps (0)	Dimensions to project

SEE ALSO

From images to scalars

 ${\tt Sum,\ Mean,\ Variance,\ Standard Deviation,\ Mean Modulus,\ Sum Modulus,\ Mean Square Modulus,\ Maximum,\ Minimum,\ Percentile}$

MedianFilter

Non-linear smoothing filter

SYNOPSIS

```
#include "dip_morphology.h"
dip_Error dip_MedianFilter ( in, out, se, boundary, param, shape )
```

DATA TYPES

integer, **float**

FUNCTION

Median filter with different filter shapes.

Only the rectangular, elliptic and diamond filter shapes are supported (DIP_FLT_SHAPE_RECTANGULAR, DIP_FLT_SHAPE_ELLIPTIC and DIP_FLT_SHAPE_DIAMOND). Other filter shapes can be implemented by setting shape to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, and passing a binary image in se. The "on" pixels define the shape of the filter window. Other values of shape are illegal.

If shape is not equal to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, se can be set to zero. When shape is set to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, param is ignored, and can be set to zero.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input image
dip_Image	out	Output image
dip_Image	se	Custom filter shape (binary)
dip_BoundaryArray	boundary	Boundary conditions
dip_FloatArray	param	Filter parameters
dip_FilterShape	shape	Filter shape

The enumerator dip_FilterShape contains the following constants:

Name	Description
DIP_FLT_SHAPE_DEFAULT	Default filter window, same as
	DIP_FLT_SHAPE_RECTANGULAR
DIP_FLT_SHAPE_RECTANGULAR	Rectangular filter window, can be even in size
DIP_FLT_SHAPE_ELLIPTIC	Elliptic filter window, always odd in size
DIP_FLT_SHAPE_DIAMOND	Diamond-shaped filter window, always odd in size
DIP_FLT_SHAPE_PARABOLIC	Parabolic filter window (morphology only)
DIP_FLT_SHAPE_DISCRETE_LINE	Rotated line structuring element (morphology only)
DIP_FLT_SHAPE_INTERPOLATED_LINE	Rotated line structuring element, through interpolation
	(morphology only)
DIP_FLT_SHAPE_PERIODIC_LINE	(not implemented)
DIP_FLT_SHAPE_STRUCTURING_ELEMENT	Use se as filter window, can be any size

SEE ALSO

PercentileFilter, Uniform, Sigma

${\tt MemoryCopy}$

Copy memory blocks

SYNOPSIS

void dip_MemoryCopy(in, out, number)

FUNCTION

Copy a memory block

ARGUMENTS

Data type	Name	Description
void *	in	pointer to memory source block
void *	out	pointer to memory destination block
dip_int	number	number of bytes to be copied

NOTE

The behaviour of this function is undefined when the in and out blocks overlap.

MemoryFree

Free a chunk of memory

SYNOPSIS

dip_Error dip_MemoryFree(pointer)

FUNCTION

Frees a chunk of memory.

ARGUMENTS

Data type	Name	Description
void *	pointer	pointer to an allocated chunk of memory

SEE ALSO

 ${\tt MemoryNew}, \, {\tt MemoryReallocate}, \, {\tt MemoryFunctionsSet}$

MemoryFunctionsSet

Sets memory allocation functions

SYNOPSIS

 $\label{thm:constraint} \begin{tabular}{ll} $\tt dip_Error \ dip_MemoryFunctionsSet(MemoryNewFunction, MemoryReallocateFunction, MemoryFreeFunction) \end{tabular}$

FUNCTION

Sets the memory allocation functions used by DIPlib.

ARGUMENTS

Data type	Name	Description
dip_MemoryNewFunction	MemoryNewFunction	pointer to a memory allocation
		function
dip_MemoryReallocateFunction	MemoryReallocateFunction	pointer to a memory
		reallocation function
dip_MemoryFreeFunction	MemoryFreeFunction	pointer to a memory freeing
		function

NOTE

The three allocation functions are defined as follows:

```
typedef void* (*dip_MemoryNewFunction)(size_t size)
```

typedef void* (*dip_MemoryReallocateFunction)(void *ptr, size_t size)

typedef void (*dip_MemoryFreeFunction)(void *ptr)

And are by default set to malloc, realloc and free.

SEE ALSO

MemoryNew, MemoryReallocate, MemoryFree

MemoryNew

Allocate and track memory

SYNOPSIS

dip_Error dip_MemoryNew(pointer, size, resources)

FUNCTION

Allocates a chunk of memory, and adds a reference to the chunk to the list of tracked resources.

ARGUMENTS

Data type	Name	Description
void **	pointer	pointer to the memory chunk pointer
size_t	size	size of the memory chunk in bytes
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

 ${\tt MemoryReallocate}, {\tt MemoryFree}, {\tt MemoryFunctionsSet}$

MemoryReallocate

Reallocate a chunk of memory

SYNOPSIS

dip_Error dip_MemoryReallocate (pointer, newsize, resources)

FUNCTION

Reallocates a chunk of memory, to change its size. resources must be the dip_Resources structure used in the call to MemoryNew when pointer was allocated.

ARGUMENTS

Data type	Name	Description
void **	pointer	pointer to the memory chunk pointer
size_t	newsize	size of the memory chunk in bytes
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

MemoryNew, MemoryFree, MemoryFunctionsSet

mErf

mathematical function

SYNOPSIS

dip_float dipm_Erf (x)

FUNCTION

Computes the error function of the input value.

ARGUMENTS

Data type	Name	Description
dip_float	х	Input value

SEE ALSO

mTruncate, mFraction, mNearestInt, mSign, mExp2, mExp10, mLog2, mSinc, mReciprocal, mBesselJ0, mBesselJ1, mBesselJN, mBesselY0, mBesselY1, mBesselYN, mLnGamma, mErfc, mGammaP, mGammaQ

mErfc

mathematical function

SYNOPSIS

dip_float dipm_Erfc (x)

FUNCTION

Computes the complementary error function of the input value.

ARGUMENTS

Data type	Name	Description
dip_float	х	Input value

SEE ALSO

mTruncate, mFraction, mNearestInt, mSign, mExp2, mExp10, mLog2, mSinc, mReciprocal, mBesselJ0, mBesselJ1, mBesselJN, mBesselY0, mBesselY1, mBesselYN, mLnGamma, mErf, mGammaP, mGammaQ

mExp10

mathematical function

SYNOPSIS

dip_float dipm_Exp10 (x)

FUNCTION

Computes the base ten exponent of the input value.

ARGUMENTS

Data type	Name	Description
dip_float	value	Value

SEE ALSO

mTruncate, mFraction, mNearestInt, mSign, mExp2, mLog2, mSinc, mReciprocal, mBesselJ0, mBesselJ1, mBesselJN, mBesselY0, mBesselY1, mBesselYN, mLnGamma, mErf, mErfc, mGammaP, mGammaQ

mExp2

mathematical function

SYNOPSIS

dip_float dipm_Exp2 (x)

FUNCTION

Computes the base two exponent of the input value.

ARGUMENTS

Data type	Name	Description
dip_float	value	Value

SEE ALSO

mTruncate, mFraction, mNearestInt, mSign, mExp10, mLog2, mSinc, mReciprocal, mBesselJ0, mBesselJ1, mBesselJN, mBesselY0, mBesselY1, mBesselYN, mLnGamma, mErf, mErfc, mGammaP, mGammaQ

mFraction

mathematical function

SYNOPSIS

dip_float dipm_Fraction (x)

FUNCTION

Computes the fraction of the input value.

ARGUMENTS

Data type	Name	Description	
dip_float	х	Input value	

SEE ALSO

mTruncate, mNearestInt, mSign, mExp2, mExp10, mLog2, mSinc, mReciprocal, mBesselJ0, mBesselJ1, mBesselJN, mBesselY0, mBesselY1, mBesselYN, mLnGamma, mErf, mErfc, mGammaP, mGammaQ

mGammaP

mathematical function

SYNOPSIS

dip_float dipm_GammaP (a, x)

FUNCTION

Computes the incomplete gamma function of the input value.

ARGUMENTS

Data type	Name	Description	
dip_float	a	A	
dip_float	х	Input value	

SEE ALSO

mTruncate, mFraction, mNearestInt, mSign, mExp2, mExp10, mLog2, mSinc, mReciprocal, mBesselJ0, mBesselJ1, mBesselJN, mBesselY0, mBesselY1, mBesselYN, mLnGamma, mErf, mErfc, mGammaQ

mGammaQ

mathematical function

SYNOPSIS

dip_float dipm_GammaQ (a, x)

FUNCTION

Computes the complementary incomplete gamma function of the input value.

ARGUMENTS

Data type	Name	Description	
dip_float	a	A	
dip_float	х	Input value	

SEE ALSO

mTruncate, mFraction, mNearestInt, mSign, mExp2, mExp10, mLog2, mSinc, mReciprocal, mBesselJ0, mBesselJ1, mBesselJN, mBesselY0, mBesselY1, mBesselYN, mLnGamma, mErf, mErfc, mGammaP

Min

arithmetic function

SYNOPSIS

dip_Error dip_Min (in1, in2, out)

DATA TYPES

binary, integer, float

FUNCTION

This function computes $\mathtt{out} = \min(\mathtt{in1}\ , \mathtt{in2})$ on a pixel by pixel basis. The data types of the $\mathtt{in1}$ and $\mathtt{in2}$ image may be of different types. See Information about dyadic operations for more information about what the type of the output will be.

ARGUMENTS

Data type	Name	Description		
$dip_{-}Image$	in1	First input		
dip_Image	in2	Second input		
dip_Image	out	Output		

SEE ALSO

Max, MaxFloat, MinFloat

MinFloat

arithmetic function

SYNOPSIS

dip_Error dip_MinFloat (in, out, constant)

DATA TYPES

binary, integer, float

FUNCTION

This function computes out = min(in, constant) on a pixel by pixel basis. The data types of the in1 image and constant may be of different types. See Information about dyadic operations for more information about what the type of the output will be.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	in	Input
dip_Image	out	Output
dip_float	constant	Constant

SEE ALSO

Max, MaxFloat, Min

Minima

Detects local minima

SYNOPSIS

```
#include "dip_filtering.h"
dip_Error dip_Minima ( in, mask, out, connectivity, booleanOutput )
```

DATA TYPES

integer, float

FUNCTION

This function detects local minima.

The algorithm finds a connected set of pixels with identical value, an no neighbours with lower value. This set is a local minimum and its pixels are set to 1 in the output image. If booleanOutput is false, the output image is a labelled image.

This function differs from LocalMinima in that it marks every minimum. LocalMinima is able to filter out unimportant minima.

For images that have large plateaus (regions of constant value) that are not local minima, this function can be quite slow. For example, an image that is one everywhere except for a small valley somewhere. For such an image it is recommended to use the mask input, for example with the output of a threshold operation.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input image
dip_Image	mask	Mask image
dip_Image	out	Binary output image
dip_int	connectivity	Connectivity
dip_Boolean	booleanOutput	Give a binary output image?

NOTE

If you are looking for the old version of Minima, it is still available through the following combination of commands:

```
dip_Erosion( in, out, se, boundary, param, shape );
dip_Equal( in, out, out );
```

SEE A	ALSO
-------	------

 ${\tt Maxima, Subpixel Minima, Local Minima, Seeded Watershed, Grow Regions}$

Minimum

statistics function

SYNOPSIS

dip_Error dip_Minimum (in, mask, out, ps)

DATA TYPES

binary, integer, float

FUNCTION

Calculates the minimum of the pixel values over all those dimensions which are specified by ps.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	mask (0)	Mask
dip_Image	out	Output
dip_BooleanArray	ps (0)	Dimensions to project

SEE ALSO

From images to scalars

 ${\tt Sum,\ Mean,\ Variance,\ Standard Deviation,\ Mean Modulus,\ Sum Modulus,\ Mean Square Modulus,\ Maximum,\ Median,\ Percentile}$

Mirror

Mirrors an image

SYNOPSIS

```
#include "dip_manipulation.h"
dip_Error dip_Mirror ( in, out, mirror )
```

DATA TYPES

binary, integer, float, complex

FUNCTION

This function mirrors the pixels in those dimensions of image as specified by mirror.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input image
dip_Image	out	Output image
dip_BooleanArray	mirror	Mirror flags

SEE ALSO

Map

mLnGamma

mathematical function

SYNOPSIS

dip_float dipm_LnGamma (x)

FUNCTION

Computes the natural logarithm of the gamma function of the input value.

ARGUMENTS

Data type	Name	Description
dip_float	value	Value

SEE ALSO

mTruncate, mFraction, mNearestInt, mSign, mExp2, mExp10, mLog2, mSinc, mReciprocal, mBesselJ0, mBesselJ1, mBesselJN, mBesselY0, mBesselY1, mBesselYN, mErf, mErfc, mGammaP, mGammaQ

mLog2

mathematical function

SYNOPSIS

dip_float dipm_Log2 (x)

FUNCTION

Computes the base two logarithm of the input value.

ARGUMENTS

Data type	Name	Description	
dip_float	value	Value	

SEE ALSO

mTruncate, mFraction, mNearestInt, mSign, mExp2, mExp10, mSinc, mReciprocal, mBesselJ0, mBesselJ1, mBesselJN, mBesselY0, mBesselY1, mBesselYN, mLnGamma, mErf, mErfc, mGammaP, mGammaQ

mNearestInt

mathematical function

SYNOPSIS

dip_float dipm_NearestInt (x)

FUNCTION

Computes the nearest int of the input value.

ARGUMENTS

Data type	Name	Description
dip_float	х	Input value

SEE ALSO

mTruncate, mFraction, mSign, mExp2, mExp10, mLog2, mSinc, mReciprocal, mBesselJ0, mBesselJ1, mBesselJN, mBesselY1, mBesselYN, mLnGamma, mErf, mErfc, mGammaP, mGammaQ

Modulo

Arithmetic function

SYNOPSIS

dip_Error dip_Modulo (in, out, period)

DATA TYPES

integer

FUNCTION

Computes the modulo of the input image values, by computing the remainder of the the division of the input image values with period.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	in	Input
$dip_{-}Image$	out	Output
dip_int	period	Period

SEE ALSO

Div, DivFloat, DivComplex, Reciprocal

Modulus

Arithmetic function

SYNOPSIS

dip_Error dip_Modulus (in, out)

DATA TYPES

binary, integer, float, $\mathbf{complex}$

FUNCTION

Computes the modulus of the input image values, and outputs a float typed image.

ARGUMENTS

Data type	Name	Description	
$dip_{-}Image$	in	Input	
dip_Image	out	Output	

SEE ALSO

Phase, Real, Imaginary

MonadicFrameWork

FrameWork for monadic operations

SYNOPSIS

dip_Error dip_MonadicFrameWork (in, out, processBoundary, processBorder, process)

FUNCTION

This function is a frontend on the SeparableFrameWork. It provides an easier interface for filters that only need to scan an image once. The dimension in which the image should be scanned can be specified or left to MonadicFrameWork by specifiying the dimension with DIP_MONADIC_OPTIMAL_DIMENSION.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
$\mathtt{dip}_{-}\mathtt{Image}$	out	Output
dip_Boundary	processBoundary	ProcessBoundary
dip_int	processBorder	ProcessBorder
dip_FrameWorkProcess	process	Process

SEE ALSO

SeparableFrameWork, SingleOutputFrameWork

${\tt MorphologicalGradientMagnitude}$

Morphological edge detector

SYNOPSIS

#include "dip_morphology.h"
dip_Error dip_MorphologicalGradientMagnitude (in, out, se, boundary, param, shape,
edgeType)

DATA TYPES

integer, float

FUNCTION

The same as MorphologicalRange.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	out	Output
$\mathtt{dip}_{-}\mathtt{Image}$	se	Custom structuring element
dip_BoundaryArray	boundary	Boundary conditions
dip_FloatArray	param	Filter parameters
dip_FilterShape	shape	Structuring element
dip_MphEdgeType	edgeType	edgeType

SEE ALSO

 ${\tt Morphological Range, Lee, MultiScale Morphological Gradient, Tophat}$

MorphologicalRange

Morphological edge detector

SYNOPSIS

#include "dip_morphology.h"
dip_Error dip_MorphologicalRange (in, out, se, boundary, param, shape, edgeType)

DATA TYPES

integer, float

FUNCTION

Implements a morphological edge detector based on the difference of two complementary morphological operations. These can be chosen through the edgeType parameter.

The rectangular, elliptic and diamond structuring elements are "flat", i.e. these structuring elements have a constant value. For these structuring elements, param determines the sizes of the structuring elements.

When shape is DIP_FLT_SHAPE_DISCRETE_LINE or DIP_FLT_SHAPE_INTERPOLATED_LINE, the structuring element is a line. param->array[0] determines the length, param->array[1] the angle. This is currently only supported for 2D images. Interpolated lines use interpolation to obtain a more accurate result, but loose the strict increasingness and extensivity (these properties are satisfied only by approximation).

When shape is set to DIP_FLT_SHAPE_PARABOLIC, params specifies the curvature of the parabola.

When shape is set to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, se is used as structuring element. It can be either a binary or a grey-value image. Its origin is the center, or one pixel to the left of the center if the size is even.

If shape is not equal to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, se can be set to zero. When shape is set to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, param is ignored, and can be set to zero.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	out	Output
dip_Image	se	Custom structuring element
dip_BoundaryArray	boundary	Boundary conditions
dip_FloatArray	param	Filter parameters
dip_FilterShape	shape	Structuring element
dip_MphEdgeType	edgeType	edgeType

The enumerator dip_FilterShape contains the following constants:

Name	Description
DIP_FLT_SHAPE_DEFAULT	Default filter window, same as
	DIP_FLT_SHAPE_RECTANGULAR
DIP_FLT_SHAPE_RECTANGULAR	Rectangular filter window, can be even in size
DIP_FLT_SHAPE_ELLIPTIC	Elliptic filter window, always odd in size
DIP_FLT_SHAPE_DIAMOND	Diamond-shaped filter window, always odd in size
DIP_FLT_SHAPE_PARABOLIC	Parabolic filter window (morphology only)
DIP_FLT_SHAPE_DISCRETE_LINE	Rotated line structuring element (morphology only)
DIP_FLT_SHAPE_INTERPOLATED_LINE	Rotated line structuring element, through interpolation
	(morphology only)
DIP_FLT_SHAPE_PERIODIC_LINE	(not implemented)
DIP_FLT_SHAPE_STRUCTURING_ELEMENT	Use se as filter window, can be any size

The enumerator dip_MphEdgeType contains the following constants:

Name	Description	
DIP_MPH_TEXTURE	Response is limited to edges in texture	
DIP_MPH_OBJECT	Response is limited to object edges	
DIP_MPH_BOTH	All edges produce equal response	

SEE ALSO

 ${\tt Morphological Gradient Magnitude, Lee, Multi Scale Morphological Gradient, Tophat}$

MorphologicalReconstruction

Morphological filter

SYNOPSIS

#include "dip_morphology.h"
dip_Error dip_MorphologicalReconstruction (marker, mask, out, connectivity)

DATA TYPES

integer, float

FUNCTION

Dilation of the image marker, constrained by the image mask. out will be smaller or equal to mask. The image is grown according to the connectivity parameter. See The connectivity parameter for more information.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{\mathtt{I}}\mathtt{Image}$	marker	Marker input image
$\mathtt{dip}_{-}\mathtt{Image}$	mask	Mask input image
$\mathtt{dip}_{-}\mathtt{Image}$	out	Output image
dip_int	connectivity	Connectivity

KNOWN LIMITATIONS

The image marker is converted to the same data type as mask. If mask is an unsigned integer, and marker has negative values, then it is possible that these negative values will be turned into large positive values, yielding an unexpected output. The solution is to make sure that mask and marker are in compatible data types.

LITERATURE

K. Robinson and P.F. Whelan, Efficient morphological reconstruction: a downhill filter, Pattern Recognition Letters 25(15):1759-1767, 2004.

SEE ALSO

Dilation, BinaryPropagation, AreaOpening

MorphologicalSmoothing

Morphological smoothing filter

SYNOPSIS

#include "dip_morphology.h"
dip_Error dip_MorphologicalSmoothing (in, out, se, boundary, param, shape, flags)

DATA TYPES

integer, float

FUNCTION

Implements a morphological smoothing based on the sequence of two complementary morphological operations. These can be chosen through the dipf_MphSmoothing parameter.

The rectangular, elliptic and diamond structuring elements are "flat", i.e. these structuring elements have a constant value. For these structuring elements, param determines the sizes of the structuring elements.

When shape is DIP_FLT_SHAPE_DISCRETE_LINE or DIP_FLT_SHAPE_INTERPOLATED_LINE, the structuring element is a line. param->array[0] determines the length, param->array[1] the angle. This is currently only supported for 2D images. Interpolated lines use interpolation to obtain a more accurate result, but loose the strict increasingness and extensivity (these properties are satisfied only by approximation).

When shape is set to DIP_FLT_SHAPE_PARABOLIC, params specifies the curvature of the parabola.

When shape is set to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, se is used as structuring element. It can be either a binary or a grey-value image. Its origin is the center, or one pixel to the left of the center if the size is even.

If shape is not equal to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, se can be set to zero. When shape is set to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, param is ignored, and can be set to zero.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	out	Output
dip_Image	se	Custom structuring element
dip_BoundaryArray	boundary	Boundary conditions
dip_FloatArray	param	Filter parameters
dip_FilterShape	shape	Structuring element
dipf_MphSmoothing	flags	flags

The enumerator dip_FilterShape contains the following constants:

Name	Description
DIP_FLT_SHAPE_DEFAULT	Default filter window, same as
	DIP_FLT_SHAPE_RECTANGULAR
DIP_FLT_SHAPE_RECTANGULAR	Rectangular filter window, can be even in size
DIP_FLT_SHAPE_ELLIPTIC	Elliptic filter window, always odd in size
DIP_FLT_SHAPE_DIAMOND	Diamond-shaped filter window, always odd in size
DIP_FLT_SHAPE_PARABOLIC	Parabolic filter window (morphology only)
DIP_FLT_SHAPE_DISCRETE_LINE	Rotated line structuring element (morphology only)
DIP_FLT_SHAPE_INTERPOLATED_LINE	Rotated line structuring element, through interpolation
	(morphology only)
DIP_FLT_SHAPE_PERIODIC_LINE	(not implemented)
DIP_FLT_SHAPE_STRUCTURING_ELEMENT	Use se as filter window, can be any size

The enumerator ${\tt dipf_MphSmoothing}$ contains the following constants:

Name	Description
DIP_MPH_OPEN_CLOSE	First the opening, then the closing
DIP_MPH_CLOSE_OPEN	First the closing, then the opening
DIP_MPH_AVERAGE	The average of the result of the two above

SEE ALSO

 ${\tt MorphologicalThreshold}, {\tt Tophat}$

MorphologicalThreshold

Morphological smoothing filter

SYNOPSIS

```
#include "dip_morphology.h"
dip_Error dip_MorphologicalThreshold ( in, out, se, boundary, param, shape, edgeType
)
```

DATA TYPES

integer, float

FUNCTION

Implements a morphological smoothing based on the average of two complementary morphological operations. These can be chosen through the edgeType parameter.

The rectangular, elliptic and diamond structuring elements are "flat", i.e. these structuring elements have a constant value. For these structuring elements, param determines the sizes of the structuring elements.

When shape is DIP_FLT_SHAPE_DISCRETE_LINE or DIP_FLT_SHAPE_INTERPOLATED_LINE, the structuring element is a line. param->array[0] determines the length, param->array[1] the angle. This is currently only supported for 2D images. Interpolated lines use interpolation to obtain a more accurate result, but loose the strict increasingness and extensivity (these properties are satisfied only by approximation).

When shape is set to DIP_FLT_SHAPE_PARABOLIC, params specifies the curvature of the parabola.

When shape is set to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, se is used as structuring element. It can be either a binary or a grey-value image. Its origin is the center, or one pixel to the left of the center if the size is even.

If shape is not equal to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, se can be set to zero. When shape is set to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, param is ignored, and can be set to zero.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
$\mathtt{dip}_{\mathtt{-}}\mathtt{Image}$	out	Output
$\mathtt{dip}_{-}\mathtt{Image}$	se	Custom structuring element
dip_BoundaryArray	boundary	Boundary conditions
dip_FloatArray	param	Filter parameters
dip_FilterShape	shape	Structuring element
dip_MphEdgeType	edgeType	edgeType

The enumerator dip_FilterShape contains the following constants:

Name	Description
DIP_FLT_SHAPE_DEFAULT	Default filter window, same as
	DIP_FLT_SHAPE_RECTANGULAR
DIP_FLT_SHAPE_RECTANGULAR	Rectangular filter window, can be even in size
DIP_FLT_SHAPE_ELLIPTIC	Elliptic filter window, always odd in size
DIP_FLT_SHAPE_DIAMOND	Diamond-shaped filter window, always odd in size
DIP_FLT_SHAPE_PARABOLIC	Parabolic filter window (morphology only)
DIP_FLT_SHAPE_DISCRETE_LINE	Rotated line structuring element (morphology only)
DIP_FLT_SHAPE_INTERPOLATED_LINE	Rotated line structuring element, through interpolation
	(morphology only)
DIP_FLT_SHAPE_PERIODIC_LINE	(not implemented)
DIP_FLT_SHAPE_STRUCTURING_ELEMENT	Use se as filter window, can be any size

The enumerator ${\tt dip_MphEdgeType}$ contains the following constants:

Name	Description
DIP_MPH_TEXTURE	Response is limited to edges in texture
DIP_MPH_OBJECT	Response is limited to object edges
DIP_MPH_BOTH	All edges produce equal response

SEE ALSO

 ${\tt MorphologicalSmoothing, Tophat}$

mReciprocal

mathematical function

SYNOPSIS

dip_float dipm_Reciprocal (x)

FUNCTION

Computes the reciprocal (1/x) of the input value.

ARGUMENTS

Data type	Name	Description
dip_float	х	Input value

SEE ALSO

mTruncate, mFraction, mNearestInt, mSign, mExp2, mExp10, mLog2, mSinc, mBesselJ0, mBesselJ1, mBesselJN, mBesselY1, mBesselYN, mLnGamma, mErf, mErfc, mGammaP, mGammaQ

mSign

mathematical function

SYNOPSIS

dip_float dipm_Sign (x)

FUNCTION

Computes the sign of the input value. The sign of zero is defined as zero.

ARGUMENTS

Data type	Name	Description
dip_float	х	Input value

SEE ALSO

mTruncate, mFraction, mNearestInt, mExp2, mExp10, mLog2, mSinc, mReciprocal, mBesselJ0, mBesselJ1, mBesselJN, mBesselY0, mBesselY1, mBesselYN, mLnGamma, mErf, mErfc, mGammaP, mGammaQ

mSinc

mathematical function

SYNOPSIS

dip_float dipm_Sinc (x)

FUNCTION

Computes the sinc $(\sin(x)/x)$ of the input value.

ARGUMENTS

Data type	Name	Description
dip_float	х	Input value

SEE ALSO

mTruncate, mFraction, mNearestInt, mSign, mExp2, mExp10, mLog2, mReciprocal, mBesselJ0, mBesselJ1, mBesselJN, mBesselY0, mBesselY1, mBesselYN, mLnGamma, mErf, mErfc, mGammaP, mGammaQ

mTruncate

mathematical function

SYNOPSIS

dip_float dipm_Truncate (x)

FUNCTION

Truncates the input value.

ARGUMENTS

Data type	Name	Description
dip_float	х	Input value

SEE ALSO

mFraction, mNearestInt, mSign, mExp2, mExp10, mLog2, mSinc, mReciprocal, mBesselJ0, mBesselJ1, mBesselJN, mBesselY0, mBesselY1, mBesselYN, mLnGamma, mErf, mErfc, mGammaP, mGammaQ

Mul

arithmetic function

SYNOPSIS

```
dip_Error dip_Mul ( in1, in2, out )
Calls Arith ( in1, in2, out, DIP_ARITHOP_MUL, DIP_DT_MINIMUM )
```

MulComplex

arithmetic function

SYNOPSIS

dip_Error dip_MulComplex (in, out, constant)

DATA TYPES

binary, integer, float, complex

FUNCTION

This function computes out = in * constant on a pixel by pixel basis. The data types of the in1 image and constant may be of different types. See Information about dyadic operations for more information about what the type of the output will be.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	in	Input
dip_Image	out	Output
dip_complex	constant	Constant

SEE ALSO

Arith, Arith_ComplexSeparated, MulInteger, MulFloat, AddComplex, SubComplex, MulConjugateComplex, DivComplex

${\tt MulConjugate}$

arithmetic function

SYNOPSIS

```
dip_Error dip_MulConjugate ( in1, in2, out )
Calls Arith ( in1, in2, out, DIP_ARITHOP_MUL_CONJUGATE, DIP_DT_MINIMUM )
```

${\tt MulConjugateComplex}$

arithmetic function

SYNOPSIS

dip_Error dip_MulConjugateComplex (in, out, constant)

DATA TYPES

binary, integer, float, complex

FUNCTION

This function computes out = in * conj(constant) on a pixel by pixel basis. The data types of the in1 image and constant may be of different types. See Information about dyadic operations for more information about what the type of the output will be.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	in	Input
$\mathtt{dip}_{-}\mathtt{Image}$	out	Output
dip_complex	constant	Constant

SEE ALSO

Arith, Arith_ComplexSeparated, AddComplex, SubComplex, MulComplex, DivComplex

MulFloat

arithmetic function

SYNOPSIS

dip_Error dip_MulFloat (in, out, constant)

DATA TYPES

binary, integer, float, complex

FUNCTION

This function computes out = in * constant on a pixel by pixel basis. The data types of the in1 image and constant may be of different types. See Information about dyadic operations for more information about what the type of the output will be.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	in	Input
dip_Image	out	Output
dip_float	constant	Constant

SEE ALSO

Arith, Arith_ComplexSeparated, MulInteger, MulComplex, AddFloat, SubFloat, DivFloat

MulInteger

arithmetic function

SYNOPSIS

dip_Error dip_MulInteger (in, out, constant)

DATA TYPES

 $binary, \, \mathbf{integer}, \, \mathbf{float}, \, \mathbf{complex}$

FUNCTION

This function computes out = in * constant on a pixel by pixel basis. The data types of the in1 image and constant may be of different types. See Information about dyadic operations for more information about what the type of the output will be.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	in	Input
dip_Image	out	Output
dip_int	constant	Constant

SEE ALSO

 ${\tt Arith, Arith_ComplexSeparated, MulFloat, MulComplex, AddInteger, SubInteger, DivInteger}$

MultiScaleMorphologicalGradient

Morphological edge detector

SYNOPSIS

#include "dip_morphology.h"
dip_Error dip_MultiScaleMorphologicalGradient (in, out, se, boundary, upperSize,
lowerSize, shape)

DATA TYPES

integer, float

FUNCTION

This function computes the average morphological gradient over a range of scales bounded by upperSize and lowerSize. The morphological gradient is computed as the difference of the dilation and erosion of the input image at a particular scale, eroded by an erosion of one size smaller. At the lowest scale, the size of the structuring element is 2 * upperSize + 1.

The rectangular, elliptic and diamond structuring elements are "flat", i.e. these structuring elements have a constant value. For these structuring elements, param determines the sizes of the structuring elements.

When shape is DIP_FLT_SHAPE_DISCRETE_LINE or DIP_FLT_SHAPE_INTERPOLATED_LINE, the structuring element is a line. param->array[0] determines the length, param->array[1] the angle. This is currently only supported for 2D images. Interpolated lines use interpolation to obtain a more accurate result, but loose the strict increasingness and extensivity (these properties are satisfied only by approximation).

When shape is set to DIP_FLT_SHAPE_PARABOLIC, params specifies the curvature of the parabola.

When shape is set to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, se is used as structuring element. It can be either a binary or a grey-value image. Its origin is the center, or one pixel to the left of the center if the size is even.

If shape is not equal to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, se can be set to zero. When shape is set to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, param is ignored, and can be set to zero.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	out	Output
dip_Image	se	Custom structuring element
dip_BoundaryArray	boundary	Boundary conditions
dip_int	upperSize	Upper size of structuring element
dip_int	lowerSize	Lower size of structuring element
dip_FilterShape	shape	Structuring element

The enumerator $\mathtt{dip_FilterShape}$ contains the following constants:

Name	Description
DIP_FLT_SHAPE_DEFAULT	Default filter window, same as
	DIP_FLT_SHAPE_RECTANGULAR
DIP_FLT_SHAPE_RECTANGULAR	Rectangular filter window, can be even in size
DIP_FLT_SHAPE_ELLIPTIC	Elliptic filter window, always odd in size
DIP_FLT_SHAPE_DIAMOND	Diamond-shaped filter window, always odd in size
DIP_FLT_SHAPE_PARABOLIC	Parabolic filter window (morphology only)
DIP_FLT_SHAPE_DISCRETE_LINE	Rotated line structuring element (morphology only)
DIP_FLT_SHAPE_INTERPOLATED_LINE	Rotated line structuring element, through interpolation
	(morphology only)
DIP_FLT_SHAPE_PERIODIC_LINE	(not implemented)
DIP_FLT_SHAPE_STRUCTURING_ELEMENT	Use se as filter window, can be any size

LITERATURE

D. Wang, Pattern Recognition, $\mathbf{30}(12)$, pp. 2043-2052, 1997

SEE ALSO

 ${\tt Lee}, {\tt Morphological Gradient Magnitude}, {\tt Morphological Range}, {\tt Tophat}$

NearestInt

Arithmetic function

SYNOPSIS

dip_Error dip_NearestInt (in, out)

DATA TYPES

binary, integer, float, $\mathbf{complex}$ binary, integer, \mathbf{float}

FUNCTION

Computes the nearest integer value of the input image values.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	out	Output

SEE ALSO

Abs, Ceil, Floor, Sign, Truncate, Fraction

NeighbourIndicesListMake

Get indices to direct neighbours

SYNOPSIS

```
#include "dip_neighbourlist.h"
dip_Error dip_NeighbourIndicesListMake ( stride, connectivity, indices, resources )
```

FUNCTION

A list indices is created with the relative indices to the direct neighbours of a pixel in an image whose strides are given by stride. How many direct neighbours are returned is controlled by connectivity, see The connectivity parameter for the available values and their meaning. indices is allocated and tracked in resources.

ARGUMENTS

Data type	Name	Description
dip_IntegerArray	stride	Stride array
dip_int	connectivity	Connectivity
dip_IntegerArray *	indices	Output neighbour indices
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

NeighbourListMake

Get list of direct neighbours

SYNOPSIS

```
#include "dip_neighbourlist.h"
dip_Error dip_NeighbourListMake ( ndims, connectivity, coords, resources )
```

FUNCTION

A list coords is created with the relative coordinates to the direct neighbours of a pixel in an ndims-dimensional image. How many direct neighbours are returned is controlled by connectivity, see The connectivity parameter for the available values and their meaning.

coords is allocated and tracked in resources.

ARGUMENTS

Data type	Name	Description
dip_int	ndims	Image dimensionality
dip_int	connectivity	Connectivity
dip_CoordinateArray *	coords	Output neighbour coordinates
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

NeighbourListMakeChamfer

Get list of neighbours based on Chamfer metric

SYNOPSIS

#include "dip_neighbourlist.h"
dip_Error dip_NeighbourListMakeChamfer (pixelsize, maxdistance, coords, distance,
resources)

FUNCTION

A list coords is created with the relative coordinates to the neighbours of a pixel in an pixelsize->size-dimensional image. Here, neighbours are all pixels within a maxdistance distance. pixelsize gives the size of a pixel, and hence controls the size of the neighbourhood with maxdistance. Anisotropic pixel grids are supported. distance contains the distance to each of the neighbours in coords.

Distances between two pixels are taken to be the Euclidean distance. There are better metrics described in the literature for small neighbourhoods, that yield a more isotropic measure when compounded over longer distances. These are not implemented in this function.

coords and distance are allocated and tracked in resources.

ARGUMENTS

Data type	Name	Description
dip_FloatArray	pixelsize	Physical dimensions of the pixels
dip_int	maxdistance	Maximum distance to which select neighbours
dip_CoordinateArray *	coords	Output neighbour coordinates
dip_FloatArray *	distance	Output distances to neighbours
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

NeighbourListMakeImage

Get list of neighbours based on metric in image

SYNOPSIS

```
#include "dip_neighbourlist.h"
dip_Error dip_NeighbourListMakeImage ( metric, coords, distance, resources )
```

FUNCTION

A list coords is created with the relative coordinates to the neighbours of a pixel in an image, with dimensionality as in metric. metric is an image that specifies the distance to each of the neighbours. This image must be odd in size, the centre pixel is the origin of the neighbourhood. Any pixel with value 0 is not considered part of the neighbourhood. distance contains the distance to each of the neighbours in coords.

Distances between two pixels are taken to be the Euclidean distance. There are better metrics described in the literature for small neighbourhoods, that yield a more isotropic measure when compounded over longer distances. These are not implemented in this function.

coords and distance are allocated and tracked in resources.

ARGUMENTS

Data type	Name	Description
dip_Image	metric	Image whose pixel values indicate the neighbour distance
dip_CoordinateArray *	coords	Output neighbour coordinates
dip_FloatArray *	distance	Output distances to neighbours
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

 $\label{lem:listMake} Neighbour List Make Chamfer, \ Neighbour List Make I mage, \\ Neighbour List To Indices, \ Neighbour Indices List Make \\$

NeighbourListToIndices

Get indices to neighbours

SYNOPSIS

```
#include "dip_neighbourlist.h"
dip_Error dip_NeighbourListToIndices ( coords, stride, indices, resources )
```

FUNCTION

This function translates the relative coordinates in coords into relative indices into an image with strides given by stride.

indices is allocated and tracked in resources.

ARGUMENTS

Data type	Name	Description
dip_CoordinateArray	coords	Input neighbour coordinates
dip_IntegerArray	stride	Stride array
dip_IntegerArray *	indices	Output neighbour indices
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

NormaliseSum

Normalise the sum of the pixel values

SYNOPSIS

```
#include "dip_manipulation.h"
dip_Error dip_NormaliseSum ( in, out, newSum )
```

DATA TYPES

binary, integer, float

FUNCTION

This function normalizes the sum of the pixel values in in to newSum, and puts the result in out.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{\mathtt{-}}\mathtt{Image}$	in	Input image
dip_Image	out	Output image
dip_float	newSum	New sum

NotEqual

Compare grey values in two images

SYNOPSIS

dip_Error dip_NotEqual (in1, in2, out)

DATA TYPES

binary, integer, float

FUNCTION

This function sets each pixel in out to "true" when corresponding pixels in in1 and in2 are different. This is the same as Compare with the DIP_SELECT_NOT_EQUAL selector flag.

in 2 can be a 0D image for comparison of pixel values with a single scalar value. This leads to the functionality of NotZero, but with more options.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{\mathtt{-}}\mathtt{Image}$	in1	First input
$\mathtt{dip}_{-}\mathtt{Image}$	in2	Second input
dip_Image	out	Output

SEE ALSO

Compare, Threshold, Equal, Greater, Lesser, NotGreater, NotLesser, SelectValue, NotZero

NotGreater

Compare grey values in two images

SYNOPSIS

dip_Error dip_NotGreater (in1, in2, out)

DATA TYPES

binary, integer, float

FUNCTION

This function sets each pixel in out to "true" when for corresponding pixels in1 <= in2. This is the same as Compare with the DIP_SELECT_LESSER_EQUAL selector flag.

in 2 can be a 0D image for comparison of pixel values with a single scalar value. This leads to a functionality similar to that of Threshold.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	in1	First input
dip_Image	in2	Second input
dip_Image	out	Output

SEE ALSO

Compare, Threshold, Equal, Greater, Lesser, NotEqual, NotLesser, SelectValue, NotZero

NotLesser

Compare grey values in two images

SYNOPSIS

dip_Error dip_NotLesser (in1, in2, out)

DATA TYPES

binary, integer, float

FUNCTION

This function sets each pixel in out to "true" when for corresponding pixels in1 >= in2. This is the same as Compare with the DIP_SELECT_GREATER_EQUAL selector flag.

in 2 can be a 0D image for comparison of pixel values with a single scalar value. This leads to a functionality similar to that of Threshold.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{\mathtt{-}}\mathtt{Image}$	in1	First input
$\mathtt{dip}_{-}\mathtt{Image}$	in2	Second input
dip_Image	out	Output

SEE ALSO

Compare, Threshold, Equal, Greater, Lesser, NotEqual, NotGreater, SelectValue, NotZero

NotZero

Point Operation

SYNOPSIS

```
#include "dip_point.h"
dip_Error dip_NotZero ( in, out )
```

DATA TYPES

integer, float

FUNCTION

This function returns a binary image with value 1 where in != 0 and value 0 elsewhere. The opposite can be accomplished with SelectValue: dip_SelectValue(in,out,0);.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	in	Input image
dip_Image	out	Output image

SEE ALSO

Threshold, SelectValue, Compare, RangeThreshold

ObjectToMeasurement

Convert object label value to measurement value

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error dip_ObjectToMeasurement ( object, intensity, out, connectivity, objectID,
featureID, measurementDim )
```

DATA TYPES

objectIm: integer intensityIm: integer, float

FUNCTION

This function produces an output image which pixel intensities are equal to the measurement value that the featureID measurement function measured on the object who label is defined by the pixel intensity of the corresponding pixel in object. This function is therefore useful to select (i.e. threshold) objects on basis of a measurement perfomed on the object. intensity provides pixel intensity information for measurements that require pixel intensity information of the objects, whose shape is defined by object.

The list of object IDs on which the measurements have to be performed is specified by objectID. If it is zero, ObjectToMeasurement will call GetObjectLabels to obtain a list of all non-zero values in objectIm.

If the featureID measurement function produces an array of measurement values, measurementDim will be used to select the desired array element.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	object	Object label image
dip_Image	intensity	Object intensity image
dip_Image	out	Output image
dip_int	connectivity	Connectivity of object's contour pixels, see The
		connectivity parameter
dip_IntegerArray	objectID	Array of Object IDs
dip_int	featureID	Measurement function ID
dip_int	measurementDim	Measurement results array index

SEE ALSO

Measure, SmallObjectsRemove, MeasurementToImage, MeasurementToHistogram

OneDimensionalSearch

Numerical algorithm

SYNOPSIS

```
#include "dip_numerical.h"
dip_Error dip_OneDimensionalSearch ( result, min, max, tol, func, dfunc, data,
searchfor )
```

FUNCTION

This function implements a numerical line-search for either the minimum or zero-crossing of a function. The obejctive is searched for in the range specified by min and man with a tolerance of tol. The search methods are based on Brent's algorithm. The dfunc parameter is preparation for support of search algorithms using derivative information. This is not supported in the current implementation, and dfunc should be set to zero.

ARGUMENTS

Data type	Name	Description
dip_float *	result	Result value
dip_float	min	Minimum value of search domain
dip_float	max	Maximum value of search domain
dip_float	tol	Tolerance
dip_OneDimensionalSearchFunction	func	Function
dip_OneDimensionalSearchFunction	dfunc	Derivative function
void *	data	User-supplied Data passed to func and
		dfunc
dipf_OneDimensionSearch	searchfor	Search for minimum of zero-crossing

Opening

Morphological opening operation

SYNOPSIS

```
#include "dip_morphology.h"
dip_Error dip_Opening ( in, out, se, boundary, param, shape )
```

DATA TYPES

integer, float, binary

FUNCTION

Grey-value opening with different structuring elements.

The rectangular, elliptic and diamond structuring elements are "flat", i.e. these structuring elements have a constant value. For these structuring elements, param determines the sizes of the structuring elements.

When shape is DIP_FLT_SHAPE_DISCRETE_LINE or DIP_FLT_SHAPE_INTERPOLATED_LINE, the structuring element is a line. param->array[0] determines the length, param->array[1] the angle. This is currently only supported for 2D images. Interpolated lines use interpolation to obtain a more accurate result, but loose the strict increasingness and extensivity (these properties are satisfied only by approximation).

When shape is set to DIP_FLT_SHAPE_PARABOLIC, params specifies the curvature of the parabola.

When shape is set to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, se is used as structuring element. It can be either a binary or a grey-value image. Its origin is the center, or one pixel to the left of the center if the size is even.

If shape is not equal to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, se can be set to zero. When shape is set to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, param is ignored, and can be set to zero.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	out	Output
dip_Image	se	Custom structuring element
dip_BoundaryArray	boundary	Boundary conditions
dip_FloatArray	param	Filter parameters
dip_FilterShape	shape	Structuring element

The enumerator dip_FilterShape contains the following constants:

Name	Description
DIP_FLT_SHAPE_DEFAULT	Default filter window, same as
	DIP_FLT_SHAPE_RECTANGULAR
DIP_FLT_SHAPE_RECTANGULAR	Rectangular filter window, can be even in size
DIP_FLT_SHAPE_ELLIPTIC	Elliptic filter window, always odd in size
DIP_FLT_SHAPE_DIAMOND	Diamond-shaped filter window, always odd in size
DIP_FLT_SHAPE_PARABOLIC	Parabolic filter window (morphology only)
DIP_FLT_SHAPE_DISCRETE_LINE	Rotated line structuring element (morphology only)
DIP_FLT_SHAPE_INTERPOLATED_LINE	Rotated line structuring element, through interpolation
	(morphology only)
DIP_FLT_SHAPE_PERIODIC_LINE	(not implemented)
DIP_FLT_SHAPE_STRUCTURING_ELEMENT	Use se as filter window, can be any size

SEE ALSO

Closing, Dilation, Erosion

0r

logic operation

SYNOPSIS

dip_Error dip_Or (in1, in2, out)

DATA TYPES

binary

FUNCTION

The function Or performs the logic OR operation between the corresponding pixels in in1 and in2, and stores the result in out.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	in1	First binary input image
dip_Image	in2	Second binary input image
dip_Image	out	Output image

SEE ALSO

Arith, And, Xor, Invert

ovl.h

Call an overloaded function

SYNOPSIS

```
dip_DataType ovlDataType = [current data type];
[#define DIP_OVL_ASSIGN [assignment target]]
#define DIP_OVL_FUNC [function base name]
#define DIP_OVL_ARGS [argument list]
#define DIP_OVL_ALLOW [type identifiers]
#include "dip_ovl.h"
```

FUNCTION

Call a type specific function based on the data type stored in the ovlDataType variable. The base name of the function is passed to dip_ovl.h by defining DIP_OVL_FUNC. The argument list is passed by defining DIP_OVL_ARGS. By defining DIP_OVL_ALLOW the list of data types for which overloading is possible can be controlled. If DIP_OVL_ALLOW is not defined, all data types are allowed. The list is specified by a logical OR of identifier and identifier group flags, see the table at DIPlib's data types. The code executed by dip_ovl.h is the following:

```
/* if ovlDataType is in the list specified by DIP_OVL_ALLOW */
DIPXJ( DIP_FUNC(DIP_OVL_FUNC,ovlDataType's extension) DIP_OVL_ARGS );
/* if ovlDataType is not in the list specified by DIP_OVL_ALLOW */
DIPSJ( DIP_E_DATA_TYPE_NOT_SUPPORTED );
```

DIP_FUNC is described in macros.h. Note that there are no brackets around DIP_OVL_ARGS, so they must be included in DIP_OVL_ARGS itself. If ovlDataType is one of the binary types, DIP_OVL_BINARY_ARGS can be defined to override DIP_OVL_ARGS.

If DIP_OVL_ASSIGN is defined, the following code will be executed by dip_ovl.h instead of the code shown above:

```
DIP_OVL_ASSIGN DIP_FUNC(DIP_OVL_FUNC, ovlDataType's extension) DIP_OVL_ARGS;
```

Note that to actually perform an assignment the "=" operator must be included in the definition of DIP_OVL_ASSIGN itself. DIP_OVL_BINARY_ASSIGN overrides DIP_OVL_ASSIGN if ovlDataType is one of the binary data types.

SEE ALSO

DIPlib's data types
DataTypeGetInfo, tpi.h

PaintBox

Paint a box

SYNOPSIS

```
#include "dip_paint.h"
dip_Error dip_PaintBox ( im, length, orign, amplitude )
```

DATA TYPES

binary, integer, float, complex

FUNCTION

Paints an box object in the image by replacing the values of the pixels in im that lie within the box (as specified by length and origin) with amplitude, and leaving the other pixel values untouched.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{\mathtt{-}}\mathtt{Image}$	im	Image
dip_FloatArray	length	Length array
dip_FloatArray	origin	Origin array
dip_float	amplitude	Pixel value of the painted ellips

SEE ALSO

PaintEllipsoid, PaintDiamond

PaintDiamond

Paint a diamond-shaped object

SYNOPSIS

```
#include "dip_paint.h"
dip_Error dip_PaintDiamond ( im, length, orign, amplitude )
```

DATA TYPES

binary, integer, float, complex

FUNCTION

Paints a diamond-shaped object in the image by replacing the values of the pixels in im that lie within the diamond (as specified by length and origin) with amplitude, and leaving the other pixel values untouched.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{\mathtt{-}}\mathtt{Image}$	im	Image
dip_FloatArray	length	Length array
dip_FloatArray	origin	Origin array
dip_float	amplitude	Pixel value of the painted ellips

SEE ALSO

PaintEllipsoid, PaintBox

PaintEllipsoid

Paint an ellipsoid

SYNOPSIS

```
#include "dip_paint.h"
dip_Error dip_PaintEllipsoid ( im, radius, orign, scale, amplitude )
```

DATA TYPES

binary, integer, float, complex

FUNCTION

Paints an elliptical object in the image by replacing the values of the pixels in im that lie within the ellips (as specified by diameter and origin) with amplitude, and leaving the other pixel values untouched.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{\mathtt{-}}\mathtt{Image}$	im	Image
dip_FloatArray	radius	Diameter array
dip_FloatArray	origin	Origin array
dip_float	amplitude	Pixel value of the painted ellips

SEE ALSO

PaintDiamond, PaintBox

PairCorrelation

Compute the pair correlation function

SYNOPSIS

```
#include "dip_analysis.h"
dip_Error dip_PairCorrelation ( object, mask, dist, probes, length, sampling,
covariance )
```

DATA TYPES

binary, integer

FUNCTION

This function computes the pair correlation function of the different phases in object. If object is a binary image, the image is a regarded as a two phase image. In case object is of the integer type, the image is regarded as a labeled image, with each integer value encoding a phase. Optionally a mask image can be provided to select which pixels in object should be used to compute the pair correlation. The probes variable specifies how many random point pairs should be drawn to compute the function. Length specifies the maximum correlation length. The correlation function can be computed using a random (DIP_CORRELATION_ESTIMATOR_RANDOM) or grid method (DIP_CORRELATION_ESTIMATOR_GRID), as specified by sampling. Finally covariance determines whether only the correlations (DIP_FALSE) or the covarianances (DIP_TRUE) have to be computed.

ARGUMENTS

Data type	Name	Description
dip_Image	object	Object image
dip_Image	mask	Mask image
dip_Distribution	dist	Ouput distribution
dip_int	probes	Number of probes
dip_int	length	Maximum correlation Length
dipf_CorrelationEstimator	sampling	Samplings method
dip_Boolean	covariance	Compute covariance

SEE ALSO

ChordLength, ProbabilisticPairCorrelation

PathOpening

Morphological filter

SYNOPSIS

```
#include "dip_morphology.h"
dip_Error PathOpening ( grey, mask, out, length, closing, constrained )
```

DATA TYPES

binary, integer, float

FUNCTION

This function applies <code>DirectedPathOpening</code> in all possible directions and takes the supremum of all results. That is, it is the supremum of all possible openings with approximately linear structuring elements of length <code>length</code>.

The number of times that DirectedPathOpening is applied is given by ((3^D)-1)/2, with D the number of image dimensions. For example, in 2D there are 4 possible values for param: [length,0], [0,length], [length,length] and [length,-length] (note that, for example, [-length,0] produces the same result as [length,0]).

See DirectedPathOpening for more information.

ARGUMENTS

Data type	Name	Description
dip_Image	grey	Grey-value input image
dip_Image	mask	Mask image for ROI processing
dip_Image	out	Output image
dip_int	length	Length of structuring element (number of pixels)
dip_Boolean	closing	DIP_FALSE for path opening, DIP_TRUE for path closing
dip_Boolean	constrained	DIP_TRUE for constrained paths, DIP_FALSE for the original path
		opening algorithm

SEE ALSO

DirectedPathOpening, Opening, Closing, AreaOpening

Percentile

statistics function

SYNOPSIS

dip_Error dip_Percentile (in, mask, out, percentile, ps)

DATA TYPES

binary, integer, float, complex

FUNCTION

Calculates the perc percentile of the pixel values over all those dimensions which are specified by ps.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
$\mathtt{dip}_{-}\mathtt{Image}$	mask (0)	Mask
$dip_{-}Image$	out	Output
dip_float	perc	Percentile
dip_BooleanArray	ps (0)	Dimensions to project

SEE ALSO

From images to scalars

 ${\tt Sum,\ Mean,\ Variance,\ Standard Deviation,\ Mean Modulus,\ Sum Modulus,\ Mean Square Modulus,\ Maximum,\ Minimum,\ Median}}$

PercentileFilter

Rank-order filter

SYNOPSIS

#include "dip_morphology.h"
dip_Error dip_PercentileFilter (in, out, se, boundary, param, shape, percentile)

DATA TYPES

integer, **float**

FUNCTION

Rank-order or percentile filter with different filter shapes.

Only the rectangular, elliptic and diamond filter shapes are supported (DIP_FLT_SHAPE_RECTANGULAR, DIP_FLT_SHAPE_ELLIPTIC and DIP_FLT_SHAPE_DIAMOND). Other filter shapes can be implemented by setting shape to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, and passing a binary image in se. The "on" pixels define the shape of the filter window. Other values of shape are illegal.

If shape is not equal to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, se can be set to zero. When shape is set to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, param is ignored, and can be set to zero.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input image
dip_Image	out	Output image
dip_Image	se	Custom filter shape (binary)
dip_BoundaryArray	boundary	Boundary conditions
dip_FloatArray	param	Filter parameters
dip_FilterShape	shape	Filter shape
dip_float	percentile	Percentile (%)

The enumerator dip_FilterShape contains the following constants:

Name	Description
DIP_FLT_SHAPE_DEFAULT	Default filter window, same as
	DIP_FLT_SHAPE_RECTANGULAR
DIP_FLT_SHAPE_RECTANGULAR	Rectangular filter window, can be even in size
DIP_FLT_SHAPE_ELLIPTIC	Elliptic filter window, always odd in size
DIP_FLT_SHAPE_DIAMOND	Diamond-shaped filter window, always odd in size
DIP_FLT_SHAPE_PARABOLIC	Parabolic filter window (morphology only)
DIP_FLT_SHAPE_DISCRETE_LINE	Rotated line structuring element (morphology only)
DIP_FLT_SHAPE_INTERPOLATED_LINE	Rotated line structuring element, through interpolation
	(morphology only)
DIP_FLT_SHAPE_PERIODIC_LINE	(not implemented)
DIP_FLT_SHAPE_STRUCTURING_ELEMENT	Use se as filter window, can be any size

SEE ALSO

 ${\tt MedianFilter}$

Phase

Arithmetic function

SYNOPSIS

dip_Error dip_Phase (in, out)

DATA TYPES

binary, integer, float, $\mathbf{complex}$

FUNCTION

Computes the phase of the input image values, and outputs a float typed image.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	in	Input
$dip_{-}Image$	out	Output

SEE ALSO

Modulus, Real, Imaginary

${\tt Physical Dimensions Copy}$

Copy a Physical Dimensions

SYNOPSIS

dip_Error dip_PhysicalDimensionsCopy (newPhysDims, src, resources)

FUNCTION

This function makes a copy of a Physical Dimensions data structure.

ARGUMENTS

Data type	Name	Description
dip_PhysicalDimensions *	newPhysDims	New Physical Dimensions data structure
dip_PhysicalDimensions	src	source data structure
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

Physical Dimensions New, Physical Dimensions Free, Physical Dimensions Is I so tropic

${\tt Physical Dimensions Free}$

Free a Physical Dimensions data structure

SYNOPSIS

dip_Error dip_PhysicalDimensionsFree (physDims)

FUNCTION

This function free the Physical Dimensions ${\tt physDims}$ structure.

ARGUMENTS

Data type	Name	Description
dip_PhysicalDimensions *	physDims	Physical Dimensions data structure

SEE ALSO

Physical Dimensions New, Physical Dimensions Copy, Physical Dimensions Is I so tropic

PhysicalDimensionsIsIsotropic

Checks if the Physical Dimensions are isotropic

SYNOPSIS

dip_Error dip_PhysicalDimensionsIsIsotropic (physDims, verdict)

FUNCTION

This function checks whether the physical dimensions physDims are isotropic by checking that all its pixel dimensions and dimension units are equal to each other. If verdict is not zero, the result (DIP_TRUE or DIP_FALSE) is stored in verdict, otherwise an error is returned in case the verification fails.

ARGUMENTS

Data type	Name	Description
dip_PhysicalDimensions	physDims	Physical Dimensions data structure
dip_Boolean *	verdict	Verdict of the test

SEE ALSO

PhysicalDimensionsNew, PhysicalDimensionsFree, PhysicalDimensionsCopy

PhysicalDimensionsNew

Allocates a new Physical Dimensions structure

SYNOPSIS

dip_Error dip_PhysicalDimensionsNew (newPhysDims, dimensionality, dims, orig, dimUnit, intensity, offset, intensUnit, resources)

FUNCTION

This function allocates a new Physical Dimensions structure.

A physical dimensions structure contains information about the physical dimensions of the data (of dimensionality dimension) in an image. It describes the position (orig) and size (dims) of a pixel in world coordinates and physical units (dimUnits), as well as the scaling (intensity) and offset (offset) of the pixel intensity in physical units (intensUnit).

Note that the initial values assigned by this function assume an isotropic pixel size. These values can be changed within the structure generated if this is not the case.

ARGUMENTS

Data type	Name	Description
dip_PhysicalDimensions *	newPhysDims	Pointer to a new Physical Dimensions data
		structure
dip_int	dimensionality	Dimensionality of the image
dip_float	dims	Initial value for dimensions along all axes
dip_float	orig	Initial value for origin along all axes
char *	dimUnit	Initial value for dimensionUnits along all axes
dip_float	intensity	Initial value for intensity
dip_float	offset	Initial value for offset
char *	intensUnit	Initial value for intensityUnit
dip_Resources	resources	Resources tracking structure. See
		ResourcesNew

The structure dip_PhysicalDimensions contains the following elements:

Data type	Name	Description
dip_FloatArray	dimensions	Dimensions of a pixel along each grid axis
dip_FloatArray	origin	Coordinates of the origin in physical units
dip_StringArray	dimensionUnits	Units for dimensions and origin
dip_float	intensity	Intensity scaling in physical units
dip_float	offset	Offset for intensity in physical units
dip_String	intensityUnit	Units for intensity and offset
dip_Resources	resources	Resource tracking; all elements within this structure are
		tracked here

SEE ALSO

Physical Dimensions Free, Physical Dimensions Copy, Physical Dimensions Is I so tropic

PixelHeapFree

Destroy heap structure

SYNOPSIS

```
#include "dip_pixelqueue.h"
dip_Error dip_PixelHeapFree ( heap )
```

FUNCTION

Frees all data associated to heap and sets heap to 0.

ARGUMENTS

Data type	Name	Description
dip_PixelHeap *	heap	The heap structure

SEE ALSO

 $\label{thm:pixelHeapNew} PixelHeapNew, PixelHeapPush, PixelHeapP$

PixelHeapIsEmpty

Query heap

SYNOPSIS

```
#include "dip_pixelqueue.h"
dip_Error dip_PixelHeapIsEmpty ( heap, result )
```

FUNCTION

Checks to see if there are any items on the heap. See PixelHeapNew for information on the heap data structure.

ARGUMENTS

Data type	Name	Description	
dip_PixelHeap	heap	The heap structure	
dip_Boolean *	result	Set to true if there are no items in the heap	

SEE ALSO

 $\label{thm:pixelHeapNew} PixelHeapNew, PixelHeapPree, PixelHeapPush, PixelHeapPop$

PixelHeapNew

Create a new heap structure

SYNOPSIS

```
#include "dip_pixelqueue.h"
dip_Error dip_PixelHeapNew ( heap, ndims, blocksize, order, resources )
```

FUNCTION

This function allocates space for a new dip_PixelHeap structure. Memory allocated is tracked in resources. The heap is dimensioned to hold pixels from an ndims-dimensional image, and initially enough space is allocated for blocksize elements. The heap will be expanded as necessary when used.

The heap stores the coordinates, the value and the pointer to a pixel in an image. Note that the value does not need to equal the data pointed to by the pointer. ndims can be set to zero, in which case no coordinates are stored; this does not affect the function of the value and the pointer.

A heap is a priority queue data structure. Just like a queue, items can be added (pushed) and subtracted (popped). However, in the priority queue the item popped is always the higherst priority one: either the one with the highest-valued item (order is DIP_GVSO_HIGH_FIRST) or lowest-valued item (order is DIP_GVSO_LOW_FIRST). However, identically-valued items are stored on the heap in unpredictable order. If this order is important (such as for the GrowRegions algorithm with integer-valued pixels, use a dip_StablePixelHeap instead. See StablePixelHeapNew for information on the stable heap structure.

ARGUMENTS

Data type	Name	Description
dip_PixelHeap *	heap	The newly allocated heap structure
dip_int	ndims	Image dimensionality
dip_int	blocksize	Size of each allocation block
dipf_GreyValueSortOrder	order	Determines the heap's sort order
dip_Resources	resources	Resources tracking structure. See ResourcesNew

The ${\tt dipf_GreyValueSortOrder}$ enumeration consists of the following values:

Name	Description	
DIP_GVSO_HIGH_FIRST	Process the pixels from high grey-value to low grey-value.	
DIP_GVSO_LOW_FIRST	Process the pixels from low grey-value to high grey-value.	

IMPLEMENTATION

When the heap grows beyond its initial size, its capacity is doubled in size by reallocating the data blocks. However, when removing pixels from the queue, the heap is not shrunk. It is assumed that

the dip_PixelHeap structure will be destroyed as soon as the algorithm using it terminates. Reducing the memory footprint of the heap therefore has no benefit.

SEE ALSO

 ${\tt StablePixelHeapNew,\ PixelQueueNew,\ PixelHeapFree,\ PixelHeapPush,\ PixelHeapPop,\ PixelHeapIsEmpty}$

PixelHeapPop

Pop item onto heap

SYNOPSIS

```
#include "dip_pixelqueue.h"
dip_Error dip_PixelHeapPop ( heap, coords, pointer, value )
```

FUNCTION

Pops the next pixel from the heap. See PixelHeapNew for information on the heap data structure. coords is a pointer to an array of dip_ints, such as that obtained with dip_IntegerArray->array. It should have as many elements as the image dimensionality. If the stack was created with ndims set to 0, the coords pointer is ignored. coords, pointer and value can be NULL if you are not interested in either those values.

ARGUMENTS

Data type	Name	Description	
dip_PixelHeap	heap	The heap structure	
dip_int *	coords	Receives the coordinates of the popped item	
void **	pointer	Receives the pointer of the popped item	
dip_sfloat *	value	Receives the value of the popped item	

SEE ALSO

PixelHeapNew, StablePixelHeapNew, PixelQueueNew, PixelHeapFree, PixelHeapPush, PixelHeapIsEmpty

PixelHeapPush

Push item onto heap

SYNOPSIS

```
#include "dip_pixelqueue.h"
dip_Error dip_PixelHeapPush ( heap, coords, pointer, value )
```

FUNCTION

Pushes a pixel onto the heap. See PixelHeapNew for information on the heap data structure. All 3 values coords, pointer and value are stored, except if the heap was created with ndims set to 0, in which case the coords pointer is ignored.

coords is a pointer to an array of dip_ints, such as that obtained with dip_IntegerArray->array. It should have as many elements as the image dimensionality. pointer is a pointer to any memory location, and value is the value to be used when sorting.

ARGUMENTS

Data type	Name	Description	
dip_PixelHeap	heap	The heap structure	
dip_int *	coords	Coordinates to be pushed	
void *	pointer	Pointer to be pushed	
dip_sfloat	value	Value to be pushed	

SEE ALSO

 $\label{thm:pixelHeapNew} PixelHeapNew, PixelQueueNew, PixelHeapFree, PixelHeapPop, PixelHeapIsEmpty$

PixelQueueFree

Destroy queue structure

SYNOPSIS

```
#include "dip_pixelqueue.h"
dip_Error dip_PixelQueueFree ( queue )
```

FUNCTION

Frees all data associated to queue and sets queue to 0.

ARGUMENTS

Data type	Name	Description	
dip_PixelQueue *	queue	The queue structure	

SEE ALSO

 $\label{thm:pixelQueue} Pixel Queue Push, Pixel Queue Push, Pixel Queue Pop, Pixel Queue Is Empty$

PixelQueueIsEmpty

Query queue

SYNOPSIS

```
#include "dip_pixelqueue.h"
dip_Error dip_PixelQueueIsEmpty ( queue, result )
```

FUNCTION

Checks to see if there are any items on the queue. See PixelQueueNew for information on the queue data structure.

ARGUMENTS

Data type	Name	Description
dip_PixelQueue	queue	The queue structure
dip_Boolean *	result	Set to true if there are no items in the queue

SEE ALSO

 ${\tt PixelQueueNew,\,PixelHeapNew,\,StablePixelHeapNew,\,PixelQueueFree,\,PixelQueuePush,\,PixelQueuePop}$

PixelQueueNew

Create a new queue structure

SYNOPSIS

```
#include "dip_pixelqueue.h"
dip_Error dip_PixelQueueNew ( queue, ndims, blocksize, resources )
```

FUNCTION

This function allocates space for a new dip_PixelQueue structure. Memory allocated is tracked in resources. The queue is dimensioned to hold pixels from an ndims-dimensional image, and initially enough space is allocated for blocksize elements. The queue will be expanded as necessary when used.

The queue stores the coordinates, and the pointer to a pixel in an image. ndims can be set to zero, in which case no coordinates are stored; this does not affect the function of the pointer.

A queue is a data structure to which items can be added (pushed) to the back, and subtracted (popped) from the front (FIFO - First In, First Out).

ARGUMENTS

Data type	Name	Description
dip_PixelQueue *	queue	The newly allocated queue structure
$\mathtt{dip}_{ ext{-}}\mathtt{int}$	ndims	Image dimensionality
dip_int	blocksize	Size of each allocation block
dip_Resources	resources	Resources tracking structure. See ResourcesNew

IMPLEMENTATION

The queue is stored in an array whose size can be incressed at will. This is accomplished by a linked list of blocks, each one holds blocksize elements. When more space is needed, a new block is simply allocated. No data needs to be moved as would be necessary when using realloc to change the size of the array. Blocks on the front of the queue that become empty are freed.

SEE ALSO

PixelHeapNew, StablePixelHeapNew, PixelQueueFree, PixelQueuePush, PixelQueuePop, PixelQueueIsEmpty

PixelQueuePop

Pop item from queue

SYNOPSIS

```
#include "dip_pixelqueue.h"
dip_Error dip_PixelQueuePop ( queue, coords, pointer, newiteration )
```

FUNCTION

Pops the next pixel from the queue. See PixelQueueNew for information on the queue data structure. coords is a pointer to an array of dip_ints, such as that obtained with dip_IntegerArray->array. It should have as many elements as the image dimensionality. If the queue was created with ndims set to 0, the coords pointer is ignored. coords and pointer can be NULL if you are not interested in either those values.

newiteration, when not NULL, will be set to DIP_TRUE if the pixel being popped is the first one in an iteration, or DIP_FALSE otherwise. When a new iteration starts, all pixels pushed onto the queue afterwards belong to the next iteration. This is useful in routines that use the queue for propagating boundaries, such as <code>GrowRegions</code>. First all boundary pixels are pushed onto the queue. The first iteration will need to process only these pixels, while at the same time push new pixels onto the queue for the second iteration. So after pushing all the initial boundary pixels onto the queue, the first iteration is started by popping the first pixel. All pixels pushed while processing this and the rest of the pixels will be pushed behind the "new iteration" marker. When the first pixel after this marker is popped, the <code>newiteration</code> boolean is set, so the program knows that the second iteration is starting. Also, the "new iteration" marker is moved to the end of the queue, so that pixels pushed subsequently will belong to iteration number 3.

ARGUMENTS

Data type	Name	Description
dip_PixelQueue	queue	The queue structure
dip_int *	coords	Receives the coordinates of the popped item
void **	pointer	Receives the pointer of the popped item
dip_Boolean *	newiteration	Set to true when the first item from an iteration is popped

SEE ALSO

PixelQueueNew, PixelHeapNew, StablePixelHeapNew, PixelQueueFree, PixelQueuePush, PixelQueueIsEmpty

PixelQueuePush

Push item onto queue

SYNOPSIS

```
#include "dip_pixelqueue.h"
dip_Error dip_PixelQueuePush ( queue, coords, pointer )
```

FUNCTION

Pushes a pixel onto the queue. See PixelQueueNew for information on the queue data structure. Both coords and pointer are stored, except if the stack was created with ndims set to 0, in which case the coords values are ignored.

coords is a pointer to an array of dip_ints, such as that obtained with dip_IntegerArray->array. It should have as many elements as the image dimensionality. pointer is a pointer to any memory location.

ARGUMENTS

Data type	Name	Description
dip_PixelQueue	queue	The queue structure
dip_int *	coords	Coordinates to be pushed
void *	pointer	Pointer to be pushed

SEE ALSO

 $\label{thm:pixelQueue} PixelQueueNew, PixelQueueFree, PixelQueuePop, PixelQueueIsEmpty$

PixelTableAddRun

Add a new run to a pixel table

SYNOPSIS

```
#include "dip_pixel_table.h"
dip_Error dip_PixelTableAddRun ( table, coordinate, length )
```

FUNCTION

Adds a new run to a pixel table. The new run is appended to the existing runs in the pixel table.

ARGUMENTS

Data type	Name	Description
dip_PixelTable	table	Pixel table
dip_IntegerArray	coordinate	Coordinate of the run
dip_int	length	Length of the run

SEE ALSO

Description of DIPlib's pixel tables

PixelTableNew, PixelTableSetRun, PixelTableGetRun, PixelTableGetRuns, PixelTableGetDimensionality, PixelTableGetDimensions, PixelTableGetOrigin, PixelTableGetSize, PixelTableGetPixelCount, PixelTableGetOffsetAndLength

PixelTableCreateFilter

Create a pixel table from a filter shape

SYNOPSIS

```
#include "dip_pixel_table.h"
dip_Error dip_PixelTableCreateFilter ( table, param, shape, se, resources )
```

FUNCTION

This function allocates and creates a new pixel table data structure. The shape and dimensionality of the pixel table is specified by the param, shape and se parameters.

Only the rectangular, elliptic and diamond filter shapes are supported (DIP_FLT_SHAPE_RECTANGULAR, DIP_FLT_SHAPE_ELLIPTIC and DIP_FLT_SHAPE_DIAMOND). Other filter shapes can be implemented by setting shape to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, and passing a binary image in se. The "on" pixels define the shape of the filter window. Other values of shape are illegal.

If shape is not equal to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, se can be set to zero. When shape is set to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, param is ignored, and can be set to zero.

ARGUMENTS

Data type	Name	Description
dip_PixelTable *	table	Pointer to a pixel table
dip_FloatArray	param	Filter size
dip_FilterShape	shape	Filter shape
dip_Image	se	Structuring element
dip_Resources	resources	Resources tracking structure. See ResourcesNew

The enumerator dip_FilterShape contains the following constants:

Name	Description
DIP_FLT_SHAPE_DEFAULT	Default filter window, same as
	DIP_FLT_SHAPE_RECTANGULAR
DIP_FLT_SHAPE_RECTANGULAR	Rectangular filter window, can be even in size
DIP_FLT_SHAPE_ELLIPTIC	Elliptic filter window, always odd in size
DIP_FLT_SHAPE_DIAMOND	Diamond-shaped filter window, always odd in size
DIP_FLT_SHAPE_PARABOLIC	Parabolic filter window (morphology only)
DIP_FLT_SHAPE_DISCRETE_LINE	Rotated line structuring element (morphology only)
DIP_FLT_SHAPE_INTERPOLATED_LINE	Rotated line structuring element, through interpolation
	(morphology only)
DIP_FLT_SHAPE_PERIODIC_LINE	(not implemented)
DIP_FLT_SHAPE_STRUCTURING_ELEMENT	Use se as filter window, can be any size

SEE ALSO

Description of DIPlib's pixel tables
BinaryImageToPixelTable, GreyValuesInPixelTable, PixelTableToBinaryImage

PixelTableFrameWork

FrameWork for PixelTable filters

SYNOPSIS

```
#include "dip_tprunlength.h"
dip_Error dip_PixelTableFrameWork ( in, out, boundary, process, table )
```

FUNCTION

This function provides a framework for filters that code the shape of their filter in a pixel table (run lengths). See SeparableFrameWork for details.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input image
dip_Image	out	Output image
dip_BoundaryArray	boundary	Boundary conditions
dip_FrameWorkProcess	process	Process
dip_PixelTable	table	Pixel table

SEE ALSO

SeparableFrameWork

PixelTableGetDimensionality

Get the dimensionality of a pixel table

SYNOPSIS

```
#include "dip_pixel_table.h"
dip_Error dip_PixelTableGetDimensionality ( table, dimension )
```

FUNCTION

Gets the dimensionality of the binary object that is encoded by the pixel table table.

ARGUMENTS

Data type	Name	Description
dip_PixelTable	table	pixel table
dip_int *	dimension	pointer to a dimensionality variable

SEE ALSO

Description of DIPlib's pixel tables

PixelTableNew, PixelTableSetRun, PixelTableGetRun, PixelTableAddRun, PixelTableGetRuns, PixelTableGetDimensions, PixelTableGetOrigin, PixelTableGetSize, PixelTableGetPixelCount, PixelTableGetOffsetAndLength

PixelTableGetDimensions

Get the dimemsions of a pixel table

SYNOPSIS

```
#include "dip_pixel_table.h"
dip_Error dip_PixelTableGetDimensions ( table, dimensions, resources )
```

FUNCTION

This functions gets the dimensions of the bounding box of the binary object that is encoded by the pixel table table.

ARGUMENTS

Data type	Name	Description
dip_PixelTable *	table	Pixel table
dip_IntegerArray *	dimensions	Pointer to a dimensions array
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

Description of DIPlib's pixel tables

PixelTableNew, PixelTableSetRun, PixelTableGetRun, PixelTableAddRun, PixelTableGetRuns, PixelTableGetDimensionality, PixelTableGetOrigin, PixelTableGetSize, PixelTableGetPixelCount, PixelTableGetOffsetAndLength

PixelTableGetOffsetAndLength

Converts the pixel table's runs

SYNOPSIS

```
#include "dip_pixel_table.h"
dip_Error dip_PixelTableGetOffsetAndLength ( table, stride, offset, length,
resources )
```

FUNCTION

This functions converts the linked-list of runs in the pixel table to two arrays of offsets and lengths. The length of these arrays equals the number of runs. The offsets are calculated by multiplying each coordinate of a run with the stride of that dimension. This function is useful when an image needs to be filtered with a filter that is encoded by a pixel table. Before processing the image. See also PixelTableFrameWork.

ARGUMENTS

Data type	Name	Description
dip_PixelTable	table	Pixel table
dip_IntegerArray	stride	Stride array
dip_IntegerArray *	offset	Pointer to offset array
dip_IntegerArray *	length	Pointer to length array
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

Description of DIPlib's pixel tables

 $\label{lem:pixelTableNew} PixelTableGetOffsetAndLength, PixelTableCreateFilter, PixelTableFrameWork$

PixelTableGetOrigin

Get the origin of the pixel table

SYNOPSIS

```
#include "dip_pixel_table.h"
dip_Error dip_PixelTableGetOrigin ( table, origin, resources )
```

FUNCTION

This function gets the origin of the pixel table table. All coordinates of the pixel table runs are defined relative to this origin. The origin is the pixel with coordinates (0,0), relative to the top left pixel.

Unless PixelTableShiftOrigin has been called, the origin is equal to the bounding box divided by 2 (integer divistion), meaning it is the middle pixel if the bounding box is odd in size, or the pixel to the right of the middle if it is even in size.

ARGUMENTS

Data type	Name	Description
dip_PixelTable *	table	Pixel table
dip_IntegerArray *	origin	Pointer to a origin array
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

Description of DIPlib's pixel tables

PixelTableNew, PixelTableSetRun, PixelTableGetRun, PixelTableAddRun, PixelTableGetRuns, PixelTableGetDimensionality, PixelTableGetDimensions, PixelTableGetSize, PixelTableGetPixelCount, PixelTableGetOffsetAndLength

PixelTableGetPixelCount

Get the number of pixels encoded in the pixel table

SYNOPSIS

```
#include "dip_pixel_table.h"
dip_Error dip_PixelTableGetPixelCount ( table, count )
```

FUNCTION

Gets the total number of pixels of the binary object that is encoded by the Pixel table table.

ARGUMENTS

Data type	Name	Description
dip_PixelTable	table	Pixel table
dip_int *	count	pointer to count

SEE ALSO

Description of DIPlib's pixel tables

PixelTableNew, PixelTableSetRun, PixelTableGetRun, PixelTableAddRun, PixelTableGetRuns, PixelTableGetDimensionality, PixelTableGetDimensions, PixelTableGetOrigin, PixelTableGetSize, PixelTableGetOffsetAndLength

PixelTableGetRun

Get the contents of a pixel table run

SYNOPSIS

```
#include "dip_pixel_table.h"
dip_Error dip_PixelTableGetRun ( table, run, coordinate, length )
```

FUNCTION

This functions get the the coordinate and length parameters of the runth run of the pixel table table.

ARGUMENTS

Data type	Name	Description
dip_PixelTable	table	Pixel table
dip_int	run	The run to be initialised
dip_IntegerArray	coordinate	Coordinate of the run
dip_int *	length	Length of the run

SEE ALSO

Description of DIPlib's pixel tables

PixelTableNew, PixelTableSetRun, PixelTableAddRun, PixelTableGetRuns, PixelTableGetDimensionality, PixelTableGetDimensions, PixelTableGetOrigin, PixelTableGetSize, PixelTableGetPixelCount, PixelTableGetOffsetAndLength

PixelTableGetRuns

Get the number of runs in a pixel table

SYNOPSIS

```
#include "dip_pixel_table.h"
dip_Error dip_PixelTableGetRuns ( table, numberOfRuns )
```

FUNCTION

Gets the number of runs in a pixel table.

ARGUMENTS

Data type	Name	Description
dip_PixelTable	table	Pixel table
dip_int *	numberOfRuns	Point to the NumberOfRuns

SEE ALSO

Description of DIPlib's pixel tables

PixelTableNew, PixelTableSetRun, PixelTableGetRun, PixelTableAddRun, PixelTableGetDimensionality, PixelTableGetDimensions, PixelTableGetOrigin, PixelTableGetSize, PixelTableGetPixelCount, PixelTableGetOffsetAndLength

PixelTableGetSize

The number of pixels in the pixel table's bounding box

SYNOPSIS

```
#include "dip_pixel_table.h"
dip_Error dip_PixelTableGetSize ( table, size )
```

FUNCTION

Gets the number of pixels in the bounding box of the pixel table table.

ARGUMENTS

Data type	Name	Description
dip_PixelTable	table	Pixel table
dip_int *	size	pointer to size

SEE ALSO

Description of DIPlib's pixel tables

PixelTableNew, PixelTableSetRun, PixelTableGetRun, PixelTableAddRun, PixelTableGetRuns, PixelTableGetDimensionality, PixelTableGetDimensions, PixelTableGetOrigin, PixelTableGetPixelCount, PixelTableGetOffsetAndLength

PixelTableNew

Allocate a new pixel table

SYNOPSIS

```
#include "dip_pixel_table.h"
dip_Error dip_PixelTableNew ( table, size, runs, resource s)
```

FUNCTION

Allocates a new pixel table table. The size array specifies the dimensionality of the coordinates in each run, and the sizes of the bounding box of the pixel table. The runs parameter specifies the number of runs in the pixel table.

ARGUMENTS

Data type	Name	Description
dip_PixelTable *	table	Pixel table
dip_IntegerArray	size	Size
dip_int	runs	Number of pixel table runs
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

Description of DIPlib's pixel tables

PixelTableSetRun, PixelTableGetRun, PixelTableAddRun, PixelTableGetRuns, PixelTableGetDimensionality, PixelTableGetDimensions, PixelTableGetOrigin, PixelTableGetSize, PixelTableGetPixelCount, PixelTableGetOffsetAndLength

PixelTableSetRun

Initialises a pixel table run

SYNOPSIS

```
#include "dip_pixel_table.h"
dip_Error dip_PixelTableSetRun ( table, run, coordinate, length )
```

FUNCTION

This function initialises the runth run of the pixel table table, by setting the run's coordinate to coordinate and its length to length. The pixel table must at least consist of run number of runs ans has to be allocated (using PixelTableNew).

ARGUMENTS

Data type	Name	Description
dip_PixelTable	table	Pixel table
dip_int	run	The run to be initialised
dip_IntegerArray	coordinate	Coordinate of the run
dip_int	length	Length of the run

SEE ALSO

Description of DIPlib's pixel tables

PixelTableNew, PixelTableGetRun, PixelTableAddRun, PixelTableGetRuns, PixelTableGetDimensionality, PixelTableGetDimensions, PixelTableGetOrigin, PixelTableGetSize, PixelTableGetPixelCount, PixelTableGetOffsetAndLength

PixelTableShiftOrigin

Changes the origin of the pixel table

SYNOPSIS

```
#include "dip_pixel_table.h"
dip_Error dip_PixelTableShiftOrigin ( table, offset )
```

FUNCTION

This function changes the origin of the pixel table table. By default, the origin is equal to the bounding box divided by 2 (integer divistion), meaning it is the middle pixel if the bounding box is odd in size, or the pixel to the right of the middle if it is even in size. After calling this function, the origin is equal to the previous origin plus the offset.

ARGUMENTS

Data type	Name	Description
dip_PixelTable	table	Pixel table
dip_IntegerArray	offset	An offset array, to be added to the origin

SEE ALSO

Description of DIPlib's pixel tables

PixelTableNew, PixelTableSetRun, PixelTableGetRun, PixelTableAddRun, PixelTableGetRuns, PixelTableGetDimensionality, PixelTableGetDimensions, PixelTableGetSize, PixelTableGetPixelCount, PixelTableGetOffsetAndLength

PixelTableToBinaryImage

Convert a pixel table to a binary image

SYNOPSIS

```
#include "dip_pixel_table.h"
dip_Error dip_PixelTableToBinaryImage ( table, im )
```

DATA TYPES

binary

FUNCTION

Converts the pixel table table to a binary image. The size of the image is set to the size of the bounding box of the pixel table.

ARGUMENTS

Data type	Name	Description
dip_PixelTable *	table	Pixel table
dip_Image	im	Binary image

SEE ALSO

Description of DIPlib's pixel tables

 ${\tt BinaryImageToPixelTable}, {\tt PixelTableCreateFilter}$

PoissonNoise

Generate an image disturbed by Poisson noise

SYNOPSIS

```
#include "dip_noise.h"
dip_Error dip_PoissonNoise ( in, out, conversion, random )
```

DATA TYPES

integer, float

FUNCTION

Generate a Poisson noise disturbed image. The intensities of the input image divided by the conversion variable are used as mean value for the Poisson distribution. The conversion factor can be used to relate the pixel values with the number of counts. For example, the simulate a photon limited image acquired by a CCD camera, the conversion factor specifies the relation between the number of photons recorded and the pixel value it is represented by.

See PoissonRandomVariable for more information on the random number generator.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	out	Output
dip_float	conversion	Conversion factor (photon/ADU)
dip_Random *	random	random

EXAMPLE

Get a Poisson disturbed image as follows:

```
dip_Image in, out;
dip_float conversion;
dip_Random random;

conversion = 2.0;
DIPXJ( dip_RandomSeed( &random, 0 ));
DIPXJ( dip_PoissonNoise( in, out, conversion, &random ));
```

SEE ALSO

 $\label{lem:poissonRandomVariable} Poisson Random Variable, Random Seed, Random Seed Vector, Uniform Noise, Gaussian Noise, Binary Noise$

PoissonRandomVariable

Poisson random variable generator

SYNOPSIS

```
#include "dip_noise.h"
dip_Error dip_PoissonRandomVariable ( random, mean, value )
```

FUNCTION

PoissonRandomVariable uses the algorithm as described in "Numerical Recipes in C, 2nd edition", section 7.3. For values of mean larger or equal to 32 the rejection method is used.

See RandomVariable for more information on the random number generator.

ARGUMENTS

Data type	Name	Description
dip_Random *	random	Pointer to a random value structure
dip_float	mean	Mean value for the distribution
dip_float *	value	Poisson distributed output value

EXAMPLE

Get a Poisson random variable as follows:

```
dip_Random random;
dip_float mean, value;

mean = 23.0;
DIPXJ( dip_RandomSeed( &random, 0 ));
DIPXJ( dip_PoissonRandomVariable( &random, mean, &value ));
```

LITERATURE

Press, W.H., Teukolsky, S.A., Vetterling, W.T., and Flannery, B.P. Numerical Recipes in C, 2nd edition, Cambridge University Press, Cambridge, 1992.

SEE ALSO

RandomVariable, RandomSeed, RandomSeedVector, UniformRandomVariable, GaussianRandomVariable, BinaryRandomVariable

ProbabilisticPairCorrelation

Compute the probabilistic pair correlation function

SYNOPSIS

#include "dip_analysis.h"
dip_Error dip_ProbabilisticPairCorrelation (phases, mask, dist, probes, length,
sampling, covariance)

DATA TYPES

float

FUNCTION

This function computes the probabilistic pair correlation function of the different phases in phases. Each image in the image array phases is treated as a separate phase. The function assumes, but does not check, that the values in these images are with the [0 1] range. Optionally a mask image can be provided to select which pixels in object should be used to compute the pair correlation. The probes variable specifies how many random point pairs should be drawn to compute the function. Length specifies the maximum correlation length. The correlation function can be computed using a random (DIP_CORRELATION_ESTIMATOR_RANDOM) or grid method (DIP_CORRELATION_ESTIMATOR_GRID), as specified by sampling. Finally covariance determines whether only the correlations (DIP_FALSE) or the covarianances (DIP_TRUE) have to be computed.

ARGUMENTS

Data type	Name	Description
dip_ImageArray	phases	Phase image array
dip_Image	mask	Mask image
dip_Distribution	dist	Ouput distribution
dip_int	probes	Number of probes
dip_int	length	Maximum chord length
dipf_CorrelationEstimator	sampling	Samplings method
dip_Boolean	covariance	Compute covariance

SEE ALSO

PairCorrelation, ChordLength

PseudoInverse

Image restoration filter

SYNOPSIS

```
#include "dip_restoration.h"
dip_Error dip_PseudoInverse ( in, psf, out, threshold, flags )
```

FUNCTION

This function performs a basic, very noise sensitive image restoration operation by inverse filtering the image with a clipped point spread function. Each frequency in the output for which the response of the PSF is smaller than threshold is set to zero.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input image
dip_Image	psf	Point spread function image
dip_Image	out	Output image
dip_float	threshold	Threshold value
dipf_Restoration	flags	Restoration flags

LITERATURE

G.M.P. van Kempen, *Image Restoration in FLuorescence Microscopy*, Ph.D. Thesis, Delft University of Technology, 1999

SEE ALSO

Wiener, TikhonovMiller

PutLine

Put a line in an image

SYNOPSIS

```
#include "dip_manipulation.h"
dip_Error dip_PutLine ( in, out, cor, dimension )
```

DATA TYPES

binary, integer, float, complex

FUNCTION

Put a line in an image. Put a line orthogonally in an image. The position of the line in the image is specified by the coordinates at which its left most pixel (cor) is be placed and on which dimension of the image, the dimension of the line maps (dimension). If in has a different type than out, it will be converted to the type of out.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	in	Input Line Image
$\mathtt{dip}_{-}\mathtt{Image}$	out	Output Image
dip_IntegerArray	cor	Coordinate in the image of the left most pixel of the line
dip_int	dimension	Dimension of the image on which the line's dimension maps

SEE ALSO

GetSlice, PutSlice, GetLine

PutSlice

Put a slice in an image

SYNOPSIS

```
#include "dip_manipulation.h"
dip_Error dip_PutSlice ( in, out, cor, dim1, dim2 )
```

DATA TYPES

binary, integer, float, complex

FUNCTION

Put a slice orthogonally in an image. The position of the slice in the image is specified by the coordinates at which its upper left corner (cor) should be placed and on which dimensions of the image, the dimensions of the slice map (dim1, dim2). If in has a different type than out, it will be converted to the type of out.

ARGUMENTS

Data type	Name	Description	
dip_Image	in	2D Input Image	
dip_Image	out	3D Output Image	
dip_IntegerArray	cor	Coordinate in out where the upper left corner of in is placed	
dip_int	dim1	Dimension of in on which the slice's first dimension maps	
dip_int	dim2	Dimension of in on which the slice's second dimensionmaps	

SEE ALSO

PutSlice, GetLine, PutLine

QuickSort

Sort a block of data

SYNOPSIS

```
#include "dip_sort.h"
dip_Error dip_QuickSort ( data, size, dataType )
```

FUNCTION

Sorts a block of data (of size size and data type dataType) using the quick sort algorithm.

ARGUMENTS

Data type	Name	Description
void *	data	Data
dip_int	size	Size
dip_DataType	dataType	Data type. See DIPlib's data types

SEE ALSO

General information about sorting

 ${\tt QuickSortIndices}, {\tt QuickSortIndices16}, {\tt Sort}, {\tt ImageSort}, {\tt SortIndices}, {\tt SortIndices16}, {\tt ImageSortIndices}$

QuickSortAnything

Sort data of any type

SYNOPSIS

```
#include "dip_sort.h"
dip_Error dip_QuickSortAnything ( data, size, compareFunction, swapFunction, tmpData
)
```

FUNCTION

Sorts a block of data (of size size) using the quick sort algorithm. This routine requires the user to write two functions in order to fully implement the sorting procedure. These are SortCompareFunction and SortSwapFunction.

ARGUMENTS

Data type	Name	Description
void *	data	Data
dip_int	size	Size
dip_SortCompareFunction	compareFunction	Function for comparing two data points
dip_SortSwapFunction	swapFunction	Function for swapping two data points, or
		copying one to the other
void *	tmpData	Pointer to a variable of the same type as the
		data, used as temporary space by some of the
		algorithms

SEE ALSO

General information about sorting

SortAnything, SortCompareFunction, SortSwapFunction

QuickSortIndices

Sort indices to a block of data

SYNOPSIS

```
#include "dip_sort.h"
dip_Error dip_QuickSortIndices ( data, indices, size, dataType )
```

FUNCTION

Sort a list of indices rather than the data itself using the quick sort algorithm.

ARGUMENTS

Data type	Name	Description
void *	data	Data
dip_sint32 *	indices	Indices
dip_int	size	Size
dip_DataType	dataType	Data type. See DIPlib's data types

SEE ALSO

General information about sorting

 ${\tt QuickSort, QuickSortIndices16, Sort, ImageSort, SortIndices, SortIndices16, ImageSortIndices}$

QuickSortIndices16

Sort indices to a block of data

SYNOPSIS

```
#include "dip_sort.h"
dip_Error dip_QuickSortIndices16 ( data, indices, size, dataType )
```

FUNCTION

Sorts a list of (16 bit) indices rather than the data itself using the quick sort algorithm.

ARGUMENTS

Data type	Name	Description
void *	data	Data
dip_sint16 *	indices	Indices
dip_int	size	Size
dip_DataType	dataType	Data type. See DIPlib's data types

SEE ALSO

General information about sorting

 ${\tt QuickSort, QuickSortIndices, Sort, ImageSort, SortIndices, SortIndices} 16, \\ {\tt ImageSortIndices}$

RadialMaximum

statistics function

SYNOPSIS

dip_Error dip_RadialMaximum (in, mask, out, ps, binSize, innerRadius, center)

DATA TYPES

binary, integer, float

FUNCTION

This function computes the radial projection of the maximum of the pixel intensities of in.

The radial projection is performed for the dimensions specified by ps. If the radial distance of a pixel to the center of the image is r, than the maximum of the intensities of all pixels with n * binSize <= r < (n + 1) * binSize is stored at position n in the radial dimension of out. The radial dimension is the first dimension to be processed (as specified by ps). If innerRadius is set to DIP_TRUE, the maximum radius that is projected is equal to the smallest dimension of the to be projected dimensions. Otherwise, the maximum radius is set equal to the diagonal length of the dimensions to be processed.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	in	Input
dip_Image	mask	Binary mask (or 0)
dip_Image	out	Output
dip_BooleanArray	ps	Dimensions to project (or 0)
dip_float	binSize	Size of radial bins
dip_Boolean	innerRadius	Maximum radius
dip_FloatArray	center	Coordinates of center (or 0)

SEE ALSO

RadialSum, RadialMean, RadialMinimum, Sum, Mean, Maximum, Minimum

RadialMean

statistics function

SYNOPSIS

dip_Error dip_RadialMean (in, mask, out, ps, binSize, innerRadius, center)

DATA TYPES

binary, integer, float

FUNCTION

This function computes the radial projection of the mean of the pixel intensities of in.

The radial projection is performed for the dimensions specified by ps. If the radial distance of a pixel to the center of the image is r, than the mean of the intensities of all pixels with n * binSize <= r < (n + 1) * binSize is stored at position <math>n in the radial dimension of out. The radial dimension is the first dimension to be processed (as specified by ps). If innerRadius is set to DIP_TRUE, the maximum radius that is projected is equal to the smallest dimension of the to be projected dimensions. Otherwise, the maximum radius is set equal to the diagonal length of the dimensions to be processed.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	mask	Binary mask (or 0)
dip_Image	out	Output
dip_BooleanArray	ps	Dimensions to project (or 0)
dip_float	binSize	Size of radial bins
dip_Boolean	innerRadius	Maximum radius
dip_FloatArray	center	Coordinates of center (or 0)

SEE ALSO

RadialSum, RadialMaximum, RadialMinimum, Sum, Mean, Maximum, Minimum

RadialMinimum

statistics function

SYNOPSIS

dip_Error dip_RadialMinimum (in, mask, out, ps, binSize, innerRadius, center)

DATA TYPES

binary, integer, float, complex

FUNCTION

This function computes the radial projection of the sum of the pixel intensities of in.

The radial projection is performed for the dimensions specified by ps. If the radial distance of a pixel to the center of the image is r, than the minimum of the intensities of all pixels with n * binSize <= r < (n + 1) * binSize is stored at position n in the radial dimension of out. The radial dimension is the first dimension to be processed (as specified by ps). If innerRadius is set to DIP_TRUE, the maximum radius that is projected is equal to the smallest dimension of the to be projected dimensions. Otherwise, the maximum radius is set equal to the diagonal length of the dimensions to be processed.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{\mathtt{-}}\mathtt{Image}$	in	Input
dip_Image	mask	Binary mask (or 0)
dip_Image	out	Output
dip_BooleanArray	ps	Dimensions to project (or 0)
dip_float	binSize	Size of radial bins
dip_Boolean	innerRadius	Maximum radius
dip_FloatArray	center	Coordinates of center (or 0)

SEE ALSO

RadialSum, RadialMean, RadialMaximum, Sum, Mean, Maximum, Minimum

RadialSum

statistics function

SYNOPSIS

dip_Error dip_RadialSum (in, mask, out, ps, binSize, innerRadius, center)

DATA TYPES

binary, integer, float

FUNCTION

This function computes the radial projection of the sum of the pixel intensities of in.

The radial projection is performed for the dimensions specified by ps. If the radial distance of a pixel to the center of the image is r, than the sum of the intensities of all pixels with n * binSize <= r < (n + 1) * binSize is stored at position <math>n in the radial dimension of out. The radial dimension is the first dimension to be processed (as specified by ps). If innerRadius is set to DIP_TRUE, the maximum radius that is projected is equal to the smallest dimension of the to be projected dimensions. Otherwise, the maximum radius is set equal to the diagonal length of the dimensions to be processed.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	mask	Binary mask (or 0)
dip_Image	out	Output
dip_BooleanArray	ps	Dimensions to project (or 0)
dip_float	binSize	Size of radial bins
dip_Boolean	innerRadius	Maximum radius
dip_FloatArray	center	Coordinates of center (or 0)

SEE ALSO

RadialMean, RadialMaximum, RadialMinimum, Sum, Mean, Maximum, Minimum

RandomSeed

Initialise random number generator

SYNOPSIS

```
#include "dip_noise.h"
dip_Error dip_RandomSeed ( random, seed )
```

FUNCTION

Initializes a dip_Random structure using a given seed value. If seed is 0, the default value of 5489 is used instead, since 0 produces a uniquely poor initialisation.

ARGUMENTS

Data type	Name	Description
dip_Random *	random	Pointer to a random value structure
dip_uint32	seed	Seed value

EXAMPLE

Initialize a dip_Random structure as follows:

```
dip_Random random;
dip_uint32 seed;
seed = 123758;
DIPXJ( dip_RandomSeed( &random, seed ));
```

SEE ALSO

 ${\tt RandomSeedVector}, \ {\tt RandomVariable}, \ {\tt UniformRandomVariable}, \ {\tt GaussianRandomVariable}, \ {\tt PoissonRandomVariable}, \ {\tt BinaryRandomVariable}$

RandomSeedVector

Initialise random number generator

SYNOPSIS

```
#include "dip_noise.h"
dip_Error RandomSeedVector ( random, seedvector )
```

FUNCTION

Initializes a dip_Random structure with a given seed value vector. The vector must have DIP_MT_STATE_SIZE (==624) values. This is an alternative to RandomSeed, which, by initialising with a 32-bit integer, only gives 4 billion different sequences. RandomSeedVector allows to initialise the whole status of the random number generator.

ARGUMENTS

Data type	Name	Description
dip_Random *	random	Pointer to a random value structure
dip_uint32[DIP_MT_STATE_SIZE]	seedvector	Seed value vector

SEE ALSO

RandomSeed, RandomVariable, UniformRandomVariable, GaussianRandomVariable, PoissonRandomVariable, BinaryRandomVariable

RandomVariable

Random number generator

SYNOPSIS

```
include "dip_noise.h"
dip_Error dip_RandomVariable ( random, value )
```

FUNCTION

Generates a random number between zero and one. The dip_Random structure must be initialized with the function RandomSeed. If the supplied dip_Random structure is not initialized, RandomVariable will initialize the dip_Random structure with the default seed. To guarantee the (psuedo) randomness between variables obtained with subsequent calls to RandomVariable (or to functions that use this function to obtain a random variable), a pointer to the same dip_Random structure has to supplied when calling RandomVariable.

The random number generator returns random numbers as 32-bit integers, which are normalised to to [0,1] range. If higher precision numbers are required, you can set random.highprecision to DIP_TRUE. This causes two random 32-bit values to be used for each floating point output value, doubling the precision of the output. There is no need to re-initialise the random structure after changing this setting.

This function is based on LGPL code by Geoff Kuenning (mtwist-0.8) implementing the Mersenne Twister pseudo-random number generator.

ARGUMENTS

Data type	Name	Description
dip_Random *	random	Pointer to a random value structure
dip_float *	value	Random value

EXAMPLE

Obtain a random number as follows:

```
dip_Random random;
dip_float val;

DIPXJ( dip_RandomSeed( &random, 0 ));
DIPXJ( dip_RandomVariable( &random, &val ));
```

LITERATURE

Makoto Matsumoto and Takuji Nishimura, Mersenne twister: a 623-dimensionally equidistributed uniform pseudo-random number generator, ACM Transactions on Modeling and Computer Simulation 8(1):3-30, 1998.

SEE ALSO

 ${\tt RandomSeed, RandomSeedVector, UniformRandomVariable, GaussianRandomVariable, PoissonRandomVariable, BinaryRandomVariable}$

Code source: mtwist-0.8 or mtwist-0.8

RangeThreshold

Point Operation

SYNOPSIS

#include "dip_point.h"
dip_Error dip_RangeThreshold (in, out, lowerBound, upperBound, foreground, background, binaryOutput)

DATA TYPES

integer, **float**

FUNCTION

out = (lowerBound <= in <= upperBound? foreground: background) If the boolean binaryOutput is true, RangeThreshold will produce a binary image. Otherwise an image of the same type as the input image is produced, with the pixels set to either foreground or background.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	in	Input image
dip_Image	out	Output image
dip_float	lowerBound	Lower bound
dip_float	upperBound	Upper bound
dip_float	foreground	Foreground value
dip_float	background	Background value
dip_Boolean	binaryOutput	Convert output image to binary

SEE ALSO

Threshold, HysteresisThreshold, IsodataThreshold

Real

Arithmetic function

SYNOPSIS

dip_Error dip_Real (in, out)

DATA TYPES

binary, integer, float, $\mathbf{complex}$

FUNCTION

Computes the real part of the input image values, and outputs a float typed image.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	in	Input
$dip_{-}Image$	out	Output

SEE ALSO

Modulus, Phase, Imaginary

Reciprocal

arithmetic function

SYNOPSIS

dip_Error dip_Reciprocal (in, out)

DATA TYPES

binary, integer, float, complex

FUNCTION

Computes the reciprocal (1/x) of the input image values. If pixel values of in are zero, the corresponding pixels in out is set to zero.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	in	Input
dip_Image	out	Output

SEE ALSO

Div, DivFloat, DivComplex, Modulo

Register

Generic registry function

SYNOPSIS

```
#include "dip_registry.h"
dip_Error dip_Register ( register )
```

FUNCTION

The Registry functions Register, Unregister, RegisterClass, RegistryList, RegistryGet RegistryValid and RegistryArrayNew are the access functions for DIPlib's generic registry framework. These functions control the access to a registry containing information of items that are registered at run-time. Each item belongs to a certain class and is identified with an ID that is unique within the item's class.

DIPlib's Registry classes are registered at run-time as well (using RegisterClass) and should be registered before an item of that class can registered.

Although the generic Registry functions can be used to register, and obtain the data of registered items of a specific clas, it is more user friendly to use class-specific Registry functions like MeasurementFeatureRegister and companions.

The dip_Register function accepts one argument, a dip_Registry structure, which contains the ID and class of the to be registered data and registry, a pointer to class-specific data. Note that this pointer, registry, is freed when the (global) registry information is freed.

The following code gives an example of a class-specific register function:

```
dip_Error dip_MeasurementFeatureRegister
(
    dip_MeasurementFeatureRegistry registry
)
{
    DIP_FN_DECLARE("dip_MeasurementFeatureRegister");
    dip_Registry globalRegistry;
    void *data;
    dip_MeasurementFeatureRegistry *reg;

    switch( registry.type )
    {
        default:
            DIPSJ( DIP_E_REGISTRY_INCOMPLETE_REGISTRY );
            break;

        case DIP_MSR_FUNCTION_LINE_BASED:
            DIPTS( ! ( registry.create &&
```

```
registry.measure.line &&
              registry.value &&
              registry.labels &&
              registry.description),
          DIP_E_REGISTRY_INCOMPLETE_REGISTRY );
         break:
      case DIP_MSR_FUNCTION_IMAGE_BASED:
         DIPTS( ! ( registry.create &&
              registry.measure.image &&
              registry.value &&
              registry.labels &&
              registry.description),
          DIP_E_REGISTRY_INCOMPLETE_REGISTRY );
         break;
      case DIP_MSR_FUNCTION_CHAINCODE_BASED:
         DIPTS( ! ( registry.create &&
              registry.measure.chaincode &&
              registry.value &&
              registry.labels &&
              registry.description),
          DIP_E_REGISTRY_INCOMPLETE_REGISTRY );
         break;
      case DIP_MSR_FUNCTION_COMPOSITE:
         DIPTS( ! ( registry.create &&
              registry.measure.composite &&
              registry.value &&
              registry.convert &&
              registry.description),
          DIP_E_REGISTRY_INCOMPLETE_REGISTRY );
         break;
   }
   /* copy the Measurement specific registry info */
   DIPXJ( dip_MemoryNew( &data, sizeof( dip_MeasurementFeatureRegistry ), 0 ));
   reg = ( dip_MeasurementFeatureRegistry * ) data;
   *reg = registry;
   globalRegistry.id
                         = registry.id.rtid;
   globalRegistry.class
                        = DIP_REGISTRY_CLASS_MEASUREMENT;
   globalRegistry.registry = reg;
   globalRegistry.free
                         = dip_MemoryFree;
   /* register this measurement registry data */
   DIPXJ( dip_Register( globalRegistry ));
dip_error:
   DIP_FN_EXIT;
```

}

ARGUMENTS

Data type	Name	Description
dip_Registry	registry	Generic registry structure

SEE ALSO

Unregister, RegisterClass, RegistryList, RegistryGet, RegistryValid, RegistryArrayNew

RegisterClass

Register a registry class

SYNOPSIS

```
#include "dip_registry.h"
dip_Error dip_RegisterClass ( class )
```

FUNCTION

This function registers a Registry class. See Register for more information about DIPlib's Registry.

ARGUMENTS

Data type	Name	Description
dip_int	class	Registry class ID

SEE ALSO

Register, Unregister, RegistryValid, RegistryList, RegistryGet, RegistryArrayNew

RegistryArrayNew

Allocate a registry array

SYNOPSIS

```
#include "dip_registry.h"
dip_Error dip_RegistryArrayNew ( array, size, resources )
```

FUNCTION

This function allocates an array of dip_Registry structures.

ARGUMENTS

Data type	Name	Description
dip_RegistryArray *	array	Pointer to the allocated array
dip_int	size	Array size
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

Register, RegistryList, RegistryGet

RegistryGet

Get a registry item

SYNOPSIS

```
#include "dip_registry.h"
dip_Error dip_RegistryGet ( id, class, registry )
```

FUNCTION

This function obtains the Registry information of the ID of the Registry class class. See Register for more information about DIPlib's Registry.

The following code gives an example of a class-specific register list function:

```
dip_Error dip_MsrRegistryGet
(
    dip_int id,
    dip_MsrRegistry *registry
)
{
    DIP_FN_DECLARE("dip_MsrRegistryGet");
    dip_MsrRegistry *reg;
    void *data;

    DIPXJ( dip_RegistryGet ( id, DIP_REGISTRY_CLASS_MEASUREMENT, &data ));
    reg = data;
    *registry = *reg;

dip_error:
    DIP_FN_EXIT;
}
```

ARGUMENTS

Data type	Name	Description
dip_int	id	Registry ID
dip_int	class	Registry class
dip_void **	registry	Pointer to registered data

SEE ALSO

Register, Unregister, RegisterClass, RegistryList, RegistryValid, RegistryArrayNew

RegistryList

Get an array of registry IDs

SYNOPSIS

```
#include "dip_registry.h"
dip_Error dip_RegistryList ( id, class, resources )
```

FUNCTION

This function obtains an array of the registered IDs of the Registry class class. See Register for more information about DIPlib's Registry.

The following code gives an example of a class-specific register list function:

ARGUMENTS

Data type	Name	Description
dip_IntegerArray *	id	Pointer to an array of Registry IDs
dip_int	class	Registry class
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

Register, Unregister, RegisterClass, RegistryGet, RegistryValid, RegistryArrayNew

RegistryValid

Validate an registry item

SYNOPSIS

```
#include "dip_registry.h"
dip_Error dip_RegistryValid ( id, class, verdict )
```

FUNCTION

This function checks whether id has been registered in the Registry in the Registry class class. If verdict is not zero, the validation information (DIP_FALSE or DIP_TRUE) is copied to verdict. Otherwise an error is returned in case the validation fails.

See Register for more information about DIPlib's Registry.

ARGUMENTS

Data type	Name	Description
dip_int	id	Registry ID
dip_int	class	Registry class
dip_Boolean *	verdict	Pointer to a boolean containing the validation data

SEE ALSO

 ${\tt Register}, {\tt Unregister}, {\tt RegisterClass}, {\tt RegistryList}, {\tt RegistryGet}, {\tt RegistryArrayNew}$

Resampling

Interpolation function

SYNOPSIS

```
#include "dip_interpolation.h"
dip_Error dip_Resampling ( in, out, zoom, shift, method )
```

DATA TYPES

binary, integer, float

FUNCTION

This function resmaples the input image in to out using various interpolation methods. Both a (subpixel) shift and a zoom factor are supported. The size of the output image is zoom times the size of in. If shift is zero, a shift of zero is assumed. If zoom is zero, a zoom of 1.0 is assumed.

ARGUMENTS

Data type	Name	Description	
dip_Image	in	Input image	
dip_Image	out	Output image	
dip_FloatArray	zoom	Zoom factor	
dip_FloatArray	shift	Shift	
dipf_Interpolation	method	Interpolation method	

The dipf_Interpolation enumaration consists of the following constants:

Name	Description
DIP_INTERPOLATION_DEFAULT	Default method, usually equivalent to
	DIP_INTERPOLATION_BSPLINE
DIP_INTERPOLATION_BSPLINE	B-Spline interpolation
DIP_INTERPOLATION_FOURTH_ORDER_CUBIC	Forth order cubic interpolation
DIP_INTERPOLATION_THIRD_ORDER_CUBIC	Third order cubic interpolation
DIP_INTERPOLATION_LINEAR	Linear interpolation
DIP_INTERPOLATION_ZERO_ORDER_HOLD	Zero order hold interpolation
DIP_INTERPOLATION_NEAREST_NEIGHBOUR	Nearest neighbour interpolation
DIP_INTERPOLATION_LANCZOS_2	Lanczos interpolation with $a=2$
DIP_INTERPOLATION_LANCZOS_3	Lanczos interpolation with $a=3$
DIP_INTERPOLATION_LANCZOS_4	Lanczos interpolation with $a=4$
DIP_INTERPOLATION_LANCZOS_6	Lanczos interpolation with $a=6$
DIP_INTERPOLATION_LANCZOS_8	Lanczos interpolation with $a=8$

All interpolation is performed separably. B-spline interpolation uses all samples on the image line.

Cubic interpolation uses a cubic spline kernel (piecewise cubic polynomial), covering 4 (third order) or 6 (fourth order) input samples. Lanczos interpolation uses a sinc function windowed by a wider sinc function, using 2a input samples. Zero order hold and nearest neighbour are the same method, but for the function Resampling, zero order hold results in a shift of half a pixel (i.e. the nearest value "to the left" is always used).

SEE ALSO

Subsampling

ResourcesFree

Free resources

SYNOPSIS

dip_Error dip_ResourcesFree(resources, flags)

FUNCTION

Free all resources registers in the resource tracking structure, as well as the resource tracking structure itself. To prevent errors, the resource tracking structure is set to 0. Passing a null pointer instead of a pointer to a dip_Resources structure is allowed and silently ignored.

ARGUMENTS

Data type	Name	Description
dip_Resources *	resources	The tracking structure which was used to register the
		resources that must be freed. Note the double pointer,
		allowing this routine to set your pointer to 0
dipf_ResourcesFree	flags	When set to DIP_RMRF_DONT_FREE, dip_ResourcesFree only
		frees the dip_Resources structure itself, not the resources
		associated with it

SEE ALSO

ResourcesNew, ResourcesMerge, ResourceSubscribe, ResourceUnsubscribe

ResourcesMerge

Add one resource list to another

SYNOPSIS

dip_Error dip_ResourcesMerge(resources, mergee)

FUNCTION

Adds one resource list to another. This function is very useful when writing functions that will support a dip_Resources parameter. Typically you want to allocate a number of resources and only add these to the user-supplied dip_Resources when all these allocations have been successful. This is where ResourcesMerge comes in. Allocate a local dip_Resources structure and register all resources with it. When no errors occured the local dip_Resources structure can be merged with the user-supplied dip_Resources structure. If an error did occur, simply free all local resources by calling ResourcesFree. In addition it is the convention that functions supporting resource tracking also accept a zero indicating that no tracking should be performed. When resources in dip_ResourcesMerge is 0, the mergee tracking structure is freed, but the resources it contains are not. In this way you get support for the "resources = 0 means no tracking" convention for free.

ARGUMENTS

Data type	Name	Description
dip_Resources	resources	The dip_Resources structure with which the additional resources
		much be merged
dip_Resources *	mergee	A pointer to the dip_Resources structure containing the
		additional resources to be merged. After the merge, mergee is
		set to 0.

SEE ALSO

ResourcesNew, ResourcesFree, ResourceSubscribe, ResourceUnsubscribe

ResourcesNew

Allocate a resource tracking structure

SYNOPSIS

dip_Error dip_ResourcesNew(resources, noItems)

FUNCTION

This function allocates a dip_Resources structure. The resource structure can be used to register various resources as they are allocated, provided that the allocating function allows you to register the resource. All resources allocated in this manner can be freed with a single call to ResourcesFree. Examples of functions supporting this registration scheme are ImageNew and special versions of the memory allocation routines.

Some operations consist of an initialization and a cleanup stage. These stages are often performed by separate routines to allow the user to execute his/her own operations in between. In DIPlib there usually is no directly callable cleanup function. Instead the initialization routine registers its cleanup routine with a dip_Resources structure provided by the user. The cleanup operation is invoked through ResourcesFree.

All functions that support the registration leave you the choice not to register the resource. This is indicated by supplying a zero instead of a resource tracking structure, unless documented otherwise. The noItems parameters can be used to give the routine a hint about the number of resources you will register. The structure grows automagically whenever more resources are registered than indicated by the hint parameter.

ARGUMENTS

Data type	Name	Description
dip_Resources *	resources	This will be used to return a dip_Resources structure
dip_int	noItems	A hint about the number of resources you are planning to
		allocate. This parameter must be ≥ 2 or 0 to indicate you
		want the default value. By the way, don't worry too much about
		this parameter, because when the structure turns out to be too
		small, it will automatically be expanded

SEE ALSO

ResourcesFree, ResourcesMerge, ResourceSubscribe, ResourceUnsubscribe, MemoryNew

ResourceSubscribe

Track a resource

SYNOPSIS

dip_Error dip_ResourceSubscribe(resource, freeResourceHandler, resources)

FUNCTION

Track a resource. The resource must be represented by a void *. A handler function to free the resource must be given. This function will be called through dip_ResourcesFree with the resource as its only parameter. If dip_ResourceSubscribe fails, the resource is not registered. It is allowed to pass a zero instead of a dip_Resources structure, in which case dip_ResourceSubscribe returns silently.

ARGUMENTS

Data type	Name	Description
void *	resource	The resource that must be registered
dip_ResourcesFreeHandler	freeResourceHandler	The handler function that will be invoked
		by dip_ResourcesFree to free the
		resource
dip_Resources	resources	dip_Resources structure to register the
		resource with

SEE ALSO

ResourcesNew, ResourcesFree, ResourcesMerge, ResourceUnsubscribe

ResourceUnsubscribe

Stop tracking a resource

SYNOPSIS

dip_Error dip_ResourceUnsubscribe(resource, resources)

FUNCTION

Stop tracking a resource. It will be removed from the dip_Resources structure. The resource itself will not be freed. If a zero is passed instead of a resource or the dip_Resources structure, dip_ResourceUnsubscribe returns silently.

ARGUMENTS

Data type	Name	Description	
void *	resource	The resource that should no longer be tracked	
dip_Resources	resources	dip_Resources structure to remove the resource from	

SEE ALSO

ResourcesNew, ResourcesFree, ResourcesMerge, ResourceSubscribe

${\tt RootMeanSquareError}$

difference measure

SYNOPSIS

dip_Error dip_RootMeanSquareError (in1, in2, mask, out)

DATA TYPES

binary, integer, float, complex

FUNCTION

Calculates the root mean square difference between each pixel value of in1 and in2. Optionally the mask image can be used to exclude pixels from the calculation by setting the value of these pixels in mask to zero.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	in1	First input
dip_Image	in2	Second input
dip_Image	mask	Mask
dip_Image	out	Output

SEE ALSO

MeanError, MeanSquareError, MeanAbsoluteError, LnNormError, IDivergence

Rotation

Interpolation function

SYNOPSIS

```
#include "dip_interpolation.h"
dip_Error dip_Rotation ( in, out, angle, method, bgval )
```

DATA TYPES

binary, integer, float

FUNCTION

This function rotates an 2-D image in over angle to out using three skews. The function implements the rotation in the mathmetical sense, but note the Y-axis is positive downwards! The rotation is over the centre of the image.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input image
dip_Image	out	Output image
dip_float	angle (radians)	Rotation angle
$dipf_Interpolation$	method	Interpolation method
dip_BackgroundValue	bgval	Background value

The dipf_Interpolation enumaration consists of the following constants:

Name	Description
DIP_INTERPOLATION_DEFAULT	Default method, usually equivalent to
	DIP_INTERPOLATION_BSPLINE
DIP_INTERPOLATION_BSPLINE	B-Spline interpolation
DIP_INTERPOLATION_FOURTH_ORDER_CUBIC	Forth order cubic interpolation
DIP_INTERPOLATION_THIRD_ORDER_CUBIC	Third order cubic interpolation
DIP_INTERPOLATION_LINEAR	Linear interpolation
DIP_INTERPOLATION_ZERO_ORDER_HOLD	Zero order hold interpolation
DIP_INTERPOLATION_NEAREST_NEIGHBOUR	Nearest neighbour interpolation
DIP_INTERPOLATION_LANCZOS_2	Lanczos interpolation with $a=2$
DIP_INTERPOLATION_LANCZOS_3	Lanczos interpolation with $a=3$
DIP_INTERPOLATION_LANCZOS_4	Lanczos interpolation with $a=4$
DIP_INTERPOLATION_LANCZOS_6	Lanczos interpolation with $a=6$
DIP_INTERPOLATION_LANCZOS_8	Lanczos interpolation with $a=8$

All interpolation is performed separably. B-spline interpolation uses all samples on the image line.

Cubic interpolation uses a cubic spline kernel (piecewise cubic polynomial), covering 4 (third order) or 6 (fourth order) input samples. Lanczos interpolation uses a sinc function windowed by a wider sinc function, using 2a input samples. Zero order hold and nearest neighbour are the same method, but for the function Resampling, zero order hold results in a shift of half a pixel (i.e. the nearest value "to the left" is always used).

The dip_BackgroundValue enumaration consists of the following flags:

Name	Description
DIP_BGV_DEFAULT	Default: fill with zeros
DIP_BGV_ZERO	Fill with zeros
DIP_BGV_MAX_VALUE	Fill with maximum value for data type
DIP_BGV_MIN_VALUE	Fill with minimum value for data type

KNOWN BUGS

This function is only implemented for 2D images, for rotating 3D images see Rotation3d.

SEE ALSO

Skewing, Rotation3d, Rotation3d_Axis

Rotation3d

Interpolation function

SYNOPSIS

```
#include "dip_interpolation.h"
dip_Error dip_Rotation3d ( in, out, alpha, beta, gamma, method, bgval )
```

DATA TYPES

binary, integer, float

FUNCTION

This function rotates an 3-D image in via the three Euler angles alpha, beta, gamma to out using nine skews. The first rotation is about alpha around the initial 3-axis. The second about beta around the intermediate 2-axis and the last about gamma around the final 3-axis. The function implements the rotation in the mathmetical sense, but note the Y-axis is positive downwards! The rotation is over the centre of the image.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	in	Input image
dip_Image	out	Output image
dip_float	alpha (radians)	Euler angle
dip_float	beta (radians)	Euler angle
dip_float	gamma (radians)	Euler angle
dipf_Interpolation	method	Interpolation method
dip_BackgroundValue	bgval	Background value

The dipf_Interpolation enumaration consists of the following constants:

Name	Description
DIP_INTERPOLATION_DEFAULT	Default method, usually equivalent to
	DIP_INTERPOLATION_BSPLINE
DIP_INTERPOLATION_BSPLINE	B-Spline interpolation
DIP_INTERPOLATION_FOURTH_ORDER_CUBIC	Forth order cubic interpolation
DIP_INTERPOLATION_THIRD_ORDER_CUBIC	Third order cubic interpolation
DIP_INTERPOLATION_LINEAR	Linear interpolation
DIP_INTERPOLATION_ZERO_ORDER_HOLD	Zero order hold interpolation
DIP_INTERPOLATION_NEAREST_NEIGHBOUR	Nearest neighbour interpolation
DIP_INTERPOLATION_LANCZOS_2	Lanczos interpolation with $a=2$
DIP_INTERPOLATION_LANCZOS_3	Lanczos interpolation with $a=3$
DIP_INTERPOLATION_LANCZOS_4	Lanczos interpolation with $a=4$
DIP_INTERPOLATION_LANCZOS_6	Lanczos interpolation with $a=6$
DIP_INTERPOLATION_LANCZOS_8	Lanczos interpolation with $a=8$

All interpolation is performed separably. B-spline interpolation uses all samples on the image line. Cubic interpolation uses a cubic spline kernel (piecewise cubic polynomial), covering 4 (third order) or 6 (fourth order) input samples. Lanczos interpolation uses a sinc function windowed by a wider sinc function, using 2a input samples. Zero order hold and nearest neighbour are the same method, but for the function Resampling, zero order hold results in a shift of half a pixel (i.e. the nearest value "to the left" is always used).

The dip_BackgroundValue enumaration consists of the following flags:

Name	Description
DIP_BGV_DEFAULT	Default: fill with zeros
DIP_BGV_ZERO	Fill with zeros
DIP_BGV_MAX_VALUE	Fill with maximum value for data type
DIP_BGV_MIN_VALUE	Fill with minimum value for data type

KNOWN BUGS

This function is only implemented for 3D images, for rotating 2D images see Rotation.

SEE ALSO

Skewing, Rotation, Rotation3d_Axis

Rotation3d_Axis

Interpolation function

SYNOPSIS

```
#include "dip_interpolation.h"
dip_Error dip_Rotation3d_Axis ( in, out, angle, axis, method, bgval )
```

DATA TYPES

binary, integer, float

FUNCTION

This function rotates an 3-D image in over angle around axis to out using three skews. The rotation axis is 0 (x), 1 (y) or 2 (z). The function implements the rotation in the mathmetical sense, but note the Y-axis is positive downwards! The rotation is over the centre of the image.

For backwards compatability, the macro Rotation3dAxis calls the function Rotation3d_Axis but uses 1, 2 and 3 to select the axis of rotation.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input image
dip_Image	out	Output image
dip_float	angle (radians)	Rotation angle
dip_int	axis	Rotation axis
${ t dipf_Interpolation}$	method	Interpolation method
dip_BackgroundValue	bgval	Background value

The dipf_Interpolation enumaration consists of the following constants:

Name	Description
DIP_INTERPOLATION_DEFAULT	Default method, usually equivalent to
	DIP_INTERPOLATION_BSPLINE
DIP_INTERPOLATION_BSPLINE	B-Spline interpolation
DIP_INTERPOLATION_FOURTH_ORDER_CUBIC	Forth order cubic interpolation
DIP_INTERPOLATION_THIRD_ORDER_CUBIC	Third order cubic interpolation
DIP_INTERPOLATION_LINEAR	Linear interpolation
DIP_INTERPOLATION_ZERO_ORDER_HOLD	Zero order hold interpolation
DIP_INTERPOLATION_NEAREST_NEIGHBOUR	Nearest neighbour interpolation
DIP_INTERPOLATION_LANCZOS_2	Lanczos interpolation with $a=2$
DIP_INTERPOLATION_LANCZOS_3	Lanczos interpolation with $a=3$
DIP_INTERPOLATION_LANCZOS_4	Lanczos interpolation with $a=4$
DIP_INTERPOLATION_LANCZOS_6	Lanczos interpolation with $a=6$
DIP_INTERPOLATION_LANCZOS_8	Lanczos interpolation with $a=8$

All interpolation is performed separably. B-spline interpolation uses all samples on the image line. Cubic interpolation uses a cubic spline kernel (piecewise cubic polynomial), covering 4 (third order) or 6 (fourth order) input samples. Lanczos interpolation uses a sinc function windowed by a wider sinc function, using 2a input samples. Zero order hold and nearest neighbour are the same method, but for the function Resampling, zero order hold results in a shift of half a pixel (i.e. the nearest value "to the left" is always used).

The dip_BackgroundValue enumaration consists of the following flags:

Name	Description
DIP_BGV_DEFAULT	Default: fill with zeros
DIP_BGV_ZERO	Fill with zeros
DIP_BGV_MAX_VALUE	Fill with maximum value for data type
DIP_BGV_MIN_VALUE	Fill with minimum value for data type

KNOWN BUGS

This function is only implemented for 3D images, for rotating 2D images see Rotation.

SEE ALSO

Skewing, Rotation, Rotation3d

${\tt ScalarImageNew}$

Allocate a scalar image

SYNOPSIS

dip_Error dip_ScalarImageNew(newImage, dataType, dimensions, resources)

FUNCTION

Allocate and forge a dip_Image structure of the DIP_IMTP_SCALAR type.

ARGUMENTS

Data type	Name	Description
dip_Image *	newImage	Used to return a pointer to your brand new dip_Image
		structure
dip_DataType	dataType	Data type. See DIPlib's data types
dip_IntegerArray	dimensions	Dimensions
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

ImageNew, ImageFree

ScanFrameWork

FrameWork for scanning multiple images

SYNOPSIS

 $\label{lem:conformed} \begin{tabular}{ll} \b$

FUNCTION

This function provides a framework for scanning nofin input images and nofout output images in one dimension of the images. The dimension in which the image should be scanned can be specified or left to ScanFrameWork by specifying the dimension with DIP_MONADIC_OPTIMAL_DIMENSION. See SeparableFrameWork for details.

ARGUMENTS

Data type	Name	Description
$ exttt{dip_ImageArray}$	in	Array of input images
dip_ImageArray	out	Array of output images
dip_FrameWorkProcess	process	Process
dip_BoundaryArray	boundary	Boundary conditions
dip_BorderArray	border	Border Array
dip_DataTypeArray	inBuffer	Array of dip_DataType of the input buffer
dip_DataTypeArray	outBuffer	Array of dip_DataType of each output buffer
dip_DataTypeArray	outImage	Array of dip_DataType of each output image

SEE ALSO

DIPlib's data types SeparableFrameWork, PixelTableFrameWork

SeededWatershed

Morphological segmentation

SYNOPSIS

```
#include "dip_morphology.h"
dip_Error dip_SeededWatershed ( seeds, in, mask, out, connectivity, order,
max_depth, max_size, binaryOutput )
```

DATA TYPES

integer, float

FUNCTION

Watershed segmentation with built-in region merging. max_depth and max_size control the merging procedure. Any region with max_size or less pixels and with max_depth grey-value difference or less will be merged to neighbouring regions when they touch (as opposed to build a watershed). max_size equal to 0 means that the size of the region is not tested when merging. To avoid merging of seeds with no grey-value difference between them, set max_size to a negative value. The regions are grown according to the connectivity parameter. See The connectivity parameter for more information. The output is either a labelled image where the pixels belonging to a catchment basin are labelled, or a binary image where the watershed pixels are 1 and the rest is 0. This is controlled by binaryOutput.

As opposed to Watershed, this function takes a seeds input image, and grows the catchment basins from there. The output image, when binaryOutput is DIP_TRUE, will have label values as given by the seed image.

If mask is not 0, only the pixels within mask will be considered. All the other pixels will be untouched.

ARGUMENTS

Data type	Name	Description
dip_Image	seeds	Binary or labelled input image
dip_Image	in	Grey-value input image
dip_Image	mask	Mask image
dip_Image	out	Output
dip_int	connectivity	Connectivity
dipf_GreyValueSortOrder	order	Whether to grow from low to high or high to low
dip_float	max_depth	Maximum depth of a region that can be merged
dip_int	max_size	Maximum size of a region that can be merged
dip_Boolean	binaryOutput	DIP_FALSE if the output should be a labelled image

The ${\tt dipf_GreyValueSortOrder}$ enumeration consists of the following values:

Name	Description
DIP_GVSO_HIGH_FIRST	Process the pixels from high grey-value to low grey-value.
DIP_GVSO_LOW_FIRST	Process the pixels from low grey-value to high grey-value.

SEE ALSO

Watershed, LocalMinima, Minima, Maxima, GrowRegions

Select

Configurable selection function

SYNOPSIS

dip_Error dip_Select (in1, in2, in3, in4, out, selector)

DATA TYPES

binary, integer, float

FUNCTION

This function can perform various pixel-by-pixel comparisons (smaller, smaller- equal, equal, not equal, greater-equal, greater) between in1 ans in2. If the result of the comparison is true, the corresponding pixel value of in3 is copied to out, otherwise it is copied from in4. In short the following operation is performed for each pixel in the five images:

```
out = in1 selector in2 ? in3 : in4
```

The images in2, in3 and in4 can be 0-D images acting as constants.

ARGUMENTS

Data type	Name	Description
dip_Image	in1	First input
dip_Image	in2	Second input
dip_Image	in3	Third input
dip_Image	in4	Fourth input
dip_Image	out	Output
dipf_Select	selector	Select flag

The dipf_Select flag can be one of:

Name	Description
DIP_SELECT_LESSER	<, Lesser than
DIP_SELECT_LESSER_EQUAL	<=, Lesser or equal
DIP_SELECT_NOT_EQUAL	!=, Unequal
DIP_SELECT_EQUAL	==, Equal
DIP_SELECT_GREATER_EQUAL	>=, Greater or equal
DIP_SELECT_GREATER	>, Greater

SEE ALSO

Compare, Max, Min

SelectValue

Point Operation

SYNOPSIS

```
#include "dip_point.h"
dip_Error dip_SelectValue ( in, out, value )
```

DATA TYPES

integer, float

FUNCTION

This function returns a binary image with value 1 where in == value and value 0 elsewhere.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	in	Input image
dip_Image	out	Output image
dip_float	value	Value to select

SEE ALSO

Threshold, NotZero, Compare, RangeThreshold

SeparableConvolution

FrameWork for separable convolution filters

SYNOPSIS

#include "dip_linear.h"
dip_Error dip_SeparableConvolution (in, out, filters, bc, process)

DATA TYPES

integer, **float**

FUNCTION

This function is a frontend to the lower level Convolve1d function. Each dimension can be processed by a different filter.

process may be zero, indicating that all dimensions should be processed.

ARGUMENTS

Data type	Name	Description
dip_Image	in Input image	
dip_Image	out	Output image
dip_SeparableConvolutionFilter *	filters	Filters
dip_BoundaryArray	bc	Boundary conditions
dip_BooleanArray	process (0)	Dimensions to process

The dip_SeparableConvolutionFilter structure contains the following elements:

Data type	Name	Description	
dip_float *	filter	Filter weights	
dip_int	filterSize	Length of filter array	
dip_int	origin	Origin of the filter, only valid in conjunction with	
		DIP_CNV_USE_ORIGIN	
dipf_Convolve	flags	Filter flags, see Convolve1d for their definitions	

SEE ALSO

General information about convolution

GeneralConvolution, ConvolveFT, SeparableFrameWork, Convolve1d

SeparableFrameWork

FrameWork for separable filters

SYNOPSIS

dip_Error dip_SeparableFrameWork (in, out, boundary, border, process)

FUNCTION

The dip_SeparableFrameWork function is a framework for separable filters. This function takes care of all the "administrative work" involved when processing a n-D DIPlib image n times with a 1-D filter. In short, using this function, one has only to create an one dimension filter function and dip_SeparableFrameWork takes care of the other stuff. The in image is filtered nrOfProcesses times using the information in each element of the process array. If nrOfProcesses is zero, only the first element of process is used to filter in all its dimensions.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input image
dip_Image	out	Output image
dip_BoundaryArray	boundary	Boundary conditions
dip_IntegerArray	border	Border array
dip_FrameWorkProcessArray	process	Array of dip_FrameWorkProcess structures

NOTE

The dip_FrameWorkProcess structure contains the following members:

Data type	Name	Description
dip_Boolean	process	flags specifying to do processing or not
dipf_FrameWorkMethod	frameWorkMethod	flags specifying the method of how
		dip_SeparableFrameWork should
		transport data from in to out
dipf_FrameWorkOperation	frameWorkOperation	flags specifying requirements of the 1-D
		filter function
dip_int	processDimension	dimension of in to be processed
dip_int	roiOrigin	ignored in current implementation
dip_int	roiSize	ignored in current implementation
dipf_FrameWorkFilter	FrameWorkFilterType	specifying the type of 1-D filter function
dip_FrameWorkFilter	FrameWorkFilter	pointer to the 1-D filter function
void *	functionParameters	parameters of the 1-D filter function, for
		all dimensions
dip_DataType	inputBufferType	data type of input buffer
dip_DataType	outputBufferType	data type of output buffer
dip_DataType	suggestedOutputType	data type of output image

The ${\tt dipf_FrameWorkMethod}$ enum contains the following elements:

Name	Description
DIP_FRAMEWORK_DEFAULT_METHOD	use dip_SeparableFrameWorks most optimal method
DIP_FRAMEWORK_CLASSICAL	use a classical method
DIP_FRAMEWORK_DOUBLE_STRIPE	use two buffers to store temporary results

It is strongly advised to use the ${\tt DIP_FRAMEWORK_DEFAULT_METHOD}$ method.

The ${\tt dipf_FrameWorkOperation}$ enum contains the following elements:

Name	Description
DIP_FRAMEWORK_DEFAULT_OPERATION	default operation
DIP_FRAMEWORK_IN_PLACE	filtering operation can be performed in-place. It is up to
	dip_SeparableFrameWork whether the actual filtering is
	done in-place
DIP_FRAMEWORK_NO_IN_BORDER	the 1-D filter does not need border extension of the input
	data
DIP_FRAMEWORK_OUT_BORDER	the 1-D filter needs border extension of the output data
DIP_FRAMEWORK_WRITE_INPUT	the 1-D filter needs to write in the input data
DIP_FRAMEWORK_USE_BUFFER_TYPES	made the input and output buffers of the
	inputBufferType and outputBufferType data type
DIP_FRAMEWORK_NO_BUFFER_STRIDE	Create input and output buffers with a stride of one
DIP_FRAMEWORK_DO_NOT_ADJUST	Do not adjust the output image, just check it
DIP_FRAMEWORK_USE_OUTPUT_TYPE	Adjust output image to the suggestedOutputType data
	type

The dipf_FrameWorkFilter enum contains the following elements:

Name	Description
DIP_FRAMEWORK_SEPARABLE_FILTER	default filter type, process one line at the time
DIP_FRAMEWORK_TWO_LINES_SEPARABLE_FILTER	process two lines in one go
DIP_FRAMEWORK_SINGLE_OUTPUT_FILTER	this filter only needs an output buffer

The union ${\tt dip_FrameWorkFunction}$ consists of the types

Name	Description
dip_SeparableFilter	one line filter function
dip_TwoLinesSeparableFilter	two lines filter function
dip_SingleOutputFilter	single output image filter

The functions have the following arguments: dip_SeparableFilter

Data type	Name	Description
void *	inData	pointer to the input data
void *	outData	pointer to the output data
dip_int	elements	number of pixels in the inData array
dip_SeparableFilterParameters	params	parameter structure for the filter function

and ${\tt dip_SeparableTwoLinesFilter}$

Data type	Name	Description
void *	inFirstLineData	pointer to the data of the
		first input line
void *	inSecondLineData	pointer to the data of the
		second input line
void *	outFirstLineData	pointer to the data of the
		first output line
void *	outSecondLineData	pointer to the data of the
		second output line
dip_int	elements	number of pixels in the
		${\tt inFirstLineData}\ {\tt array}$
$ ext{dip}_{-} ext{TwoLinesSeparableFilterParameters}$	params	parameter structure for the
		two lines filter function

The inData, inFirstLineData and inSecondLineData will always point to the first pixel of the line of in that is processed. Therefore, the 1-D filter can access pixels with indices ranging from -border[dimension] up to elements + border[dimension]. If the flag DIP_FRAMEWORK_OUT_BOUNDARY is specified, the same holds for the output data pointers.

The structure dip_SeparableFilterParameters contains the following elements:

Data type	Name	Description	
void *	functionParameters	parameters of the 1-D filter function per dimension	
dip_int	dimension	the dimension in which direction the input buffer is	
		taken from the input image	
dip_int	processNumber	number of times dip_SeparableFrameWork has	
		already filtered in with an 1-D filter including	
		current filtering	
dip_DataType	inType	dip_DataType of the input buffer	
dip_int	inStride	stride of the elements in the input array	
dip_int	inPlane	plane number in case in is a binary image	
dip_DataType	outType	dip_DataType of the output buffer	
dip_int	outStride	stride of the elements in the output array	
dip_int	outPlane	plane number in case out is a binary image	
dip_int	outDimension	the dimension in which direction the output buffer is	
		taken from the output image	
dip_IntegerArray	position	coordinate of the first pixel of the input buffer in the	
		input image	

The structure $\mathtt{dip_TwoLinesSeparableFilterParameters}$ contains the following elements:

Data type	Name	Description
void *	functionParameters	parameters of the 1-D filter function per dimension
dip_int	dimension	the dimension in which direction the input buffer is
		taken from the input image
dip_int	processNumber	number of times dip_SeparableFrameWork has
		already filtered in with an 1-D filter including
		current filtering
dip_DataType	inType	dip_DataType of the input buffer
dip_int	inStride	stride of the elements in the input array
dip_int	inPlane	plane number in case in is a binary image
dip_DataType	outType	dip_DataType of the output buffer
dip_int	outStride	stride of the elements in the output array
dip_int	outPlane	plane number in case out is a binary image
dip_int	outDimension	the dimension in which direction the output buffer is
		taken from the output image
dip_IntegerArray	position	coordinate of the first pixel of the input buffer in the
		input image

The structure ${\tt dip_SingleOutputFilterParameters}$ contains the following elements:

Data type	Name	Description
void *	functionParameters	parameters of the 1-D filter function per dimension
dip_int	dimension	the dimension in which direction the input buffer is
		taken from the input image
$\mathtt{dip}_{-}\mathtt{int}$	processNumber	number of times dip_SeparableFrameWork has
		already filtered in with an 1-D filter including
		current filtering
dip_DataType	type	dip_DataType of the input buffer
dip_int	stride	stride of the elements in the input array
$\mathtt{dip}_{-}\mathtt{int}$	plane	plane number in case in is a binary image
dip_IntegerArray	position	coordinate of the first pixel of the input buffer in the
		input image

SEE ALSO

 $DIPlib's\ data\ types\ {\tt SeparableConvolution},\ {\tt MonadicFrameWork},\ {\tt SingleOutputFrameWork},\ {\tt PixelTableFrameWork},\ {\tt ScanFrameWork}$

Set

the value of a pixel

SYNOPSIS

```
#include "dip_manipulation.h"
dip_Error dip_Set ( out, const, cor, adjust )
```

FUNCTION

This function set a value of a pixel at position cor in the image out to the value const.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	out	Output image
$\mathtt{dip}_{-}\mathtt{Image}$	constImage	0-D image
dip_IntegerArray	cor	Pixel coordinate
dip_Boolean	adjust	Adjust data type of output image

SEE ALSO

 ${\tt SetInteger, SetFloat, SetComplex, dip_PixelSetInteger, dip_PixelSetFloat, GetEnteger, dip_PixelSetFloat, GetEnteger, dip_PixelSetFloat, GetEnteger, dip_PixelSetFloat, GetEnteger, dip_PixelSetFloat, GetEnteger, dip_PixelSetFloat, GetEnteger, di$

${\tt SetComplex}$

Set a pixel value

SYNOPSIS

```
#include "dip_manipulation.h"
dip_Error dip_SetComplex ( out, constant, cor, adjust )
```

FUNCTION

This function set a value of a pixel at position cor in the image out to the value constant.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	out	Output
dip_complex	constant	Constant
dip_IntegerArray	cor	Pixel coordinate
dip_Boolean	adjust	Adjust data type of output image

SEE ALSO

 ${\tt Set}, {\tt SetInteger}, {\tt SetFloat}, {\tt dip_PixelSetInteger}, {\tt dip_PixelSetFloat}, {\tt Get}$

SetFloat

Set a pixel value

SYNOPSIS

```
#include "dip_manipulation.h"
dip_Error dip_SetFloat ( out, constant, cor, adjust )
```

FUNCTION

This function set a value of a pixel at position cor in the image out to the value constant.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	out	Output
dip_float	constant	Constant
dip_IntegerArray	cor	Pixel coordinate
dip_Boolean	adjust	Adjust data type of output image

SEE ALSO

 ${\tt Set}, {\tt SetInteger}, {\tt SetComplex}, {\tt dip_PixelSetInteger}, {\tt dip_PixelSetFloat}, {\tt GetInteger}, {\tt GetInteger}, {\tt GitInteger}, {$

SetInteger

Set a pixel value

SYNOPSIS

```
#include "dip_manipulation.h"
dip_Error dip_SetInteger ( out, constant, cor, adjust )
```

FUNCTION

This function set a value of a pixel at position cor in the image out to the value constant.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	out	Output
dip_int	constant	Constant
dip_IntegerArray	cor	Pixel coordinate
dip_Boolean	adjust	Adjust data type of output image

SEE ALSO

 ${\tt Set}, {\tt SetFloat}, {\tt SetComplex}, {\tt dip_PixelSetInteger}, {\tt dip_PixelSetFloat}, {\tt Get}$

Sharpen

Enhance an image

SYNOPSIS

```
#include "dip_derivatives.h"
dip_Error dip_Sharpen ( in, out, weight, bc, ps, sigmas, tc, flavour )
```

DATA TYPES

See Laplace

FUNCTION

This function enhances the high frequencies ("sharpens") of the input image in by subtracting a Laplace filtered version of in from it. The weight parameter determines by which amount the laplace information is subtracted from the original: output = input - weight * laplace(input) The sigmas are the Gaussian smoothing parameters of the Laplace operation, and determine how strongly the high-frequency noise in in is suppressed.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input image
dip_Image	out	Output image
dip_float	weight	Laplacian weight
dip_BoundaryArray	boundary	Boundary conditions
dip_BooleanArray	process (0)	Dimensions to process
dip_FloatArray	sigmas	Sigma of Gaussian
dip_float	truncation (<0)	Truncation of Gaussian, see
		GlobalGaussianTruncationGet
dip_DerivativeFlavour	flavour	Derivative Flavour

SEE ALSO

Laplace

Shift

an image manipulation function

SYNOPSIS

```
#include "dip_manipulation.h"
dip_Error dip_Shift ( in, out, shift, killNy )
```

DATA TYPES

binary, integer, float, complex

FUNCTION

This function shifts an image in the Fourier Domain. All frequiencies larger than the Nyquist frequency are set to zero if killNy is true. It performs:

out = Real(InverseFourierTransform(GeneratePhase(shift) * FourierTransform(in))

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	in	Input image
dip_Image	out	Output image
dip_FloatArray	shift	Shift array
dip_Boolean	killNy	set frequencies > Nyquist to zero?

SEE ALSO

Crop, Wrap, FourierTransform, Real

Sigma

Adaptive uniform smoothing filter

SYNOPSIS

```
#include "dip_filtering.h"
dip_Error dip_Sigma ( in, out, se, boundary, param, shape, sigma, outputCount )
```

DATA TYPES

integer, float

FUNCTION

The Sigma filter is an adaptive Uniform smoothing filter. The value of the pixel underinvestigation is replaced by the average of the pixelvalues in the filter region (as specified by param, shape and se) which lie in the interval +/- 2 sigma from the value of the pixel that is filtered. If outputCount is DIP_TRUE, the output values represent the number of pixels over which the average has been calculated. When threshold is DIP_TRUE, the pixel intensities are thresholded at +/- 2 sigma, when it is set to DIP_FALSE, the intensities are weighted with the Gaussian difference with the intensity of the center pixel.

Only the rectangular, elliptic and diamond filter shapes are supported (DIP_FLT_SHAPE_RECTANGULAR, DIP_FLT_SHAPE_ELLIPTIC and DIP_FLT_SHAPE_DIAMOND). Other filter shapes can be implemented by setting shape to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, and passing a binary image in se. The "on" pixels define the shape of the filter window. Other values of shape are illegal.

If shape is not equal to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, se can be set to zero. When shape is set to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, param is ignored, and can be set to zero.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	in	Input image
dip_Image	out	Output image
dip_Image	se	Custom filter window (binary)
dip_BoundaryArray	boundary	Boundary conditions
dip_FloatArray	param	Filter sizes
dip_FilterShape	shape	Filter shape
dip_float	sigma	Sigma
dip_Boolean	outputCount	Output the Count

The enumerator dip_FilterShape contains the following constants:

Name	Description
DIP_FLT_SHAPE_DEFAULT	Default filter window, same as
	DIP_FLT_SHAPE_RECTANGULAR
DIP_FLT_SHAPE_RECTANGULAR	Rectangular filter window, can be even in size
DIP_FLT_SHAPE_ELLIPTIC	Elliptic filter window, always odd in size
DIP_FLT_SHAPE_DIAMOND	Diamond-shaped filter window, always odd in size
DIP_FLT_SHAPE_PARABOLIC	Parabolic filter window (morphology only)
DIP_FLT_SHAPE_DISCRETE_LINE	Rotated line structuring element (morphology only)
DIP_FLT_SHAPE_INTERPOLATED_LINE	Rotated line structuring element, through interpolation
	(morphology only)
DIP_FLT_SHAPE_PERIODIC_LINE	(not implemented)
DIP_FLT_SHAPE_STRUCTURING_ELEMENT	Use se as filter window, can be any size

LITERATURE

John-Sen Lee, $Digital\ Image\ Smoothing\ and\ the\ Sigma\ Filter,$ Computer Vision, Graphics and Image Processing, 24, 255-269, 1983

SEE ALSO

 ${\tt BiasedSigma,\,GaussianSigma,\,Uniform}$

Sign

Arithmetic function

SYNOPSIS

dip_Error dip_Sign (in, out)

DATA TYPES

binary, integer, float

FUNCTION

Computes the sign of the input image values, and outputs a signed integer typed image. The sign of zero is defined as zero.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	in	Input
dip_Image	out	Output

SEE ALSO

Abs, Ceil, Floor, Truncate, Fraction, NearestInt

SimulatedAttenuation

Simulation of the attenuation process

SYNOPSIS

```
#include "dip_microscopy.h"
dip_Error dip_SimulatedAttenuation ( in, out, fAttenuation, bAttenuation, NA,
refIndex, zxratio, oversample, rayStep )
```

DATA TYPES

binary, integer, float

FUNCTION

This function simulates an attenuation based on the model of a CSLM, using a ray tracing method.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input image
$dip_{-}Image$	out	Output image
dip_float	fAttenuation	Forward attenuation factor
dip_float	bAttenuation	Backward attenuation factor
dip_float	NA	Numerical aperture
dip_float	refIndex	Refractive index
dip_float	zxratio	Z/X sampling ratio
dip_int	oversample	Ray casting oversampling
dip_float	rayStep	Ray step

LITERATURE

K.C. Strasters, H.T.M. van der Voort, J.M. Geusebroek, and A.W.M. Smeulders, *Fast attenuation correction in fluorescence confocal imaging: a recursive approach*, BioImaging, vol. 2, no. 2, 1994, 78-92.

AUTHOR

Karel Strasters, adapted to DIPlib by Geert van Kempen.

SEE A	LSO
-------	------------

 ${\tt AttenuationCorrection}, {\tt ExponentialFitCorrection}$

Sin

trigonometric function

SYNOPSIS

dip_Error dip_Sin (in, out)

DATA TYPES

binary, integer, float, complex

FUNCTION

Computes the sine of the input image values.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	out	Output

SEE ALSO

Cos, Tan, Asin, Acos, Atan, Sinh, Cosh, Tanh

Sinc

mathematical function

SYNOPSIS

dip_Error dip_Sinc (in, out)

DATA TYPES

binary, integer, ${f float}$

FUNCTION

Computes the sinc $(\sin(x)/x)$ of the input image values.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	out	Output

SEE ALSO

BesselJO, BesselJI, BesselJN, BesselYO, BesselYI, BesselYN, LnGamma, Erf, Erfc

SingleOutputFrameWork

FrameWork for generation functions

SYNOPSIS

dip_Error dip_SingleOutputFrameWork (out, processBoundary, processBorder, process)

FUNCTION

This function is a frontend on the SeparableFrameWork. It provides an easier interface for filters that only need to scan an single output image. The dimension in which the image should be scanned can be specified or left to SingleOutputFrameWork by specifying the dimension with DIP_MONADIC_OPTIMAL_DIMENSION.

ARGUMENTS

Data type	Name	Description
dip_Image	out	Output
dip_Boundary	processBoundary	ProcessBoundary
dip_int	processBorder	ProcessBorder
dip_FrameWorkProcess	process	Process

SEE ALSO

SeparableFrameWork, MonadicFrameWork, PixelTableFrameWork, ScanFrameWork

${\tt SingularValueDecomposition}$

Singular value decomposition

SYNOPSIS

dip_Error dip_SingularValueDecomposition (in, sz, u, s, v)

DATA TYPES

float

FUNCTION

Computes the SVD of the ImageArray in, such that in = u * s * transpose(v), with s being diagonal. The size of the in matrix is passed to the function via the integer array sz. If the input is of size MxN, then the outputs must be u: nxM, s: NxN, and v: NxN.

Optionally, set u and v to NULL, and let s have N elements, it will contain only the singular values.

ARGUMENTS

Data type	Name	Description
dip_ImageArray	in	Input
dip_IntegerArray	sz	Matrix size of Input
dip_ImageArray	u	Output
dip_ImageArray	s	Output
dip_ImageArray	v	Output

Sinh

trigonometric function

SYNOPSIS

dip_Error dip_Sinh (in, out)

DATA TYPES

binary, integer, ${f float}$

FUNCTION

Computes the hyperbolic sine of the input image values.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{\mathtt{-}}\mathtt{Image}$	in	Input
$\mathtt{dip}_{-}\mathtt{Image}$	out	Output

SEE ALSO

Sin, Cos, Tan, Asin, Acos, Atan, Atan2, Cosh, Tanh

Skewing

Interpolation function

SYNOPSIS

```
#include "dip_interpolation.h"
dip_Error dip_Skewing ( in, out, shear, skew, axis, method, bgval, periodicSkew )
```

DATA TYPES

binary, integer, float

FUNCTION

This function skews the axis axis of in over an angle angle to out using the interpolation method method. The skew is over the centre of the image. If periodicSkew is set to DIP_TRUE, the output image will be of the same size as the input image, and its pixels in the skew dimension wrapped around the image boundaries. bgval is not used in this case.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input image
dip_Image	out	Output image
dip_float	shear (radians)	Shear angle
dip_int	skew	Skew dimension
dip_int	axis	Skew axis
dipf_Interpolation	method	Interpolation method
dip_BackgroundValue	bgval	Background value
dip_Boolean	periodicSkew	Skew using periodic image boundaries

The ${\tt dipf_Interpolation}$ enumaration consists of the following constants:

Name	Description
DIP_INTERPOLATION_DEFAULT	Default method, usually equivalent to
	DIP_INTERPOLATION_BSPLINE
DIP_INTERPOLATION_BSPLINE	B-Spline interpolation
DIP_INTERPOLATION_FOURTH_ORDER_CUBIC	Forth order cubic interpolation
DIP_INTERPOLATION_THIRD_ORDER_CUBIC	Third order cubic interpolation
DIP_INTERPOLATION_LINEAR	Linear interpolation
DIP_INTERPOLATION_ZERO_ORDER_HOLD	Zero order hold interpolation
DIP_INTERPOLATION_NEAREST_NEIGHBOUR	Nearest neighbour interpolation
DIP_INTERPOLATION_LANCZOS_2	Lanczos interpolation with $a=2$
DIP_INTERPOLATION_LANCZOS_3	Lanczos interpolation with $a=3$
DIP_INTERPOLATION_LANCZOS_4	Lanczos interpolation with $a=4$
DIP_INTERPOLATION_LANCZOS_6	Lanczos interpolation with $a=6$
DIP_INTERPOLATION_LANCZOS_8	Lanczos interpolation with $a=8$

All interpolation is performed separably. B-spline interpolation uses all samples on the image line. Cubic interpolation uses a cubic spline kernel (piecewise cubic polynomial), covering 4 (third order) or 6 (fourth order) input samples. Lanczos interpolation uses a sinc function windowed by a wider sinc function, using 2a input samples. Zero order hold and nearest neighbour are the same method, but for the function Resampling, zero order hold results in a shift of half a pixel (i.e. the nearest value "to the left" is always used).

The dip_BackgroundValue enumaration consists of the following flags:

Name	Description	
DIP_BGV_DEFAULT	ULT Default: fill with zeros	
DIP_BGV_ZERO	Fill with zeros	
DIP_BGV_MAX_VALUE	Fill with maximum value for data type	
DIP_BGV_MIN_VALUE	Fill with minimum value for data type	

SEE ALSO

Rotation, Rotation3d, Rotation3d_Axis

SmallObjectsRemove

Remove small objects from an image

SYNOPSIS

```
#include "dip_measurement.h"
dip_Error dip_SmallObjectsRemove ( in, out, threshold )
```

DATA TYPES

integer

FUNCTION

This function removes from the labeled image in those objects whose size (measured in the number of pixels) is smaller than threshold.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	in	Input image
dip_Image	out	Output image
dip_int	threshold	Minimum object size

SEE ALSO

Measure, ObjectToMeasurement, Label

SobelGradient

A linear gradient filter

SYNOPSIS

```
#include "dip_linear.h"
dip_Error dip_SobelGradient ( in, out, boundary, processDim )
```

DATA TYPES

binary, integer, float, complex

FUNCTION

The SobelGradient filter computes a finite difference gradient $(1\ 0\ -1)/2$ in the processDim, and performs a local $(1\ 2\ 1)/4$ smoothing in the other dimensions. Note that in 2D, this differs by a multiplication factor of 1/8 to the original definition by Sobel.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{\mathtt{-}}\mathtt{Image}$	in	Input
dip_Image	out	Output
dip_BoundaryArray	boundary	Boundary conditions
dip_int	processDim	ProcessDim

SEE ALSO

General information about convolution

FiniteDifference, Uniform, Gauss, SeparableConvolution, Convolve1d, Derivative

Sort

Sort a block of data

SYNOPSIS

```
#include "dip_sort.h"
dip_Error dip_Sort ( data, size, algorithm, dataType )
```

FUNCTION

Sorts a block of data (of size $\tt size$ and data type $\tt dataType$) using the algorithm specified by algorithm.

ARGUMENTS

Data type	Name	Description
void *	data	Data
dip_int	size	Size
dip_Sort	algorithm	Sort algorithm
dip_DataType	dataType	Data type. See DIPlib's data types

The sortType parameter is one of:

Name	Description
DIP_SORT_DEFAULT	Default sort algorithm
DIP_SORT_QUICK_SORT	Quick sort
DIP_SORT_DISTRIBUTION_SORT	Distribution sort
DIP_SORT_INSERTION_SORT	Insertion sort

SEE ALSO

General information about sorting

 ${\tt DistributionSort, InsertionSort, QuickSort, ImageSort, SortIndices, SortIndices16, ImageSortIndices}$

SortAnything

Sort data of any type

SYNOPSIS

```
#include "dip_sort.h"
dip_Error dip_SortAnything ( data, size, compareFunction, swapFunction, tmpData,
algorithm )
```

FUNCTION

Sorts a block of data (of size size) using the algorithm specified by algorithm. This routine requires the user to write two functions in order to fully implement the sorting procedure. These are SortCompareFunction and SortSwapFunction.

ARGUMENTS

Data type	Name	Description
void *	data	Data
dip_int	size	Size
dip_SortCompareFunction	compareFunction	Function for comparing two data points
dip_SortSwapFunction	swapFunction	Function for swapping two data points, or
		copying one to the other
void *	tmpData	Pointer to a variable of the same type as the
		data, used as temporary space by some of the
		algorithms
dip_Sort	algorithm	Sort algorithm

SEE ALSO

General information about sorting

 ${\tt QuickSortAnything, SortCompareFunction, SortSwapFunction}$

SortCompareFunction

Typedef for comparison function (sorting)

SYNOPSIS

```
#include "dip_sort.h"
dip_Boolean (*dip_SortCompareFunction) ( data1, index1, data2, index2 )
```

FUNCTION

A function of this type must be supplied to the sorting algorithms for data of arbitrary type. It should return DIP_TRUE if data1[index1] > data2[index2].

Example:

```
dip_Boolean MyComplexCompare( void *data1, dip_int index1, void *data2, dip_int index2 )
   dip_complex *cmplx1, *cmplx2;
   dip_float magnitude1, magnitude2;
   cmplx1 = data1;
   cmplx2 = data2;
   cmplx1 += index1;
   cmplx2 += index2;
   magnitude1 = sqrt( cmplx1->re * cmplx1->re + cmplx1->im * cmplx1->im );
   magnitude2 = sqrt( cmplx2->re * cmplx2->re + cmplx2->im * cmplx2->im );
   if ( magnitude1 > magnitude2 )
   {
      return( DIP_TRUE );
   }
   else
      return( DIP_FALSE );
   }
}
```

ARGUMENTS

Data type	Name	Description
void *	data1	Pointer to first data array
$\mathtt{dip}_{-}\mathtt{int}$	index1	Index to element in first data array
void *	data2	Pointer to second data array
dip_int	index2	Index to element in second data array

SEE ALSO

 $\label{thm:conting} General \ information \ about \ sorting \\ SortAnything, \ QuickSortAnything, \ SortSwapFunction \\$

SortIndices

Sort indices to a block of data

SYNOPSIS

```
#include "dip_sort.h"
dip_Error dip_SortIndices ( data, indices, size, algorithm, dataType, indexType )
```

FUNCTION

Sorts a list of indices rather than the data itself using the algorithm specified by algorithm.

ARGUMENTS

Data type	Name	Description
void *	data	Data
void *	indices	Indices
dip_int	size	Size
dip_Sort	algorithm	Sort algorithm
dip_DataType	dataType	Data type. See DIPlib's data types
dip_DataType	indexType	Data type of the index array. Must be either DIP_DT_SINT32 or
		DIP_DT_SINT16.

The sortType parameter is one of:

Name	Description
DIP_SORT_DEFAULT	Default sort algorithm
DIP_SORT_QUICK_SORT	Quick sort
DIP_SORT_DISTRIBUTION_SORT	Distribution sort
DIP_SORT_INSERTION_SORT	Insertion sort

SEE ALSO

General information about sorting

 ${\tt DistributionSort, InsertionSort, QuickSort, Sort, ImageSort, SortIndices 16, ImageSortIndices}$

SortIndices16

Sort indices to a block of data

SYNOPSIS

```
#include "dip_sort.h"
dip_Error dip_SortIndices16 ( data, indices, size, algorithm, dataType )
```

FUNCTION

Sorts a list of (16 bit) indices rather than the data itself using the algorithm specified by algorithm.

ARGUMENTS

Data type	Name	Description
void *	data	Data
dip_sint16 *	indices	Indices
dip_int	size	Size
dip_Sort	algorithm	Sort algorithm
dip_DataType	dataType	Data type. See DIPlib's data types

The sortType parameter is one of:

Name	Description
DIP_SORT_DEFAULT	Default sort algorithm
DIP_SORT_QUICK_SORT	Quick sort
DIP_SORT_DISTRIBUTION_SORT	Distribution sort
DIP_SORT_INSERTION_SORT	Insertion sort

SEE ALSO

General information about sorting

 ${\tt DistributionSort, InsertionSort, QuickSort, Sort, ImageSort, SortIndices, ImageSortIndices}$

SortSwapFunction

Typedef for swap and copy function (sorting)

SYNOPSIS

```
#include "dip_sort.h"
void (*dip_SortSwapFunction) ( data1, index1, data2, index2, copy )
```

FUNCTION

A function of this type must be supplied to the sorting algorithms for data of arbitrary type. It should swap data1[index1] and data2[index2] if copy = DIP_FALSE, and copy data1[index1] to data2[index2] if copy = DIP_TRUE.

Example:

}

```
void dip_MyComplexSwap( void *data1, dip_int index1, void *data2, dip_int index2, dip_Boolean
{
   dip_complex *cmplx1, *cmplx2, tmpValue;
   cmplx1 = data1;
   cmplx2 = data2;
   cmplx1 += index1;
   cmplx2 += index2;
   if ( copy == DIP_TRUE )
      cmplx2->re = cmplx1->re;
      cmplx2->im = cmplx1->im;
   }
   else
      tmpValue.re = cmplx2->re;
      tmpValue.im = cmplx2->im;
      cmplx2->re = cmplx1->re;
      cmplx2->im = cmplx1->im;
      cmplx1->re = tmpValue.re;
      cmplx1->im = tmpValue.im;
   }
   return;
```

ARGUMENTS

Data type	Name	Description
void *	data1	Pointer to first data array
dip_int	index1	Index to element in first data array
void *	data2	Pointer to second data array
dip_int	index2	Index to element in second data array
dip_Boolean	сору	if DIP_FALSE, swap data. if DIP_TRUE copy data from data1 to data2

SEE ALSO

General information about sorting

SortAnything, QuickSortAnything, SortCompareFunction

Sqrt

arithmetic function

SYNOPSIS

dip_Error dip_Sqrt (in, out)

DATA TYPES

binary, integer, ${f float}$

FUNCTION

Computes the square root of the input image values.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	in	Input
$dip_{-}Image$	out	Output

SEE ALSO

Exp, Exp2, Exp10, Ln, Log2, Log10

${\tt StablePixelHeapFree}$

Destroy heap structure

SYNOPSIS

```
#include "dip_pixelqueue.h"
dip_Error dip_StablePixelHeapFree ( heap )
```

FUNCTION

Frees all data associated to heap and sets heap to 0.

ARGUMENTS

Data type	Name	Description
dip_StablePixelHeap *	heap	The heap structure

SEE ALSO

StablePixelHeapNew, PixelHeapNew, PixelQueueNew, StablePixelHeapPush, StablePixelHeapPop, StablePixelHeapIsEmpty

StablePixelHeapIsEmpty

Query heap

SYNOPSIS

```
#include "dip_pixelqueue.h"
dip_Error dip_StablePixelHeapIsEmpty ( heap, result )
```

FUNCTION

Checks to see if there are any items on the heap. See StablePixelHeapNew for information on the heap data structure.

ARGUMENTS

Data type	Name	Description
dip_StablePixelHeap	heap	The heap structure
dip_Boolean *	result	Set to true if there are no items in the heap

SEE ALSO

StablePixelHeapNew, PixelHeapNew, PixelQueueNew, StablePixelHeapFree, StablePixelHeapPush, StablePixelHeapPop

StablePixelHeapNew

Create a new heap structure

SYNOPSIS

```
#include "dip_pixelqueue.h"
dip_Error dip_StablePixelHeapNew ( heap, ndims, blocksize, order, resources )
```

FUNCTION

This function allocates space for a new dip_StablePixelHeap structure. Memory allocated is tracked in resources. The heap is dimensioned to hold pixels from an ndims-dimensional image, and initially enough space is allocated for blocksize elements. The heap will be expanded as necessary when used.

The heap stores the coordinates, the value and the pointer to a pixel in an image. Note that the value does not need to equal the data pointed to by the pointer. ndims can be set to zero, in which case no coordinates are stored; this does not affect the function of the value and the pointer.

A heap is a priority queue data structure. Just like a queue, items can be added (pushed) and subtracted (popped). However, in the priority queue the item popped is always the higherst priority one: either the one with the highest-valued item (order is DIP_GVSO_HIGH_FIRST) or lowest-valued item (order is DIP_GVSO_LOW_FIRST). When various identically-valued items are stored on the heap, they will be extracted in the same order as they were insterted (FIFO - first in, first out). If this order is unimportant (such as for the GrowRegionsWeighted algorithm, use the more efficient dip_PixelHeap instead. See PixelHeapNew for information on the unstable heap structure.

ARGUMENTS

Data type	Name	Description
dip_StablePixelHeap *	heap	The newly allocated heap structure
dip_int	ndims	Image dimensionality
dip_int	blocksize	Size of each allocation block
dipf_GreyValueSortOrder	order	Determines the heap's sort order
dip_Resources	resources	Resources tracking structure. See ResourcesNew

The dipf_GreyValueSortOrder enumeration consists of the following values:

Name	Description
DIP_GVSO_HIGH_FIRST	Process the pixels from high grey-value to low grey-value.
DIP_GVSO_LOW_FIRST	Process the pixels from low grey-value to high grey-value.

IMPLEMENTATION

This data structure is implemented identically to PixelHeapNew (see that function's description for details), but an insertion order value is attached to each pixel pushed onto the heap. This is used to

maintain stability.

SEE ALSO

 $\label{lem:pixelHeapNew} PixelQueueNew, StablePixelHeapFree, StablePixelHeapPush, StablePixelHeapIop, StablePixelHeapIsEmpty$

StablePixelHeapPop

Pop item onto heap

SYNOPSIS

```
#include "dip_pixelqueue.h"
dip_Error dip_StablePixelHeapPop ( heap, coords, pointer, value )
```

FUNCTION

Pops the next pixel from the heap. See StablePixelHeapNew for information on the heap data structure. coords is a pointer to an array of dip_ints, such as that obtained with dip_IntegerArray->array. It should have as many elements as the image dimensionality. If the stack was created with ndims set to 0, the coords pointer is ignored. coords, pointer and value can be NULL if you are not interested in either those values.

ARGUMENTS

Data type	Name	Description
dip_StablePixelHeap	heap	The heap structure
dip_int *	coords	Receives the coordinates of the popped item
void **	pointer	Receives the pointer of the popped item
dip_sfloat *	value	Receives the value of the popped item

SEE ALSO

StablePixelHeapNew, PixelHeapNew, PixelQueueNew, StablePixelHeapFree, StablePixelHeapPush, StablePixelHeapIsEmpty

StablePixelHeapPush

Push item onto heap

SYNOPSIS

```
#include "dip_pixelqueue.h"
dip_Error dip_StablePixelHeapPush ( heap, coords, pointer, value )
```

FUNCTION

Pushes a pixel onto the heap. See StablePixelHeapNew for information on the heap data structure. All 3 values coords, pointer and value are stored, except if the heap was created with ndims set to 0, in which case the coords pointer is ignored.

coords is a pointer to an array of dip_ints, such as that obtained with dip_IntegerArray->array. It should have as many elements as the image dimensionality. pointer is a pointer to any memory location, and value is the value to be used when sorting.

ARGUMENTS

Data type	Name	Description
dip_StablePixelHeap	heap	The heap structure
dip_int *	coords	Coordinates to be pushed
void *	pointer	Pointer to be pushed
dip_sfloat	value	Value to be pushed

SEE ALSO

StablePixelHeapNew, PixelHeapNew, PixelQueueNew, StablePixelHeapFree, StablePixelHeapPop, StablePixelHeapIsEmpty

StandardDeviation

statistics function

SYNOPSIS

dip_Error dip_StandardDeviation (in, mask, out, ps)

DATA TYPES

binary, integer, float

FUNCTION

Calculates the standard deviation of the pixel values over all those dimensions which are specified by ps.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	in	Input
$\mathtt{dip}_{-}\mathtt{Image}$	mask (0)	Mask
dip_Image	out	Output
dip_BooleanArray	ps (0)	Dimensions to project

SEE ALSO

From images to scalars

Sum, Mean, Variance, MeanModulus, SumModulus, MeanSquareModulus, Maximum, Minimum, Median, Percentile

StringAppend

Append a string to another

SYNOPSIS

```
#include "dip_string.h"
dip_Error dip_StringAppend ( str1, str2, cstr )
```

FUNCTION

Concatenates str1 and str2 and puts the result in str1, which is increased in size if necessary. If str2 is 0, cstr is used instead.

ARGUMENTS

Data type	Name	Description
dip_String	str1	First string
dip_String	str2	Second string
char *	cstr	Second string

SEE ALSO

 ${\tt StringCat}, {\tt StringCompare}, {\tt StringCompareCaseInsensitive}, {\tt StringCopy}, {\tt StringCopy}, {\tt StringReplace}, {\tt UnderscoreSpaces}$

StringArrayCopy

Copy a string array

SYNOPSIS

```
#include "dip_string.h"
dip_Error dip_StringArrayCopy ( new, src, resources )
```

FUNCTION

This function copies the complete src string array to new.

ARGUMENTS

Data type	Name	Description
dip_StringArray *	new	Pointer to the destination dip_StringArray structure
dip_StringArray	src	Source string array
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

StringArrayNew, StringArrayFree

IntegerArrayCopy, FloatArrayCopy, ComplexArrayCopy, DataTypeArrayCopy, BooleanArrayCopy,
VoidPointerArrayCopy, StringArrayCopy

StringArrayFree

Array free function

SYNOPSIS

dip_Error dip_StringArrayFree (array)

FUNCTION

This function frees *array, and sets array to zero.

ARGUMENTS

Data type	Name	Description
dip_ImageArray *	array	string array

SEE ALSO

StringArrayNew, StringArrayCopy

ArrayFree, IntegerArrayFree, FloatArrayFree, ComplexArrayFree, BoundaryArrayFree, FrameWorkProcessArrayFree, DataTypeArrayFree, ImageArrayFree, BooleanArrayFree, VoidPointerArrayFree, StringArrayFree, CoordinateArrayFree

StringArrayNew

Allocate an array of strings

SYNOPSIS

```
#include "dip_string.h"
dip_Error dip_StringArrayNew ( array, size, stringSize, init, resources )
```

FUNCTION

This function allocates an array of strings. size specifies the size of the array, stringSize the size of the individual strings, which are allocated too. If StringSize is zero, and init is not, the strings in the array are initialised with init (in size and content). If both stringSize and init are zero, strings of zero length are created. (see StringNew).

ARGUMENTS

Data type	Name	Description
dip_StringArray *	array	Pointer to the array
dip_int	size	Size of the array
dip_int	stringSize	Size of the strings
char *	init	Initialisation string
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

 ${\tt StringArrayFree}, {\tt StringArrayCopy}$

ArrayNew, IntegerArrayNew, FloatArrayNew, ComplexArrayNew, BoundaryArrayNew, FrameWorkProcessArrayNew, DataTypeArrayNew, ImageArrayNew, BooleanArrayNew, VoidPointerArrayNew, StringArrayNew, CoordinateArrayNew

StringCat

Concatenate two strings

SYNOPSIS

```
#include "dip_string.h"
dip_Error dip_StringCat ( newStr, str1, str2, cstr, resources )
```

FUNCTION

Concatenates str1 and str2 and puts the result in newStr, which is allocated. If str2 is 0, cstr is used instead.

ARGUMENTS

Data type	Name	Description
dip_String *	newStr	Destination
dip_String	str1	First string
dip_String	str2	Second string
char *	cstr	Second string
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

 ${\tt StringAppend}, {\tt StringCompare}, {\tt StringCompareCaseInsensitive}, {\tt StringCopy}, {\tt StringCopy}, {\tt StringReplace}, {\tt UnderscoreSpaces}$

StringCompare

Compare two strings

SYNOPSIS

```
#include "dip_string.h"
dip_Error dip_StringCompare ( orig, copy, verdict )
```

FUNCTION

This function uses the strcmp function to compare orig and copy. If the strings are different, an error is generated, or verdict obtains the value DIP_FALSE, if it is not zero.

ARGUMENTS

Data type	Name	Description
dip_String	orig	The original string
dip_String	сору	The fake (or not) string
dip_Boolean *	verdict	Verdict of the comparison

SEE ALSO

 ${\tt StringAppend, StringCat, StringCompareCaseInsensitive, StringCopy, StringCrop, StringNew, StringReplace, UnderscoreSpaces}$

${\tt StringCompareCaseInsensitive}$

Compare two strings without minding case

SYNOPSIS

```
#include "dip_string.h"
dip_Error dip_StringCompareCaseInsensitive ( orig, copy, verdict )
```

FUNCTION

This function uses the strcasecmp (or stricmp) function to compare orig and copy. If the strings are different, an error is generated, or verdict obtains the value DIP_FALSE, if it is not zero.

ARGUMENTS

Data type	Name	Description	
dip_String	orig	The original string	
dip_String	сору	The fake (or not) string	
dip_Boolean *	verdict	Verdict of the comparison	

SEE ALSO

 ${\tt StringAppend}, {\tt StringCompare}, {\tt StringCopy}, {\tt StringCrop}, {\tt StringNew}, {\tt StringReplace}, {\tt UnderscoreSpaces}$

StringCopy

Copy a String

SYNOPSIS

```
#include "dip_string.h"
dip_Error dip_StringCopy ( new, src, resources )
```

FUNCTION

Thsi function copies string src to new.

ARGUMENTS

Data type	Name	Description	
dip_String *	new	Pointer to a destination dip_String strcture	
dip_String	src	Source string	
dip_Resources	resources	Resources tracking structure. See ResourcesNew	

SEE ALSO

 ${\tt StringAppend}, {\tt StringCat}, {\tt StringCompare}, {\tt StringCompareCaseInsensitive}, {\tt StringCrop}, {\tt StringNew}, {\tt StringReplace}, {\tt UnderscoreSpaces}$

StringCrop

Crop a string

SYNOPSIS

```
#include "dip_string.h"
dip_Error dip_StringCrop ( str, length )
```

FUNCTION

Crops str to length characters.

ARGUMENTS

Data type	Name	Description
dip_String	str	String to be cropped
dip_int	length	New string length

SEE ALSO

 ${\tt StringAppend}, {\tt StringCat}, {\tt StringCompare}, {\tt StringCompareCaseInsensitive}, {\tt StringCopy}, {\tt StringNew}, {\tt StringReplace}, {\tt UnderscoreSpaces}$

StringFree

Free a string

SYNOPSIS

```
#include "dip_string.h"
dip_Error dip_StringFree ( void )
```

FUNCTION

This function frees a string data structure that has been allocated using StringNew.

ARGUMENTS

Data type	Name	Description	
dip_String *	string	Pointer to the string to be freed	

SEE ALSO

StringNew

StringNew

Allocate a string

SYNOPSIS

```
#include "dip_string.h"
dip_Error dip_StringNew ( string, size, init, resources )
```

FUNCTION

This function allocates a string of size size. If init is not zero, its contents is copied into the new string. If size is zero, and init is not, the size of string is made equal to init plus one.

ARGUMENTS

Data type	Name	Description	
dip_String *	string	Pointer to the new string	
dip_int	size	Size of the string	
char *	init	Initialisation string	
dip_Resources	resources	Resources tracking structure. See ResourcesNew	

SEE ALSO

 ${\tt StringArrayNew,\,StringAppend,\,StringCat,\,StringCompare,\,StringCompareCaseInsensitive,\,StringCopy,\,StringCrop,\,StringReplace,\,UnderscoreSpaces}$

StringReplace

Replace the contents of one string with that of another

SYNOPSIS

```
#include "dip_string.h"
dip_Error dip_StringReplace ( str1, str2, cstr )
```

FUNCTION

Replaces the content of str1 with str2. str1 is increased in size if necessary. If str2 is 0, cstr is used instead.

ARGUMENTS

Data type	Name	Description
dip_String	str1	Destination string
dip_String	str2	Source string
char *	cstr	Source string

SEE ALSO

 ${\tt StringAppend}, {\tt StringCompare}, {\tt StringCompareCaseInsensitive}, {\tt StringCopy}, {\tt StringCrop}, {\tt StringNew}, {\tt UnderscoreSpaces}$

StructureTensor2D

Two dimensional Structure Tensor

SYNOPSIS

#include "dip_structure.h"

dip_Error dip_StructureTensor2D(in, mask, orientation, energy, 11, 12, anisotropy1,
anisotropy2, curvature, boundary, gradSpec, gradSigmas, tensorSpec, tensorSigmas,
curvatureSpec, curvatureSigmas)

DATA TYPES

integer, float

FUNCTION

This function computes the Structure Tensor (ST) at each point in the image. For a description of this technique see the references. There are two stages in the computation. The first stage computes the gradient vector at each point, using Derivative with parameters gradSpec and gradSigmas. The second stage, the tensor smoothing, is also performed using Derivative (with order = 0). The parameters used are tensorSpec and tensorSigmas.

If a mask image is given, a technique called normalised convolution (see references) is used to "fill in" the missing data.

The routine has a number of output images. Each of these can be set to zero. If set to zero, the corresponding result will not be computed. The following quantities are computed by this routine:

orientation	Orientation. Lies in the interval (-pi/2,pi/2).
energy	Sum of the two eigenvalues 11 and 12.
11	The largest eigenvalue.
12	The smallest eigenvalue.
anisotropy1	Measure for local anisotropy: (11 - 12) /
	(11 + 12).
anisotropy2	Measure for local anisotropy: 1 - 12 / 11.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	mask	Mask image (0=missing data)
dip_Image	orientation	Orientation
dip_Image	energy	Energy (11+12)
dip_Image	11	Largest eigenvalue
dip_Image	12	Smallest eigenvalue
dip_Image	anisotropy1	Local anisotropy: (11-12)/(11+12)
dip_Image	anisotropy2	Local anisotropy: 1-12/11
dip_Image	curvature	undocumented, set to 0
dip_BoundaryArray	boundary	Boundary conditions
dip_DerivativeSpec	gradSpec	Parameters for derivative to compute gradient (see
		DerivativeSpec data structure)
dip_FloatArray	gradSigmas	Sigmas of derivative to compute gradient
dip_DerivativeSpec	tensorSpec	Parameters for Gaussian for tensor smoothing (see
		DerivativeSpec data structure)
dip_FloatArray	tensorSigmas	Sigmas of Gaussian for tensor smoothing
dip_DerivativeSpec	curvatureSpec	undocumented, set to 0
dip_FloatArray	curvatureSigmas	undocumented, set to 0

LITERATURE

Bernd Jahne, Practical Handbook on Image Processing for Scientific Applications, chapter 13, CRC Press, 1997

L.J. van Vliet and P.W. Verbeek, *Estimators for Orientation and Anisotropy in Digitized Images*, in: J. van Katwijk, J.J. Gerbrands, M.R. van Steen, J.F.M. Tonino (eds.), ASCI'95, Proc. First Annual Conference of the Advanced School for Computing and Imaging (Heijen, NL, May 16-18), ASCI, Delft, 1995, pp. 442-450.

C.F. Westin, A $Tensor\ Framework\ for\ Multidimensional\ Signal\ Processing,$ PhD thesis, Linkoping University, Sweden, 1994

SEE ALSO

Derivative

Sub

arithmetic function

SYNOPSIS

```
dip_Error dip_Sub ( in1, in2, out )
Calls Arith ( in1, in2, out, DIP_ARITHOP_SUB, DIP_DT_MINIMUM )
```

SubComplex

arithmetic function

SYNOPSIS

dip_Error dip_SubComplex (in, out, constant)

DATA TYPES

binary, integer, float, complex

FUNCTION

This function computes out = in - constant on a pixel by pixel basis. The data types of the in1 image and constant may be of different types. See Information about dyadic operations for more information about what the type of the output will be.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	out	Output
dip_complex	constant	Constant

SEE ALSO

SubFloat

arithmetic function

SYNOPSIS

dip_Error dip_SubFloat (in, out, constant)

DATA TYPES

binary, integer, float, complex

FUNCTION

This function computes out = in - constant on a pixel by pixel basis. The data types of the in1 image and constant may be of different types. See Information about dyadic operations for more information about what the type of the output will be.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	in	Input
dip_Image	out	Output
dip_float	constant	Constant

SEE ALSO

Arith, Arith_ComplexSeparated, SubInteger, SubComplex, AddFloat, MulFloat, DivFloat

SubInteger

arithmetic function

SYNOPSIS

dip_Error dip_SubInteger (in, out, constant)

DATA TYPES

 $binary, \, \mathbf{integer}, \, \mathbf{float}, \, \mathbf{complex}$

FUNCTION

This function computes out = in - constant on a pixel by pixel basis. The data types of the in1 image and constant may be of different types. See Information about dyadic operations for more information about what the type of the output will be.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	in	Input
dip_Image	out	Output
dip_int	constant	Constant

SEE ALSO

Arith, Arith_ComplexSeparated, SubFloat, SubComplex, AddInteger, MulInteger, DivInteger

SubpixelLocation

Gets coordinates of an extremum with sub-pixel precision

SYNOPSIS

```
#include "dip_analysis.h"
dip_Error dip_SubpixelLocation ( in, pos, coords, val, method, polarity )
```

DATA TYPES

integer, float

FUNCTION

Determines the sub-pixel location of a local maximum or minimum close to pos. pos should point to a pixel that is larger than its direct neighbours (if polarity is DIP_SEP_MAXIMUM) or smaller than its direct neighbours (polarity is DIP_SEP_MINIMUM). coords will contain the the sub-pixel location of this local extremum. val will contain the interpolated grey value at the location of the extremum. method determines which method is used.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input grayscale image
dip_IntegerArray	pos	Input coordinates
dip_FloatArray	coords	Output coordinates
dip_float*	val	Output grey value
dipf_SubpixelExtremumMethod	method	Sub-pixel detection method
dipf_SubpixelExtremumPolarity	pol	Maximum or minimum?

The dipf_SubpixelExtremumMethod flag can be any of these values:

Name	Description
DIP_SEM_DEFAULT	Same as DIP_SEM_PARABOLIC_SEPARABLE.
DIP_SEM_LINEAR	Computes the center of gravity of 3 pixels around the
	extremum, in each dimension independently. The val returned
	is that of the pixel at pos.
DIP_SEM_PARABOLIC_SEPARABLE	Fits a parabola to 3 pixels around the extremum, for each
	dimension independently. The val returned is the maximum
	(minimum) value of these 1D extrema, and thus not equivalent
	to the grey value obtained by true interpolation.
DIP_SEM_PARABOLIC	Fits a parabolic patch to a region 3x3 or 3x3x3 pixels around
	the extremum (only for 2D or 3D images).
DIP_SEM_GAUSSIAN_SEPARABLE	Same as DIP_SEM_PARABOLIC_SEPARABLE, but using the log of
	the pixel values, very accurate if peak is a Gaussian.
DIP_SEM_GAUSSIAN	Same as DIP_SEM_PARABOLIC, but using the log of the pixel
	values (only for 2D or 3D images).
DIP_SEM_BSPLINE	Fits a B-spline to 11 pixels around the extremum, in each
	dimension independently. The val returned is the maximum
	(minimum) value of these 1D extrema, and thus not equivalent
	to the grey value obtained by true interpolation

SEE ALSO

SubpixelMaxima, SubpixelMinima

SubpixelMaxima

Gets coordinates of local maxima with sub-pixel precision

SYNOPSIS

```
#include "dip_analysis.h"
dip_Error dip_SubpixelMaxima ( in, mask, out_coord, out_val, method )
```

DATA TYPES

integer, float

FUNCTION

Detects local maxima in the image, and returns their coordinates, with sub-pixel precision, in the output image out_coord. Only pixels where mask is on will be examined. Local maxima are detected using Maxima, then their position is determined accurately using SubpixelLocation. out_coord will have ndims pixels along the first dimension (ndims being the number of dimensions in in), and as many pixels along the second dimension as there are local maxima in in. Thus, each row of the image out_coord contains the coordinates of one local maximum. out_coord is always dip_float. out_val, when not 0, will contain the interpolated values of the image at the local maxima. out_val will have the same size and type as out_coord, except only one pixel along the first dimension.

A local maximum can not touch the edge of the image. That is, its integer location must be one pixel away from the edge.

See SubpixelLocation for the definition of the method parameter.

ARGUMENTS

Data type	Name	Description
$ exttt{dip_Image}$	in	Input grayscale image
$ extstyle{dip_Image}$	mask	Binary mask for ROI processing
dip_Image	out_coord	Output coordinates, image will be N_dims x
		N_maxima
dip_Image	out_val	Output values, image will be 1 x N_maxima
dipf_SubpixelExtremumMethod	method	Sub-pixel detection method

SEE ALSO

SubpixelMinima, SubpixelLocation, Maxima

SubpixelMinima

Gets coordinates of local minima with sub-pixel precision

SYNOPSIS

```
#include "dip_analysis.h"
dip_Error dip_SubpixelMinima ( in, mask, out_coord, out_val, method )
```

DATA TYPES

integer, float

FUNCTION

Detects local minima in the image, and returns their coordinates, with sub-pixel precision, in the output image out_coord. Only pixels where mask is on will be examined. Local minima are detected using Minima, then their position is determined accurately using SubpixelLocation. out_coord will have ndims pixels along the first dimension (ndims being the number of dimensions in in), and as many pixels along the second dimension as there are local minima in in. Thus, each row of the image out_coord contains the coordinates of one local minimum. out_coord is always dip_float. out_val, when not 0, will contain the interpolated values of the image at the local minima. out_val will have the same size and type as out_coord, except only one pixel along the first dimension.

A local minimum can not touch the edge of the image. That is, its integer location must be one pixel away from the edge.

See SubpixelLocation for the definition of the method parameter.

ARGUMENTS

Data type	Name	Description
$ exttt{dip_Image}$	in	Input grayscale image
$ extstyle{dip_Image}$	mask	Binary mask for ROI processing
dip_Image	out_coord	Output coordinates, image will be N_dims x
		N_minima
dip_Image	out_val	Output values, image will be 1 x N_minima
dipf_SubpixelExtremumMethod	method	Sub-pixel detection method

SEE ALSO

SubpixelMaxima, SubpixelLocation, Minima

Subsampling

Interpolation function

SYNOPSIS

```
#include "dip_interpolation.h"
dip_Error dip_Subsampling ( in, out, sample )
```

DATA TYPES

binary, integer, float, complex

FUNCTION

This function subsamples in by copying each sampleth pixel to out.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input image
dip_Image	out	Output image
dip_IntegerArray	sample	Sample spacing

SEE ALSO

Resampling

Sum

statistics function

SYNOPSIS

dip_Error dip_Sum (in, mask, out, ps)

DATA TYPES

binary, integer, float, complex

FUNCTION

Calculates the sum of the pixel values over all those dimensions which are specified by ps.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	mask (0)	Mask
$dip_{-}Image$	out	Output
dip_BooleanArray	ps (0)	Dimensions to project

SEE ALSO

From images to scalars

 ${\tt Mean, Variance, Standard Deviation, Mean Modulus, Sum Modulus, Mean Square Modulus, Maximum, Minimum, Median, Percentile}$

SumModulus

statistics function

SYNOPSIS

dip_Error dip_SumModulus (in, mask, out, ps)

DATA TYPES

binary, integer, float, complex

FUNCTION

Calculates the sum of the modulus the pixel values over all those dimensions which are specified by ps.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	mask (0)	Mask
dip_Image	out	Output
dip_BooleanArray	ps (0)	Dimensions to project

SEE ALSO

From images to scalars

Sum, Mean, Variance, StandardDeviation, MeanModulus, MeanSquareModulus, Maximum, Minimum, Median, Percentile

Tan

trigonometric function

SYNOPSIS

dip_Error dip_Tan (in, out)

DATA TYPES

binary, integer, float, complex

FUNCTION

Computes the tangent of the input image values.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	in	Input
dip_Image	out	Output

SEE ALSO

Sin, Cos, Tan, Asin, Acos, Atan, Atan2, Sinh, Cosh, Tanh

Tanh

trigonometric function

SYNOPSIS

dip_Error dip_Tanh (in, out)

DATA TYPES

binary, integer, ${f float}$

FUNCTION

Computes the hyperbolic tangent of the input image values.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	out	Output

SEE ALSO

Sin, Cos, Tan, Asin, Acos, Atan, Atan2, Sinh, Cosh

${\tt TensorImageInverse}$

Invert tensor image

SYNOPSIS

dip_Error dip_TensorImageInverse (in, out)

DATA TYPES

float

FUNCTION

Inverts the NxN tensor image in (stored as an array with N*N elements) using LU decomposition.

ARGUMENTS

Data type	Name	Description
dip_ImageArray	in	Input
dip_ImageArray	out	Output

TestObjectAddNoise

TestObject generation function

SYNOPSIS

```
#include "dip_generation.h"
dip_Error dip_TestObjectAddNoise ( object, noisy, background, backvalue,
gaussianNoise, poissonNoise, snr, conversion, variance, random )
```

DATA TYPES

binary, integer, float

FUNCTION

This function adds a mixture of Gaussian and Poisson noise to a testobject at a specified signal-to-noise ratio. The SNR is defined as the average object energy divided by the average noise power.

ARGUMENTS

Data type	Name	Description
dip_Image	object	Input Object Image
$\mathtt{dip}_{-}\mathtt{Image}$	noisy	Output Image
dip_Image	background	Background Image
dip_float	backvalue	Constant Background Value
dip_float	gaussianNoise	Relative Amount of Gaussian Noise
dip_float	poissonNoise	Relative Amount of Poisson Noise
dip_float	snr	Signal to Noise Ratio
dip_float *	conversion (0)	Pointer to the Poisson Conversion Factor
dip_float *	variance (0)	Pointer to the Gaussian Variance
dip_Random *	random	Pointer to a random value structure

SEE ALSO

TestObjectCreate, TestObjectModulate, TestObjectBlur

TestObjectBlur

TestObject generation function

SYNOPSIS

```
#include "dip_generation.h"
dip_Error dip_TestObjectBlur ( object, psf, convolved, xNyquist, testPSF )
```

DATA TYPES

binary, integer, float

FUNCTION

This function blurs a testobject with a Gaussian psf, with a two dimensional in focus diffraction limited incoherent PSF or with an user-supplied PSF. The xNyquist parameter specifies the oversampling factor of the incoherent PSF and Gaussian PSF. The sigma of the Gaussian PSF is equal to 0.9 * xNyquist.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	object	Input Object Image
$\mathtt{dip}_{-}\mathtt{Image}$	psf	User supplied PSF
dip_Image	convolved	Output Image
dip_float	xNyquist	Oversampling Factor
dipf_TestPSF	testPSF	TestPSF

The dipf_TestPSF enumaration consists of the following flags:

Name	Description
DIP_TEST_PSF_GAUSSIAN	Gaussian PSF
DIP_TEST_PSF_INCOHERENT_OTF	in-focus, diffraction limited, incoherent PSF
DIP_TEST_PSF_USER_SUPPLIED	User supplied PSF with the psf image
DIP_TEST_PSF_NONE	no blurring

SEE ALSO

TestObjectCreate, TestObjectModulate, TestObjectAddNoise

TestObjectCreate

TestObject generation function

SYNOPSIS

#include "dip_generation.h"

dip_Error dip_TestObjectCreate (object, testObject, objectHeight, objectRadius, scale, scaleRadius, scaleAmplitude, objSigma, position, random)

DATA TYPES

Output: sfloat

FUNCTION

This function can generate an aliasing free object (ellips, box, ellipsoid shell, box shell) or uses an user-supplied object. The generated objects have their origin at the center in the image, but can be generated with a sub-pixel random shift around the center, to average out dicretization effects over several instances of the same generated object. Optinally the generated object can be convolved with an isotropic Gaussian with a width specified by objSigma. Elliptical objects are only supported for images with a dimsnionality equal or less than three. The position boolean variable specifies whether a subpixel random shift should be applied to the object. This can be used to average out digitisation error over a repetition of the generation of the same object.

ARGUMENTS

Data type	Name	Description
dip_Image	object	Output Object Image
dipf_TestObject	testObject	Type of Test Object
dip_float	objectHeight	Object Height
dip_float	objectRadius	Object Radius
dip_FloatArray	scale	Relative Radii for each dimension
dip_float	scaleRadius	ScaleRadius
dip_float	scaleAmplitude	ScaleAmplitude
dip_float	objSigma	Sigma of Gaussian Object Blur
dip_Boolean	position	Random Subpixel Position Shift
dip_Random *	random	Pointer to a random value structure

SEE ALSO

TestObjectModulate, TestObjectBlur, TestObjectAddNoise

TestObjectModulate

TestObject generation function

SYNOPSIS

```
#include "dip_generation.h"
dip_Error dip_TestObjectModulate ( in, out, modulation, modulationDepth )
```

DATA TYPES

Output: sfloat

FUNCTION

This function adds a sine modulation to a test object, with modulation the modulation frequency and modulationDepth the modulation depth.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{\mathtt{-}}\mathtt{Image}$	in	Input
dip_Image	out	Output
dip_FloatArray	modulation	Modulation Frequency
dip_float	modulationDepth	ModulationDepth

SEE ALSO

TestObjectCreate, TestObjectBlur, TestObjectAddNoise

Threshold

Point Operation

SYNOPSIS

#include "dip_point.h"
dip_Error dip_Threshold (in, out, threshold, foreground, background, binaryOutput)

DATA TYPES

integer, **float**

FUNCTION

This function thresholds an image at the threshold value. If the boolean binaryOutput is true, Threshold will produce a binary image. Otherwise an image of the same type as the input image is produced, with the pixels set to either foreground or background. In other words: out = (in >= threshold? foreground: background)

ARGUMENTS

Data type	Name	Description	
$dip_{-}Image$	in	Input image	
$\mathtt{dip}_{-}\mathtt{Image}$	out	Output image	
dip_float	threshold	Threshold value	
dip_float	foreground	Foreground value	
dip_float	background	Background value	
dip_Boolean	binaryOutput	Convert output image to binary	

SEE ALSO

See section 10.3, "Segmentation", in Fundamentals of Image Processing.

 ${\tt RangeThreshold}, \, {\tt SelectValue}, \, {\tt NotZero}, \, {\tt Compare}, \, {\tt HysteresisThreshold}, \, {\tt IsodataThreshold}, \, {\tt Clip}$

TikhonovMiller

Image restoration filter

SYNOPSIS

#include "dip_restoration.h" #include "dip_transform.h"
dip_Error dip_TikhonovMiller (in, psf, out, reg, background, method, var, lambda,
flags)

FUNCTION

The TikhonovMiller restoration filter is a linear least squares restoration algorithm.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input image
dip_Image	psf	Point spread function image
dip_Image	out	Output image
dip_Image	reg	Regularisation filter image
dip_Image	background (0)	Background image
dipf_RegularizationParameter	method	Method used to determine the
		regularisation parameter
dip_float	var	Noise variance
dip_float *	lambda	Regularisation parameter
${ t dipf_ImageRestoration}$	flags	Restoration flags

LITERATURE

G.M.P. van Kempen, *Image Restoration in FLuorescence Microscopy*, Ph.D. Thesis, Delft University of Technology, 1999

SEE ALSO

 ${\tt Wiener,\,TikhonovRegularizationParameter}$

TikhonovRegularizationParameter

Determine the value of the regularisation parameter

SYNOPSIS

#include "dip_restoration.h"
dip_Error dip_TikhonovRegularizationParameter (in, psf, reg, background, max, min,
lambda, method, var, flags)

FUNCTION

This function implements different methods to estimate the value of the regularistion parameter lambda of the TikhonovMiller restoration filter.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input image
dip_Image	psf	Point spread function image
dip_Image	reg	Regularisation filter rimage
dip_Image	background (0)	Background image
dip_float	max	Maximum value of lambda
dip_float	min	Minimum value of lambda
dip_float *	lambda	pointer to the regularisation parameter
dipf_RegularizationParameter	method	Method used to determine lambda
dip_float	var	Noise variance
dipf_ImageRestoration	flags	Restoration flags

LITERATURE

G.M.P. van Kempen, *Image Restoration in FLuorescence Microscopy*, Ph.D. Thesis, Delft University of Technology, 1999

SEE ALSO

 ${\tt TikhonovRegularizationParameter}$

TimerGet

Timing functions

SYNOPSIS

```
#include "dip_timer.h"
dip_Error dip_TimerGet ( timer )
```

FUNCTION

This function gets three timer values elapsed since the last call to TimerSet.

The dip_Timer struct contains the following values:

Data type	Name	Description	
$\mathtt{dip}_{-}\mathtt{int}$	setTime	Time stamp when TimerSet was called.	
$\mathtt{dip}_{\mathtt{-}}\mathtt{int}$	getTime	Time stamp when TimerGet was called.	
dip_float	getClockTime	Amount of CPU time (in seconds) between TimerSet and	
		TimerGet.	
dip_float	getSystemTime	Amount of CPU time (in seconds) executing system calls in the	
		process.	
dip_float	getUserTime	Amount of CPU time (in seconds) executing instructions in the	
		process.	
dip_float	setClockTime	Set by TimerSet, just ignore!	
dip_float	setSystemTime	Set by TimerSet, just ignore!	
dip_float	setUserTime	Set by TimerSet, just ignore!	

setTime and getTime give the time, in seconds, elapsed since the Epoch (00:00:00 UTC, January 1, 1970). The C function ctime can convert this value into a date string.

getClockTime gives the CPU time, in seconds, between the call to TimerSet and TimerGet. The number of significant digits depends on your system. getUserTime contains the portion of this time that was used by the CPU to process instructions for the current process. getSystemTime contains the portion of time spent in the system while executing tasks on behalf of the current process (e.g. doing file I/O). getUserTime and getSystemTime do not necessarily add up to getClockTime if there are other processes running on the same processor.

ARGUMENTS

Data type	Name	Description	
dip_Timer *	timer	Pointer to a dip_Timer struct	

NOTES

Note that getClockTime, getUserTime and getSystemTime can wrap around. The system returns these values as a clock_t value. If this is a 32-bit integer, these timers wrap around after only 72

minutes.

 ${\tt getUserTime}$ and ${\tt getSystemTime}$ are not supported on some systems, it is possible that these values are always 0.

SEE ALSO

 ${\tt TimerSet}$

TimerSet

Timing functions

SYNOPSIS

```
#include "dip_timer.h"
dip_Error dip_TimerSet ( timer )
```

FUNCTION

This function resets three timers that can be obtained by TimerGet.

ARGUMENTS

Data type	Name	Description	
dip_Timer *	timer	pointer to a dip_Timer structure	

SEE ALSO

TimerGet

Tophat

Morphological high-pass filter

SYNOPSIS

```
#include "dip_morphology.h"
dip_Error dip_Tophat ( in, out, se, boundary, param, shape, edgeType, polarity )
```

DATA TYPES

integer, float

FUNCTION

The top-hat is the difference between a morphological operation and the original image, comparable to a high-pass filter. Which operation is used can be chosen through the dip_MphEdgeType and dip_MphTophatPolarity parameters.

The rectangular, elliptic and diamond structuring elements are "flat", i.e. these structuring elements have a constant value. For these structuring elements, param determines the sizes of the structuring elements.

When shape is DIP_FLT_SHAPE_DISCRETE_LINE or DIP_FLT_SHAPE_INTERPOLATED_LINE, the structuring element is a line. param->array[0] determines the length, param->array[1] the angle. This is currently only supported for 2D images. Interpolated lines use interpolation to obtain a more accurate result, but loose the strict increasingness and extensivity (these properties are satisfied only by approximation).

When shape is set to DIP_FLT_SHAPE_PARABOLIC, params specifies the curvature of the parabola.

When shape is set to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, se is used as structuring element. It can be either a binary or a grey-value image. Its origin is the center, or one pixel to the left of the center if the size is even.

If shape is not equal to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, se can be set to zero. When shape is set to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, param is ignored, and can be set to zero.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	out	Output
dip_Image	se	Custom structuring element
dip_BoundaryArray	boundary	Boundary conditions
dip_FloatArray	param	Filter parameters
dip_FilterShape	shape	Structuring element
dip_MphEdgeType	edgeType	edgeType
dip_MphTophatPolarity	polarity	polarity

The enumerator dip_FilterShape contains the following constants:

Name	Description
DIP_FLT_SHAPE_DEFAULT	Default filter window, same as
	DIP_FLT_SHAPE_RECTANGULAR
DIP_FLT_SHAPE_RECTANGULAR	Rectangular filter window, can be even in size
DIP_FLT_SHAPE_ELLIPTIC	Elliptic filter window, always odd in size
DIP_FLT_SHAPE_DIAMOND	Diamond-shaped filter window, always odd in size
DIP_FLT_SHAPE_PARABOLIC	Parabolic filter window (morphology only)
DIP_FLT_SHAPE_DISCRETE_LINE	Rotated line structuring element (morphology only)
DIP_FLT_SHAPE_INTERPOLATED_LINE	Rotated line structuring element, through interpolation
	(morphology only)
DIP_FLT_SHAPE_PERIODIC_LINE	(not implemented)
DIP_FLT_SHAPE_STRUCTURING_ELEMENT	Use se as filter window, can be any size

The enumerator ${\tt dip_MphEdgeType}$ contains the following constants:

Name	Description	
DIP_MPH_TEXTURE	Response is limited to edges in texture	
DIP_MPH_OBJECT	Response is limited to object edges	
DIP_MPH_BOTH	All edges produce equal response	

The enumerator dip_MphTophatPolarity contains the following constants:

Name	Description
DIP_MPH_TEXTURE	Response is limited to edges in texture
DIP_MPH_OBJECT	Response is limited to object edges
DIP_MPH_BOTH	All edges produce equal response

SEE ALSO

Lee, MorphologicalGradientMagnitude, MorphologicalRange, MultiScaleMorphologicalGradient, MorphologicalSmoothing, MorphologicalThreshold

tpi.h

Type iterator

SYNOPSIS

#include "dip_tpi.h"

FUNCTION

Type iterator. For each data type specified by the define DIP_TPI_ALLOW, dip_tpi.h will include the file specified by the define DIP_TPI_FILE. If DIP_TPI_ALLOW is not defined the file will be included for all data types. DIP_TPI_TYPES must be defined as a logical OR of identifier flags and identifier group flags, as given in DIPlib's data types and the table below. During each "iteration" the main symbols defined by dip_tpi.h are DIP_TPI_DATA_TYPE, DIP_TPI_IDENTIFIER and DIP_TPI_EXTENSION. The following table shows how these are defined for each data type:

DID EDI	DID MDI DAMA MUDD	DID MDI IDDUMININ	DID EDI DVEDNOTON
DIP_TPI	DIP_TPI_DATA_TYPE	DIP_TPI_IDENTIFIER	DIP_TPI_EXTENSION
dip_bin8	DIP_DT_BIN8	DIP_DTID_BIN8	_b8
dip_bin16	DIP_DT_BIN16	DIP_DTID_BIN16	_b16
dip_bin32	DIP_DT_BIN32	DIP_DTID_BIN32	_b32
dip_uint8	DIP_DT_UINT8	DIP_DTID_UINT8	_u8
dip_uint16	DIP_DT_UINT16	DIP_DTID_UINT16	_u16
dip_uint32	DIP_DT_UINT32	DIP_DTID_UINT32	_u32
dip_sint8	DIP_DT_SINT8	DIP_DTID_SINT8	_s8
dip_sint16	DIP_DT_SINT16	DIP_DTID_SINT16	_s16
dip_sint32	DIP_DT_SINT32	DIP_DTID_SINT32	_s32
dip_sfloat	DIP_DT_SFLOAT	DIP_DTID_SFLOAT	_sfl
dip_dfloat	DIP_DT_DFLOAT	DIP_DTID_DFLOAT	_dfl
dip_scomplex	DIP_DT_SCOMPLEX	DIP_DTID_SCOMPLEX	_SCX
dip_dcomplex	DIP_DT_DCOMPLEX	DIP_DTID_DCOMPLEX	_dcx

Using this include file it is possible to compile source code for different data types. We recommend that instead of splitting your code into two files, one for generic code and one for type specific code, that you use dip_tpi.h to let the source file include itself. This also prevents dependency problems with makefiles. A source file that includes itself through dip_tpi.h should have the following format:

```
contents of example.c:
#ifndef DIP_TPI

#include "diplib.h"

#define DIP_TPI_FILE "example.c"
#include "dip_tpi.h"

/* This is where the generic code should be */
```

#else

/* This is where the type specific code should be */

#endif

In addition to the main defines as described above, there are a number of macro's that are defined by dip_tpi.h:

DIP_TPI_FUNC (function name)	attaches the current type suffix to the function
	name.
DIP_TPI_DEFINE (function name)	equivalent to: dip_Error DIP_TPI_FUNC(
	function name) useful for function
	definitions.
DIP_TPI_DECLARE (function name)	equivalent to: dip_Error DIP_TPI_FUNC(
	function name) useful for function
	declarations. Don't forget the trailing ";".
DIP_TPI_NAME (function name)	attaches the current type suffix to the function
	name and puts double quotes around the result,
	thus creating a string.

There are also a couple of defines that are only available for some of the data types:

When DIP_TPI is	
dip_sfloat	DIP_TPI_CAST_R2C is defined as dip_scomplex
dip_dfloat	DIP_TPI_CAST_R2C is defined as dip_dcomplex
dip_scomplex	DIP_TPI_CAST_C2R is defined as dip_sfloat
dip_dcomplex	DIP_TPI_CAST_C2R is defined as dip_dfloat

Other type iterators may be created by making a copy of the dip_tpi.h file and replacing DIP_TPI throughout the file by a different name for the new type iterator.

ARGUMENTS

Name	Description
DIP_TPI_ALLOW	logical OR of data type identifier and identifier group flags to indicate for which
	data types the file should be included
DIP_TPI_FILE	Name of the file to be included by dip_tpi.h

SEE ALSO

DIPlib's data types

DataTypeGetInfo, ovl.h

Truncate

Arithmetic function

SYNOPSIS

dip_Error dip_Truncate (in, out)

DATA TYPES

binary, integer, float

FUNCTION

Computes the truncation of the input image values, and outputs a signed integer typed image.

ARGUMENTS

Data type	Name	Description	
$dip_{-}Image$	in	Input	
$dip_{-}Image$	out	Output	

SEE ALSO

Abs, Ceil, Floor, Sign, Fraction, NearestInt

UnderscoreSpaces

Replace spaces with underscores

SYNOPSIS

```
#include "dip_string.h"
dip_Error dip_UnderscoreSpaces ( string )
```

FUNCTION

This function replaces spaces in string with underscores. This function works in-place.

ARGUMENTS

Data type	Name	Description
dip_String	string	String to be examined

SEE ALSO

StringAppend, StringCat, StringCompare, StringCompareCaseInsensitive, StringCopy, StringCrop, StringNew, StringReplace

Uniform

Uniform filter

SYNOPSIS

```
#include "dip_linear.h"
dip_Error dip_Uniform ( in, out, se, boundary, param, shape )
```

DATA TYPES

binary, integer, float, complex

FUNCTION

This functions implements an uniform convolution filter with support for various filter shapes.

Only the rectangular, elliptic and diamond filter shapes are supported (DIP_FLT_SHAPE_RECTANGULAR, DIP_FLT_SHAPE_ELLIPTIC and DIP_FLT_SHAPE_DIAMOND). Other filter shapes can be implemented by setting shape to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, and passing a binary image in se. The "on" pixels define the shape of the filter window. Other values of shape are illegal.

If shape is not equal to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, se can be set to zero. When shape is set to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, param is ignored, and can be set to zero.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
$\mathtt{dip}_{-}\mathtt{Image}$	out	Output
dip_Image	se	Custom filter shape (binary)
dip_BoundaryArray	boundary	Boundary conditions
dip_FloatArray	param	Filter parameters
dip_FilterShape	shape	Filter shape

The enumerator dip_FilterShape contains the following constants:

Name	Description
DIP_FLT_SHAPE_DEFAULT	Default filter window, same as
	DIP_FLT_SHAPE_RECTANGULAR
DIP_FLT_SHAPE_RECTANGULAR	Rectangular filter window, can be even in size
DIP_FLT_SHAPE_ELLIPTIC	Elliptic filter window, always odd in size
DIP_FLT_SHAPE_DIAMOND	Diamond-shaped filter window, always odd in size
DIP_FLT_SHAPE_PARABOLIC	Parabolic filter window (morphology only)
DIP_FLT_SHAPE_DISCRETE_LINE	Rotated line structuring element (morphology only)
DIP_FLT_SHAPE_INTERPOLATED_LINE	Rotated line structuring element, through interpolation
	(morphology only)
DIP_FLT_SHAPE_PERIODIC_LINE	(not implemented)
DIP_FLT_SHAPE_STRUCTURING_ELEMENT	Use se as filter window, can be any size

SEE ALSO

 $\label{thm:convolution} General information about convolution \\ \textbf{Gauss, GeneralConvolution}$

UniformNoise

Generate an image disturbed by uniform noise

SYNOPSIS

```
#include "dip_noise.h"
dip_Error dip_UniformNoise ( in, out, lowerBound, upperBound, random )
```

DATA TYPES

integer, **float**

FUNCTION

Generate an image disturbed by additive uniform noise. See UniformRandomVariable for more information on the random number generator.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	out	Output
dip_float	lowerBound	Lower bound of the uniform distribution the noise is drawn from
dip_float	upperBound	Upper bound of the uniform distribution the noise is drawn from
dip_Random *	random	Pointer to a random value structure

EXAMPLE

Get a image with additive uniform noise as follows:

```
dip_Image in, out;
dip_float lower, upper;
dip_Random random;

lower = 1.0;
upper = 10.0;
DIPXJ( dip_RandomSeed( &random, 0 ));
DIPXJ( dip_UniformNoise( in, out, lower, upper, &random ));
```

SEE ALSO

 $\label{thm:comparison} Uniform Random Variable, Random Seed, Random Seed Vector, Gaussian Noise, Poisson Noise, Binary Noise$

UniformRandomVariable

Uniform random variable generator

SYNOPSIS

```
#include "dip_noise.h"
dip_Error dip_UniformRandomVariable ( random, lowerBound, upperBound, output)
```

FUNCTION

Generate an uniform distributed random variable. See RandomVariable for more information on the random number generator.

ARGUMENTS

Data type	Name	Description
dip_Random *	random	Pointer to a random value structure
dip_float	lowerBound	Lower bound of the uniform distribution the variable is drawn from
dip_float	upperBound	Upper bound of the uniform distribution the variable is drawn
		from
dip_float*	output	output

EXAMPLE

Get a uniform random variable as follows:

```
dip_Random random;
dip_float lower, upper, value;

lower = -1.0;
upper = 1.0;
DIPXJ( dip_RandomSeed( &random, 0 ));
DIPXJ( dip_UniformRandomVariable( &random, lower, upper, &value ));
```

SEE ALSO

RandomVariable, RandomSeed, RandomSeedVector, GaussianRandomVariable, PoissonRandomVariable, BinaryRandomVariable

Unregister

Remove a registry item

SYNOPSIS

```
#include "dip_registry.h"
dip_Error dip_Unregister ( id, class )
```

FUNCTION

This function removes the Registry information of the ID of the Registry class class. See Register for more information about DIPlib's Registry.

ARGUMENTS

Data type	Name	Description	
dip_int	id	Registry ID	
dip_int	class	Registry class	

SEE ALSO

 ${\tt Register}, {\tt RegistryList}, {\tt RegistryGet}, {\tt RegistryArrayNew}$

UpperEnvelope

Upper envelope transform (a flooding and an algebraic closing)

SYNOPSIS

```
#include "dip_morphology.h"
dip_Error dip_UpperEnvelope ( in, out, bottom, labels, connectivity, max_depth,
max_size )
```

DATA TYPES

integer, float

FUNCTION

The Upper envelope transform produces a flooding of the input image (which is an algebraic closing). See any article by F. Meyer for further explanations.

The Upper envelope is based on the watershed transform, each region being filled up to the level where it meets a neighbouring region. See Watershed for information on the parameters.

The bottom image is a second output image that contains the whole watershed region painted with the lowest value in it. It is useful for stretching the input image: (out - in) / (in - bottom). labels returns the label image used during region growing.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	out	Output
$\mathtt{dip}_{-}\mathtt{Image}$	bottom	Optional output
$\mathtt{dip}_{-}\mathtt{Image}$	labels	Optional output
dip_int	connectivity	Connectivity
dip_float	max_depth	Maximum depth of a region that can be merged
dip_int	max_size	Maximum size of a region that can be merged

SEE ALSO

Watershed, LocalMinima

Variance

statistics function

SYNOPSIS

dip_Error dip_Variance (in, mask, out, ps)

DATA TYPES

binary, integer, float

FUNCTION

Calculates the variance of the pixel values over all those dimensions which are specified by ps.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	mask (0)	Mask
$\mathtt{dip}_{-}\mathtt{Image}$	out	Output
dip_BooleanArray	ps (0)	Dimensions to project

SEE ALSO

From images to scalars

 ${\tt Sum}, \, {\tt Mean}, \, {\tt StandardDeviation}, \, {\tt MeanModulus}, \, {\tt SumModulus}, \, {\tt MeanSquareModulus}, \, {\tt Maximum}, \, {\tt Minimum}, \, {\tt Median}, \, {\tt Percentile}$

VarianceFilter

Sample Variance Filter

SYNOPSIS

```
#include "dip_filtering.h"
dip_Error dip_VarianceFilter ( in, out, se, boundary, param, shape )
```

DATA TYPES

binary, integer, float

FUNCTION

This function calculates for every pixel the sample variance of the pixels in the filter window (its size specified by param).

Only the rectangular, elliptic and diamond filter shapes are supported (DIP_FLT_SHAPE_RECTANGULAR, DIP_FLT_SHAPE_ELLIPTIC and DIP_FLT_SHAPE_DIAMOND). Other filter shapes can be implemented by setting shape to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, and passing a binary image in se. The "on" pixels define the shape of the filter window. Other values of shape are illegal.

If shape is not equal to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, se can be set to zero. When shape is set to DIP_FLT_SHAPE_STRUCTURING_ELEMENT, param is ignored, and can be set to zero.

ARGUMENTS

Data type	Name	Description
dip_Image	in	Input
dip_Image	out	Output
dip_Image	se	Custom filter window (binary)
dip_BoundaryArray	boundary	Boundary conditions
dip_FloatArray	param	Filter sizes
dip_FilterShape	shape	Filter shape

The enumerator dip_FilterShape contains the following constants:

Name	Description
DIP_FLT_SHAPE_DEFAULT	Default filter window, same as
	DIP_FLT_SHAPE_RECTANGULAR
DIP_FLT_SHAPE_RECTANGULAR	Rectangular filter window, can be even in size
DIP_FLT_SHAPE_ELLIPTIC	Elliptic filter window, always odd in size
DIP_FLT_SHAPE_DIAMOND	Diamond-shaped filter window, always odd in size
DIP_FLT_SHAPE_PARABOLIC	Parabolic filter window (morphology only)
DIP_FLT_SHAPE_DISCRETE_LINE	Rotated line structuring element (morphology only)
DIP_FLT_SHAPE_INTERPOLATED_LINE	Rotated line structuring element, through interpolation
	(morphology only)
DIP_FLT_SHAPE_PERIODIC_LINE	(not implemented)
DIP_FLT_SHAPE_STRUCTURING_ELEMENT	Use se as filter window, can be any size

SEE ALSO

Kuwahara

VectorDistanceTransform

Euclidean vector distance transform

SYNOPSIS

#include "dip_distance.h"
dip_Error dip_VectorDistanceTransform (in, outx, outy, outz, distance, border,
method)

DATA TYPES

binary

FUNCTION

This function produces the vector components of the Euclidean distance transform. These are stored in the output images, one for each dimension of the input image. See the EuclideanDistanceTransform for detailed information about the parameters.

To compute the Euclidean distance from the vector components produced by this function, one needs to multiply each component with the sampling distance, square the result, sum the results for all components and take the square root of the sum.

ARGUMENTS

Data type	Name	Description
$ exttt{dip_Image}$	in	Input image
dip_ImageArray	out	Output images
dip_FloatArray	distance	Sampling distances
dip_Boolean	border	Image border type
dipf_DistanceTransform	method	Transform method

dipf_DistanceTransform defines the following distance transform types:

Name	Description
DIP_EDT_FAST	fastest, but most errors
DIP_EDT_TIES	slower, but fewer errors
DIP_EDT_TRUE	slow, uses lots of memory, but is "error free"
DIP_EDT_BRUTE_FORCE	gives a result from which errors are calculated for the other methods. This
	method is extremly slow and should only be used for testing purposes.

LITERATURE

See EuclideanDistanceTransform

KNOWN BUGS

 $See \ {\tt Euclidean Distance Transform}$

AUTHOR

James C. Mullikin, adapted to DIPlib by Geert M.P. van Kempen

SEE ALSO

 ${\tt Euclidean Distance Transform, Grey Weighted Distance Transform}$

VoidPointerArrayCopy

Copy an array

SYNOPSIS

dip_Error dip_VoidPointerArrayCopy (dest, src, resources)

FUNCTION

This function copies the void pointer array src to dest. The array dest is created by this function as well.

ARGUMENTS

Data type	Name	Description
dip_IntegerArray *	dest	Destination array
dip_IntegerArray	src	Source array
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

VoidPointerArrayNew, VoidPointerArrayFree, VoidPointerArrayCopy, VoidPointerArrayFind IntegerArrayCopy, FloatArrayCopy, ComplexArrayCopy, DataTypeArrayCopy, BooleanArrayCopy, VoidPointerArrayCopy, StringArrayCopy

VoidPointerArrayFind

Find value in array

SYNOPSIS

dip_Error dip_VoidPointerArrayFind (array, value, index, found)

FUNCTION

Finds a value in an array and "returns" its index in the array. If found is zero, VoidPointerArrayFind will produce an error if value is not found, otherwise found obtains the search result (DIP_FALSE if value is not found).

ARGUMENTS

Data type	Name	Description
dip_VoidPointerArray	array	Array to find value in
void *	value	Value to find
dip_int *	index	Index of the found value
dip_VoidPointer *	found	Value found or not

SEE ALSO

 $\label{thm:condition} VoidPointerArrayFree, VoidPointerArrayCopy, VoidPointerArrayFind\\ IntegerArrayFind, FloatArrayFind, ComplexArrayFind, DataTypeArrayFind, BooleanArrayFind, VoidPointerArrayFind\\$

VoidPointerArrayFree

Array free function

SYNOPSIS

dip_Error dip_VoidPointerArrayFree (array)

FUNCTION

This function frees *array, and sets array to zero.

ARGUMENTS

Data type	Name	Description
dip_VoidPointerArray *	array	Array

SEE ALSO

BooleanArrayNew, BooleanArrayFree, BooleanArrayCopy, BooleanArrayFind

ArrayFree, IntegerArrayFree, FloatArrayFree, ComplexArrayFree, BoundaryArrayFree, FrameWorkProcessArrayFree, DataTypeArrayFree, ImageArrayFree, BooleanArrayFree, VoidPointerArrayFree, StringArrayFree, CoordinateArrayFree

VoidPointerArrayNew

Array allocation function

SYNOPSIS

dip_Error dip_VoidPointerArrayNew (array, size, resources)

FUNCTION

This function allocates the size elements of a dip_VoidPointerArrayNew and sets the size of the array to size.

ARGUMENTS

Data type	Name	Description
dip_VoidPointerArray *	array	Array
dip_int	size	Size
dip_Resources	resources	Resources tracking structure. See ResourcesNew

SEE ALSO

VoidPointerArrayNew, VoidPointerArrayFree, VoidPointerArrayCopy, VoidPointerArrayFind ArrayNew, IntegerArrayNew, FloatArrayNew, ComplexArrayNew, BoundaryArrayNew, FrameWorkProcessArrayNew, DataTypeArrayNew, ImageArrayNew, BooleanArrayNew, VoidPointerArrayNew, StringArrayNew, CoordinateArrayNew

Watershed

Morphological segmentation

SYNOPSIS

```
#include "dip_morphology.h"
dip_Error dip_Watershed ( in, mask, out, connectivity, max_depth, max_size, binaryOutput )
```

DATA TYPES

integer, float

FUNCTION

Watershed segmentation with built-in region merging. max_depth and max_size control the merging procedure. Any region with max_size or less pixels and with max_depth grey-value difference or less will be merged to neighbouring regions when they touch (as opposed to build a watershed). max_size equal to 0 means that the size of the region is not tested when merging. The regions are grown according to the connectivity parameter. See The connectivity parameter for more information. The output is either a labelled image where the pixels belonging to a catchment basin are labelled, or a binary image where the watershed pixels are 1 and the rest is 0. This is controlled by binaryOutput.

If mask is not 0, only the pixels within mask will be considered. All the other pixels will be marked as watershed pixels.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	in	Input
$dip_{-}Image$	mask	Mask
dip_Image	out	Output
dip_int	connectivity	Connectivity
dip_float	max_depth	Maximum depth of a region that can be merged
dip_int	max_size	Maximum size of a region that can be merged
dip_Boolean	binaryOutput	DIP_FALSE if the output should be a labelled image

SEE ALSO

SeededWatershed, UpperEnvelope, LocalMinima, GrowRegions

WeightedAdd

arithmetic function

SYNOPSIS

dip_Error dip_WeightedAdd (in1, in2, out, weight)

DATA TYPES

binary, integer, float, complex

FUNCTION

This function calculates out = in1 + weight * in2;

ARGUMENTS

Data type	Name	Description
dip_Image	in1	First input
dip_Image	in2	Second input
dip_Image	out	Output
dip_float	weight	Weight

SEE ALSO

WeightedMul, WeightedSub, WeightedDiv

WeightedDiv

arithmetic function

SYNOPSIS

dip_Error dip_WeightedDiv (in1, in2, out, weight)

DATA TYPES

binary, integer, float, complex

FUNCTION

This function calculates $\mathtt{out} = \mathtt{in1} \ / \ \mathtt{weight} \ ^* \ \mathtt{in2}; \ \mathrm{If} \ (\mathtt{weight} \ ^* \ \mathtt{in2}) \ \mathrm{is} \ \mathrm{zero}, \ \mathtt{out} \ \mathrm{will} \ \mathrm{be} \ \mathrm{set} \ \mathrm{to} \ \mathrm{zero} \ \mathrm{as} \ \mathrm{well}.$

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	in1	First input
dip_Image	in2	Second input
dip_Image	out	Output
dip_float	weight	Weight

SEE ALSO

WeightedAdd, WeightedMul, WeightedSub

WeightedMul

arithmetic function

SYNOPSIS

dip_Error dip_WeightedMul (in1, in2, out, weight)

DATA TYPES

binary, integer, float, complex

FUNCTION

This function calculates out = in1 * weight * in2;

ARGUMENTS

Data type	Name	Description
dip_Image	in1	First input
dip_Image	in2	Second input
$dip_{-}Image$	out	Output
dip_float	weight	Weight

SEE ALSO

WeightedAdd, WeightedSub, WeightedDiv

WeightedSub

arithmetic function

SYNOPSIS

dip_Error dip_WeightedSub (in1, in2, out, weight)

DATA TYPES

binary, integer, float, complex

FUNCTION

This function calculates out = in1 - weight * in2;

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{\mathtt{-}}\mathtt{Image}$	in1	First input
dip_Image	in2	Second input
dip_Image	out	Output
dip_float	weight	Weight

SEE ALSO

WeightedAdd, WeightedMul, WeightedDiv

Wiener

Image Restoration Filter

SYNOPSIS

```
#include "dip_restoration.h"
dip_Error dip_Wiener ( in, psf, signalPower, noisePower, out, flags )
```

FUNCTION

This function performs an image restoration using the Wiener filter. The Wiener filter is the linear restoration filter that is optimal in mean square error sense.

ARGUMENTS

Data type	Name	Description
$\mathtt{dip}_{-}\mathtt{Image}$	in	Input image
dip_Image	psf	Point spread function image
dip_Image	signalPower	SignalPower image
dip_Image	noisePower	NoisePower image
dip_Image	out	Output image
dipf_Restoration	flags	Restoration flags

LITERATURE

 ${\rm G.M.P.}$ van Kempen, $Image\ Restoration\ in\ FLuorescence\ Microscopy,\ Ph.D.\ Thesis,\ Delft\ University\ of\ Technology,\ 1999$

SEE ALSO

PseudoInverse, TikhonovMiller

Wrap

Wrap an image

SYNOPSIS

```
#include "dip_manipulation.h"
dip_Error dip_Wrap ( in, out, wrap )
```

DATA TYPES

binary, integer, float, complex

FUNCTION

This function wraps the in image around its image borders. wrap specifies the number of pixels over which the image has to wrapped in each dimension.

ARGUMENTS

Data type	Name	Description
${\tt dip_Image}$	in	Input image
dip_Image	out	Output image
dip_IntegerArray	wrap	Wrap parametrs

SEE ALSO

Wrap, Crop, Shift

Xor

logic operation

SYNOPSIS

dip_Error dip_Xor (in1, in2, out)

DATA TYPES

binary

FUNCTION

The function Xor performs the logic XOR operation between the corresponding pixels in in1 and in2, and stores the result in out.

ARGUMENTS

Data type	Name	Description
$dip_{-}Image$	in1	First binary input image
dip_Image	in2	Second binary input image
dip_Image	out	Output image

SEE ALSO

 ${\tt Arith,\,And,\,Or,\,Invert}$

Chapter 3

Assorted topics

3.1 Boundary conditions

Neighbourhood operations pose a problem. What happens when the neighbourhood operator operates near the border of the image and needs data from the area outside the image? The usual solution, also adopted by DIPlib, is to silently extend the image. There are various ways of extending the boundary. Below is a list of the possible methods. More details can be found in the user guide. Note that not all functions support all of these.

Name	Description
DIP_BC_SYM_MIRROR	Symmetric mirroring
DIP_BC_ASYM_MIRROR	Asymmetric mirroring
DIP_BC_PERIODIC	Periodic copying
DIP_BC_ASYM_PERIODIC	Asymmetric periodic copying
DIP_BC_ADD_ZEROS	Extending the image with zeros
DIP_BC_ADD_MAX_VALUE	Extending the image with +infinity
DIP_BC_ADD_MIN_VALUE	Extending the image with -infinity

SEE ALSO

BoundaryArrayNew, BoundaryArrayFree FillBoundaryArray, SeparableFrameWork

3.2 Compression methods for image files

The dipio_Compression structure

The structure dipio_Compression specifies the compression method to use when writing an image file, and contains the following elements:

Data type	Name	Description
dipio_CompressionMethod	method	Compression method
dip_int	level	Compression parameter, dependent on method

dipio_CompressionMethod is an enum with the known compression methods. File formats typically only support one or a few of these, and most of these methods do not have a parameter to set, in which case level is ignored. If an unsupported compression method is selected, no compression is done. The dipio_CompressionMethod has the following values:

Name	Description
DIPIO_CMP_DEFAULT	Default compression method for the file format
DIPIO_CMP_NONE	No compression
DIPIO_CMP_GZIP	ZIP compression, using zlib. The level parameter is between 0 and 9,
	1 being the faster, lesser compression and 9 being the slower, higher
	compression. 0 indicates no compression.
DIPIO_CMP_DEFLATE	Deflate (same as DIPIO_CMP_GZIP)
DIPIO_CMP_COMPRESS	Using UNIX's "compress" utility, which uses the LZW algorithm
DIPIO_CMP_LZW	LZW compression (same as DIPIO_CMP_COMPRESS)
DIPIO_CMP_JPEG	Lossy JPEG compression. The level parameter is between 1 and 100,
	higher numbers giving better quality output but larger files.
DIPIO_CMP_PACKBITS	PackBits
DIPIO_CMP_THUNDERSCAN	ThunderScan
DIPIO_CMP_NEXT	NeXT
DIPIO_CMP_CCITTRLE	CCITT RLE
DIPIO_CMP_CCITTRLEW	CCITT RLE/W
DIPIO_CMP_CCITTFAX3	CCITT Group 3
DIPIO_CMP_CCITTFAX4	CCITT Group 4

Thus only DIPIO_CMP_GZIP and DIPIO_CMP_JPEG currently have a level to set.

Supported compression methods for the various file formats

The TIFF file writer understand the methods DIPIO_CMP_NONE, DIPIO_CMP_DEFLATE, DIPIO_CMP_LZW, DIPIO_CMP_JPEG, DIPIO_CMP_PACKBITS, DIPIO_CMP_THUNDERSCAN, DIPIO_CMP_NEXT, DIPIO_CMP_CCITTRLE, DIPIO_CMP_CCITTRLEW, DIPIO_CMP_CCITTFAX3 and DIPIO_CMP_CCITTFAX4. It defaults to DIPIO_CMP_DEFLATE. The level parameter is currently not used.

The ICS file writer understands DIPIO_CMP_NONE, DIPIO_CMP_GZIP and DIPIO_CMP_COMPRESS, although DIPIO_CMP_COMPRESS is currently not implemented. It defaults to DIPIO_CMP_GZIP.

The GIF file writer only understands <code>DIPIO_CMP_LZW</code>. The compression method selected is simply ignored.

The JPEG file writer only understands $\mathtt{DIPIO_CMP_JPEG}$. The compression method selected is simply ignored.

The PNG file writer only understands <code>DIPIO_CMP_GZIP</code>. The compression method selected is simply ignored. [The PNG writer is not yet implemented!]

All other file writers do not compress, and simply ignore the compression method requested.

3.3 DerivativeSpec data structure

STRUCTURE

This structure is an aggregate of common parameters for derivative operators. Its current definition is:

```
typedef struct
{
    dip_DerivativeFlavour flavour;
    dip_float truncation;
} dip_DerivativeSpec;
```

The enumerator flavour parameter is one of:

Name	Description
DIP_DF_DEFAULT	Default derivative flavour (==DIP_DF_FIRGAUSS)
DIP_DF_FIRGAUSS	Gaussian family, FIR implementation, Gauss
DIP_DF_IIRGAUSS	Gaussian family, IIR implementation, GaussIIR
DIP_DF_FTGAUSS	Gaussian family, FT implementation, GaussFT
DIP_DF_FINITEDIFF	Finite difference implementation, FiniteDifferenceEx

SEE ALSO

StructureTensor2D, Derivative

3.4 DIPlib's data types

Pixel values are represented by	different types,	called data	types.	DIPlib	supports the	data types
given in the following table:						

data type	dip_DataType	data type identifier	suffix
dip_bin8	DIP_DT_BIN8	DIP_DTID_BIN8	_b8
dip_bin16	DIP_DT_BIN16	DIP_DTID_BIN16	_b16
dip_bin32	DIP_DT_BIN32	DIP_DTID_BIN32	_b32
dip_uint8	DIP_DT_UINT8	DIP_DTID_UINT8	_u8
dip_uint16	DIP_DT_UINT16	DIP_DTID_UINT16	_u16
dip_uint32	DIP_DT_UINT32	DIP_DTID_UINT32	_u32
dip_sint8	DIP_DT_SINT8	DIP_DTID_SINT8	_s8
dip_sint16	DIP_DT_SINT16	DIP_DTID_SINT16	_s16
dip_sint32	DIP_DT_SINT32	DIP_DTID_SINT32	_s32
dip_sfloat	DIP_DT_SFLOAT	DIP_DTID_SFLOAT	_sfl
dip_dfloat	DIP_DT_DFLOAT	DIP_DTID_DFLOAT	_dfl
dip_scomplex	DIP_DT_SCOMPLEX	DIP_DTID_SCOMPLEX	_SCX
dip_dcomplex	DIP_DT_DCOMPLEX	DIP_DTID_DCOMPLEX	_dcx

The data types can be divided into five classes: the binary, unsigned integer, signed integer, floating point and complex classes. Different data types in the same class (e.g. dip_uint8 and dip_uint16) provide a different range of values they can represent.

The complex data types are defines as follows:

```
typedef struct
{
    dip_sfloat re;
    dip_sfloat im;
} dip_scomplex;

typedef struct
{
    dip_dfloat re;
    dip_dfloat im;
} dip_dcomplex;
```

The binary data types are simply aliases for a set of corresponding unsigned integer types. The reason for having a separate typedef for the binary types is that they are not used like ordinary integers. Each bit of the integer can store one binary value. When manipulating binary data, care must be taken not to change any of the other bits of the integer used for storing it.

The dip_DataType enumeration is used to represent data types symbolically. It is used in dip_Image's to indicate what the data type of the image is. Data type identifiers are used by the type iterator (see tpi.h) and overload schemes (see ovl.h and overload.h). Type suffixes are used to give type specific routines a unique name. Using a standard set of suffixes enables the type iterator and overload schemes to deal with these type specific routines. The dip_DataType enumeration, data type identifiers and suffixes can be found in the table above.

In addition to the data type identifiers for individual data types, there are also defines to represent an entire group. These are given in the following table:

Data type identifier group	data types
DIP_DTGID_UINT	unsigned integer
DIP_DTGID_UNSIGNED	unsigned integer
DIP_DTGID_SINT	signed integer
DIP_DTGID_INT	signed and unsigned integer
DIP_DTGID_INTEGER	signed and unsigned intege
DIP_DTGID_FLOAT	floating-point
DIP_DTGID_REAL	integer and floating-point
DIP_DTGID_COMPLEX	complex floating-point
DIP_DTGID_SIGNED	signed integer, floating-point and complex
DIP_DTGID_BINARY	binary
DIP_DTGID_ALL	all

SEE ALSO

 ${\tt DataTypeGetInfo}$

3.5 Description of DIPlib's pixel tables

Pixel tables provide an efficient way to encode a multi-dimensional binary object. DIPlib's dip_PixelTable implements this using runlength encoding (in 2-D this coding scheme is known as pxy-tables).

A DIPlib pixel table is a structure (defined in dip_pixel_table.h) that incorporates a link-list of runlengths. Each run-length consists of a n-D coordinate (integer array) and the length of the run along the X dimension. All the runlengths in total encode the binary object.

LITERATURE

See section 3.6, "Contour representations", in Fundamentals of Image Processing.

I.T. Young, R.L. Peverini, P.W. Verbeek and P.J. van Otterloo, *A New Implementation for Binary and Minkowski Operators*, Computer Graphics and Image Processing, Volume 17, No. 3, 189-210, 1981

3.6 File formats recognized by dipIO

The Registry

A number of file reading and writing functions are included in dipIO. These are registered in the ImageReadRegistry and the ImageWriteRegistry. Through this registry, ImageRead and ImageWrite are able to read from and write to any registered file format. You can add your own functions to these (the interface functions for this are not documented yet), thereby increasing the possibilities of ImageRead and ImageWrite.

Below you can find a list of currently supported file formats for both reading and writing. To obtain the format ID from the registry, you need to include the specified file and call the specified function.

Reading

These are the file formats currently supported for reading:

format	include file	registry ID retrieval function	dimension- ality	colour	data types
ICS (Image Cytometry	dipio_ics.h	dipio_ReadICS	any ID	yes	any
Standard) TIFF (Tagged Image File Format)	dipio_tiff.h	dipio_ReadTIF	2D FID	yes	any
PNG (Portable Network Graphics) [not yet implemented!]	dipio_png.h	dipio_ReadPNG	2D ID	yes	uint8 and uint16
JPEG (JPEG File Interchange Format)	dipio_jpeg.h	dipio_ReadJPE	2D GID	yes	uint8
GIF (Graphics Interchange Format)	dipio_gif.h	dipio_ReadGIF	2D ID	yes	uint8
LSM (Zeiss LSM file format)	dipio_lsm.h	dipio_ReadLSM	1D - 4D ID	no	uint8, uint16 and sfloat
PIC (BioRad PIC file format)	dipio_pic.h	dipio_ReadPIC	2D and 3D ID	no	uint8
CVS (Comma Separated Values)	dipio_csv.h	dipio_ReadCSV	2D ID	no	sfloat

Writing

These are the file formats currently supported for writing:

format	include file	registry ID retrieval function	dimension- ality	colour	data types
ICS v1 (Image Cytometry Standard)	dipio_ics.h	dipio_WriteIC	any Sv1ID	yes	any
ICS v2 (Image Cytometry Standard)	dipio_ics.h	dipio_WriteIC	any Sv2ID	yes	any
TIFF (Tagged Image File Format)	dipio_tiff.h	dipio_WriteTI	2D FFID	yes	any in grey-value, uint8 in colour
PNG (Portable Network Graphics) [not yet implemented!]	dipio_png.h	dipio_WritePN	2D GID	yes	uint8 and uint16
JPEG (JPEG File Interchange Format)	dipio_jpeg.h	dipio_WriteJP	2D EGID	yes	uint8
GIF (Graphics Interchange Format)	dipio_gif.h	dipio_WriteGI	2D FID	no	uint8
CVS (Comma Separated Values)	dipio_csv.h	dipio_WriteCS	2D VID	no	any except complex
FLD (AVS field file)	dipio_fld.h	dipio_WriteFL	any DID	no	any
PS (PostScript)	dipio_ps.h	dipio_WritePS	2D ID	yes	uint8, others automatically converted
EPS (Encapulated PostScript)	dipio_ps.h	dipio_WriteEP	2D SID	yes	uint8, others automatically converted

3.7 From images to scalars

Within DIPlib all data, i.e. multi-dimensional data, such as images, and scalar, are all represented by the same object: the image. Scalars are stored as zero dimensional images. Examine, for example, the following code to compute the sum over all the grey values:

```
dip_Image img;
dip_Image value;
...
dip_Sum ( img, 0, value, 0 );
```

Which stores the sum over all the pixel values of img in the 0-D image value. We often want to directly manipulate scalars, in which case we need to extract the value. This can be accomplished easily with the GetInteger, GetFloat or the GetComplex functions:

```
dip_Image img;
dip_Image valueimg;
dip_float value;
...
dip_Sum ( img, 0, valueimg, 0 );
dip_GetFloat ( valueimg, &value, 0 );
printf ( "The sum is: %f\n", value );
```

3.8 General information about convolution

Convolution can be explained in just a few words: it is a local weighted average (the weights can be negative). This of course does not explain how to use it or what its properties are. For this we refer to the following sources:

Ian T. Young, Jan J. Gerbrands and Lucas J. van Vliet, Fundamentals of Image Processing.

Alan V. Oppenheim, Alan S. Willsky and I.T. Young, "Signals and Systems", Prentice-Hall, 1983.

Anil K. Jain, "Fundamentals of Digital Image Processing", Prentice-Hall, 1989.

"The Digital Signal Processing Handbook", Vijay K. Madisetti and Douglas B. Williams (eds), CRC Press + IEEE Press, 1998.

Kenneth R. Castleman, "Digital Image Processing", Prentice-Hall, 1996.

3.9 General information about sorting

There are two kinds of sorting routines in DIPlib. The first sorts a one-dimensional array of data, the second sorts a set of indices to a one-dimensional array of data. The result of the sort routines can be summarised as follows:

```
Sort: data[ i ] <= data[ i + 1 ]
Sort indices: data[ indices[ i ] ] <= data[ indices[ i + 1 ] ]</pre>
```

Note that the number of indices does not have to be equal to the amount of pixels in the image, it may be either smaller or larger. The indices themselves should of course "point" to a valid pixel.

The sorting algorithms are described in the following reference:

Donald E. Knuth, "The Art of Computer Programming, volume 3: Sorting and Searching", second edition, Addison-Wesley, 1998.

3.10 Information about dyadic operations

There are two types of dyadic operations. First there are operations such as Add, Sub, etc... which take two input images. The second category consists of functions such as AddFloat, AddComplex etc... The data type of the output image given the data types of the input images is given by the following table:

	dcomplex	scomplex	dfloat	sfloat	sint32	sint16
dcomplex						
scomplex	dcomplex	scomplex	dcomplex	scomplex	scomplex	scomplex
dfloat	dcomplex	dcomplex	dfloat	dfloat	dfloat	dfloat
sfloat	dcomplex	scomplex	dfloat	sfloat	sfloat	sfloat
sint32	dcomplex	scomplex	dfloat	sfloat	sint32	sint32
sint16	dcomplex	scomplex	dfloat	sfloat	sint32	sint16
sint8	dcomplex	scomplex	dfloat	sfloat	sint32	sint16
uint32	dcomplex	scomplex	dfloat	sfloat	sint32	sint32
uint16	dcomplex	scomplex	dfloat	sfloat	sint32	sint16
uint8	dcomplex	scomplex	dfloat	sfloat	sint32	sint16
binary	dcomplex	scomplex	dfloat	sfloat	sint32	sint16
	sint8	uint32	uint16	uint8	binary	
dcomplex	dcomplex	dcomplex	dcomplex	dcomplex	dcomplex	
scomplex	scomplex	scomplex	scomplex	scomplex	scomplex	
dfloat	dfloat	dfloat	dfloat	dfloat	dfloat	
sfloat	sfloat	sfloat	sfloat	sfloat	sfloat	
sint32	sint32	sint32	sint32	sint32	sint32	
sint16	sint16	sint32	sint16	sint16	sint16	
sint8	sint8	sint32	sint16	sint8	sint8	
uint32	sint32	uint32	uint32	uint32	uint32	
uint16	sint16	uint32	uint16	uint16	uint16	
uint8	sint8	uint32	uint16	uint8	uint8	
binary	sint8	uint32	uint16	uint8	sint8	
		·				

The output data type of an operation involving an image and a constant of one of the types: dip_complex, dip_float, dip_int, is given by the following table:

	dip_complex	dip_float	dip_int
dcomplex	dcomplex	dcomplex	dcomplex
scomplex	scomplex	scomplex	scomplex
dfloat	dcomplex	dfloat	dfloat
sfloat	scomplex	sfloat	sfloat
sint32	scomplex	sint32	sint32
sint16	scomplex	sint16	sint16
sint8	scomplex	sint8	sint8
uint32	scomplex	uint32	uint32
uint16	scomplex	uint16	uint16
uint8	scomplex	uint8	uint8
binary	scomplex	sint8	sint8

3.11 The image structure

DESCRIPTION

dip_Image is the structure that is used to store images in DIPlib. It contains a number of fields that are used to describe an image. The type field stores the type of the image using a dip_ImageType enumeration. Currently scalar images are the only supported type (DIP_IMTP_SCALAR). The DIP_IMTP_ALIEN type is used internally by DIPlib for creating interfaces to other packages. Whether the other fields in the dip_Image are meaningful depends on the image type. A dip_Image may contain fields specific to the current image type. These will be discussed on the pages pertaining to the type in question. The standard fields that are always present are:

field type	short description	access functions
dip_ImageType	The image type	ImageGetType, ImageSetType
$ ext{dip_ImageState}$	The image state	(none)
dip_DataType	Data type used to store pixel	ImageGetDataType,
	values	ImageSetDataType
dip_IntegerArray	Dimensions of the image	ImageGetDimensions,
		ImageGetDimensionality,
		ImageSetDimensions
void *	Pointer to the pixel data	ImageGetData
dip_int	Plane number, for binary	ImageGetPlane
	images	
dip_IntegerArray	Stride array (see below)	ImageGetStride

Pixel values are stored in the data type specified by the data type field. For a list of possible data types see DIPlib's data types.

The dimensionality of the image and the size of each individual dimension is stored in the dimensions Array.

The data pointer points to the pixel at the origin of the image. For each dimension the stride array holds the interleave between two neighbouring pixels in memory. The following equation may be used to compute the address of a pixel at a coordinate specified by an array called cor[]:

A dip_Image structure does not necessarily have pixel data associated with it. When a dip_Image does not contain pixel data, it is said to be in the "raw" state. A dip_Image that does contain data, is said to be "forged". For binary images the plane field holds the number of the bit in which the binary data is stored. Access to the fields of a dip_Image is restricted to a number of functions, which are given in the table above. The "set" functions can only be used on "raw" images.

SEE ALSO

DIPlib's data types

ImageNew

3.12 The connectivity parameter

DIPlib uses a different name for the various possible connectivites than you might be used to. This is to generalize this parameter to images of any dimensionality. It is defined as follows: if connectivity is 1 all pixels for which only one coordinate differs from the pixel's coordinates by maximally 1 are considered neighbours; if it is 2, all pixels for which one or two coordinates differ maximally 1 are considered neighbours. The connectivity can never be larger than the image dimensionality.

In terms of the obsolete connectivity definitions we have:

In 2-D	this connectivity	corresponds to	and forms this
			structuring element
	1	4 connectivity	diamond
	2	8 connectivity	square
	-1	4-8 connectivity	octagon
	-2	8-4 connectivity	octagon
In 3-D	this connectivity	corresponds to	and forms this
			structuring element
	1	6 connectivity	octahedron
	2	18 connectivity	cuboctahedron
	3	26 connectivity	cube
	-1	6-26 connectivity	small
			rhombicuboctahedron
	-3	26-6 connectivity	small
			rhombicuboctahedron

The negative connectivities are only defined for the functions in binary morphology such as **BinaryDilation** and **BinaryErosion**. These alternate steps with different connectivity to produce a better approximation to an isotropic structuring element.